



The Verifit Directional Mic Test
Evaluating Modern Directional Microphone Technologies

Presenter — David J. Smriga, M.A.
Senior Audiology Consultant
Audiocan Education


audioscan
Hearing Instruments Fitting Systems

www.audioscan.com USA 800-365-2093 519-368-3313 info@audioscan.com
Copyright 2015 Etymonic Design Incorporated. All rights reserved.

IMPORTANT NOTICE

- Unless otherwise noted, the contents of this workshop/presentation are copyright of Audioscan, a Division of Etymonic Design Inc. and are intended for use by workshop attendees only
- Any recording, presentation, adaptation, publishing or sharing of any kind is prohibited without the express written permission of Audioscan

Audioscan® Education
Copyright 2015 Etymonic Design Incorporated. All rights reserved.




Directional Microphone Technology

- Directional technology has been available in hearing aids since 1969
- Remains to this day the single most effective “head-worn” approach to improving signal-to-noise ratio for certain listening conditions

Lybarger, S.F. & Lybarger, E.H. (2000). A historical overview. In R. Sandlin (Ed.), *Textbook of hearing aid amplification: Technical and clinical considerations, 2nd edition*. San Diego, California: Singular Thomson Learning.

Audioscan® Education
Copyright 2015 Etymonic Design Incorporated. All rights reserved.



SNR Impact on Intelligibility

BKB Sentences: Moore, Johnson, Clark & Pluvinaige, 1992

From Dillon, H., "Hearing Aid Technology & The Future", ASA Presentation, April 2006

Audioscan® Education
Copyright 2015 Etymonic Design Incorporated. All rights reserved.

Directional Microphones Make A Difference On Key Satisfaction Factors

Customer satisfaction variable	Multiple Microphone Digital (%)	Single Microphone Digital (%)
Overall satisfaction	17	3
Comfort with loud sounds	17	11
Use in noisy situations	16	7
Sound clarity	16	4
Natural sounding	16	9
Will repurchase brand	15	18
Outdoor situations	14	5
Hours worn	14	18
Battery life	14	11
Whistling/feedback	13	11
Workplace	12	4
Quality of life	12	7
Recommend dispenser	12	9
On-going expense	10	8

% difference in customer satisfaction ratings

Kochkin S. On the issue of value: Hearing aid benefit, price, satisfaction, and brand repurchase rates. Hearing Review. 2003;10(2):12-26

Audioscan® Education
Copyright 2015 Etymonic Design Incorporated. All rights reserved.

Directional Microphone History

circa. 1970's to app. 1995

- The pressure-gradient directional microphone
 - Two mic ports
 - One "front-facing"
 - One "rear-facing"
 - Ports directed to opposite sides of mic diaphragm
 - Port distance = external delay
 - Baffle in RP = internal delay
 - Polar pattern is function of ratio between external and internal delays


Described in Knowles Electronics Technical Bulletin 21

Audioscan® Education
Copyright 2015 Etymonic Design Incorporated. All rights reserved.

Directional Microphone History

App. 1995 to present

- Two omni-directional microphones
- Internal delay is part of DSP algorithm
 - Allows for different polar patterns
 - Allows directional mic to be "On" or "Off"
 - Facilitates automated directional mic activation



H. Dillon; NAL, CRC for CI and HAI

Audioscan® Education
Copyright 2015 Etymonic Design Incorporated. All rights reserved.



Modern Directional Mic Design Possibilities

		<u>Polar Pattern</u>	
		Fixed	Adaptive*
<u>Activation</u>	Manual	Fixed	Manual/Adaptive
	Automatic	Automatic	Adaptive

* Includes broad-band and multi-band adaptive directionality

Fabry, D., "Adaptive Directional Microphone Technology and Hearing Aids: Theoretical and Clinical Implications", The Hearing Review, April 2, 2005

Audioscan® Education
Copyright 2015 Etymonic Design Incorporated. All rights reserved.



Remote Microphone




www.oticon.ca

Audioscan® Education
Copyright 2015 Etymonic Design Incorporated. All rights reserved.




Standard Specification of Directionality

Polar Plots




Audioscan® Education
Copyright 2015 Etymotic Design Incorporated. All rights reserved.





Audioscan® Education
Copyright 2015 Etymotic Design Incorporated. All rights reserved.







Fabry, D., "Adaptive Directional Microphone Technology and Hearing Aids: Theoretical and Clinical Implications", The Hearing Review, April 2, 2005

Audioscan® Education
Copyright 2015 Etymotic Design Incorporated. All rights reserved.







**Considerations Regarding
Conventional Polar Plot Testing**

- Polar plots are obtained in the presence of a single pure tone frequency
- Polar plots do not measure in the presence of multiple input sources
- Polar plots can be less robust in the presence of non-linear (compression) amplification


Audioscan® Education
Copyright 2015 Etymotic Design Incorporated. All rights reserved.

Methods of Quantifying Directional Advantage

- Research Lab
 - The Directivity Index (DI)
 - Compares output for sounds originating from the front to output for sounds originating from all other locations
 - Expressed in dB
 - Omni directional aid would have a 0dB DI
 - Two-mic array devices can have as much as a 6dB DI
 - The AI-DI*
 - An AI (or SII) weighted version of the DI, taking into account the effects of the DI on certain speech frequencies

*Killion, M., et. al., "Real World Performance of an ITE Directional Microphone", The Hearing Journal 51(4), 1998 pp 24-38


Audioscan® Education
Copyright 2015 Etymonic Design Incorporated. All rights reserved.



Methods of Quantifying Directional Advantage

- Clinical
 - The Front-to-Back Ratio (FBR)
 - Comparison of output curves for the same broadband stimulus obtained with the hearing instrument at 0° and 180° orientation to the stimulus
 - Speech In Noise Tests
 - Fixed SNR approach: Difference between omni-directional and directional % correct = directional advantage
 - Variable SNR approach: Difference in SNR required to reach 50% correct = directional advantage
 - SIN, QuickSIN, HINT tests (and others) can be used for this purpose

Audioscan® Education
Copyright 2015 Etymonic Design Incorporated. All rights reserved.




Effects of Compression on FBR

TABLE 2. Single-source traditional and two-source modified front-to-back ratio values averaged within two hearing aid brands and across the four test frequencies of 500, 1000, 2000, and 4000 Hz.

	Traditional	Modified
Type B (compression off)	7.7 dB	7.7 dB
Type B (compression on)	3.9 dB	7.7 dB
Type C (compression off)	11.2 dB	11.2 dB
Type C (compression on)	3.7 dB	11.1 dB

Ricketts, T., "Directivity Quantification in Hearing Aids: Fitting & Measurement Effects" Ear & Hearing, 2000; 21: 45-58

Audioscan® Education
Copyright 2015 Etymonic Design Incorporated. All rights reserved.



Education
The Science of Hearing Aid Fitting


A Unique Way of Measuring Directional Mic Performance

Multiple Input Locations
FFT Analysis

audioscan
Hearing Instrument Fitting Systems

www.audioscan.com USA 800-265-2093 519-368-3313 info@audioscan.com
Copyright 2015 Etymonic Design Incorporated. All rights reserved.


Recommended Positioning - Verifit



A photograph showing a Verifit hearing aid connected to a test rig. Two large speakers are positioned above the hearing aid, creating a sound field. The hearing aid is placed on a dark surface with a microphone and a cable connected to it.

Audioscan® Education
Copyright 2015 Etymonic Design Incorporated. All rights reserved.


Recommended Positioning - Verifit



A diagram showing the Verifit hearing aid on a test rig. The hearing aid is positioned between two speakers. Dotted lines and arrows indicate the directional microphone's field of view at three angles: 0 deg, 135 deg, and 180 deg.

Audioscan® Education
Copyright 2015 Etymonic Design Incorporated. All rights reserved.

Verifit 2 Binaural Coupler




HA stabilizer can be "rocked" until front and rear mic ports of each instrument are aligned parallel to the floorboard of the text box

Putty-less TRIC adapter couples RIC and Thin Tube products to wideband .4cc coupler on both right and left coupler mics

Audioscan® Education
Copyright 2015 Etymotic Design Incorporated. All rights reserved.

Recommended Positioning – Verifit2



Directional Test Speaker

Directional Test Speaker

Front Speaker

Audioscan® Education
Copyright 2015 Etymotic Design Incorporated. All rights reserved.

Directional in Test-box - BTE

Front speaker

Rear speaker



Axis of microphone openings

Audioscan® Education
Copyright 2015 Etymotic Design Incorporated. All rights reserved.

Directional in Test-box - ITE

Front speaker
Rear speaker


Audioscan® Education
Copyright 2015 Etymotic Design Incorporated. All rights reserved.

Directional testing

- Stimulus presented simultaneously from front and rear to condition hearing instrument directional operation.
- Periodically, short burst from only front speaker and only rear speaker used for measurement.
- Duration of measurement bursts short enough to not affect HI output.
- Two versions – Noise + Noise, Noise + Speech

Audioscan® Education
Copyright 2015 Etymotic Design Incorporated. All rights reserved.


Noise + Noise Directional Test



NOISE

FRONT SPEAKER

NOISE




REAR SPEAKER

- Simultaneous composite noise signals
- Alternates for measurement
- Rear stimulus presented longer to condition adaptive directional instruments.

Audioscan® Education
Copyright 2015 Etymotic Design Incorporated. All rights reserved.

Noise + Speech Directional Test



- Speech from front speaker, noise from rear to condition speech activated directionals
- Alternating short noise bursts during measurement intervals

Audioscan® Education
Copyright 2015 Eymonic Design Incorporated. All rights reserved.



DEMONSTRATION ONE:

OMNI DIRECTIONAL SWITCHING TO DIRECTIONAL

Audioscan® Education
Copyright 2015 Eymonic Design Incorporated. All rights reserved.



DEMONSTRATION TWO:


DIRECTIONAL RESULTS AS A FUNCTION OF LEVEL

Audioscan® Education
Copyright 2015 Eymonic Design Incorporated. All rights reserved.




DEMONSTRATION THREE:
DIRECTIONAL RESULTS AS A FUNCTION OF SIGNAL TO NOISE RATIO

Audioscan® Education
Copyright 2015 Elymnic Design Incorporated. All rights reserved.




DEMONSTRATION FOUR:
ADAPTIVE DIRECTIONAL INTERACTION

Audioscan® Education
Copyright 2015 Elymnic Design Incorporated. All rights reserved.



DEMONSTRATION FIVE:
TESTING A DIRECTIONAL MIC IMBEDDED IN A REMOTE MIC

Audioscan® Education
Copyright 2015 Elymnic Design Incorporated. All rights reserved.



Remote Mic Considerations

- Remote microphone options (often connected to paired hearing instruments via blue tooth) can be either hand held (table top) or talker worn
- It can be a separate device or integrated into a wireless remote control
- If the remote microphone contains a directional mic, it too can be tested with the Verifit D-Mic test construct

Audioscan® Education
Copyright 2015 Etymonic Design Incorporated. All rights reserved.



Positioning Remote Mic In Test Box for Directional Testing of Right HA



Audioscan® Education
Copyright 2015 Etymonic Design Incorporated. All rights reserved.



Placement of Coupler Mic Assembly Outside of the Test Box



Audioscan® Education
Copyright 2015 Etymonic Design Incorporated. All rights reserved.



Omni-directional Remote Mic Setting



Audioscan® Education
Copyright 2015 Etymotic Design Incorporated. All rights reserved.



Directional Mic Test Result



Audioscan® Education
Copyright 2015 Etymotic Design Incorporated. All rights reserved.




Directional Remote Mic Setting



Audioscan® Education
Copyright 2015 Etymotic Design Incorporated. All rights reserved.





DIRECTIONAL MEASUREMENT TROUBLESHOOTING

Audioscan® Education
Copyright 2015 Etymonic Design Incorporated. All rights reserved.

- ### No Curve Separation
- Possible Causes:
 - Intensity level may be inadequate
 - SNR may need to be decreased
 - HA placement may be incorrect
 - The directional mic feature of the aid has not been activated
 - One or more mic ports are plugged with debris
 - The directional mic is not working
- Audioscan® Education
Copyright 2015 Etymonic Design Incorporated. All rights reserved.

Curves Are Reversed (when not expected)

- Possible Causes:
 - HA placement may be incorrect (backward)
 - Directional mic is incorrectly wired

Audioscan® Education
Copyright 2015 Etymotic Design Incorporated. All rights reserved.

Poorer Curve Separation Than Initially Benchmarked

- Possible Causes:
 - HA placement may be incorrect or different
 - Current coupling \neq original coupling
 - Slit leak
 - One or more mic ports are compromised by debris
 - Microphone drift

Audioscan® Education
Copyright 2015 Etymotic Design Incorporated. All rights reserved.

Directionality Test (REM)

Rear speaker Aided ear with probe tube positioned Verifit in On-ear directional mode

Equal distance

Audioscan® Education
Copyright 2015 Etymotic Design Incorporated. All rights reserved.

CONSIDERATIONS REGARDING ON-EAR DIRECTIONAL MIC TESTING

Audioscan® Education
Copyright 2015 Etymotic Design Incorporated. All rights reserved.




Venting

- Venting does effect directivity, especially at the low frequencies
 - Best directivity with no vent
- Low frequency sounds pass through the vent w/out attenuation

Audioscan® Education
Copyright 2015 Etymotic Design Incorporated. All rights reserved.





RIC w Power Dome (On-Ear)



Audioscan® Education
Copyright 2015 Etymotic Design Incorporated. All rights reserved.







What is the impact of Venting?

- Data suggest that largest impact of vent in low frequencies (500 & 1000 Hz)
- Largest effect seen for open mold (OM)
- OM condition still resulted in significantly greater directivity than omnidirectional hearing aid
- Listeners receive significant directivity from many directional aids, regardless of venting
 - unless aid only provides directivity mainly in low frequencies

Audioscan® Education
Copyright 2015 Etymotic Design Incorporated. All rights reserved.

Mic Port Orientation

- Custom Products:
 - Modeled (or marked) representation of the horizontal plane for mic port placement on faceplate
 - May be difficult to achieve in some ears
- BTE Products:
 - Horizontal positioning will be influence by:
 - Earmold tubing length (BTE)
 - Receiver wire length (RIC)
 - Thin tube length (Thin tube)

Audioscan® Education
Copyright 2015 Etymotic Design Incorporated. All rights reserved.

Quantifying Mic Port Orientation

Horizontal
-22 degrees

Ricketts, T., "Directivity Quantification in Hearing Aids: Fitting & Measurement Effects"
Ear & Hearing, 2000; 21; 45-58


Audioscan® Education
Copyright 2015 Etymotic Design Incorporated. All rights reserved.


Mic Orientation Effects on DI

Relative to optimal (horizontal) mic port orientation

Ricketts, T., "Directivity Quantification in Hearing Aids: Fitting & Measurement Effects"
Ear & Hearing, 2000; 21; 45-58

Audioscan® Education
Copyright 2015 Etymotic Design Incorporated. All rights reserved.





Summary

- D-mics remain the most effective approach to “headworn” SNR improvement
- Through Audioscan’s unique real-time independent analysis test, a great deal of practical clinical information can be obtained about D-mic functionality
- Test-box D-mic testing is both stable and robust
- On-Ear D-mic testing can be done with consideration of additional factors

Audioscan® Education
Copyright 2015 Etyonic Design Incorporated. All rights reserved.



Education
The Science of Hearing Aid Fitting

Thank You!

daves@audioscan.com



audioscan
Hearing Instrument Fitting Systems

www.audioscan.com USA 800-265-2093 519-368-3313 info@audioscan.com
Copyright 2015 Elymoric Design Incorporated. All rights reserved.
