

**GSI Corti: New Features to Enhance Workflow**

Karen Morris, MS, CCC/A



---

---

---

---

---

---

---

---

**Introduction**

- Experience as a Clinical Audiologist in hospital and neurotology settings for many years.
- Past 10 years worked for audiology equipment manufacturing in education and support role



---

---

---

---

---

---

---

---

**Agenda**

1. Introduction and Review
2. DPOAE Normative Data Review
3. Interpretation of DPOAE Test Results
4. Accessing OAE Collection Parameters and Settings
5. Data Management: Reporting and Printing Options



---

---

---

---

---

---

---

---

### Course Objectives

- Participants will be able to use published DPOAE normative data as a guide when interpreting DPOAE test results
- Participants will be able to list 3 ways to electronically print and save OAE test results
- Participants will be able to access 3 advanced OAE test options included with the GSI Corti



---

---

---

---

---

---

---

---

### OAE Review

- Described by Kemp in 1978
- Acoustic sounds generated by movement of the hair cells in the cochlea
- Emissions are usually present when 'hearing' is normal
- Pre-neural and represent cochlear function



---

---

---

---

---

---

---

---

### OAE Review

- Ear, site and frequency specific
- Very sensitive test of outer hair cell function: OHC damage is closely linked to hearing loss
- Objective in that patient response is not required
- Pediatric and Adult applications



---

---

---

---

---

---

---

---

### 4 Types of OAE

1. **Stimulus Frequency** OAE (SFOAE): elicited by one continuous, low-level sinusoidal signal.
2. **Spontaneous** OAE (SOAE): recorded in the ear canal with no stimulus (not evoked)
3. **Transient** or TEOAE: Stimulus is a broadband click. Activates a wide frequency region. Level is 80 dB SPL



---

---

---

---

---

---

---

---

### 4 Types of OAE

4. **Distortion Product DPOAE**  
Elicited by 2 pure tones (F1 and F2)  
(F2/F1 ratio 1.22 example: F2=2000 and F1=1639Hz)  
At intensity levels (L1 and L2) 65/55  
DP measured at: 2F1-F2 (largest distortion)  
F2 is cochlear place being stimulated  
Use different pure tone pairs to elicit OAE's over a wide range of frequencies



---

---

---

---

---

---

---

---

### How is it measured?

- Probe contains speakers and a microphone
- Speakers deliver the evoking stimuli
- Microphone picks up the sound energy in the ear canal which includes the emission and noise
- Measure the sound pressure level (SPL) present in the *sealed* ear canal



---

---

---

---

---

---

---

---

### How is it analyzed?

- OAE system analyses the spectrum of the sound and separates the “noise” and the “emission”
- An emission usually should be about 4-6 dB higher than the noise floor (NF)
- Can use additional criteria for screening when analyzing data that requires a minimum OAE amplitude or compare to normative data

(Example NF = -20, OAE = -14)



---

---

---

---

---

---

---

---

### Pediatric Applications

- Newborn Hearing Screening
- Mandatory component of pediatric test battery (AAP)
- Dx of ANSD (OAE present/ABR absent-ABN)
- Ototoxicity Monitoring
- Pre-school and school screenings
- Malingering

**“Cross-check Principle”**



---

---

---

---

---

---

---

---

### Adult Applications

- Dx cochlear vs retro cochlear dysfunction
- Malingering
- Ototoxicity Monitoring
- Noise Exposure
- Tinnitus Assessment
- Obscure Auditory Dysfunction (normal audiogram with complaints of hearing loss)



---

---

---

---

---

---

---

---

## Normative Data



---

---

---

---

---

---

---

---

## Normative Data Reference

---

Ear & Hearing  
December 1997 - Volume 18 - Issue 6 - pp 440-455  
Articles

**From Laboratory to Clinic: A Large Scale Study of Distortion Product Otoacoustic Emissions in Ears with Normal Hearing and Ears with Hearing Loss**

Gorga, Michael P.; Neely, Stephen T.; Ohlrich, Brenda; Hoover, Brenda; Redner, Joelle; Peters, Jo



---

---

---

---

---

---

---

---

## A Few Study Goals

- Describe DP's in subjects of all hearing levels (normal and HI) using a large subject sample
- Develop a way to determine the probability that a response was coming from a normal or impaired ear
- Develop a clinical tool



---

---

---

---

---

---

---

---

### Details

- 1267 ears
- 806 individual subjects
- Ranging in age from 1-96 yrs
- No middle ear dysfunction determined by tympanograms performed at 226 Hz and sometimes otoscopic exam by ENT
- Patients with A-B gap of 10 dB at any single frequency were not included



---

---

---

---

---

---

---

---

### Details

- Subjects audiogram was considered the gold standard
- All DPOAE test results were compared to was the subjects audiogram (750-8000 Hz)
- Measurements conducted in a quiet room or "typical clinical condition"



---

---

---

---

---

---

---

---

### DP Stimuli

- F2 (test) frequencies ranged from 750-8000Hz
- F2/F1 ratio of 1.22 (this ratio produces the largest DP)
- Most prominent DP occurs at 2f1-f2
- L1=65, L2=55 (10dB separation, under ~70dB, f2 area of cochlea stimulated/most accurate separation of normal and HI)



---

---

---

---

---

---

---

---

## DP Stimuli

- 500Hz was not included
- Not reliable in predicting auditory status
- $F2/F1 = 1.22$

$F2=500\text{Hz}$      $F1=410\text{Hz}$

$2F1-F2 = 820 - 500 = 320\text{Hz}$




---

---

---

---

---

---

---

---

---

---

## Results

- DP amplitudes were largest for ears with normal hearing
- DP's obtained from hearing impaired ears were not completely separated from DPs from normal hearing ears—**overlap**
- When comparing to audiogram, DPs at 1500-6000 Hz best
- Best correlation when normal hearing is defined as 20-30 dB HL
- Based on test data, were able to develop a way to determine the probability that a response was coming from a normal or impaired ear that could be used clinically




---

---

---

---

---

---

---

---

---

---

### Fig 8: Amplitude/F2




Figure 8. DPOAE amplitude (dB SPL) and DPOAE/noise (dB) as a function of  $F_2$  frequency. The top row represents DPOAE amplitude measurements while the bottom row represents DPOAE/noise estimates. The left column represents the case when no effort was made to account for any floor effects. The right column shows results when DPOAE and noise amplitudes less than -20 dB SPL were set equal to -20 dB SPL. The solid line limiting the top of the hashed areas represents the 95th percentile from the distributions of responses derived from impaired ears. The solid line limiting the bottom of the hashed areas represents the 5th percentile from the distributions of responses from ears with normal hearing. The dashed line limiting the top of the shaded areas represents the 90th percentile from the impaired distributions. The dashed line limiting the bottom of the shaded areas represents the 10th percentile from the normal distributions. Normal hearing was defined as audiometric thresholds of 20 dB HL or less. The values represented in this figure are provided in an Appendix.

---

---

---

---

---

---


---

---

---

---

## Clinical Tool: BT DP Report



Plotting DP amplitude. SNR = 6dB




---

---

---

---

---

---

---

---

---

---

---

---

### HISTORY/REFERRAL REASON:

Referring doctor: \_\_\_\_\_  
 Birth Hospital: \_\_\_\_\_  
 Results from Newborn Hearing Screen: 


Right	Pass	Refer
Left		

Y    N  
 Hearing concerns \_\_\_\_\_  
 Family Hx of hearing loss \_\_\_\_\_  
 Birth complications \_\_\_\_\_  
 History of otitis media \_\_\_\_\_

### RESULTS:

LEFT EAR			RIGHT EAR		
DP level	Noise level	Result	DP level	Noise level	Result
6000 Hz			6000 Hz		
4000 Hz			4000 Hz		
2000 Hz			2000 Hz		
1000 Hz			1000 Hz		

Left Tympanogram: Normal/Abnormal/Tube \_\_\_\_\_  
 Right Tympanogram: Normal/Abnormal/Tube \_\_\_\_\_




---

---

---

---

---

---

---

---

---

---

---

---

## Norms in Data Manager

**GSI Data Manager**

For eligible DPOAE results, the program will display the Expanded Boys Town Norms template as an overlay on the DP-Gram. The light shaded area at the top of the normative curve represents the 50th to 95th percentile of DP amplitudes from the hearing impaired population. DP amplitudes within or above this range indicate a high probability of normal hearing. The light shaded area at the bottom of the normative curve represents the 5th to 10th percentile of DP amplitudes from hearing individuals. DP amplitudes within or below this range indicate a high probability of hearing loss. The dark shaded area in between represents a range of uncertainty where the normal hearing and hearing impaired populations overlap.

The values used to create the template are as shown in table A2 from Gorga, M.P., Neely, S.T., Ohnicks, B., Hoover, B., Redner, J., and Peters, I. (1997). "From laboratory to clinic: a large scale study of distortion product otoacoustic emissions in ears with normal hearing and ears with hearing loss." *Ear & Hearing*, 18, 440-455.

Note: The template will only show if target L1 and L2 levels are 65 and 55 dB SPL respectively, and frequency ratio is between 1.20 and 1.22.

See Operating Manual for more information.

OK




---

---

---

---

---

---

---

---

---

---

---

---

Coming Soon 

Same normative data will be available on the Corti



---

---

---

---

---

---

---

---

Test Interpretation



---

---

---

---

---

---

---

---

US CPT Codes

- **92558**: DP 4S/2S or TE64/TE32 Screening
- **92587**: DP1.5-6.0. TE.7-4.0
  - DPOAE, limited 3-6 frequencies with interpretation and reporting in medical record
  - TEOAE, with interpretation and reporting in medical record
- **92588**: DP1.6-8.0, DP1.5-12 DPOAE, comprehensive, minimum 12 frequencies, with interpretation and reporting in medical record



---

---

---

---

---

---

---

---

## Screening Test: 92558

- Automated pass/refer screening at fixed number of frequencies at a single intensity level
- Test interpretation is not required as all tests will end with a Pass or Refer result
- Used in NHS, nurseries, preschools for this reason




---

---

---

---

---

---

---

---

## Screening Results: DM




---

---

---

---

---

---

---

---

## US CPT Codes

- **92558:** DP 4S/2S or TE64/TE32 Screening
- **92587:** DP1.5-6.0. TE.7-4.0
  - DPOAE, limited 3-6 frequencies with interpretation and reporting in medical record
  - TEOAE, with interpretation and reporting in medical record
  - No Automated pass/refer by default
- **92588:** DP1.6-8.0, DP1.5-12
  - DPOAE, comprehensive, minimum 12 frequencies, with interpretation and reporting in medical record
  - No Automated pass/refer




---

---

---

---

---


---

---

---



### Left Ear 1.5-6.0




---

---

---

---


---

---

---

---

### Right Ear 1.6-8.0




---

---

---

---

---


---

---

---

### Reporting

Normative data will appear in reports




---

---

---

---

---

---

---

---

**Accessing OAE Collection Parameters and Settings**

NOTE 1: Have your device in hand  
NOTE 2: Some of these are coming soon



---

---

---

---

---

---

---

---

**Corti Default Settings**

- L/R Mode
- **Min Value OFF**
- **Graph display: SNR bar graph**
- No **Normative data**
- Diagnostic protocols do not have pass/refer requirements



---

---

---

---

---

---

---

---

**2 'Hidden' Sub-Menus**

- 1. Overall System Options Menu**  
(changes are global)
- 2. Advanced TE and DP OAE (diagnostic) Settings Menu**  
(changes will affect selected protocol)



---

---

---

---

---

---

---

---

### How to Get to Overall System Menu?

From the Corti main screen, press the down button until you reach the time and date screen:



---

---

---

---

---


---

---

---

### From the Time and Date screen:

2. Press and **hold the down** button for 3 seconds or until the green "READY" light turns off. You should see:



New BT is the first item in General system Sub-menu



---

---

---

---

---

---

---

---

### From the Time and Date screen:

2. Press and **hold the down** button for 3 seconds or until the green "READY" light turns off. You should see:



New BT is the first item in General system Sub-menu



---

---

---

---

---

---

---

---

### Overall Systems Menu

9-OCT 02:36 PM

**NEW BT DEVICE**

+ DISCOVER +

+ NEXT

9-OCT 02:36 PM

**POWER OFF**

1 minute

+ CHANGE +

+ NEXT

9-OCT 02:36 PM

**SAVE L/R TESTS**

+ CHANGE +

+ NEXT

09/1/02 15-DEC 07:43 AM

**CLOCK MODE**

12 Hour

+ CHANGE +

+ NEXT

9-OCT 02:36 PM

**1 TEST RESULTS**

+ CLEAR +

+ NEXT

9-OCT 02:37 PM

**MIN VALUE OFF**

+ CHANGE +

+ NEXT

9-OCT 02:37 PM

**RESET TO DEFAULT**

+ RESET +

+ NEXT

9-OCT 02:37 PM

**LANGUAGE**

English

+ CHANGE +

+ NEXT

9-OCT 02:37 PM

**GRAPH SNR**

+ CHANGE +

+ NEXT

---

---

---

---

---

---

---

---

---

---

---

---

### Making a Change: Screening

9-OCT 02:37 PM

**MIN VALUE OFF**

+ CHANGE +

+ NEXT

**DPOAE:** Min Amp ON enables -5dB as an additional requirement to PASS

**TEOAE:** Min Amp ON enables -5 OR -10dB as an additional requirement to PASS

000/002 15-DEC 07:43 AM

**MIN VALUE ON**

+ CHANGE +

+ NEXT

---

---

---

---

---

---

---

---

---

---

---

---

### Making a change: SNR or Value Graph?

9-OCT 02:37 PM

**GRAPH SNR**

+ CHANGE +

+ NEXT

000/002 15-DEC 07:43 AM

**GRAPH VALUE**

+ CHANGE +

+ NEXT

---

---

---

---

---

---

---

---

---

---

---

---

### Making a Change: Value Graph with Norms



Set to ON: Boys Town Norms will appear on Corti



---

---

---

---

---

---

---

---

### Why Norms on Corti?

- Norms will display for all diagnostic DP protocols
- Can be used as a guide when interpreting results
- Provides instant interpretation assistance before results are imported to the Data Manager



---

---

---

---


---

---

---

---

### What Do Norms Look Like?



---

---

---

---


---

---

---

---

### What Do Norms Look Like?



---

---

---

---

---


---

---

---

### Making a Change: Min Amp

000/002 15-DEC 07:48 AM  
MIN VALUE ON  
+ CHANGE +  
+ NEXT



---

---

---

---


---

---

---

---

### Min Amp + Norms



---

---

---

---

---

---

---

---

## Advanced Options Diagnostic



---

---

---

---

---

---

---

---

## Advanced DP Menu

Diagnostic Protocols

- DP 1.5 to 6.0 (6 freq)
- DP 1.6 to 8.0 (12 freq)
- DP 1.5 to 12.0 (12 freq)

Diagnostic protocols do not have pass/refer results. Audiologist is tasked with test interpretation.



---

---

---

---



---

---

---

---

## DPOAE: Advanced Parameters

- L1/L2 Intensity 
- SNR 
- Averaging Time per Frequency
- # Freq to Pass



---

---

---

---

---

---

---

---

### How to Get to DP Advanced Options Menu?

1. Select a diagnostic DP protocol
2. From the Corti main screen, press the down button until you reach the time and date screen:



---

---

---

---

---

---

---

---

### From the Date screen:

2. Press and **hold the down** button for 3 seconds or until the green "READY" light turns off. You should see:



---

---

---

---

---

---

---

---

### From the New BT Device:

- Press and **hold the down** button for 3 seconds or until the green "READY" light turns off. You should see:



---

---

---

---

---

---

---

---

### DP Advance Options:



Norms will not appear if L1 or L1 levels are changed.




---

---

---

---

---

---

---

---

---

---

### 6 Freq Protocols TE or DP



Allows user to enable Pass criteria for the specified number of frequencies

SNR/Min Amp/# of Freq

If implemented a Pass/Refer will appear at the end of the test




---

---

---

---

---

---

---

---

---

---

### 12 Freq Protocols DP



Allows user to enable Pass criteria for the specified number of frequencies

SNR/Min Amp/# of Freq

P will appear on report/DM per freq. Overall pass/refer will not appear




---

---

---

---

---

---

---

---

---

---

### How to Save Changes?

0 FREQS FOR PASS  
+ CHANGE  
+ NEXT

RESET PROTOCOL?  
+ RESET  
+ NEXT

SAVE PROTOCOL  
+ SAVE  
+ DONE

CUSTOM PROTOCOL SAVED  
+ DONE

---

---

---

---

---

---

---

---

### Modified Protocol “\*”

000/002 13-JAN 12:53

DP\* 1.5-6.0  
+L TEST R  
+ CHANGE

---

---

---

---

---

---

---

---

### TE Advanced Options

AVERAGING TIME = 64  
+ CHANGE  
+ NEXT

PASS SNR = 4  
+ CHANGE  
+ NEXT

0 FREQS FOR PASS  
+ CHANGE  
+ NEXT

---

---

---

---

---

---

---

---



New Feature: Printing Feature



---

---

---

---

---

---

---

---

### Printing Options

- Print to small thermal printer
- Print full page results using Data Manager
- NEW: Auto Print



---

---

---

---

---

---

---

---

### Auto Print

- Tool where user can by-pass data management functionality while maintaining the ability to quickly print OAE test results.
- For sites that want to print test results but use a different system for reporting test results.
- For customers who do not need a database or do not want to transfer names and data back and forth.



---

---

---

---

---

---

---

---

## NEW: Auto Print



PDF, any PC based printer OR PC based label printer



---

---

---

---

---

---

---

---

## Report Options: Norms, Name, Logo



Name can be added if in L/R mode.



---

---

---

---

---

---

---

---

## How does it work?

- Set Up printer and report settings
- Minimize Auto Print software
- As soon as Corti is connected to PC, results will be accessed and printed
- Can connect via micro USB cable OR while in the Corti cradle



---

---

---

---

---

---

---

---

### Print to PDF

- When in L/R mode, enter patient name
- Report will print to any location on the PC
- Report can be uploaded into EMR system or attached to audiometric report



---

---

---

---

---

---

---

---

### Print to PC Based Label Printer

- Can rapidly print any number of results to high quality labels
- When in L/R mode, can add name to label



---

---

---

---

---

---

---

---



### New Feature: Import Names



---

---

---

---

---

---

---

---

## Import Patient Names

- Quick way to add names to the Data Manager
- Can import an unlimited number of patient names from a spreadsheet




---

---

---

---

---

---

---

---

## How does it work?

Import patient list from excel spreadsheet




---

---

---

---

---

---

---

---

## How?



- Match GSI Database field name to the spreadsheet field name.
- Click OK!
- Names will appear in the Patients Tab




---

---

---

---

---

---

---

---

## Conclusion

- All new features will be available with the next release of software
- Software is in final testing



---

---

---

---

---

---

---

---

## Have Questions?

Sherrie Weller  
[srw@grason-stadler.com](mailto:srw@grason-stadler.com)

Karen Morris  
[Karen@grason-stadler.com](mailto:Karen@grason-stadler.com)

Laura Prigge  
[lcp@grason-stadler.com](mailto:lcp@grason-stadler.com)

GSI Marketing  
[marketing@grason-stadler.com](mailto:marketing@grason-stadler.com)



---

---

---

---

---

---

---

---

## Thank You!

Visit the Grason-Stadler Website for additional CEU opportunities  
[www.grason-stadler.com/e-learning](http://www.grason-stadler.com/e-learning)

Sign-up for E-Learning Alerts  
<http://www.grason-stadler.com/e-learning/e-learning-alerts>



---

---

---

---

---

---

---

---