
 **Education**
The Science of Hearing Aid Fitting

On-Ear Verification of Open Fittings

Presenter — David J. Smriga, M.A.
Senior Audiology Consultant
Audioscan Education

 **audioscan**[®]
Hearing Instrument Fitting Systems

www.audioscan.com USA 800-265-2093 519-268-3313 info@audioscan.com
Copyright 2016 Etymonic Design Incorporated. All rights reserved.

IMPORTANT NOTICE

- Unless otherwise noted, the contents of this webinar are copyright of Audioscan, a Division of Etymonic Design Inc. and are intended for use by webinar attendees only
- The name 'Audioscan' and 'Education' as used in these materials are trademarks of Etymonic Design Inc.
- Any subsequent recording, presentation, adaptation, publishing or sharing of any kind of this material is prohibited without first securing the express written permission of Audioscan.



Learning Objectives

- Correctly identify and objectively verify an open fitting condition with any patient.
- Correctly identify open fit candidates using real ear measures.
- Correctly fit open fit candidates using real ear measures.


Audioscan® Education

Copyright 2016 Etymonic Design Incorporated. All rights reserved.



What Is An Open Fitting?

- *Any style of hearing instrument (and its associated “plumbing”) that creates minimal occlusion of the ear canal when worn.*



Staab, W., Open Earmold Hearing Aid Fittings”, HHTM January 2016

McCabe, E., Galster, J. “Reducing Occlusion In CIC’s With Custom Venting” online pdf

Audioscan® Education

Copyright 2016 Etymonic Design Incorporated. All rights reserved.



Why An Open Fitting?

- Primary value:
 - Addresses occlusion and own voice complaints
 - Particularly for patients with little hearing loss in the lower frequencies
 - Provides a more natural sound quality
- Primary limitation:
 - Increased feedback potential will limit available gain
 - Best suited for high frequency loss

Audioscan® Education

Copyright 2016 Etymonic Design Incorporated. All rights reserved.



Digital/Thin-Tube/RIC

- Contributions:
 - Digital feedback management technology has facilitated greater access to more useable gain in an open ear canal environment
 - Thin-tube and RIC styles have advanced the cosmetic acceptance of the BTE form
 - Comfortable, non-occluding ear tips make these products easier to fit, adjust to and wear



Audioscan® Education

Copyright 2016 Etymonic Design Incorporated. All rights reserved.



As A Result. . .

- In the first quarter of 2016:
 - RIC products accounted for 64.5% of sales
 - BTE products accounted for 80.7% of sales




Audioscan® Education

Copyright 2016 Etymonic Design Incorporated. All rights reserved.



Not All RIC/Thin-Tubes Are Open



The six possibilities of BTE couplings used in today's fittings

Kuk, F., Baekgaard, L., "Hearing Aid Selection and BTE's: Choosing Among Various Open Ear and Receiver-in-Canal Options", Hearing Review, March 2008

Audioscan® Education

Copyright 2016 Etymonic Design Incorporated. All rights reserved.



What Is An “Open” Instrument?

- “Open” should NOT be defined by hearing aid design
 - Thin Tube
 - RIC
- “Open” defined by hearing aid coupling to the ear:
 - If coupling does not occlude the ear canal, it is an “Open” fitting

Audioscan® Education

Copyright 2016 Etymonic Design Incorporated. All rights reserved.



METHODS FOR CONFIRMING AN OPEN FITTING

- 1) The Real Ear Unaided Response (REUR) vs. the Rear Ear Occluded Response (REOR)

Probe tube placement

- Begin with otoscopic inspection of ear canal
- Place probe tip within 5mm of the eardrum
- Placement methods include:
 - Visually-assisted positioning
 - Use anatomy of external ear to guide positioning



Audioscan® Education

Copyright 2016 Etymonic Design Incorporated. All rights reserved.



Interested in learning more?

Visit the education and article section of Audioscan's expo page on AudiologyOnline

www.audiologyonline.com/expo/audioscan

Video tab: "Probe tube Placement"



Probe Tube Placement

September 18, 2015

Tips on proper probe tube placement for Speechmap and other real ear measures on Audioscan hearing instrument fitting systems.


Watch Video

Audioscan® Education

Copyright 2016 Etymonic Design Incorporated. All rights reserved.



Method For Determining an “Open” Fit



Once the probe tube has been properly inserted into the ear canal, select “Speechmap” in the On-ear column of the Verifit test selection menu.

Audioscan® Education

Copyright 2016 Etymonic Design Incorporated. All rights reserved.



REUR




On-Ear result obtained when only the probe tube is in the ear canal

Audioscan® Education

Copyright 2016 Etymonic Design Incorporated. All rights reserved.



REOR (Open dome)




On-Ear result obtained when the probe tube and a **muted open-dome** RIC instrument is in the ear canal

Audioscan® Education

Copyright 2016 Etymonic Design Incorporated. All rights reserved.



REUR & REOR Superimposed




Note that the two measurements are virtually identical, indicating that the physical presence of the hearing instrument/plumbing does not change ear acoustics

Audioscan® Education

Copyright 2016 Etymonic Design Incorporated. All rights reserved.



REOR (Closed Dome)




On-Ear result obtained when the probe tube and a muted **single-wall closed-dome** RIC instrument is in the ear canal

Audioscan® Education

Copyright 2016 Etymonic Design Incorporated. All rights reserved.



REUR & REOR Superimposed




Now, the two measurements are not identical, indicating that the presence of the hearing instrument/plumbing is indeed changing ear acoustics

Audioscan® Education

Copyright 2016 Etymonic Design Incorporated. All rights reserved.



REOR (Ventless Earmold)




On-Ear result obtained when the probe tube and a muted BTE with a **Double-walled power dome** is in the ear canal

Audioscan® Education

Copyright 2016 Etymonic Design Incorporated. All rights reserved.



REUR & REOR Superimposed



Note the significant difference between the two measures when this magnitude of occlusion is present

Audioscan® Education

Copyright 2016 Etymonic Design Incorporated. All rights reserved.



METHODS FOR DETERMINING AN OPEN FITTING

2) OCCLUSION EFFECT TEST




Occlusion Effect Test

The screenshot displays the Verifit software interface, which is used for hearing aid verification. The interface is divided into several sections:

- Verifit Logo:** Located at the top left, with the tagline "Total Control & Complete Verification".
- Test box:** A central area with two columns of settings. The left column is labeled "On-ear" and includes:
 - Speechmap
 - Directional
 - Noise reduction
 - Manual control
 - Calibration
 The right column is labeled "Test box" and includes:
 - WRECD
 - ANSI / IEC
 - Distortion
 - Input/Output
 - Battery drain
 - Multicurve
 - Telecoil
- Bottom Bar:** Contains icons for "Left/Right", "Print", "Session", "Help", and "Setup".
- Right Panel:** Displays the "Occlusion Effect Test" results. It includes:
 - Max TM SPL:** Set to 120.
 - Stimulus Filter:** Set to "Off" and "C Weighted".
 - Start test:** A button to initiate the test.
 - Probe microphone:** A graph showing the frequency response of the probe microphone, with a scale from 0 to 130 dB.
 - Reference microphone:** A graph showing the frequency response of the reference microphone, with a scale from 0 to 130 dB.
 - Occlusion:** A graph showing the frequency response of the occlusion effect, with a scale from 0 to 30 dB.



Occlusion Test Result: Open Dome




Audioscan® Education

Copyright 2016 Etymonic Design Incorporated. All rights reserved.



Occlusion Test Result: Closed Dome




Audioscan® Education

Copyright 2016 Etymonic Design Incorporated. All rights reserved.



Occlusion Test: Power Dome



Audioscan® Education


Copyright 2016 Etymonic Design Incorporated. All rights reserved.



OPEN FIT ACOUSTIC CONSIDERATIONS

1) VENTING

Comparative Vent Effects



Effect of vent diameter and open fitting on hearing aid output


Kuk, F., Baekgaard, L., "Hearing Aid Selection and BTE's: Choosing Among Various Open Ear and Receiver-in-Canal Options", Hearing Review, March 2008

Audioscan® Education

Copyright 2016 Etymonic Design Incorporated. All rights reserved.



NAL-NL2



Audioscan® Education

Copyright 2016 Etymonic Design Incorporated. All rights reserved.




OPEN FIT ACOUSTIC CONSIDERATIONS

2) EAR CANAL RESONANCE



Comparative Vent Effects




Mueller, H.G., Ricketts, T.A., "Open-canal Fittings: Ten Take-home Tips: Hearing Journal, November, 2006

Audioscan® Education

Copyright 2016 Etymonic Design Incorporated. All rights reserved.



Ear Canal Resonance



Difference in probe measures obtained in a Zwislocki coupler driven with steady state noise between a response obtained with the probe sealed at the coupler opening (zero line) and a response obtained with the probe simply placed in the coupler (measured line).

Mueller, H.G., Ricketts, T.A., "Open-canal Fittings: Ten Take-home Tips: Hearing Journal, November, 2006

Audioscan® Education

Copyright 2016 Etymonic Design Incorporated. All rights reserved.



OPEN FIT ACOUSTIC CONSIDERATIONS 3) ACOUSTIC PATHWAYS

In an Open Fit, There Are Several Pathways to Get SPL to the Eardrum


- Direct pathway
 - Sounds that reach the eardrum unprocessed by the hearing aid
- Amplification (Indirect) pathway
 - Sounds that reach the eardrum via the output signal of the hearing aid

Audioscan® Education

Copyright 2016 Etymonic Design Incorporated. All rights reserved.



Acoustic Pathways



Audioscan® Education

Copyright 2016 Etymonic Design Incorporated. All rights reserved.



The Advantage of REAR Measurement

- An REAR displays the combined acoustic pathway result
 - It measures and displays the aided SPL at the probe tube tip regardless of how that SPL was delivered
- An REIG displays only the hearing aid pathway

Audioscan® Education

Copyright 2016 Etymonic Design Incorporated. All rights reserved.



OPEN FIT VERIFICATION CONSIDERATIONS

 **Education**
The Science of Hearing Aid Fitting

If an Open Fitting Is Present. . .

- Pre-fit (in the test box) is off the table
 - There is no such thing as an “Open” coupler
- All fitting measures need to be done on the ear
- “Concurrent” calibration (equalization) should not be used

Audioscan® Education


Copyright 2016 Etymonic Design Incorporated. All rights reserved.



Impact of sound calibration method

For open fitting:

- Outflow from ear canal received by reference mic.
 - Lowers speaker output
- Resulting measured output will be reduced




Audioscan® Education

Copyright 2016 Etymonic Design Incorporated. All rights reserved.



Open Fitting Recommendation

- Use 'stored' equalization instead of 'concurrent' equalization
 - Presents the calibration signal independent from the speech signal
 - Hearing aid is on the ear, but muted during the calibration signal event
 - Thus, outflow from ear canal can not affect speaker level and your measurement



Audioscan® Education

Copyright 2016 Etymonic Design Incorporated. All rights reserved.



Effect of equalization approach




Figure 11. REAR findings for an OC fitting using concurrent equalization (reference microphone on) versus stored equalization (reference mic off)

Mueller, H.G., Ricketts, T.A., "Open Canal Fittings" Ten Take Home Tips", Hearing Journal, v9:11, Nov. 2006

Audioscan® Education

Copyright 2016 Etymonic Design Incorporated. All rights reserved.



OPEN FIT VERIFICATION PROCEDURE


1) VERIFYING PATIENT CANDIDACY



Patient Candidacy


- Select “Speechmap” from the On-Ear menu
- Enter audiogram into verification system
- Select a fitting target (i.e., NAL-NL2)
- Place the probe tube only in the patient’s ear(s)
- Run and record the REUR for 65dB speech
- Compare the displayed LTASS to the displayed target





Audioscan® Education



Copyright 2016 Etymonic Design Incorporated. All rights reserved.



Audioscan® Education

Copyright 2016 Etymonic Design Incorporated. All rights reserved.






OPEN FIT VERIFICATION PROCEDURE

2) VERIFYING AIDED PERFORMANCE (ON-EAR)




Open-Fit On-Ear Verification Procedure



Select "Speechmap" in the On-ear column of the Verifit test selection menu.



Audiogram Being Used In This Example




Audioscan® Education

Copyright 2016 Etymonic Design Incorporated. All rights reserved.



1) Select "Open" in Instrument menu



Audioscan® Education

Copyright 2016 Etymonic Design Incorporated. All rights reserved.



2) Place probe mic and HI on ear per usual




Turn OFF or Mute the Hearing Instrument

Audioscan® Education

Copyright 2016 Etymonic Design Incorporated. All rights reserved.



3) Select Play in Test 1 and store equalization when prompted



Once Equalization is completed, turn ON the hearing instrument

Audioscan® Education


Copyright 2016 Etymonic Design Incorporated. All rights reserved.




Re-equalizing sound field

If sound field changes at *any time* during testing, including if patient or clinician moves, the sound field eq must be repeated.

Click the “EQ” button to interrupt the running test to allow a new equalization to be stored.



Verifit2




Axiom

Audioscan® Education

Copyright 2016 Etymonic Design Incorporated. All rights reserved.



4) Run Test 1: Calibrated Soft Speech




Adjust HI gain until best target match has been achieved

Audioscan® Education

Copyright 2016 Etymonic Design Incorporated. All rights reserved.



5) Run Test 2: Calibrated Average Speech




Audioscan® Education

Copyright 2016 Etymonic Design Incorporated. All rights reserved.



5) Run Test 2: Calibrated Average Speech




“Showing” unaided speech banana (grey) helps define HI difference

Audioscan® Education

Copyright 2016 Etymonic Design Incorporated. All rights reserved.



6) Run Test 3: MPO




Adjust hearing aid output to approximate MPO target on Speechmap Screen

Audioscan® Education

Copyright 2016 Etymonic Design Incorporated. All rights reserved.



Optional: Run REUR in Test 4




Actual unaided speech in ear canal (gold) vs. predicted unaided speech (grey)

Audioscan® Education

Copyright 2016 Etymonic Design Incorporated. All rights reserved.



Optional: Run REUR In Test 4



Actual unaided (gold) and aided (pink) normal conversational speech in ear canal

Audioscan® Education

Copyright 2016 Etymonic Design Incorporated. All rights reserved.



AudiologyOnline

Home Contact Help Downloads 188-755-2945

Continuing Education Career Center Journal AudiologExpo Resources

The new

Verifit 2
Total Control & Complete Verification

- Because people have 2 ears
- Because 2 aids interact
- Because it's 2x faster
- Because it measures 2x the bandwidth
- Because Verifit2 integrates seamlessly

Home / Expo Home / Audioscan



Audioscan

Hearing Instrument Fitting Systems



Audioscan® is the #1 hearing instrument fitting system manufacturer in North America with more systems in use than all other manufacturers combined! Audioscan products are renowned world-wide because their users experience tangible direct savings, scientific verification of audiology, reduced returns, and more satisfied patients. Audioscan has recently introduced its most anticipated product in its 25 year history. The new Verifit2 is the culmination

Contact Information

Audioscan
A Division of Etymonic Design Inc.
271 Loring St.
Dorchester, ON N9L 1G4
CA
Phone: 519-268-7513 CANADA
800-265-2929 USA

General: info@audioscan.com

Support: support@audioscan.com

Location: Find a nearby distributor

Visit us: audioscan.com

Sign up for software update notifications!

Click Here

Interested in learning more?

Visit the education and article section of Audioscan's expo page on AudiologyOnline

www.audiologyonline.com/expo/audioscan

Education tab: "Speechmapping As A Fitting and Counseling Tool"

Articles tab: "Speechmapping As A Fitting and Counseling Tool"

Audioscan® Education

Copyright 2016 Etymonic Design Incorporated. All rights reserved.



Directional Mic Test

OPEN FITTING CONSIDERATIONS
ON EAR AND TEST BOX MODES




Directionality Test (ON EAR)


Rear speaker

Aided ear with
probe tube positioned

Verifit in On-ear
directional mode



On Ear Directional Mic Test Result




Audioscan® Education

Copyright 2016 Etymonic Design Incorporated. All rights reserved.



Aided On Ear REAR for Speech Compared to On Ear Directional Mic Test



Frequency Range of Aided Region

Frequency Range of Aided Region

Audioscan® Education

Copyright 2016 Etymonic Design Incorporated. All rights reserved.



On Ear Directional Mic Measurement Considerations


- Open venting will likely mask directional mic functionality in the lower frequencies
- Directional mic port angle can effect the magnitude of curve separation
 - The closer to horizontal, the better the curve separation
- Head and pinna effects play a role in the directional mic result obtained.

Audioscan® Education

Copyright 2016 Etymonic Design Incorporated. All rights reserved.



Directionality Test Test Box



Audioscan® Education

Copyright 2016 Etymonic Design Incorporated. All rights reserved.



Directional in Test-box - RIC

Front speaker

Rear speaker



Audioscan® Education

Copyright 2016 Etyronic Design Incorporated. All rights reserved.



Recommended Positioning – Verifit2




Audioscan® Education

Copyright 2016 Etyronic Design Incorporated. All rights reserved.



Directional Test Box Result




Audioscan® Education

Copyright 2016 Etymonic Design Incorporated. All rights reserved.




Aided On Ear REAR for Speech Compared to Test Box Directional Mic Test



Audioscan® Education

Copyright 2016 Etymonic Design Incorporated. All rights reserved.





AudiologyOnline

Continuing Education Career Center Journal AudiologyExpo Resources

The new Verifit 2
Real-World & Complete Verification

- Because people have 2 ears
- Because 2 aids interact
- Because it's 2x faster
- Because it measures 2x the bandwidth
- Because Verifit2 integrates seamlessly

Audioscan
Hearing Instrument Fitting Systems

Audioscan® is the #1 hearing instrument fitting system manufacturer in North America with more systems in use than all other manufacturers combined. Audioscan products are renowned world-wide because their owners experience lengthy three average scientific verification of audibility, reduced returns, and more satisfied patients. Audioscan has recently introduced its most anticipated product in its 25 year history. The new Verifit2 is the culmination of years of research and development.

Contact Information
Audioscan
A Division of Tympanic Design Inc.
20 Luning Dr.
Dorchester, ON N2L 1G4
CA
Phone: 519 265-2213 CANADA
800 365 2861 USA

General: info@audioscan.com
Support: support@audioscan.com
Location: Find a nearby distributor
Visit us: audioscan.com

Sign up for software update notifications!

Audioscan® Education
Copyright 2016 Tympanic Design Incorporated. All rights reserved.

Noise Reduction Test

On Ear and Test box

Digital Noise Reduction Properties


- Digital algorithm programmed to recognize “non-speech” elements of incoming stimulus
 - Operates independently in bands
 - Analyzes incoming signal modulation
- Can vary in terms of time constants
 - Typically, slow attack, fast release

Audioscan® Education

Copyright 2016 Etymonic Design Incorporated. All rights reserved.



On Ear Noise Reduction Test Result




Audioscan® Education

Copyright 2016 Etymonic Design Incorporated. All rights reserved.




Aided On Ear REAR for Speech Compared to On Ear Noise Reduction Test



Audioscan® Education
Copyright 2016 Etymonic Design Incorporated. All rights reserved.




Test Box Noise Reduction Test Result



Audioscan® Education
Copyright 2016 Etymonic Design Incorporated. All rights reserved.



Aided On Ear REAR for Speech Compared to Test Box Directional Mic Test



Audioscan® Education

Copyright 2016 Etymonic Design Incorporated. All rights reserved.



Education
The Science of Hearing Aid Fitting

Summary

audioscan
Hearing Instrument Fitting Systems

www.audioscan.com USA 800-265-2093 519-268-3313 info@audioscan.com
Copyright 2016 Etymonic Design Incorporated. All rights reserved.

Key Open Fit Verification Considerations

- An open fitting exists when any hearing aid is coupled to the ear in a way that does not occlude the ear
 - Defined more by plumbing than HA style
 - An occlusion-free ear canal can be verified
 - REUR/REOR test
 - Occlusion effect test

Audioscan® Education

Copyright 2016 Etymonic Design Incorporated. All rights reserved.



Key Open Fit Verification Considerations

- There are several important acoustic properties associated with open fittings
 - Vent effects
 - Ear canal resonance
 - Combined acoustic pathways

Audioscan® Education

Copyright 2016 Etymonic Design Incorporated. All rights reserved.



Key Open Fit Verification Considerations

- Regarding open fit verification:
 - Pre-fitting (test-box Speechmap) is not an option
 - Concurrent microphone equalization (in REM) is replaced with stored microphone equalization
 - All fitting/fine-tuning steps will need to be completed On-Ear

Audioscan® Education

Copyright 2016 Etymonic Design Incorporated. All rights reserved.



Key Open Fit Verification Considerations

- Advantages associated with the REAR:
 - Can be used to “verify” open-fit candidacy
 - Displays an aided result that accounts for the combined acoustic signal conditions unique to an open fitting

Audioscan® Education

Copyright 2016 Etymonic Design Incorporated. All rights reserved.



Key Open Fit Verification Considerations

- Directional mic can be verified On-Ear with an auxiliary speaker
 - Demonstrates the effects of open venting on directional performance
 - Shows the value of D-Mic under these conditions
- Noise reduction can be verified On-Ear
 - Demonstrates the effects of open venting on this feature
 - Shows the value of noise reduction under these conditions.

Audioscan® Education

Copyright 2016 Etymonic Design Incorporated. All rights reserved.



Key Open Fit Verification Considerations

- The test box can still be a valuable tool in open fit verification
 - Test box directional mic test results can be compared to REAR results to identify directional mic advantage potential
 - Test box noise reduction test results can be compared to REAR results to identify noise reduction advantage potential

Audioscan® Education

Copyright 2016 Etymonic Design Incorporated. All rights reserved.





 **Education**
The Science of Hearing Aid Fitting

Thank You!

info@audioscan.com

 **audioscan**[®]
Hearing Instrument Fitting Systems

www.audioscan.com USA 800-265-2093 519-268-3313 info@audioscan.com
Copyright 2016 Etymonic Design Incorporated. All rights reserved.