



## uLearn series: Verification of Unitron's signature features with Otometrics Aurical Verification System

Presenters:

Josie Helmbrecht AuD

**Unitron**

Sarah Carls AuD

**Unitron**

Wendy Switalski MBA AuD

**Audiology Systems**

**unitron.** Hearing matters

## Agenda

Introduction

Advantages and application of advanced features

Unitron fitting philosophy

Noise reduction verification

Speech Enhancement verification

Directionality verification

Conclusion

**unitron.**

## Unitron's philosophy


### Industry Standard Fitting Formulas

NAL and DSL are based on a combination of theory and empirical evidence and have been clinical standards for several decades.

The benefit of an evidence based approach is that it is proven to provide maximum speech intelligibility.

unitron.

## TrueFit 3.3 Default fitting formula



We default to NAL-NL2.  
Fitting formula defaults can  
be changed under  
preferences in the software.

We offer DSL adult and  
pediatric, NAL-NL1, and NAL-  
NL2 for tonal languages as  
other options.

## Unitron's Feature philosophy

Based on a proven strategy  
for providing the greatest improvement  
in speech intelligibility

Balancing listening  
goals across  
7 environments

# Sound Conductor

Sound Conductor knows the listening **goals** across SoundNav's 7 environments.

## Conversations

- Conversation in quiet
- Conversation in a small group
- Conversation in a crowd
- Conversation in noise

Sound Conductor **goal** is understanding conversations

## No conversation

- Quiet
- Noise

Sound Conductor **goal** is comfort with awareness

## Music

- Music

Sound Conductor **goal** is beautiful sound quality with enhanced lyrics


unitron.

## Goals:

Understanding conversations

Comfort with awareness

Natural Sound quality



Sound Conductor

unitron.

Three features have  
been optimized to work  
together

SC28

Verifying Advanced Features –  
Otometrics Aurical  
FreeFit PMM and HIT systems



## Slide 12

---

### SC28 Wendy Content

Advanced Features verification

-interactive- turn off certain features

Carls, Sarah, 1/30/2017

## Verification using Aurical

Aurical PMM and Aurical HIT can be used to confirm that advanced features are working and also to quantify the how the HI reacts in response to various signals.

Due to the interactive nature of adaptive features, it is helpful to turn off all features except for the one being evaluated.

Duration, level, and type of signals needed will vary. Keep this in mind especially if making comparisons between instruments.

While advanced features can be evaluated with both On-Ear (PMM) and Coupler (HIT) measurements, using HIT provides the following advantages:

- More consistent and acoustically isolated environment

- Multiple speaker sources for providing input from both front and rear directions

unitron.

## Verification using Aurical

RECD Unaided Response Occluded Response Aided Response Noise Reduction FreeStyle

Aurical is designed for advanced feature verification using both the Noise Reduction tab for NR, and the FreeStyle tab for all other features

- Offers access to 60+ input signals, including speech, noise, and common environmental sounds, including Speech + Noise pre-mixed options

- Also allows for mixing of multiple signals at varied SNR to replicate challenging listening environments

The customizable 'User Tests' are a valuable tool when performing these measurements. For copies of the User Tests used throughout the presentation today, please contact us.


unitron.

SC5  
SC29

# 1. Noise reduction

SC6

## Noise reduction (NR)



## Slide 15

---

**SC5** Josie explains noise reduction

Carls, Sarah, 1/5/2017

**SC29** Sarah

Carls, Sarah, 1/30/2017

## Slide 16


---

**SC6** Josie

Carls, Sarah, 1/5/2017


SC7

## Noise reduction LD



SC8

## Noise reduction LD



## Slide 17

---

**SC7** Josie

Carls, Sarah, 1/5/2017

## Slide 18


---

**SC8** Josie

Carls, Sarah, 1/5/2017

SC30

# Noise reduction LD



Unitron TrueFit 3.3.0 Verification - AG

File Edit Window Help

Unitron TrueFit File Forms Hearing Instrument Options Help

Client Instruments Fitting End Fitting Technology Change Acoustics Accessories Pre-Fitting In Situ

Quick Fit Conversation in quiet Detect X L Flextrial Moxi Fit Pro North xS RIC Moxi Fit Pro North xS RIC Connection +

Detect iCube II

Simulate

Stride Pro 800 700 600 500

Moxi Pro 800 700 600 500

Moxi<sup>2</sup> Pro 20 16 10 E

Quantum<sup>2</sup> Pro 20 16 10 E

Max 20 6 E

Moxi

Current Product in this Session

Moxi Fit Pro North xS RIC (Flextrial)

Moxi Fit Pro North xS RIC

Unitron TrueFit

Connection successful.  
 Moxi Fit Pro North RIC (Flextrial)  
 SN: 1509X085P

No hearing instrument detected.

Moxi Fit Pro North

- SoundNav and Sound Conductor work together in 7 environments
- SpeechZone 2 to focus on speech coming from any direction even in the toughest listening situations
- MyMusic with automatic binaural treatment
- Unwanted sounds are faded into the background
- Hear phone in both ears with Binaural Phone
- Restores ear's natural acoustics with Pinna Effect
- Wireless streaming accessories
- 20 Channels



312 battery  
 113/47 dB (2cc)  
 Push Button  
 Telecoil  
 Easy-1


## Slide 19

---

**SC9** Josie  
Carls, Sarah, 1/5/2017

**SC30** Josie  
Carls, Sarah, 1/30/2017





# Noise Reduction Video Capture

Wendy narration


## 2. Speech enhancement

SC4

### Speech enhancement

- Speech enhancement provides additional gain to speech sounds detected in each band
- Graph explanation:
  - SE is applied starting at the expansion kneepoint
  - SE applies the same amount of gain to soft and average inputs
  - SE gain is rolled off for loud inputs as it approaches saturation to avoid distortion and maintain excellent sound quality

(The Gain function shown here is linear for ease of understanding, however the gain function often also has compression.)




unitron.

**SC4** Josie explains speech enhancement

Carls, Sarah, 1/5/2017

The screenshot shows the UniTron TrueFit software interface. The main window is titled 'Configure Features' and is set to the 'Input' configuration. On the left, there is a sidebar with categories: 'Universal' (Sound Conductor, Input, Wind Control, AntiShock, Phase Canceller) and 'Specific' (Binaural Phone, MyMusic, Automatic Controls). The main area is divided into sections: 'SoundNav', 'Manual Programs', and 'Wireless Programs'. The 'SoundNav' section is currently active and shows settings for 'Conversations' and 'Music'. The 'Conversations' section has four rows: 'In quiet', 'In a small group', 'In a crowd', and 'In noise'. Each row has two columns for 'R' (Right) and 'L' (Left) settings, with dropdown menus for 'Pinna Effect', 'Fixed Directional', and 'Fixed Wide with SpeechZone 2'. The 'Music' section has two rows: 'Music' and 'Music', with dropdown menus for 'Omni (Mic)'. The 'Manual Programs' section shows 'Telephone' and the 'Wireless Programs' section shows 'uPhone', 'uAudio', and 'uFM'. The top of the window has a menu bar (File, Forms, Hearing Instrument, Options, Help) and a toolbar with icons for Client Instruments, Fitting, End Fitting, Feedback Optimization, Program Manager, Tuning, and Data Logging. The bottom of the window shows a status bar with 'Ready' and 'NUM'.



**SC26** Wendy

Carls, Sarah, 1/19/2017

# Speech Enhancement Video Capture

Wendy narration

unitron.

SC11

## 3. Directional strength

**SC11** Josie explains directional microphones


Carls, Sarah, 1/5/2017

SC12

# Directional input options

Available input options:

- Omni
- Pinna Effect
- Fixed Wide directional
- Fixed directional
- Adaptive directional



SC19 trueFit 3.10 Verification AG

Unitron TrueFit File Forms Hearing Instrument Options Help

Client Instruments Fitting End Fitting

Feedback Optimization Program Manager Tuning Data Logging

Quick Fit

Conversation in a crowd

R Detect X L

Flextrial Maxi Fit Pro North X S RIC Maxi Fit Pro North X S RIC

Connection +

Configure Features

Universal

Sound Conductor

Input

Wind Control

AntiShock

Phase Canceller

Specific

Binaural Phone

MyMusic

Automatic Controls

Select a sound file

Input

SpeechZone 2

Conversations R L

In quiet	Pinna Effect	Pinna Effect
In a small group	Fixed Directional	Fixed Directional
In a crowd	Adaptive with SpeechZone 2	Adaptive with SpeechZone 2
In noise	Fixed Wide with SpeechZone 2	Fixed Wide with SpeechZone 2
No conversations		
Quiet	Pinna Effect	Pinna Effect
Noise	Omni (Mic)	Omni (Mic)
Music		
Music	Omni (Mic)	Omni (Mic)

+ Manual Programs Telephone

+ Wireless Programs uPhone uAudio uFM

## Slide 31

---

### SC12 Josie



Carls, Sarah, 1/5/2017

## Slide 32

---

### SC19 Wendy reviews video capture

Carls, Sarah, 1/5/2017



### Slide 33

---

**SC20** Josie show's set up through screenshots or video capture



Carls, Sarah, 1/5/2017

### Slide 34

---

**SC21** Wendy shows functionality via video capture

Carls, Sarah, 1/5/2017



# Directionality Video Capture

Wendy narration

unitron.

SC15

# Conversations


## Slide 38

---

**SC15** Not sure if we want to use the following slides, but I left them in bc they show how much NR and SE we use.

Carls, Sarah, 1/5/2017

SC16




SC17

# Conversations environments

## Sound Conductor

### Maximum strengths:

- Speech enhancement **4dB**
- Noise reduction **-7dB**
- Directionality independently configurable



unitron

**Slide 39**

---

**SC16 ?**

Carls, Sarah, 1/5/2017

**Slide 40**

---

**SC17** Carls, Sarah, 1/5/2017

# Comfort with awareness: Quiet and Noise

Sound Conductor

**Maximum strengths:**

- Speech enhancement strength **0**
- Noise reduction strength **-7dB**
- Directionality independently configurable



unitron.

# Music

Sound Conductor

**Maximum strengths:**

- Speech enhancement strength **4dB**
- Noise reduction strength **0dB**
- Directionality independently configurable



unitron.

## Noise Reduction with Speech Added Video Capture

Wendy narration of noise reduction releasing when speech is present


Good place to answer question / address why is it preferred to verify features in isolation versus a more real world combine approach

unitron.

## Personalization


# Personalization Sound Conductor

Sound Conductor can be adjusted independently in all programs



# Personalization Sound Conductor

Sound Conductor can be linked/unlinked in all programs



# Personalization of Sound Conductor

Adjusting Sound Conductor will influence the strength of speech enhancement, noise reduction, as well as the strength of the directionality for **Fixed Directional and Adaptive Directional inputs** (Fixed directional shown here)

**Feature Details**

To increase speech understanding move the slider to the right

Speech Enhancement  
Noise Reduction

**SoundNav**

- Conversation in quiet
- Conversation in a small group
- Conversation in a crowd
- Conversation in noise
- Quiet
- Noise
- Music

unitron.

# Personalization of Sound Conductor

Directionality is independently configurable in all environments

**SoundNav** SpeechZone 2

Conversations		R	L
In quiet		Pinna	Pinna
In a small group	◀	Fixed Directional	Fixed Directional
In a crowd		Fixed Wide with SpeechZone 2	Fixed Wide with SpeechZone 2
In noise		Fixed Wide with SpeechZone 2	Fixed Wide with SpeechZone 2
No conversations			
Quiet		Pinna	Pinna
Noise		Fixed Wide Directional	Fixed Wide Directional
Music			
Music		Pinna	Pinna

unitron.

# Conclusion

Questions?

unitron.