

Product Feature	Implementation	Laboratory Evidence	Real-World Evidence	Informing the Patient ("When I adjust this feature, it will allow the hearing to...")
Multiple Channels of Compression	Repackage sounds into residual dynamic range of the user, better able to shape gain within specific frequency response	Data shows speech intelligibility is improved with up to 5 channels in quiet (Woods et al, 2006) and up to 8 to 16 in noise (Yund & Buckles, 1995)	None	Effectively place more of the speech sounds into your "restricted" hearing range. Maximizing audibility while maintaining comfort will result in improved speech intelligibility.
Multi-channel AGC-O	Keep loud sounds under the user's threshold of discomfort while maintaining headroom across the entire frequency response	None	None	Keep virtually all loud sounds from becoming uncomfortably loud or intolerable.
Expansion	Minimize the negative effects of soft ambient noise, like microphone noise	Even when intelligibility for soft speech is negatively affected, patients still prefer expansion (Plyler, et al, 2005, 2006, 2007)	Limited evidence suggests "expansion on" is preferred to "expansion off" in everyday listening conditions	Minimize low levels of room noise that can be extremely annoying in quiet listening situations.
Automatic Feedback Reduction	Reduce acoustic feedback without sacrificing gain and headroom	Data (Johnson et al, 2007) shows that phase cancellation automatic feedback reduction systems do not negatively effect sound quality	None	Automatically detect whistling and eliminate or reduce it almost instantaneously. There are large individual differences in performance of this feature.
DNR - Modulation detection	Reduce gain in each channel where the signal is classified as noise	In paired comparison studies, patients prefer "DNR on" when listening in both low and high noise conditions (E.g., Ricketts and Hornsby, 2005)	Studies published in non-peer reviewed journals indicate that DNR is advantageous in everyday listening situations (Numerous studies - see Bentler, 2006 for review)	Improve listening comfort in noise but will not directly improve intelligibility in noise. DNR does no harm to your understanding of speech.
DNR - Impulse detection	Reduces gain of impulse sounds	Paired comparison study indicates that "impulse noise reduction on" is preferred over "impulse noise reduction off" (Keldser, et al, 2007)	Non-peer reviewed published study indicates this type of DNR makes impulse sounds less annoying. (Chalupper & Powers, 2007)	Makes intense impulse sounds less annoying.
DNR - Spectral subtraction or Wiener	Filter out the noise that occurs between gaps of syllables	Limited data suggests that Wiener filtering may improve perceived annoyance of noise (Mueller, et al, 2006)	Non-peer reviewed published study indicates this type of "DNR on" is preferred in everyday listening situations (Powers, et al 2006)	Improves listening comfort in noise, does no harm in quiet.

Audibility & Comfort

Comfort in Noise

Speech Intelligence in Noise	Convenience & Value-Added Extras	DNR - "reverberation tail reduction"	Reduces energy of sound envelope based on time of arrival and amplitude of signal	None	None	Makes highly reverberant places (e.g. gym) more tolerable
Fixed Directional Microphones	Improves speech understanding in noise, assuming you can spatially separate the speech & noise, and remember to push the memory button on the hearing aid.	When distance and reverberation effects are controlled, DMT technology improve signal to noise ratio 3 to 6dB (See Ricketts, 2005 for review)	Improves signal to noise ratio of listening area when speech and noise are spatially separated	30 to 50% of users prefer DMT on in noisy situations (Walden et al, 2004)	None	Improve speech understanding in noise, assuming you can spatially separate the speech & noise, and remember to push the memory button on the hearing aid.
Adaptive Directional Technology	Adaptive DMT does not outperform conventional fixed DMT, but it may be easier to use for many patients	In paired comparison studies, no significant differences between fixed and adaptive DMT benefit (Bentler, 2004, Ricketts, 2005)	The null of the polar pattern automatically changes as the primary noise source moves in space	There is little basis for customizing automatic directionality feature based on user preference data (Walden et al, 2007)	None	Adaptive DMT does not outperform conventional fixed DMT, but it may be easier to use for many patients
Multi-Channel Adaptive DMT	Adaptive DMT does not outperform conventional fixed DMT, but it may be easier to use for many patients	None	Polar patterns tuned to multiple frequencies automatically change depending on the input levels	None	None	Adaptive DMT does not outperform conventional fixed DMT, but it may be easier to use for many patients
Input Signal Classification	The hearing aid analyzes your listening situations and classifies the signal based on loudness duration and spectrum. "The brain of the hearing aid."	None	Uses modulation, intensity and time characteristics to classify sound and change electroacoustic parameters, if necessary	None	None	The hearing aid analyzes your listening situations and classifies the signal based on loudness duration and spectrum. "The brain of the hearing aid."
Trainable Volume Control	The hearing aid calculates volume control setting, etc, and automatically determines where you prefer it, and remembers the setting. For many people this feature makes the hearing aid easier to use	None	Hearing aid automatically sets the gain and/or compression parameters based on real time usage	Majority of participants preferred "trained" compression and noise suppression parameters to "untrained" condition (Zakis, et al, 2007)	None	The hearing aid calculates volume control setting, etc, and automatically determines where you prefer it, and remembers the setting. For many people this feature makes the hearing aid easier to use
Automatic Program Switching	Automatic program switching may make the hearing aid easier to use because it makes the decisions as to what the optimal program should be for you in various listening situations.	None	On-board classification system allows the hearing aid to change programs automatically depending in the characteristics of the incoming signal (automatic directional technology)	None	None	Automatic program switching may make the hearing aid easier to use because it makes the decisions as to what the optimal program should be for you in various listening situations.
Datalogging	Datalogging allows me to do a better job making precise adjustments to your hearing aids any time you are in the office for an appointment	None	Objectively records use of hearing aids, programmed used, etc.	Experts agree that data logging helps them do a better job of tailoring the acoustic parameters to the hearing loss (Mueller, 2007)	None	Datalogging allows me to do a better job making precise adjustments to your hearing aids any time you are in the office for an appointment
Data was compiled April, 2008						