

# Optimizing Your Patient Outcomes Using a Complete Solutions Approach

Phonak Naída V & Roger



**PHONAK**  
life is on

## Disclosures

- **Cheri Hebeisen, Au.D., CCC-A**
- Cheri Hebeisen, has been an audiologist for over 15 years. She earned her Doctor of Audiology degree in 2003 from the University of Florida and her Masters of Communication Disorders in 2000 from Louisiana State University Health Science Center in New Orleans. Prior to working as a clinical trainer for Phonak, her clinical experience included diagnostics for pediatrics through geriatrics, hearing aid and FM fittings and evaluations, as well as, auditory processing disorder evaluations.
- *Financial-Phonak employee who receives a salary for employment*
- *Nonfinancial-No relevant nonfinancial relationships exists*

**PHONAK**  
life is on

---


## Learner objectives

- Participants will be able to describe the solutions based approach for a patient with severe to profound hearing loss.
- Participants will be able to describe the new features and benefits of Phonak Naida V and Roger technology.
- Participants will be able to outline characteristics of the patient profile that can be fit with Naida V as well as the benefits associated with wearing Naida V.


---

**PHONAK**  
life is on

## Significant variability in performance in noise and in quiet



Speech in noise (SNR) vs pure tone average



Speech in quiet (% correct) vs pure tone average

- Mild to Moderate
- Moderately Severe to Severe


Souza, P. (2009). Severe Hearing Loss Recommendations for Fitting Amplification. *Audiology Online*, January 19

**PHONAK**  
life is on

## Defining the power patient

By the patient need?

- **Comfort and durability**
- **Natural sound quality and clarity**
- Hearing in **quiet**, near field is a challenge – in **noise** and over far field seem impossible without technology
- A hearing aid that can **adapt** to changing environments



**PHONAK**  
life is on

## Seeing the whole patient

A holistic approach:

1. Hearing and listening needs assessment, including the situations most important to the patient
2. Full audiological assessment
3. Recommend a complete solution that aims to meet the listening and communication needs of the patient and their communication partners



A complete solution often includes elements beyond the hearing aids. A solution that not only solves face to face communication but helps patients navigate hearing at home, at work, with family and in their social life.

**PHONAK**  
life is on

1994 | **AudioZoom**: Pioneering multi-microphone technology

2003 | **Multi-frequency FM**: First frequency FM system

A history of  
POWER INNOVATION

2002 | **Supero**: First super high power digital product

2007 | **Dynamic FM**: First adaptive FM

2008 | **Naída**: Introduced SoundRecover



2013 | **Roger**: First to use adaptive wireless transmission at 2.4 GHz



2013 | **Naída Q**



2011 | **Naída S**



2016 | **Naída V**

**PHONAK**  
life is on

## Phonak Portfolio of Solutions

**Bolero V**  
3 new designs



**Venture Chip**



**Music Performance**



**Audéo V**  
4 designs



**Virto V**



**AutoSense OS**



**NEW!**

**Naída V**



**30% less battery consumption**



**Roger EasyPen**



**CROS II**  
2 new designs



**PHONAK**  
life is on

## Phonak Power Portfolio: Naída V

### Naída V



Naída V-UP

Naída V-SP

Naída V-RIC

### Naída V + Roger



Naída V-UP  
Roger 19

Naída V-SP  
Roger 18

Naída V-RIC  
Roger 18

All with Volume Control, Push Button, Telecoil and EasyPhone

**PHONAK**  
life is on

## Colors for Naída V and Roger receivers

### Hair & Skin



P1

Sand Beige

P2

Amber Beige

P3

Sandalwood

P4

Chestnut

P5

Champagne

P6

Silver Gray

P7

Graphite Gray

P8

Velvet Black

### Fashion



P9

Ruby

Q1

Petrol

### Traditional




O1

Beige

Same color palette as Audéo V and Bolero V

**PHONAK**  
life is on

## Naída V fitting ranges



## Acoustic coupling for Naída V-RIC

Fitting Range	xReceiver	Hearing Instrument	Coupling Options
Mild to moderately-severe hearing loss	Standard (xS) 45/112 (2cc)	Phonak Audio V 10 312 312T 13 Phonak Naída V RIC	Open Dome Closed Dome Power Dome SlimTip Hard SlimTip Soft cShell (Hard or Soft)
Mild to severe hearing loss	Power (xP) 55/126 (2cc)	Phonak Audio V 10 312 312T 13 Phonak Naída V RIC	Open Dome Closed Dome Power Dome SlimTip Hard cShell (Hard or Soft)
Moderate to severe hearing loss	SuperPower plus (xSP plus) 65/133 (2cc)	Phonak Audio V 312 312T 13	cShell (Hard only)
Moderate to profound hearing loss	UltraPower (xUP) 66/136 (2cc)	Phonak Naída V RIC	cShell (Hard only)

**PHONAK**  
life is on

## Simply ingenious solutions you already know



### AutoSense OS

Improved speech understanding by 20%.



### Music

Fuller, richer, more natural listening experience.



### Better Battery Life

Reduced power consumption by up to 30% while streaming.



### Tinnitus Balance Noise Generator

Flexible tinnitus solution in all performance levels.



### Speech in Loud Noise

1.5 dB SNR improvement over static StereoZoom.



### CROS II

The smart solution for single sided hearing loss.

**PHONAK**  
life is on

## 4 performance levels

### 90 – Premium

- Speech in Car\*\*
- Speech in Loud Noise\*\*
- Comfort in Echo\*\*
- Music\*\*
- Comfort in Noise\*\*
- Speech in Noise\*\*
- Calm Situation\*\*
- Speech in Wind\*
- Speech in 360°\*
- StereoZoom\*
- EchoBlock
- WindBlock
- SoundRelax
- SNR-Boost
- FlexControl
- Flex Volume
- DuoPhone\*
- WhistleBlock
- NoiseBlock
- QuickSync
- SoundRecover
- Tinnitus Balance
- auto Acclimatization
- UltraZoom Premium
- User Preference Tuning
- 20 channels
- WCP compatible
- CROS II
- AOV

### 70 – Advanced

- Music\*\*
- Comfort in Noise\*\*
- Speech in Noise\*\*
- Calm Situation\*\*
- Speech in 360°\*
- StereoZoom\*
- WindBlock
- SoundRelax
- SNR-Boost
- FlexControl
- Flex Volume
- DuoPhone\*
- WhistleBlock
- NoiseBlock
- QuickSync
- SoundRecover
- Tinnitus Balance
- auto Acclimatization
- UltraZoom Advanced
- User Preference Tuning
- 16 channels
- WCP compatible
- CROS II
- AOV

### 50 – Standard

- Comfort in Noise\*\*
- Speech in Noise\*\*
- Calm Situation\*\*
- SNR-Boost
- FlexControl
- Flex Volume
- DuoPhone\*
- WhistleBlock
- NoiseBlock
- QuickSync
- SoundRecover
- Tinnitus Balance
- auto Acclimatization
- UltraZoom Standard
- User Preference Tuning
- 12 channels
- WCP compatible
- CROS II
- AOV

### 30 – Essential

- Speech in Noise\*\*
- Calm Situation\*\*
- WhistleBlock
- NoiseBlock
- QuickSync
- SoundRecover
- Tinnitus Balance
- auto Acclimatization
- UltraZoom Essential
- User Preference Tuning
- 8 channels
- WCP compatible
- CROS II
- AOV



\*\*

AutoSense OS



\*

Binaural  
VoiceStream  
Technology™

**PHONAK**  
life is on

16



## Naída V Benefits

### Enhanced Hearing Performance

- Improvement in high frequency audibility (with SoundRecover2)
- Increased audibility:
  - up to 5 dB more output with BroadbandBooster - up to 6 dB more output with xUP receiver

### Increased confidence and reliability

- 60% more robust with high tech composite material
- IP68 – water and dust resistant

### Smaller devices for more people

- More power allows for broader fitting ranges in smaller devices (xUP receiver and BroadbandBooster)
- 25% slimmer hearing aids

**PHONAK**  
life is on

## Goals of Non-Linear Frequency Compression

- Increase audibility for high frequency phonemes and sounds, specifically 's' and 'sh'
- Improve detection, recognition and perception of high frequencies
- Overcome high frequency receiver output limitations



**PHONAK**  
life is on

## Results of Non-Linear Frequency Compression

“Clinical field studies have shown that the use of SoundRecover restores high frequency audibility and results in high spontaneous user acceptance, quick acclimatization, and improvements in the wearer’s own voice quality.”<sup>1,2,3</sup>

“SoundRecover increased the pleasantness of sound and of the subject’s own voice, resulting in a highly satisfactory overall impression.”<sup>4</sup>

1 Simpson A, Hersbach AA, McDermott HJ: Improvements in speech perception with an experimental nonlinear frequency-compression hearing device. IJA 2005;44(5):281–292.

2 Scollie S, Glista D, Seewald R: Speech quality ratings of nonlinear frequency compressed speech by normal and hearing-impaired listeners. Ear Hear (submitted).

3 Bagatto M, Scollie S, Glista D, Pasa V, et al.: Case study outcomes of hearing-impaired listeners using nonlinear frequency compression technology. Audiol Online, 2008.

4 Nyffeler, M. (2008). Study finds that non-linear frequency compression boosts speech intelligibility. The hearing journal, 61(12), 22-24.



**PHONAK**  
life is on

## SoundRecover2

Restores access to high frequency information while preserving high frequency discrimination and sound quality

Retains the essence of our SoundRecover strategy, while allowing lower cut-off frequencies with weaker compression ratios

Maintains the familiarity of low and mid frequency sounds



Extends the benefits of SoundRecover to:

- more profound losses
- left corner audiograms
- ski-slope audiograms

**PHONAK**  
life is on

## A specific look at those with more profound hearing loss


- Audiogram with no response about 2000 Hz
- Even with SoundRecover (original) enabled, high frequency audibility is not possible due to the lowest possible cut off frequency of 1500 Hz



**PHONAK**  
life is on


## SoundRecover rules – severe to profound hearing loss

Increase the compression ratio to try to make more sounds audible in a smaller bandwidth?



Not the solution: could create a negative impact on sound quality


Lower the cut-off (CT) required to improve audibility?



Lowering the CT below 1.5 kHz could result in distortion of vowels

**PHONAK**  
life is on


## The solution is SoundRecover2 - adaptive frequency lowering algorithm




**PHONAK**  
life is on

## Adaptive nature of SoundRecover2

Input with significant **high frequency** energy present




Input with significant **low/mid frequency** energy present



25


**PHONAK**  
life is on

## Original sound spectrogram



**PHONAK**  
life is on

## SoundRecover vs SoundRecover2



**PHONAK**  
life is on

## Adjusting SoundRecover2 in Target 4.3




**PHONAK**  
life is on

## SR2 Controls

- CT1 @ 1.3 kHz
- CT2 @ 2.7 kHz
- Maximum output frequency (fmax) @ 4 kHz

- ← Adjustment toward “Audibility”
- Lowering of CT1, CT2 and fmax
  - “Stronger” frequency lowering
- Adjustment toward “Distinction”
- Raising of CT1, CT2 and fmax
  - “Weaker” frequency lowering

Adaptive frequency compression




29

**PHONAK**  
life is on

## SR2 Clarity/Comfort

- CT1 @ 1.5 kHz
- CT2 @ 3.2 kHz
- fmax @ 4 kHz

Adaptive frequency compression




30

**PHONAK**  
life is on

## SR2 Clarity/Comfort

- CT1 @ 1.6 kHz
- CT2 @ 3.5 kHz
- fmax @ 4 kHz

Adaptive frequency compression



31


**PHONAK**  
life is on

## SR2 Clarity/Comfort

- CT1 @ 1.8 kHz
- CT2 @ 4 kHz
- fmax @ 4 kHz

Total area of frequency compression is in adaptive region between CT1 and CT2

Adaptive frequency compression



32

**PHONAK**  
life is on






Figure 4: Phoneme perception test: median detection thresholds for profound hearing losses, Naída Q SR (original SoundRecover)<sup>2</sup> vs Naída V SR2 (SoundRecover2). Detection of 3 out of 4 stimuli tested (sh5, sh6, sh9) was significantly better with SoundRecover2.

- » 8 male hearing impaired adults
- » Profound, symmetrical, sensorineural or mixed loss
  - » Average loss >90dB 250Hz-8kHz
- » Phoneme Perception Test
  - » Significantly better detection for 3 out of 4 tested stimuli

**PHONAK**  
life is on

Baumann, Silvia; Jha, Siddhartha; Rehmann, Julia. Phonak Insight 2016. "SoundRecover2 – the first adaptive frequency compression algorithm."

## Benefits of SoundRecover2 pre-calculation

**Maximize stimulation of the hearing nerve** by setting the **Maximum Output Frequency** at the upper limit of the individual's audible bandwidth.

No risk of deprivation

**Protect the low and mid frequencies** by setting the upper cut-off (CT2) high enough so that audible speech (without compression) is not affected.

No risk for distortion of low and mid frequencies

**Extend the benefits to severe-profound losses** by combining a lower starting point for compression (CT1 or lower cut-off) with a weaker compression ratio.

Overcomes the constraints of SoundRecover and provides benefit for more users

**PHONAK**  
life is on

## BroadbandBooster: What is it?

When activated, BroadbandBooster increases the overall level of the broadband output of the hearing instrument by up to 6 dB.

- The gain-frequency response and MPO are not changed.
- When activated, the BroadbandBooster setting may vary between 0 to 10, where 0 equals no boost and 10 equals the largest boost in output.


When do you use BroadbandBooster?

1. You need the output of an xP receiver but only an xS will fit in a narrow ear canal.
2. You have already maxed-out a Naída V-UP BTE or a RIC with xUP receiver and still need more broadband output.

Naída V	DEFAULT setting
RIC xS	Off
RIC xP	Off
RIC xUP	On
SP BTE	On
UP BTE	On

**PHONAK**  
life is on


## BroadbandBooster in the test-box



Speech (male voice)


boost on (green)

boost off (purple)



**PHONAK**  
life is on

## BroadbandBooster in Target 4.3




**PHONAK**  
life is on

## xUP UltraPower receiver – backvent

Venting from the back cavity into the shell

Reduces the stiffness of the membrane

Increases the amplitude of the membrane vibration without feedback




\*Requires a custom cShell

- Up to 6 dB more output in the low frequencies
- Significantly less battery consumption by the receiver compared to non-backvented receivers

**PHONAK**  
life is on

## xUP UltraPower receiver: The most powerful RIC receiver from Phonak



**PHONAK**  
life is on

## Naída V Benefits

### Enhanced Hearing Performance

- Improvement in high frequency audibility (with SoundRecover2)
- Increased audibility:
  - up to 5 dB more output with BroadbandBooster
  - up to 6 dB more output with xUP receiver

### Increased confidence and reliability

- 60% more robust with high tech composite material
- IP68 – water and dust resistant

### Smaller devices for more people

- More power allows for broader fitting ranges in smaller devices (xUP receiver and BroadbandBooster)
- 25% slimmer hearing aids

**PHONAK**  
life is on

Naída V: 60% more robust and IP68



### **More Robust**


- Naída V features a new composite housing material.
- Pressure test validation found that 60% more pressure was needed to damage the housing of the instrument compared to its predecessor.

### **IP68**


- Naída V achieves an ingress protection (IP) rating of 68, the highest possible rating for hearing instruments.
- IP68 rating is maintained when Roger 18 and Roger 19 receivers are attached.

**PHONAK**  
life is on

Naída V tested in the water




Naída V survives continuous immersion in 3 feet of water for 60 minutes



**PHONAK**  
life is on

## Naída V tested in the dust chamber



After 8 hours in a dust chamber  
no traces of dust were evident  
within the housing



**PHONAK**  
life is on

## Naída V Benefits

### Enhanced Hearing Performance

- Improvement in high frequency audibility (with SoundRecover2)
- Increased audibility:
  - up to 5dB more output with BroadbandBooster
  - up to 6dB more output with xUP receiver

### Increased confidence and reliability


- 60% more robust with high tech composite material
- IP68 – water and dust resistant

### Smaller devices for more people

- More power allows for broader fitting ranges in smaller devices (xUP receiver and BroadbandBooster)
- 25% slimmer hearing aids

**PHONAK**  
life is on

Severe to profound clients want small hearing instruments too!




- Naída V-SP is 25% thinner than Naída Q-SP and 19% smaller in volume.
- Naída V-UP is also smaller!

**PHONAK**  
life is on




## A history of innovation in wireless for adults



**PHONAK**  
life is on

## Why Roger?



**PHONAK**  
life is on

## What is a Roger system?

- A Roger system consists of **one or more Roger microphones** with accompanying **Roger receivers**
- The Roger system functions by **picking up the voice of a speaker and wirelessly transmitting it to the listener** while reducing background noise



**PHONAK**  
life is on

## Roger Microphones



Roger Pen	Roger EasyPen	Roger Clip-On Mic
A discreet, wireless Roger microphone that delivers superior speech-in-noise and over distance performance. It's the ideal solution for an all-inclusive listening experience including Bluetooth™ connectivity.	An easy-to-use wireless microphone for one-on-one or group conversations, featuring full Roger speech-in-noise and over distance performance.	A discreet and lightweight stand-alone microphone for a conversation partner, with industry leading speech-in-noise and over distance performance.

**PHONAK**  
life is on

## Roger Table Mic

- Roger Table Mic is a **dedicated** microphone for meetings
- Designed purely for use on a **table**
- Extending the existing Roger for adult portfolio including Roger Pen, Roger EasyPen and Roger Clip-On Mic
- **Simple** - only two buttons
- Combining multiple Roger Table Mics will make it possible to cover large meeting rooms with many participants



**PHONAK**  
life is on

## Roger Receivers



Roger Design Integrated




Roger X



Roger MyLink

**PHONAK**  
life is on

## Roger relative to traditional FM



## Roger benefits



### Performance

- Sound Quality
- Speech understanding in noise and over distance
- Adaptive behavior
- **Roger Directional** NEW!

### Ease of Use

- Simple - No programming required
- Automatic mic mode
- **RogerReady** NEW!
- No frequency planning - 2.4 GHz

### Multi-functionality

- One device for everything
- TV/multimedia connectivity
- Wideband Bluetooth (Roger Pen)
- Expandable multi-talker network

### Design

- Discreet
- Cool design
- Appealing colors

**PHONAK**  
life is on

# RogerReady

## RogerReady



The Roger/DAI + Mic program is always in the program structure


No additional programming required

Simple 2 step process: Attach the Roger receiver(s) and switch on the Roger microphone

Ensures every power user's Naída V hearing aids are ready for Roger

**PHONAK**  
life is on

## RogerReady - Automatic activation of Roger



Field Study News, RogerReady: Automating the Roger Experience, April 2016

**PHONAK**  
life is on

## RogerReady

- Roger/DAI+Mic program is on by default and cannot be deleted from the program structure
- RogerReady is checked as default in the (Program manager)
  - Replaces EasyRoger
  - Can be unchecked – warning message if no manual program
- Hearing instrument microphone is always enabled by default
  - Can be disabled or attenuated



**PHONAK**  
life is on

## The best of both worlds: Roger + Directional



NEW!


- Now Naída V has directional microphones active when using a Roger microphone
- Patients wearing Naída V with Roger + Directional will be able to enjoy the benefits of the hearing instrument directionality in the near field and the Roger microphone in the far field
- Roger + Directional setting offers 3 hearing instrument microphone modes:
  - **Omni:** always omnidirectional microphone
  - **Real Ear Sound:** always Real Ear Sound
  - **Fixed directional:** adaptive activation of fixed beamformer

---

**PHONAK**  
life is on

## Program options in Target 4.3

- Roger/DAL+Mic program has three microphone mode settings
  - Omni-directional
  - Real Ear Sound
  - Fixed directional



**PHONAK**  
life is on

**The Phonak power hearing solution covers all listening situations and communication needs in the near and far field**

**PHONAK**  
life is on



## A solutions approach

- A complete solution that is aimed to meet all the listening and communication needs.
- For power patients, a single product solution is unlikely to meet all of their needs.
- Support encompasses *more* than just hearing aids and involves *more* than just the person with hearing loss.

**Timing is everything: To drive adoption, the complete solution must be recommended at the same time.**



**PHONAK**  
life is on

## Subjective Needs Assessment

**Hearing Solutions**

Name: \_\_\_\_\_ Date: \_\_\_\_\_

Needs Assessment

Which hearing problem do you wish to improve?

Hearing problem	How often	How much
Understanding speech in noisy situations	Always	Very much
Understanding speech in quiet situations	Always	Very much
Understanding speech in noisy situations	Often	Very much
Understanding speech in quiet situations	Often	Very much
Understanding speech in noisy situations	Sometimes	Very much
Understanding speech in quiet situations	Sometimes	Very much
Understanding speech in noisy situations	Rarely	Very much
Understanding speech in quiet situations	Rarely	Very much
Understanding speech in noisy situations	Never	Very much
Understanding speech in quiet situations	Never	Very much

Primary Hearing Goals

PHONAK

## Objective evaluation




## Recommendation



**PHONAK**  
life is on

## Indication for Roger based on need

- Patients who express improvements in any of these situations:
  - Hearing in noise
  - Hearing in quiet or noise at a distance
  - Listening to a passenger in a car, train, plane or bus
  - TV / Multimedia (computer, music player, etc.)
  - Cell phone
  - Landline phone
  - Business meetings
  - Conferences/lectures



Naida V + Roger is the best solution to address these situational needs.  
In quiet and in noise, near and far.

**PHONAK**  
life is on



Thank you

**PHONAK**  
life is on