



**More Than Just Sound – It's About Audibility:
The Oticon Medical Power
Solution**

Rebecca Cihocki AuD
May 10, 2016

Because sound matters

oticon
MEDICAL

Learning Objectives

- Explain to patients the benefits of Power bone anchored devices in everyday life.
- Describe the new surgical technique MIPS (Minimally Invasive Ponto Surgery).
- Explain how to connect wireless accessories to both Oticon Medical Ponto as well as Oticon Inc hearing aids.


Because sound matters


oticon
MEDICAL


**Introduction to Bone
Anchored Hearing Solutions**

Because sound matters

oticon
MEDICAL









Audiometry, Measuring AC and BC

Measurement of Bone-Conduction thresholds represent loss due to inner ear



Because sound matters


oticon
MEDICAL

Traditional Earphones/Headphones



Because sound matters


oticon
MEDICAL



The Batband uses three transducers to relay sound waves to the wearer's inner ear via the bones of the skull
(Image Gallery (4 Images))

Because sound matters

oticon
MEDICAL



Google Glass to use bone vibration instead of traditional headphones

TECHNOLOGY / 04 FEBRUARY 13 / by PHILIPPA WARR

Here are a few tips to help improve Glass' sound quality:


- Press the BCT against your head so that it is flush against your skull. This improves the conduction of sound through your bones and to your inner ear.
- Cup your right ear to reduce ambient sounds that might interfere with Glass' sound quality.
- Sometimes you may feel a tickling vibration from the BCT. That's normal, but if it's bugging you, turn down the volume.

Because sound matters

oticon MEDICAL

July 12, 2013 12:55am EDT

Bone conduction: the new front in guerilla advertising



Bone conduction technology enables commuters to hear advertisements via their skull – but how? *M. F. Fox's Travel Photos*

The audio information of these ads is transmitted via the window of the train. When you rest your head against the glass the signal is heard inside your head, utilising a phenomenon known as bone conduction. See some reactions to the technology in



Because sound matters



oticon MEDICAL



So What's Wrong With Traditional Hearing Aids?

Nothing - but for some people they just don't work.

1. Conductive/mix loss due to chronic otitis
2. Conductive/mix loss due to atresia
3. Unilateral (single sided) loss

Because sound matters

oticon MEDICAL

What Exactly Are We Making?




The Ponto is an auditory osseointegrated bone conduction implant system

Because sound matters

oticon MEDICAL

Osseointegration




Living Bone Tissue (Osseo)Integrating with Titanium

Because sound matters

oticon MEDICAL

Osseointegration and Direct Bone Conduction Simulation




Sound processor coupling ISQ screw Abutment Ti-Implant

Because sound matters

oticon MEDICAL

Ponto Plus Sound Processors



Ponto Plus
Up to 45 dB BHL

Ponto Plus Power
Up to 55 dB BHL


Because sound matters

oticon
MEDICAL

Normal Hearing

Hearing Threshold Level (HL)

Uncomfortable Level (UCL)




Auditory Dynamic Range (ADR)

Because sound matters

oticon
MEDICAL


Type of Hearing Losses

Sensorineural
BC ~ AC




ADR is reduced

Conductive
BC ~ 0 dB HL



ADR is not reduced

Mixed
Sensori-neural and Conductive




ADR is somewhat reduced

Because sound matters

oticon
MEDICAL

What Amount of Amplification is Needed?


1. Sounds should be amplified so they are within the Dynamic Range.



Because sound matters

oticon
MEDICAL

Hearing Loss and Traditional Hearing Aids



Because sound matters

oticon
MEDICAL

What Amount of Amplification is Needed?

- People with sensorineural, conductive or mixed HL have different ADR; therefore they need different amount of amplification.
- Fitting rules (=rationale) tell how much amplification is needed to compensate for the HL. The needed gain is often expressed in Insertion Gain (dB IG).

Because sound matters

oticon
MEDICAL

Need For Gain in Traditional HI

- is determined by AC and BC thresholds

Sensori-neural

AC = 40 @2k Hz

1/3 - 1/2 of the hearing loss,
~ 16 dB IG

Conductive

AC = 40 @2k Hz

1:1 of the hearing loss,
40 dB IG

Mixed

Combination of what's is
needed for sensori-neural
and conductive part of the
HL

Because sound matters

oticon
MEDICAL

Need For Gain in Bone Anchored Solutions

- is determined by BC thresholds

Sensori-neural

AC=BC = 40 @2k Hz

~ 16 dB IG

Conductive

BC = 0 @2k Hz

0 dB IG

Mixed

Determined by BC

Because sound matters


oticon
MEDICAL


Bone Anchored Hearing Aids


How does it work?

Because sound matters

oticon
MEDICAL








Power Across the Frequency Spectrum

- BAHs are inherently weak because the vibrator has to vibrate the whole skull
- More power means more audibility of everyday sounds
- Lower power means:
 - Reduced audibility
 - Poorer reproduction of own voice because of saturation issues
 - Tinny sound quality
 - Poorer result for patients with mixed hearing losses – the 45 dB HL fitting range may certainly be a stretch
 - Poorer general acceptance by patients

Because
sound matters

oticon
MEDICAL

Ponto Plus Power and Single-Sided Deafness



- Single-sided deaf patients with slight cochlea loss will get more HF amplification when fitted with Ponto Plus Power

Because
sound matters

oticon
MEDICAL

Known Consequences of Low MFO


- BAHs are inherently weak because the vibrator has to vibrate the whole skull – focus ought to be on making them stronger – not weaker
- What are the potential issues?
 - Poor reproduction of own voice because of saturation issues
 - Reduction of patient's ability to hear speech at a distance
 - Tinny sound quality
 - Poor general acceptance by patients – especially upgrades
 - Poor result with patients with mixed hearing losses – the 45 dB HL fitting range may certainly be a stretch

Because
sound matters

oticon
MEDICAL

Powerful Transducer = More of Everyday Sounds

Sound levels in different everyday situations*



*Based on recordings of actual everyday situations and a typical Ponto user. The exact effect for a user will depend on the situation, hearing loss and choice of sound processor (see Wegmann et al., 2006).

Because sound matters

oticon
MEDICAL

Known Consequences of Reduced Bandwidth

- Good high frequency performance is very important as patients with BAHs normally have good cochlear function and therefore can benefit from hearing the high frequencies
- What are the potential issues?
 - Reduced general sound quality
 - For patients with good cochlear function a reduction in speech understanding in certain background noise situations
 - A harder and less balanced sound experience as the reduced bandwidth might force audiologists to provide more gain in the mid frequencies


Because sound matters

oticon
MEDICAL

Ponto Plus – Provides More Output

BAHs are inherently weak because the transducer has to vibrate the whole skull.
Make sure never to compromise on sound!

- Increased low frequencies helps to provide better loudness sensation and voice quality.
- More high frequency gain delivers important speech cues.
- Wider bandwidth helps with speech perception, sound quality and understanding speech in noise.




Because sound matters

oticon
MEDICAL

Considerations For Mixed HL



Mixed Hearing losses




- More gain and headroom with Ponto Plus Power – benefits a wide range of mixed HL



Mixed Hearing Losses

- The need for amplification is high because:
 - Compensation for the sensori-neural hearing loss
 - Additional amplification is needed in a softband / skin drive solution to compensate for the skin attenuation
- Only a **Power** bone anchored hearing device can provide the needed gain




Mixed Hearing Losses

- The power sound processor needs to have an effective feedback management system, to ensure that:

- The prescribed gain will not be limited by the feedback limit
- The patients are not bothered by feedback

- A power sound processor is also needed in order to provide the patient with the possibility to turn up the volume.




Because
sound matters

oticon
MEDICAL

Bone Conduction Direct Drive vs Skin Drive systems

- There are two types of bone conduction devices:
- Direct Drive** devices send vibrations via direct route to bone
- Skin Drive** devices send vibrations through the skin to bone
- Softband and magnet solutions provide similar performance¹



Because
sound matters

¹ Briggs (2015) Clinical Performance of a New Magnetic Bone Conduction Hearing Implant System: Results from a Prospective, Multicenter, Clinical Investigation. *Otol Neurotol*, 2015 Jan;36(1):88-95.

oticon
MEDICAL

Avoiding Poor Results For SSD Patients

- Only direct drive solutions can amplify the needed frequencies.
- Fitting a skin drive solution to an SSD patient may lead to poor results.
- Thanks to the Feedback Shield, Ponto Plus provides more high frequency amplification than other bone conduction devices, and without feedback.




Because
sound matters

oticon
MEDICAL

On The Importance of Bandwidth For Children

- Stelmakowich et al. (2001)
- Recognition of high-frequency consonants
- Difference between speakers



Because
sound matters

oticon
MEDICAL

Conclusion

Medwetsky, 1992; Boothroyd, Erickson, and Medwetsky, 1994; Jongman, Wayland, and Wong, 2000), these results suggest that, when fitting hearing aids to young children, audible energy through 5 kHz is needed to maximize perception of /s/ produced by a male talker and that a 9-kHz upper limit is needed for /s/ produced by a female. From a clinical perspective, it is important to distinguish between cases where performance *cannot be improved* by high-frequency


Because
sound matters

oticon
MEDICAL

The Importance of High Frequencies For Language Development in Children

Clinical study

- Purpose: Determine learning rate for words presented with 4 kHz and 9 kHz bandwidths
- Method: Five nonsense words, paired with 5 novel pictures
- Outcome: Children need **3 times as many trials with limited bandwidth** to learn new words as compared to the children who listen to extended bandwidth.⁷



⁷ Pittman (2008) Short-term word learning rate in children with normal hearing and children with hearing loss in limited and extended high-frequency bandwidths. *Journal of Speech, Language and Hearing Research*, 51(5), 1265-1271

Because
sound matters

oticon
MEDICAL

The Importance of /S/

- 3rd or 4th most frequently occurring phoneme in the English language
- Makes a word plural
- Shows tense
- Shows possession
- Is that Mary? vs Is that Mary's?
- She put it on vs. she puts it on
- Keep Vs. Keeps
- Bring the Chair Vs. Bring the chairs

Because
sound matters

oticon
MEDICAL


Daily Usage And Comfort

Ponto

- Minimal physical sensation of the device
- Device stays attached to abutment
- Nothing prevents the user from wearing it during all waking hours (> 8 hours)

Magnetic solution

- Low avg. usage: 7.0 h/day
- 22 % of patients had numbness 9 months post-surgery.




² Dettl et al. Day to day use and service related issues with bone anchored hearing aids: The tiffic medical system questionnaire. *Journal of Laryngology and Otolaryngology* 2002; 122: 20-28
⁴ Ridge et al (2015) Clinical performance of a new magnetic bone conduction hearing implant system: Results from a prospective, multicenter, clinical investigation. *Otolaryngology & Neurotology*. Jan;36(5):838-42

Because
sound matters

oticon
MEDICAL


Surgical Considerations Pertinent to Patients with SSD Due to Acoustic Neuroma

- MRI
- A big artifact causes shadowing. Nothing in the shadow is visible
- MRI is widely used. There is a high likelihood patients will undergo MRI examination during their lifetime



<http://www.magnetic-resonance.org/021-01.html>

- Acoustic neuroma patients undergo a number of MRI's post surgery in follow up program



¹⁾ Oticon Medical MRI Safety card
²⁾ Imaging guide, Oticon Medical's Magnetic System, Surgical procedure


Because
sound matters

oticon
MEDICAL

Many BAHS Options Available

Some important considerations:

- Does the solution provide enough amplification / output?
- Is it comfortable enough to wear for a full day?
- Cosmetically, how does it compare?



Because sound matters

oticon MEDICAL

Wireless




Because sound matters

oticon MEDICAL

Wireless Connectivity

•Ponto Steamer



Because sound matters

oticon MEDICAL

Bimodal Use Of Oticon Medical Streamer



Because sound matters


oticon
MEDICAL

Which Oticon Hearing Aids Can Be Used in a Bimodal Setup?

• Currently, there are six models that can be used bimodally with OM Streamer*:


• **Dynamo:**

- Dynamo SP4
- Dynamo SP6
- Dynamo SP8
- Dynamo SP10



• **Sensei SP:**

- Sensei SP
- Sensei SP Pro



* Find the full list of Oticon Medical Streamer-compatible hearing aids at www.oticonmedical.com/OMStreamerHA

Because sound matters

oticon
MEDICAL

How Do We Recommend That Oticon Medical Streamer is Used in Bimodal Fitting Cases?

• Leave Oticon Medical Streamer **OPEN**.


• The wireless Ponto sound processor and the Oticon hearing aid can both be used with the Oticon Medical Streamer without linking, right out of the box, in the **OPEN** condition



Because sound matters

oticon
MEDICAL

Wireless Possibilities Should Never Compromise on the Power Which Users Need!




Because sound matters

*The MaxGain (max audiology gain) setting was used for all devices during testing - with and without iPhone audio streaming link active - The MFO MaxGain curves shows results without actively streaming a sound file through the wireless link.

oticon MEDICAL

Wireless Possibilities Should Never Compromise on the Power Which Users Need!




Because sound matters

*The MaxGain (max audiology gain) setting was used for all devices during testing - with and without iPhone audio streaming link active - The MFO MaxGain curves shows results without actively streaming a sound file through the wireless link. The Streaming Silence curves show the results when streaming an empty (zero) sound file from the iPhone during the MFO measurement. This shows how the wireless link itself influences the MFO.

oticon MEDICAL


Ponto Plus - Wireless Possibilities Does Not Affect the Power Users Need!




Because sound matters

*The MaxGain (max audiology gain) setting was used for all devices during testing - with and without iPhone audio streaming link active - The MFO MaxGain curves shows results without actively streaming a sound file through the wireless link.

oticon MEDICAL








MIPS Procedure, Summary

1. Decide position and measure skin thickness
2. Punch
3. Insert Cannula
4. Drill. Guide- and widening drilling
5. Insert implant
6. Fit Soft Healing cap



Because sound matters

oticon MEDICAL

A Truly New Perspective on Tissue Preservation




Because sound matters

oticon MEDICAL

MIPS Procedure, Summary


1. Decide position and measure skin thickness
2. Punch
3. Insert Cannula
4. Drill. Guide- and widening drilling
5. Insert implant
6. Fit Soft Healing cap



Because sound matters

oticon MEDICAL

Minimizing Post-op Complications



Tissue preservation surgery
3 months post-op


MIPS surgery 1 week post-op

*Johansson M, Holmberg M, Hakkarinen M, "Bone anchored hearing
implant surgery with tissue preservation – A systematic literature review,"
Otitis Media white paper, MS2017, 2014

Because
sound matters

oticon
MEDICAL

Patient Outcomes 7 Days Post-op



Because
sound matters

oticon
MEDICAL

Abutment Extension



Because
sound matters

oticon
MEDICAL

Small Change. Big Difference

- No surgery
 - The Abutment extension is connected on the existing abutment without any surgery involved.
 - It is delivered non-sterile.
 - Placement can take place outside the surgical theatre.
 - The Abutment extension does not work with BAHA processors. It must be removed if a BAHA is fitted.



Because
sound matters

oticon
MEDICAL

Small Change. Big Difference

- Seamless fit on existing abutment
 - The Abutment extension has full compatibility to most bone anchored hearing aid abutments available on the market*.




Because
sound matters

*Please read the Abutment extension IFU for further instructions

oticon
MEDICAL

Is the Sound Transmission Affected by the Abutment Extension?

- No
 - The graph shows the output as measured with accelerometer on counterweight mounted on patient's abutment; with and without abutment extension.
 - The curves show that they are on par with no major difference.




Because
sound matters

*When normal Gaze Medical fitting flow is followed

oticon
MEDICAL

Abutment Extension




Because sound matters

oticon
MEDICAL

Future Technologies

- More powerful sound processor*



- Requires FDA clearance to commercially market in the United States
- Not currently available.

Because sound matters

oticon
MEDICAL
