

SRII


Taking the high road to high frequency access



Christine Jones, AuD

PARC
Phonak Audiology
Research Center

PHONAK
life is on



PHONAK
life is on

Slide 3

A1 Add screen shot and referencec
Author, 3/2/2016

Searching for the optimal solution...



5 Phonak Audiology
Research Center

AudOnline Aug 2016

PHONAK
life is on

SoundRecover evidence

Peer reviewed

Glista D, Scollie S, Sulkers J Perceptual Acclimatization Post Nonlinear Frequency Compression Hearing Aid Fitting in Older Children. J Speech Lang Hear Res. 2012 May 21. [Epub ahead of print]

Wolfe J, John A, Schafer E, Nyffeler M, Boretzki M, Caraway T, Hudson M. Long-term effects of non-linear frequency compression for children with moderate hearing loss. Int J Audiol. 2011 Jun;50(6):396-404

Wolfe J, John A, Schafer E, Nyffeler M, Boretzki M, Caraway T. Evaluation of nonlinear frequency compression for school-age children with moderate to moderately severe hearing loss. J Am Acad Audiol. 2010 Nov-Dec;21(10):618-28

Bohnert A, Nyffeler M, Keilmann A. Advantages of a non-linear frequency compression algorithm in noise. Eur Arch Otorhinolaryngol. 2010 Jul;267(7):1045-53

Glista D, Scollie S, Bagatto M, Seewald R, Parsa V, Johnson A. Evaluation of nonlinear frequency compression: clinical outcomes. Int J Audiol. 2009;48(9):632-44

Simpson A. Frequency-lowering devices for managing high-frequency hearing loss: a review. Trends Amplif. 2009;13:87-106.

Stelmachowicz, P., Pittman, A., Hoover, B., Lewis, D. (2002). Aided perception of the /s/ and /z/ by hearing-impaired children. Ear and Hearing, 23 (4), 316-324

Non-peer reviewed

McDermott HJ. A technical comparison of digital frequency-lowering algorithms available in two current hearing aids. PLoS One. 2011;6(7):e22358


Timmer B. 2010: Neue Ansätze bei direktionalen multi-mikrofonsystemen. Hörakustik 11/2010:12-16

Nyffeler M. 2010: Geschaffen für erstklassigen Hörgenuss, tatsächlich binaural. AudiolInfos 10:78-82

Wolfe, J., Caraway, T., John, A., Schafer, E. C., & Nyffeler, M. (2009). Study suggests that non-linear frequency

PHONAK
life is on

Word learning improves as a function of bandwidth




PARC

7 Phonak Audiology Research Center


AudOnline Aug 2016

» Pittman, 2010

PHONAK
life is on



SR activated commonly for all ages



A3

Clinician recognized benefits

Better detection,
distinction
recognition of high
frequency
consonants

Improved access to
environmental
sounds

Improved vocal
quality

Reduced feedback

Slide 10

A3 Include picture icons here
Author, 3/2/2016

Why SoundRecover is disabled

Of those surveyed:

- ☐ 40% routinely disable SoundRecover for mild losses
- ☐ 85% report that age does not seem to be a factor

In open question, they report

- ☐ Disabled for rising, flat or mild hearing losses
- ☐ Disabled when Targets could be matched without it

Concerns mentioned:


- ☐ Always enabled in the fitting software
- ☐ Default settings seem too aggressive
- ☐ Patients report unnatural sound quality



Challenges from our own experience


We know for more profound hearing losses and for left corner and ski slope hearing losses

Trade-off between audiological benefit and sound quality




Trade-off – audiological benefit and sound quality

- Severe to profound hearing loss with no responses beyond 2 kHz
- SoundRecover enabled at maximum




- Minimal benefit – only loud inputs
- Stronger setting needed




But such settings introduce unacceptable sound quality

Phonak has led the way in modern day frequency lowering and we won't rest on our success.

We will continue to innovate and develop this technology so that more people can experience the benefit of high frequency audibility and improved sound quality.



SoundRecover2 - intelligent adaptive frequency lowering algorithm




Parameter explanation

- **CT1**: below this lower cut-off frequency the output is always linear, no frequency compression applied
- **CT2**: beyond this upper cut-off, compression is always active, as with the original algorithm
 - between the cut-off frequencies, compression is applied only when significant high frequency energy is present thereby protecting and maintaining the familiarity of mid and low frequencies sounds
- **fOutMax**: maximum output frequency after frequency lowering - depends on audiogram
- **Compression Ratio (CR)**: constant, fixed across the compressed range. SR2 allows an overall lower compression ratio.


The philosophy of SoundRecover2

- Maximize stimulation of the hearing nerve
- Protect the mid and low frequencies.
- Extend the benefits for more severe-profound hearing losses




Philosophy of the SoundRecover2

- **Maximize stimulation of the hearing nerve** by setting the initial **fOutMax** at the upper limit of audible bandwidth.



Philosophy of the SoundRecover2 pre-calculation

- **Protect the mid-frequencies** by setting **(CT2 or upper cut-off)** high enough so that audible speech (without compression) is not affected.



Benefits of SoundRecover2 pre-calculation

Maximize stimulation of the hearing nerve by setting the **fOutMax** at the upper limit of the individual's audible bandwidth.

No risk of deprivation


Protect the mid and low frequencies by setting **(CT2 or upper cut-off)** high enough so that audible speech (without compression) is not affected.

No risk for distortion of mid and low frequencies

Extend the benefits of SoundRecover2 for severe-profound losses by combining a lower starting point for compression (CT1 or lower cut-off) with a weaker compression ratio

Overcomes the constraints of SoundRecover and provides benefit for more users

Original sound spectrogram



My name is Asa


21

AudOnline Aug
2016

PHONAK
life is on


SoundRecover vs SoundRecover2

SoundRecover at strongest settings



My name is Asa

SoundRecover2




My name is Asa

PARC
Phonetic Auditory Research Center
22 PhonAudOnline/Aug
Rese2016center

PHONAK
life is on


Protecting the vowel formants at various settings




23 AudOnline Aug 2016

PHONAK
life is on

Phoneme perception test



PPT: Detection




25

AudOnline Aug
2016

PHONAK
life is on

PPT: Recognition



26

AudOnline Aug
2016

PHONAK
life is on

Scientific evidence – in progress

- Studies ongoing with
 - Susan Scollie and Danielle Glista at Western University (UWO)
 - Jace Wolfe at Hearts for Hearing, OK
 - Andrea Bohnert at Mainz University, Germany
- IJA in review: Jace Wolfe
- Publications:
 - Field Study News SoundRecover2
 - SoundRecover Insight
 - Best Practice Protocol Verification of SoundRecover2 for Pediatrics

SoundRecover2 and Phonak Target 4.3

SoundRecover2 will be available in all new pediatric and power form factors and technology levels.

- Naida V SP and RIC
- Preorders for UP

- Enabled by default for moderate to profound hearing losses.
- Can be enabled for mild losses

Parameter limits

	SoundRecover2
CT1	$\geq 800\text{Hz}$
CT2	$\geq 1.6\text{kHz}$
CR	1:1.05 to 1:4
fOutMax	$\geq 2.7\text{kHz}$

SoundRecover2 – overcomes the constraints

Restores access to high frequency information while preserving high frequency discrimination and sound quality

Maintains the familiarity of low and mid frequency sounds

Retains the essence of our SoundRecover strategy, while allowing lower cut-off frequencies with weaker compression ratios

Extends the benefits of SoundRecover to those with

- more profound hearing loss
- left corner audiogram
- ski slope audiogram

