


Technology in 10: Benefits of Binaural Beamforming for Individuals with Severe Hearing Loss



Hearing
Systems

SIEMENS

Benefits of Binaural Beamforming for Individuals with Severe Hearing Loss

Presented by:
Dianna Theokas, AuD

Transmitted to you in cooperation with our colleagues at AudiologyOnline.
If you are experiencing audio or visual difficulties,
please contact AudiologyOnline (800-753-2160)

Hearing
Systems

SIEMENS

<http://www.hearingreview.com/2016/04/benefits-binaural-beamforming-individuals-severe-hearing-loss/>



TECH TOPIC

Benefits of Binaural Beamforming for Individuals with Severe Hearing Loss

Published on April 29, 2016

Tech Topic | May 2016 *Hearing Review*

By Thomas A. Powers, PhD, and Veronika Littmann, PhD


Historically, one of the most effective methods to improve speech understanding in background noise has been the use of directional technology. In recent years, through the use of wireless full-audio transmission between a bilateral pair of hearing aids, binaural beamforming technology

New research shows it is reasonable to expect the same speech recognition benefit for individuals with severe hearing impairment when using binaural beamforming technology as previously reported for those with mild-to-moderate losses.

has taken the benefit of directional technology to a new level.^{1,2} Clinical trials with these new products have shown a substantial benefit over omnidirectional, and in fact two independent studies have shown that hearing-impaired individuals fitted with these instruments experience significantly better speech recognition in background noise than their normal-hearing counterparts.^{3,4}

- Through the use of wireless full-audio transmission between a bilateral pair of hearing aids, binaural beamforming technology has taken the benefit of directional technology to a new level
- Clinical trials with these new products have shown substantial benefit over omnidirectional
- Research has focused on hearing-impaired with mild-to-moderate hearing losses, since this is the profile of the most common hearing aid user
- Can we expect the same speech recognition benefit from those with severe hearing losses?

- The research with these beamforming instruments has focused on hearing-impaired individuals with mild-to-moderate hearing losses
 - The primary goal is to restore lost audibility due to outer hair cell damage
- However, a significant minority of people who use hearing aids have a severe impairment
 - When a severe hearing loss is present, we expect there will also be inner hair cell damage
 1. This is associated with broadened auditory filters causing loss of frequency selectivity
 2. Reducing the ability to separate two signals in the frequency domain
 3. Broad filters allow more noise to pass



TruEar (Mild Directional)


Simulates the acoustics of the pinna to allow for better front/back localization of BTE/RIC fittings

Standard Directionality (Adaptive directional)

Attenuates noise from the sides to behind the listener, but not in the listener's angle of focus

Narrow Directionality (Frontal beamforming)

Narrows directional microphone focus and results in better-than-normal hearing in demanding environments



Subjects


- Experienced hearing aid wearers
- Bilateral symmetrical severe SNHL
- Age 61-78 years (mean age = 70.3)

Instruments:

- Binaural fitting of Pure 7 bx with DoubleDomes

Programmed for each individual based on their hearing thresholds:

- First Fit using binaxFit
- Mild directional (TruEar)
- Adaptive directional
- Frontal beamforming (ND)



Test set-up:

- The array for the presentation of the target and competing speech material consisted of eight loudspeakers surrounding the participant
- Sentences of the English version of the Oldenburger Satztest (OLSA) presented from speaker at 0 degrees
- Competing OLSA sentences presented at 68 dBA through seven, equally spaced loudspeakers surrounding subject
- The competing signal remained constant for all testing
- Following a practice list, one list of 20 sentences of the OLSA was administered and scored
- The threshold in noise signal-to-noise ratio (SNR) was automatically calculated

Clinical Evaluation

Binaural beamforming with severe hearing loss

Results

- Adaptive directional ~ 3dB better than mild directional
- Binaural beamforming added another 2 dB adaptive directional benefit
- Same speech recognition benefit as reported for those with mild-to-moderate hearing loss as reported in separate clinical study (Froehlich, et al 2015)
Hearing Review 2015

- Purpose: Compare the benefit observed for this severe hearing loss group to previous research for individuals with mild-to-moderate hearing loss
- As reported in the earlier research (Froehlich, et al 2015) and using the same experimental design, the benefit of beamforming, compared to mild directional, was an identical 5.0 dB for both test sites (using the German OLSA and the English HINT material)
- In the present research (as shown in Figure 2), the mild directional SRT-in-noise was -3.0 dB, and mean performance for the beamforming technology was -7.9 dB minus a 4.9 dB advantage
- Results: When using the binax beamforming technology, it is reasonable to expect the same speech recognition benefit for individuals with severe hearing impairment as previously reported for those with mild-to-moderate losses

signia

Life sounds brilliant.



Hearing
Systems

SIEMENS

Upcoming Technology in 10 Webinars



July 27, 2016 at 1 PM EDT

- A Comparison of Prescriptive Fitting Algorithms for Tonal Languages

August 2, 2016 at 1 PM EDT

- Objective Measurement of Listening Effort

August 9, 2016 at 12 PM EDT

- eWindScreen binaural

August 16, 2016 at 12 PM EDT

- Signia wireless CROS and BiCROS

August 24, 2016 at 11 AM EDT

- Directional microphones

August 30, 2016 at 4 PM EDT

- Spatial SpeechFocus

September 6, 2016 at 4 PM EDT

- EchoShield

signia

Life sounds brilliant.

primax

Redefining the ease of listening



(800) 766-4500
www.usa.siemens.com/hearing

Hearing
Systems

SIEMENS