

signia
Life sounds brilliant.

Technology in 10: Directional Microphones



Hearing
Systems
SIEMENS

signia
Life sounds brilliant.

Directional Microphones

Presented by:
Jennifer Gehlen, AuD


Transmitted to you in cooperation with our colleagues at AudiologyOnline.
If you are experiencing audio or visual difficulties,
please contact AudiologyOnline (800-753-2160)

Hearing
Systems
SIEMENS

Directional microphone technology Background

signia
Life sounds brilliant.

- Available in hearing aids in the U.S. since the 1970s
- Improves signal-to-noise ratio by reducing sound outside of the angle of focus
- Utilizing hearing aid processing alone, directional microphone technology is the only proven way to increase signal-to-noise ratio
- Directional microphones have been shown to reduce listening effort in background noise in older adults with sensorineural hearing loss¹



¹ The Effects of Hearing Aid Directional Microphone and Noise Reduction Processing on Listening Effort in Older Adults with Hearing Loss, Desjardins, J Am Acad Audiol. 2016 Jan;27(1):29-41

Advanced directionality


signia
Life sounds brilliant.



Bilateral beamforming

e2e wireless 3.0 Binaural audio processing

signia
Life sounds brilliant.



Hearing
Systems
SIEMENS

Signiausa.com Whitepaper

signia
Life sounds brilliant.

www.siemens.com/hearing

e2e 3.0 – energy efficient inter-aural audio transmission

Thomas Lotter, Ph.D., Ulrich Schätzle, Dipl.-Ing., Thomas Fischer, Ph.D.


Lotter, T., Schätzle, U., Fischer, T. (2014). e2e 3.0 – energy efficient inter-aural audio transmission. Siemens Technology Backgrounder

Hearing
Systems
SIEMENS

Directionality algorithms Automatic in Universal program


signia
Life sounds brilliant.

Narrow directionality




e2e wireless 3.0

Binaural OneMic directionality




e2e wireless 3.0

Spatial SpeechFocus



e2e wireless 3.0


Directional speech enhancement



Narrow Directionality Automatic in Universal program

signia
Life sounds brilliant.

Narrow directionality




e2e wireless 3.0

Hearing
Systems
SIEMENS

Narrow Directionality Automatic in Universal program

signia
Life sounds brilliant.

NOTE: Our automatic ND algorithm is only engaged in extremely noisy environments and only in the channels where the SNR is negative




Narrow Directionality Wearer Benefits


signia
Life sounds brilliant.

1 dB SRT corresponds to 11% improvement in speech understanding for the HINT and 17% improvement for the OLSA

Univ. Northern Colorado



Hörzentrum Oldenburg



*Two clinical studies have shown that binax provides better than normal hearing in certain demanding environments (University of Northern Colorado, 2014; Oldenburg Hörzentrum, 2013). Speech Reception Thresholds (SRT) in cocktail-party situations improved up to 2.9 dB for wearers with mild to moderate hearing loss using Carat binax or Pure binax hearing aids with narrow directionality, compared to people with normal hearing.

*Studies conducted at University of Northern Colorado (2014) and Oldenburg Hörzentrum (2013) showed that Speech Reception Thresholds (SRT) in cocktail-party situations improved up to 2.9 dB for wearers with mild to moderate hearing loss using binax with Narrow Directionality, compared to people with normal hearing. This corresponds to over 25% improvement in speech understanding.

Published Research Article

signia
Life sounds brilliant.

May 2015 *Hearing Review* by Rebecca Herbig, AuD and Matthias Froehlich, PhD

Binaural Beamforming: The Natural Evolution

<http://www.hearingreview.com/2015/04/binaural-beamforming-natural-evolution/>

TECH TOPIC | READER CORNER

Binaural Beamforming: The Natural Evolution

BY REBECCA HERBIG, AuD, AND MATTHIAS FROELICH, PhD

This paper reports on a new measurement protocol based on the Hagmann-Ohlsson Method that yields a superior directivity index (DI), and demonstrates that the Signia binaural beamformer Narrow Directionality in a new 2-microphone device is as effective as the RIC form factor—the latter of which has been shown in clinical studies to result in better SII benefits for patients in certain noisy environments than those with normal hearing. Additionally, a new binaural beamformer designed for single-microphone devices, like CICs, is shown to have significantly higher SII/DI values than relying on omnidirectional processing and the natural pinna effect alone.

With the introduction of binaural use in ear-wearable audio transmission, such hearing devices in a bilateral pair can now operate with input not only from microphones on its own hearing, but it also receives acoustic signals directed and processed by the other hearing device. This has opened up new possibilities for advancement in many aspects of hearing treatment processing, most notably binaural beamforming technology.

For ear-wearable audio transmission, there is an algorithm that can provide a spatial focus on speech signals originating from the side or, in the back, to the more common frontal pattern. As illustrated in Figure 1, we have been able to create a binaural beamforming algorithm, which can narrow the directional focus of the beam even further so that essentially everything outside of what is immediately in

front of the wearer is attenuated (for a review, see Karkar, Park, Tietzer, and Aubertle). This technology, commonly marketed as Signia narrow directionality, has been shown in two clinical studies to allow hearing-impaired wearers to have a better speech reception threshold (SRT) in noise than those with normal hearing in certain noisy listening situations (see Figure 2 from Froehlich and Powers).

The benefit we have observed in clinical studies, however, is not reflected in conventional lab measurements used to quantify directivity, such as the directivity index (DI). This is because standard measurement techniques do not reflect the real-world benefits of modern dynamic beamforming technology. The primary goal of conventional measurement has been to optimize directivity to a controlled and repeatable way, accepting data from real-world application where noise

very following the rationale, directivity has often been measured in anechoic chambers, without sharing any formal target signal, or using pure-tone signals as the reference—all of which do not represent real-world applications.

The Hagmann-Ohlsson Method Using Speech Signals

While clinical study results are highly relevant, it is not efficient or feasible to use this efficacy measure to assess the effectiveness of every new directional microphone technology for every hearing aid model. Therefore, a lab measurement utilizing measurement conditions close to real-world environments is necessary. This is the only means by which we can quickly assess and quantify directional performance for new beamforming techniques as they are applied in various hearing aid form factors. This was the motivation behind the development of a new protocol for assessment of directivity in modern hearing aids, which has the benefit of simulating both a target source and an interfering source at the same time using speech signals. This protocol uses the established Hagmann-Ohlsson Method, but approximates real-world situations, such as noisy restaurants or cocktail parties where the wearer wants to listen to a specific speech source, but the unwanted interfering noise is other people talking in the background. This new measurement method results in a superior directivity index or SII (Auditory Noise Index) as the same "target" to reflect the Hagmann-Ohlsson approach, where signals only allow desired components as measured sequentially in order to obtain the final result.



Rebecca Herbig, AuD, is Manager & Editor of Hearing Services, and Matthias Froehlich, PhD, is Head of Product Management, Audiology for Research Center, Falmouth, Germany.




Figure 1. Compared with omnidirectional microphones and relying on the natural pinna effect, the Signia narrow directionality technology can focus on the target source more effectively.

24 | HEARINGREVIEW.COM | MAY 2015

Published Research Article

signia
Life sounds brilliant.

May 2016 *Hearing Review* by Thomas A. Powers, PhD and Veronika Littmann, PhD

Benefits of Binaural Beamforming for Individuals with Severe Hearing Loss

<http://www.hearingreview.com/2016/04/benefits-binaural-beamforming-individuals-severe-hearing-loss/>

hearingreview.com

<http://www.hearingreview.com/2016/04/benefits-binaural-beamforming-individuals-severe-hearing-loss/>

Benefits of Binaural Beamforming for Individuals with Severe Hearing Loss

By Thomas A. Powers, PhD, and Veronika Littmann, PhD

Tech Topic | May 2016 *Hearing Review*

New research shows it is reasonable to expect the same speech recognition benefit for individuals with severe hearing impairment when using binaural beamforming technology as previously reported for those with mild-to-moderate losses.

Historically, one of the most effective methods to improve speech understanding in background noise has been the use of directional technology. In recent years, through the use of wireless full-audio transmission between a bilateral pair of hearing aids, binaural beamforming technology has taken the benefit of directional technology to a new level.^{1,2} Clinical trials with these new products have shown a substantial benefit over omnidirectional, and in fact two independent studies have shown that hearing-impaired individuals fitted with these instruments experience significantly better speech recognition in background noise than their normal-hearing counterparts.^{3,4}

In general, the research with these beamforming instruments has focused on hearing-impaired individuals with mild-to-moderate hearing losses, as this is the audiometric profile of the most common hearing aid user. However, a significant minority of people who use hearing aids have a severe impairment. Can we expect the same speech recognition benefit for these hearing aid users?

When fitting hearing aids to individuals with mild-to-moderate hearing impairment, the primary goal is to restore the lost audibility due to their outer hair cell damage. However, when a severe hearing loss is present, we expect that there will also be inner hair cell damage. This is associated with broadened auditory filters causing loss of frequency selectivity, reducing the ability to separate two signals in the frequency domain. Broad filters allow more noise to pass, which then makes it even more difficult to detect the desired target speech in a competing background noise environment. Moreover, when signals are made audible, there may be a greater susceptibility to distortions.⁵

A related issue when fitting individuals with severe hearing losses is the consideration of "effective audibility." As reviewed by Ching et al.,⁶ when hearing loss becomes more severe, the ability of the ear to extract useful information from the signal is reduced, and there comes a point where more audibility will not provide more benefit, and in fact may even be detrimental, particularly if cochlear dead regions are present. Therefore, the amplified signal must be audible, but effective, which can be a small residual dynamic range window for many patients with severe hearing loss.

Simply said, there are many more factors present that can negatively impact the benefit obtained from hearing aids for the severe hearing loss patient than for those with mild-to-moderate losses. Limited research in this area suggests that people with severe losses will obtain at least some additional benefit when fitted with directional technology,⁷ although this has not been studied with modern beamforming processing. Therefore, the current study was designed to determine if the large benefit obtained in previous research for individuals with mild-to-moderate losses using binaural beamforming technology would also be present when the participants had severe hearing impairments.


Signia research results

Narrow Directionality—Severe Hearing Loss

signia
Life sounds brilliant.

Results

- Adaptive directional ~ 3dB better than mild directional
- Binaural beamforming ~ additional 2 dB adaptive directional benefit
- Same benefit as reported for those with mild-to-moderate hearing loss as reported in the separate clinical study (Herbig & Froehlich, et al 2015)



Displayed are the mean SRT-In-Noise findings and error bars (English OLSA) for three different microphone conditions: Mild Directional (omnidirectional with pinna compensation), Adaptive Directional, and Beamforming

Narrow Directionality

Wearer control via TouchControl App

signia
Life sounds brilliant.



If the wearer wants the greatest impact of bilateral beamforming in very demanding environments they can select the narrowest beam via the TouchControl App (Narrow Directionality is now fully engaged in all channels)


Available on the
App Store

ANDROID APP ON
Google play

Hearing
Systems
SIEMENS

Spatial SpeechFocus Automatic in the car

signia
Life sounds brilliant.



Hearing
Systems
SIEMENS

Published Research Article

signia
Life sounds brilliant.

September 2014 Directional Steering for Special Listening Situations: Benefit Supported by Research Evidence

The Siemens SpeechFocus algorithm has been available since 2010, and research with this product is starting to emerge in referred publications. In this paper, we review three of these articles.

Eric Branda, AuD, Joel Beilin, MScEE, Thomas A. Powers, PhD

<http://www.audiologyonline.com/articles/directional-steering-for-special-listening-12974/>

Directional Steering for Special Listening Situations: Benefit Supported by Research Evidence

Eric Branda, AuD, Joel Beilin, MScEE, Thomas A. Powers, PhD

September 23, 2014

Would you like to earn CEEs for this article?



Like?

Print

Since the introduction of the directional hearing aid over 40 years ago, this technology has proven to be one of the most successful solutions for improving speech understanding in background noise. Well-performing bilaterally-fitted directional instruments can improve the signal-to-noise ratio (SNR) by 6 dB or more (Powers and Beilin, 2013). This degree of improvement is substantial, and often determines whether a patient can follow a conversation. Moreover, the benefit is present for individuals with a wide range of hearing loss, different audiometric configurations and listening needs.


Historically, directional amplification has been designed so that maximum amplification is applied to signals from the front of the user, with maximum attenuation applied to signals originating from the back. This of course makes sense, as in most listening situations the wearer faces the speaker of interest. Even when a conversation originates from a location other than the front, most listeners will turn their head to face the speaker. In this position, the fundamental design of directional technology is working at or near its maximum efficiency.

So while we would expect an SNR improvement using traditional directional microphone technology for most listening situations, there are some cases when the wearer is unable to face the speaker. One such commonplace situation is when the hearing aid user is driving a car and a passenger in the back seat is speaking. Here, the noise comes from the sides and front while, speech comes from the back—yet the hearing aid user cannot turn to face the speaker. In fact, in this situation we would predict that traditional directional technology would actually be worse than omnidirectional. Other similar communication situations might occur in group events, large gatherings, when the wearer is in a wheelchair, or when having a conversation while walking.

Spatial SpeechFocus

Manual control via TouchControl App

signia
Life sounds brilliant.



Summary

signia
Life sounds brilliant.

Benefits of directionality and binaural beamforming technology:

- Automatic SNR benefit in difficult environments with Narrow Directionality, including listeners with severe hearing loss
- Automatic SNR benefit for CICs using directionality with binaural beamforming technology
- Automatic directional benefit toward targeted speech source in the car
- Manual access to directional preference via the discrete app available for Apple® and Android® phones



Life sounds brilliant.



Hearing
Systems

SIEMENS

2016 Expert Series



Life sounds brilliant.

2016 Expert Series.

Join us all summer long in a series of free online seminars. All sessions are recorded and will be available at [AudiologyOnline.com](http://www.audiologyonline.com/c/signa-siemens).

To register for these courses and many more, visit <http://www.audiologyonline.com/c/signa-siemens>



Life sounds brilliant.

CEUs Available



June 27, 2016 at 12 PM ET

Driving a Hearing Aid on Highway 80

Presented by Yu-Hsiang Wu, MD, Ph.D.

It is important to know whether new interventions, such as hearing aid technologies or fitting strategies, deliver greater benefit to listeners with hearing impairment than older interventions. Measuring the intervention benefits - or the outcomes - is not as straightforward because many factors, such as the acoustic characteristics of listening environments, can



June 27, 2016 at 12 PM ET

Driving a Hearing Aid on Highway 80

Presented by Yu-Hsiang Wu, MD, Ph.D.

It is important to know whether new interventions, such as hearing aid technologies or fitting strategies, deliver greater benefit to listeners with hearing impairment than older interventions. Measuring the intervention benefits - or the outcomes - is not as straightforward because many factors, such as the acoustic characteristics of listening environments, can

effect the outcomes. In this talk, a series of studies will be presented to illustrate the effect of hearing aid directional microphone technology on speech understanding and listening effort in a laboratory setting and in an automobile. The use of dual task paradigms to measure listening effort will also be discussed.



July 8, 2016 at 12 PM ET

Speech-in-noise testing for selection and fitting of hearing aids: Worth the effort?

Presented by H. Guenter Mueller, Ph.D.

Hearing aid fitting protocols must be efficient, and some have questioned if there is value in including speech-in-noise testing. In this course, we'll review evidence showing that indeed, the findings from this testing can be used for the selection and adjustment of special features, selection of accessories, and most importantly, patient counseling. We'll also discuss what specific speech tests easily can be implemented into a busy schedule.



July 26, 2016 at 12 PM ET

Impact of auditory access on speech and language development

Presented by Mary Spafford, Au.D.

Recent results from the Outcomes of Children with Hearing Loss (OCHL) study indicate that auditory access - aided auditory levels (i.e., speech intelligibility index) and amount and duration of hearing aid use - contributes to speech and language development for children who are hard of hearing. Longitudinal trends in auditory and hearing aid use show some children to be at constant risk for decreased auditory access, thus also at risk for speech and language delay. Auditory access profiles for individual children will be discussed relative to speech and language outcomes and recommendations for audiological/ amplification management.

Hearing
Systems

SIEMENS

Copyright © 2016 Signia GmbH. All rights reserved. Siemens AG is a trademark of Siemens AG. 016-517658-13

signia
Life sounds brilliant.

primax
Clinically proven better than normal
hearing with less listening effort

Thank you!

Hearing
Systems
SIEMENS

