

## ReSound Sound Shaper

Understanding, fitting and verifying for 2016 and beyond

Jenn Schumacher, AuD • GN ReSound

GN Making Life Sound Better

## Today's Agenda

Introduction to frequency compression

Description of Sound Shaper

Update on current research

Candidacy for Sound Shaper

Fitting and verifying Sound Shaper

Case studies

Q & A

GN Making Life Sound Better

ReSound GN


## Learning Objectives

After this course learners will be able to describe how the Sound Shaper feature works.

After this course learners will be able to define candidacy criteria for fitting Sound Shaper.

After this course learners will be able to explain how to apply tips, current research and case studies in this presentation to fit and verify Sound Shaper.

## Starts with a problem



## Introduction to frequency compression

### **Audibility of high frequency sounds cannot be adequately restored with conventional amplification**

Consonants- high frequency and/or softer intensity

More difficulties in background noise

Environmental sounds and localization cues

Speech production

## Introduction to frequency compression

### **Dead regions**

Areas of dead inner cells in cochlea

Amplifying these frequencies can decrease speech understanding

Listeners with dead regions may still benefit from wider bandwidth amplification

### **Issues with physical/acoustic limitations**

High-powered receivers

Must have tight seal to ensure high frequency amplification is reaching TM and to prevent feedback

## The Solution: ReSound Sound Shaper

### **Provide audibility for high frequency sounds not aidable via conventional amplification**

High frequency sounds are moved to a lower frequency region where there is aidable hearing


### **Preserve sound quality**

Attempts to minimize perceptual distortion by using linear (proportional) frequency compression


Frequencies below cutoff are not altered

Implementation of Sound Shaper has been simplified

## Frequency compression



## Introduction to Sound Shaper



GN Making Life Sound Better

ReSound GN

## Introduction to Sound Shaper

### Nonlinear (non-proportional) frequency compression

Most common type used in today's hearing aids

Frequencies closer to start frequency are shifted less than frequencies at the higher range of the source region

When activated, it is always operating

### Proportional frequency compression

Frequencies that are lowered maintain a "constant" relationship with target frequencies

*Source frequency \* constant = target frequency*

$$6000 \text{ Hz} * 0.6 = 3600 \text{ Hz}$$

$$5000 \text{ Hz} * 0.6 = 3000 \text{ Hz}$$

$$4000 \text{ Hz} * 0.6 = 2400 \text{ Hz}$$


Harmonic relationships are maintained

When activated, algorithm only lowers signal when speech is detected

GN Making Life Sound Better

ReSound GN

## Introduction to Sound Shaper



Unprocessed  
Pure tone sweep


Non-proportional

Proportional

GN Making Life Sound Better

ReSound GN

## Introduction to Sound Shaper



GN Making Life Sound Better

ReSound GN

## What outcomes do we see with frequency compression?

### Data reviews from:


Simpson (2009) 4 studies  
Alexander (2013) 5 studies

Adult hearing impaired listeners

Using simulated and commercially-available frequency compression techniques

Speech in quiet

### When using frequency compression over conventional amplification, is speech understanding:



GN Making Life Sound Better

ReSound GN

## What other outcomes have been measured?

### Speech understanding in noise

Benefit from FC in noise (McCreery et al, 2014)

No benefit from FC in noise (Ellis & Munro, 2013; Millet et al, 2016)

Some benefit in spatially-separated noise for those with poor SNR abilities (Shehorn et al, 2013)

### Sound quality

Cutoff frequency and CR affect sound quality for HI listeners (Souza et al, 2013; Parsa et al, 2013; Johnson & Light, 2015)

Sound quality is poorer than conventional, even at low CF/CR for speech and music (Souza et al, 2013; Mussoi & Bentler, 2015)

A balance between audibility and sound quality can often be reached (Johnson & Light, 2015)

### Binaural cues

FC could reduce perception of binaural cues (Brown et al, 2016)

GN Making Life Sound Better

ReSound GN

What's the conclusion on frequency lowering?

**It depends!**

What are your goals?

Who is your patient?

How are you setting the feature?

Research can also guide our decisions on candidacy and choosing settings...

What factors may affect candidacy?

**Hearing loss**

Still the most powerful predictor of frequency compression benefit (Souza et al, 2013; Ellis & Munro, 2015)

**Age**

Older listeners may benefit more from FC than younger listeners (Kokx-Ryan et al, 2015)

**Cognitive abilities**

Reduced FC benefit seen in listeners with poorer working memory (Arehart et al, 2013)

WM + hearing loss appears to interact (Souza et al, 2015)

Impact of cognitive status may be minimized by using conservative settings that are individualized for patients (Ellis & Munro, 2015)



## Sound Shaper Candidacy

### Which patients are candidates for frequency compression?

Sound Shaper is designed for “typical” patient profile:

Severe to profound high frequency sensorineural hearing loss

Steeply sloping thresholds

Not so “typical” but may still benefit:

Mild-moderate hearing loss

Do see some benefits for high frequency consonant identification for this group (Alexander et al, 2014)

Flat severe to profound hearing loss

### If you are concerned about distortion

Consider measure of speech in noise, such as QuickSIN

Unaided QSIN performance best predictor of speech abilities using FC (Kates et al, 2013)

### Patients should be evaluated on an individual basis

GN Making Life Sound Better



## Sound Shaper Candidacy

### Which patients may Sound Shaper NOT be recommend for?

Recommended fitting software setting is Off

Patient has known cognitive issues, especially if they are elderly

GN Making Life Sound Better



## Sound Shaper Candidacy

**If you suspect a patient is not or will not receive audibility of high frequency sounds, measure it!**

## Fitting Sound Shaper

**If you feel your patient would benefit from Sound Shaper, how do you determine settings?**

Activate most conservative setting that provides high frequency audibility your patient is missing with conventional amplification

As FC settings become more aggressive, speech intelligibility performance can decrease (e.g., Souza et al, 2013)

## Fitting Sound Shaper

### Three settings for Sound Shaper in fitting software:

Mild (4 kHz)

Moderate (3.5 kHz)

Strong (2.5 kHz)

**Note that the default setting is Off**

Sound Shaper Setting	FC Knee Point	FC Ratio
Mild	4000 Hz	1.33:1
Moderate	3500 Hz	2:1
Strong	2500 Hz	2:1

## Fitting Sound Shaper

### “Mild”

Recommended if the audiogram has a slope  $\geq 10$  dB per octave frequency, and the slope begins  $\geq 4000$  Hz

### “Moderate”


Recommended if the audiogram has a slope  $\geq 10$  dB per octave frequency, and the slope begins at 2000 Hz

### “Strong”

Recommended if the audiogram has a slope  $\geq 10$  dB per octave frequency, and the slope ends  $\geq 2000$  Hz

## Sound Shaper settings

Setting ID	Cut-off [Hz]	CR
"Off"	n/a	n/a
1 (weakest)	5000	1.33
2	4000	1.33
3	4000	2
4	3500	2
5	3000	2
6	2500	2
7	2250	2
8 (strongest)	2000	2



GN Making Life Sound Better

ReSound GN

**Aventa 3.6**

File Edit View Instrument Fitting Help

**Start** Pre-Fit Fit Summary

Right Ear: UN961-DRW HP RIE Left Ear: UN961-DRW HP RIE

Gain Curves: [x] 50 [x] 65 [x] 80 Target Curves: [x] 50 [x] 65 [x] 80

dB

125 250 500 1K 2K 4K 8K 16K Hz

P1: Binaural Directionality P2: Restaurant P3: Acoustic Phone P4: None PhoneNow: None Streamers: TV / Mini Mic Phone Clip: Phone Clip

Remove Reorder... Create Comparison Recalculate Autorelate Copy Paste... Manage Program Names

**Features**

Binaural Directionality: Automatic

Directional Mix: Very Low

DPS Ultra II: Off

Auto DFS: On

Expansion: Off

Sound Chimes: Off

**NoiseTracker II:** Off (4K) Moderate (3.5K) **Strong (2.5K)**

Windguard: Off

Select Sound: 00:00/00:00

Calibrate DFS

Right side: Automatic Very Low Off On Off Off Off Off Off


Left side: Automatic Very Low Off On Off Off Off Off Off

Off Off Mild (4K) Moderate (3.5K) **Strong (2.5K)**

ReSound GN

GN Making Life Sound Better

ReSound GN



## Fitting Sound Shaper

### Asymmetric hearing loss

- Sound Shaper can be set differently for each ear in the fitting software
- Consider choosing one Sound Shaper setting and using it on both sides
- May want to use more conservative setting in both ears
- Different Sound Shaper settings may be warranted in some cases

## Verifying Sound Shaper

### **Number one goal is to ensure there is audibility using real ear measurements**

Sound Shaper is providing adequate audibility of high frequencies

Also, that there is not audibility without Sound Shaper (if you are not sure)

### **Real ear verification will not be able to establish possible speech understanding improvements**

Nonsense syllable/phoneme comparison tests are best for showing benefits from FC algorithm

## Verifying Sound Shaper

### **Recommended real ear stimuli**

Long term average speech spectrum (LTASS) at various input levels


Ling /s/ in Aurical FreeFit

Filtered speech in Audioscan Verifit

## Verifying Sound Shaper

### 1) Measure LTASS at various input levels without and with Sound Shaper

Ensure that there is a need for Sound Shaper, and that Sound Shaper does not reduce high-frequency bandwidth




GN Making Life Sound Better

ReSound GN

## Verifying Sound Shaper

### 2) Use high-frequency stimuli to verify audibility




Verifit  
filtered  
speech

GN Making Life Sound Better

ReSound GN

## Verifying Sound Shaper

### 2) Use high-frequency stimuli to verify audibility



Aurical  
FreeFit  
Ling  
sounds

GN Making Life Sound Better

ReSound GN

## Verifying Sound Shaper

### 3) Determine other effects Sound Shaper may have on signal

Measure separation of /s/ and /sh/ phonemes

GN Making Life Sound Better

ReSound GN



## Case study 1- Mr. L

### 70 year old male, mild to profound noise-induced sensorineural hearing loss

Experienced hearing aid user (10 years)

Reported extreme distortion when listening to speech, that his own hearing aids were not very helpful

Confirmed dead region bilaterally using TEN test (Moore et al, 2000)

Importance of wider bandwidth for amplification


Tested Mr. L using varying degrees of high frequency bandwidth- He performed best with greater high-frequency bandwidth



GN Making Life Sound Better

ReSound GN

## Case study 1- Mr. L



### Hearing aid fitting:

RIE with Ultra power receiver

High frequency gain was maximized before activating Sound Shaper, but ran into feedback issues

Sound Shaper = Strong (2.5k) bilaterally

### Mr. L's experience:

Mr. L required several weeks of wear time with Sound Shaper to begin demonstrating benefits for phoneme perception

Benefits were still small, even after experience with Sound Shaper. But there was improvement compared to own hearing aids

GN Making Life Sound Better

ReSound GN


## Case study 2- Mr. J

### 40 year old male, severe congenital sensorineural hearing loss

Lifelong hearing aid user- first time frequency compression user

Father of two young kids


Avid music lover



GN Making Life Sound Better

ReSound GN

## Case study 2- Mr. J



GN Making Life Sound Better

### Hearing aid fitting:

RIE with Ultra power receiver

High frequency gain was maximized before activating Sound Shaper

Sound Shaper = Strong (2.5k) bilaterally

It was helpful to give Mr. J 2 take-home programs when fitted: Program 1 used Sound Shaper and Program 2 used conventional amplification

### Mr. J's experience:

Better audibility of /s/, though some confusions with /s/ and /sh/

Reported better environmental awareness, noticed audibility of high frequency consonants for the first time in other people's speech

Music was more difficult to accept- He could pick up on changes in pitch with Sound Shaper compared to conventional amplification

ReSound GN


### Case study 3- Ms. R

#### 30 year old female, mild to moderate sensorineural hearing loss

Experienced hearing aid user (10 years), desired more benefit

Based on her hearing loss, Ms. R is not a "typical" FC candidate


She was willing to try something novel for more audibility of high frequencies



GN Making Life Sound Better

ReSound GN

### Case study 3- Ms. R



#### Hearing aid fitting:

Open fit BTE with power dome

High frequency gain was maximized before activating Sound Shaper

Sound Shaper = Mild (4k) bilaterally

#### Ms. R's experience:


Subjectively preferred FC over conventional amplification

She appreciated the new stimuli she was hearing

Noticed the change in sound quality- but liked it as it meant she was hearing more high frequency sounds

GN Making Life Sound Better

ReSound GN



Thank you!

Any questions?

[jschumacher@gnresound.com](mailto:jschumacher@gnresound.com)

GN Making Life Sound Better

ReSound GN