



**Hybrid L24 Study Update:
Where Are They Now?**

Ginger Grant, AuD
Audiology Online: January 18, 2017

Hear now. And always



Learning Objectives



1. After the completion of this course participants will be able to list the benefits of acoustic hearing.
2. After the completion of this course participants will be able to describe the expanded indications for the Nucleus® Hybrid™ system.
3. After the completion of this course participants will be able to outline the steps to programming a Nucleus Hybrid device.



Why Hybrid?



Audibility 125 → 8000Hz



High frequency electric hearing is critical for improved speech understanding

Benefits of Acoustic Hearing

- ✓ Enhanced sound quality
- ✓ Improved understanding in background noise
- ✓ Fundamental frequency cues for both pitch and vowel discrimination
- ✓ Interaural timing difference cues

Good frequency resolution

Restores HF Sensitivity



Pivotal Clinical Trial

Hear now. And always.



Background



More than 15 years of research

- 8 year FDA trial of 50 adult subjects

FDA-approved March 2014

- 1st Hybrid implant system (U.S.) April 16, 2014
- 1st cochlear implant labeled for potential preservation of residual hearing
- 1st Indication using CNC words
- 1st Integrated Electric + Acoustic Sound Processor



Study Design



- Multicenter pivotal trial
- 14 surgeons at 10 US-based implant centers
- 50 adult subjects implanted
- Primary endpoint was 6 months postactivation




Pivotal IDE Study



Subject Demographics	N=50
Gender	
Male	25 (50%)
Female	25 (50%)
Mean Age at Implantation (yrs)	64.1 [37-86]
Mean Duration Hearing Loss (yrs)	28.1 [3-74]
Mean Duration Sev/Prof Loss (yrs)	13.1 [2-30]
Mean Preoperative CNC Score CI Ear	28.4% [9-64%]
Mean Preoperative AzBio Sentence Score CI Ear	16.3% [0-64%]

Hybrid L24: An Expanded Indication



Ear to be Implanted

- Aided CNC word score between 10% and 60% correct, inclusively

Contralateral Ear

- Aided CNC word score better than ear to be implanted but less than 80% correct

Audiometric

- Severe to profound HF SNHL bilaterally

Adults

- aged **18 years** and older

Key Findings of Pivotal Study




United States Multicenter Clinical Trial of the Cochlear Nucleus Hybrid Implant System

J. Thomas Roland Jr, MD; Bruce J. Gantz, MD; Susan B. Waltzman, PhD; Aaron J. Parkinson, PhD;
The Multicenter Clinical Trial Group

Published results demonstrate:

- Significant improvement in speech understanding in quiet and noise
- Expanded indications for individuals with severe high frequency hearing loss
- Increased satisfaction, consistent with speech perception results

Laryngoscope_2016 Jan; 126(1): 175-81.



Recent Updates

Hear now. And always.

Reimbursement

Hear now. And always.



Reimbursement Success



Since August 2016...

142,000,000

Covered Lives Now Have Access To Hybrid Coverage



BlueCross
BlueShield
62 Million

Anthem

40 Million



23 Million



9 Million

FEHB
Federal Employees
Health Benefits
Program

8 Million

Acoustic Component

Hear now. And always.



Cochlear™ Hybrid Hearing™



- Amplify low frequency hearing and provide access to high frequencies via electrical stimulation
- The Acoustic Component is compatible with all N6 Sound Processors (CP910 and CP920)
- No processor exchange needed
- Available for use by all recipients, pediatric and adult, with **all Nucleus implants**



The Acoustic Component should only be used when behavioral audiometric thresholds can be obtained and the recipient can provide feedback regarding sound quality. Not available for the Nucleus 22 series of implants.

Hybrid™ Hearing Availability Today:



Nucleus® Hybrid™ Hearing Now FDA
Approved for **ALL** Nucleus Implants



Cochlear™ Nucleus® Profile Slim Modiolar Electrode



Cochlear™ Nucleus® Contour Advance®



Cochlear™ Nucleus® Slim Straight



Cochlear™ Nucleus® Full-Band Straight



Cochlear™ Nucleus® Hybrid™ L24



*The Acoustic Component should only be used when behavioral audiometric thresholds can be obtained and the recipient can provide feedback regarding sound quality.

Hybrid L24 Clinical Trial Outcomes Update

Hear now. And always.



Background:



- All Hybrid implant subjects upgraded to the N6 sound processor
 - 39/50 (78%) consented to long term follow-up and upgrade study
- Questions:
 1. What does long-term hearing preservation look like?
 2. What does Hybrid performance look like?
 - > Electric Alone?
 - > Electric + Acoustic?

Hearing Preservation



Established Definitions:

1. Preservation of measurable^{1,2} acoustic hearing
 - Definition used by industry and in most publications
 - Measure of detection/audibility but not usability
2. Preservation of functional^{1,2} acoustic hearing
 - A **severe or better** degree of hearing loss post-implantation that may provide a level of acoustic hearing to enhance electrical stimulation

1. Lenarz et. AL (2013) European Multi-center Study of the Nucleus Hybrid L24 cochlear implant. International Journal of Audiology. 1-11
 2. Roland, et AL (2015) United States Multicenter Clinical Trial Cochlear Nucleus Hybrid Implant System. Laryngoscope, 2-7

Evolution of Sentence Materials



- CUNY Sentences were used in our earliest clinical trials (1999) with Hybrid implants but many recipients quickly topped out
- Members of AAA and AAO-HNS established use of HINT sentences in noise as part of the Minimum Speech Test Battery developed in 1996
- AzBio Sentences were developed by Spahr and Dorman in 2004
- Cochlear's Hybrid L24 trial was implemented with use of AzBio Sentences in noise in an attempt to avoid ceiling effects
- The recommended MSTB protocol was updated in 2011 to replace HINT sentences with AzBio Sentences

Power of the Nucleus® 6:



ACOUSTIC

- ASC
- **SCAN**
- **Wind Noise Reduction**

ELECTRIC

- ADRO®
- ASC
- **SCAN**
- **Wind Noise Reduction**
- **Background Noise Reduction**
- **Whisper™**




TRUE WIRELESS

Speech Perception in Noise



AzBio Sentences in Noise (+5dB SNR)



N=39

Next Steps



- Long term follow up post approval study
 - Original pivotal Hybrid L24 subjects followed to 5 years postactivation
 - The majority of subjects have completed all requirements
 - Publication planned
- Newly implanted post approval study
 - Cochlear's initiative to gather even more 5 year data
 - 19 implant centers currently enrolling subjects
 - Data collection and analysis ongoing




Programming Considerations

Hear now. And always




Smart Hearing Alliance




Cochlear®

Industry's only
TRUE WIRELESS
bimodal solution

ReSound



What have we learned about programming?



Hybrid Mode: Cut-off Frequency



What's changed in Custom Sound® 4.4?

Now each stimulation mode has its own cut-off frequency. We no longer have what we called a 'cross-over frequency'.

- **Acoustic** cut-off at ≤ 90 dB HL up to 2 kHz *remains unchanged*
- **Electric** cut-off has been *lowered* to >70 dB HL up to 2 kHz

This results in more overlap for some audiograms between acoustic and electric stimulation, especially for shallow sloped audiograms.


NOTE: This change *only applies when creating new MAPs* (for existing or new recipients). The software uses the currently entered audiogram and applies the above calculations.

Acoustic Component Specifications



Measurements According to ANSI S3.22-2009


Output Sound Pressure Level (OSPL 90)



SPA-OSPL90 = 116.5 dB SPL

OSPL90 max = 117.5 dB SPL

Full-on Gain and Reference Test Gain



Full-on gain
Maximum = 54.7 dB

SPA-full-on gain = 48.8 dB

Reference test gain
SPA-Reference test
gain = 39.9 dB


Frequency range:
Between 100 Hz and 2200 Hz

Hybrid Mode: Frequency Bands



≤90 dB HL
cut-off
ACOUSTIC
bands enabled
up to 2 kHz

Acoustic
Unchanged






Electric
Revised

>70 dB HL
cut-off
ELECTRIC
Bands
enabled

Robust Custom Sound algorithm:

- Checks the Acoustic cut-off frequency at 90 dB HL, up to 2 kHz
- Checks the Electric cut-off frequency at 70 dB HL, up to 2 kHz

Electro-Acoustic Cross-over Ranges




Results in electro-acoustic overlap between 563 Hz and 685 Hz

Steps in Programming



1. Create a new recipient



2. Measure Impedances




3. Create a New MAP (will default to 37 PW) & check the box for Hybrid

Create a New MAP

Processor	Strategy	Mode	Rate	Maxima	Pulse Width	Hybrid Mode
CP900 Series sound processor	ACE / ACE (RE)	MP1+2	900	8	37	<input checked="" type="checkbox"/>

Steps in Programming



4. Remove the coil from the head



- Open the Acoustics screen

Set Prescription

Method: **NAL-RP**

☒ WDRC ☐ Linear


Mould Type: **Custom**

Vent Size (mm): **0**

Recalculate

- Select the fitting prescription, mould type, compression and vent size.

Steps in Programming



5. Go Live with acoustic stimulation only

- Complete Real Ear measurements
- Adjust gain and MPO if needed
- Take into account the subjective feedback of the recipient.


If the desired gain cannot be met acoustically:

- enable adjacent channel to obtain additional gain
- disable acoustic channels not meeting target
- the lower CI frequency boundary may need to be lowered further in the next step

Steps in Programming



6. Measure T levels using the streamlined programming method. Go live to obtain C levels so sound is comfortably loud.



Coil must be on the head


7. Loudness
Balance C levels

Acoustics Panel Options




What's new in the Acoustics panel in CS 4.4?


It is now possible to make very quick loudness comparisons between Hybrid mode and full electrical mode.



Hybrid MAP Checkbox




Mute Acoustic Checkbox



Summary of Quick Loudness Checks



CI only
(full frequency allocation)





Listen in CI mode

Finalize Programming



8. Check battery suitability and personalize MAPs by adding SmartSound® iQ
 - Program 1 primary MAP with SCAN
 - Program 2 default
 - Program 3 and 4 can be used to try other cut offs with SCAN.
9. Write MAPs to the processor and finalize
10. Verification of audibility




Summary



- Hybrid indication is continuing to gain momentum
- Programming of both the electric and acoustic portions can be completed with Custom Sound® Software
 - Overlap between acoustic and electric stimulation will vary by patient
- Hybrid L24 + Nucleus 6
 - 95% of subjects perform the same or better on AzBio +5 dB SNR

Gilden et al (2015). Improved hearing in noise using new signal processing algorithms with the Cochlear Nucleus 6 sound processor. Jour of Otol 10: 51-56.

Summary



References



1. Dorman MF & Spahr, A.J. Speech perception by adults with multichannel cochlear implants. In Cochlear Implants, S.B. Waltzman & J.T. Roland (Eds.) Thieme Medical Publishers, NY, Stuttgart, 2006.
2. Ebrahimi-Madiseh et al (2016). Speech perception scores in cochlear implant recipients: An analysis of ceiling effects in the CUNY sentence test (Quiet) in post-lingually deafened cochlear implant recipients. Cochlear Implants Int 2016 17(2): 75-80.
3. Gifford RH and Dorman MF. The Psychophysics of Low-Frequency Acoustic Hearing in Electric and Acoustic Stimulation (EAS) and Bimodal Patients. J Hear Sci 2012 May 1; 2(2): 33-44.
4. Gilden et al (2015). Improved hearing in noise using new signal processing algorithms with the Cochlear Nucleus 6 sound processor. Jour of Otol 10: 51-56.
5. Lenarz et. Al (2013) European Multi-center Study of the Nucleus Hybrid L24 cochlear implant. International Journal of Audiology. 1-11
6. Roland, et al (2016) United States Multicenter Clinical Trial Cochlear Nucleus Hybrid Implant System. Laryngoscope, 126(1): 175-81.
7. Sheffield SW, Jahn K, Gifford RH. Preserved Acoustic Hearing in Cochlear Implantation Improves Speech Perception. J Am Acad Audiol 2015 Feb; 26(2): 145-154.



Hear now. And always



Like us on Facebook: Cochlear Americas



Follow us on LinkedIn: Cochlear/medical devices



Follow us on Twitter: @CochlearUS

Visit our website: www.Cochlear.com/us



Hear now. And always

©Cochlear Limited 2017. All rights reserved. Hear now. And always and other trademarks and registered trademarks are the property of Cochlear Limited.

The Acoustic Component should only be used when behavioral audiometric thresholds can be obtained and the recipient can provide feedback regarding sound quality; The Hybrid L24 Implant is approved in the US for adults ages 18 and older; SNR-WR, WNR and SCAN are approved for use with any recipient ages 6 years and older, who is able to: 1) complete objective speech perception testing in quiet and in noise in order to determine and document performance; and 2) report a preference for different program settings.

Cochlear Americas is pleased to provide this guidance on insurance coverage about Cochlear's products and technology, but please realize that this information is provided as guidance only. Cochlear Americas makes no representation, warranty or guarantee regarding the fitness for your particular situation of the guidance provided or that following these guidelines will result in any form of coverage or reimbursement. Please be sure to seek legal advice and consult with your own reimbursement advisors to confirm the applicability of our guidance to your particular needs. And remember that all services or products billed to an insurance company or payer must be medically necessary, actually performed or purchased, and appropriately documented.