



Connecting Patients to Their Favorite Activities: With ReSound's Multi and Micro Mics

Neil A. Wright, AuD

GN Making Life Sound Better

Welcome!

Quick Bio:

- **Neil Wright, AuD**
- **Audiologist at GN ReSound in Glenview, IL**
- **Received AuD from Northwestern University**
- **Financial Disclosure: I am an employee of GN ReSound**

GN Making Life Sound Better

ReSound GN

Outline

- **(0-5) – Introductions**
- **(5-10) – Introduction to remote microphone technology**
- **(10-20) – Restaurants: Improving SNR**
- **(20-30) – Cars: Compensating for lack of visual cues**
- **(30-40) – Classrooms and lectures: Listening at a distance**
- **(40-50) – Meetings and Large groups: When its all working against you**
- **(50-55) – Other functions and scenarios of benefit**
- **(55-60) – Conclusions and Q&As**

Learning Goals

- **To quantify the benefits of Remote Microphone Technology**
- **To identify environmental factors that negatively impact hearing aid user performance**
- **To demonstrate programming strategies to maximize hearing aid and remote microphone user performance in difficult listening environments**



ReSound Micro Mic



ReSound Multi Mic



Remote Microphone Technology

GN Making Life Sound Better

The Current State of Hearing Aid Users

- **Performance of hearing aid technology has perceptibly improved overall**
 - Marketrak IX¹: 80% satisfaction rate across all environments for devices <5 years old
- **BUT, many environments still cause hearing aid users problems**
 - Marketrak VIII² (2010): Satisfaction by environment →
 - Marketrak IX (2015): Listening in noise: lowest satisfaction
- **The Usual Suspects**
 - Noise
 - Distance
 - Reverberation
 - Limited Visual Cues / Multiple Talkers
 - Often in combo

GN Making Life Sound Better




Figure 9. Consumer satisfaction with hearing aids in various listening situations where hearing aids are <=4 years of age.


Kochkin S. MarkeTrak VIII: Consumer satisfaction with hearing aids is slowly increasing. *Hear Jour.* 2010;63(1):19-27

ReSound GN

How do we address these environmental struggles?

HA Programming Options


- Directional Microphones³
 - Digital Noise reduction
 - Gain Reduction
 - Environment-specific programs
 - **Measureable benefits, but have limitations**
- When these aren't enough, what now?**
- **Remote Microphones!**



Kochkin S. MarkeTrak VIII: Consumer satisfaction with hearing aids is slowly increasing. *Hear Jour.* 2010;63(1):19-27

Remote Microphone Technology: At a Glance

- **Worn by conversation partner**
 - Improves SNR in noise, speech perception and understanding in noise and at a distance
- **History: a wired connection**
- **Key piece of FM systems**
 - First wireless remote microphone
 - Significant benefits in noise⁴
 - Universal but more expensive
- **Digital Wireless has made them more affordable, practical, & accessible**



1st Generation: The Mini Mic

- **Released: Spring 2012**
- **Omni directional remote mic**
- **Line-in functionality**
- **Truly Wireless connectivity**
 - 2.4 GHz wireless connection
 - No intermediate device
 - 30 ft. range
- **Significant speech perception benefit in noise^{5,6}**
- **Significant benefit at a distance^{5,6}**



GN Making Life Sound Better

ReSound GN

Next Generation: Multi Mic and Micro Mic


- **Multi and Micro Mic advantages over 1st Gen Mini Mic**
 - Directional mic towards the speaker's mouth
 - Longer distance of transmission: 90 feet – clear line of sight
 - Drop Detection
- **Multi Mic Advantages**
 - Omni Directional Table Mode
 - Built-in Telecoil
 - DAI compatible
 - Line in functionality
- **Micro Mic**
 - small, discrete size



ReSound Micro Mic



ReSound Multi Mic



GN Making Life Sound Better

ReSound GN

Streaming Program Default Settings

- **Hearing aid mics always in Omni**
- **Mic Balance function: defaults to +6 towards the Remote Mic**
- **Wireless-specific programming: Gain settings mimic All-Around program**
 - Independent Gain Settings
- **Most advanced features default “OFF”**
- **Noise reduction defaults to “Mild”**

Manage Programs Create Comparison Manage Program Names

Features

	Right side	Left side
Time Constants:	Syllabic	Syllabic
Streamer BassBoost:	Mild	Mild
DFS Ultra II:	Off	Off
Auto DFS:	On	On
Expansion:	Off	Off
Sound Shaper:	Off	Off
NoiseTracker II:	Mild	Mild
HI Mic / TV:	3	3
HI Mic / Micro-Multi Mic:	6	6

GN M2



Listening in Restaurants:

Improving the Signal to Noise Ratio

Restaurant: Challenges

Biggest Challenges


- **Signal to Noise Ratio (SNR)**
- **Background noise comes from multiple directions**
- **Poor room acoustics**
- **Poor lighting**

Why They Matter

- **Hard to pick out target speech from background**
- **Background noise typically speech**
- **Fixed directionality may not be enough**
- **Sound rebounding, reverberation adds to noise**
- **Can impact visual cues**


Remote Mic Benefits in Noise

- **Data from Rodemerk and Galster⁽²⁰¹⁵⁾⁶**
- **HINT measurement tested across different types of remote microphones**
 - 6' from speaker and remote microphone
 - Hearing aid mics muted
- **Remote microphone technology provides significant SNR and speech perception benefit over HAs alone**
- **HA mics were set to Omni**
 - What about directionality?



Rodemerk K, Galster J. The benefit of remote microphones using wireless protocols. *J Acad Audiol*. 2015;26:724-731.

Improving SNR: Improvement over Mini Mic



- The new ReSound Multi Mic on average provides 15dB of SRT benefit above what directionality provides.
- It provides a significantly better SRT improvement (Tukey's Honest Significant Difference) than the existing ReSound Mini Mic and an available smartphone solution.
- The new Multi Mic on average provides 5 dB more SRT benefit compared to the first generation Mini Mic.

Jespersen & Kirkwood, 2016

GN Making Life Sound Better



Programming Recommendations for Restaurants

- **Low Frequency Gain Settings**
 - Decrease overall gain at 750 Hz & below
- **Microphone Balance**
 - Slide balance towards the remote microphone
 - Will increase input from remote mic relative to hearing aid microphones
- **Mid & High Frequency Gain Settings**
 - Increase mid level gain in mid & high frequencies to increase target speaker's voice clarity
 - 1000-6000 Hz, less increase at 6 kHz
 - Similar gains to restaurant program, but in mid-level gain instead of soft gain




GN Making Life Sound Better




Programming Recommendations for Restaurants: Continued

- **Increase Noise Tracker**
 - Some built in noise reduction in Multi Mic
- **Utilize Expansion**
- **Vertical Orientation key for lanyard (clipped on) usage**
 - Directional Microphone will activate in noise for Multi Mic
 - Micro Mic always in directional



GN Making Life Sound Better

ReSound GN



Listening in the Car:

Compensating for a Lack of Visual Cues

GN Making Life Sound Better

Listening in the Car: Challenges

Biggest Challenges

- **Signal to Noise Ratio (SNR)**
 - High levels of background noise⁷
- **Driving demands high level of attention**
- **Limited Visual Cues**
- **Indirect Speech – coming from behind or side of patient**

Why They Matter

- **Where the rubber meets the road: makes for a high noise level**
- **Driver cannot turn to speaker, and/or driver should not turn to listener**
- **Visual Cues contribute to speech understanding for hearing impaired**
- **Fixed directionality – a conventional approach to noise reduction – is detrimental to speech understanding from back/side**

GN Making Life Sound Better

ReSound GN

Importance of Visual Cues

- **Visual cues critical for speech understanding**
 - Lipreading and Speechreading
- **Visual cues can contribute up to +4-6 dB SNR benefit⁸**
 - Often missing in car listening environment
- **Face-to-face conversations uncommon in the car**
 - Negative impact on speech perception
- **Speakers often to the side, in front, or behind the listener**
 - Fixed directionality will decrease speech understanding in this scenario



GN Making Life Sound Better

ReSound GN

Programming Considerations for the Car

- **Reducing road noise input**
 - Decreasing wideband loud input gain
 - Decrease low frequency (below 1 kHz) gain
- **Noise Reduction**
 - Some built-in noise reduction in Multi Mic
 - Can increase Noise Reduction setting in fitting software
- **Microphone balance**
 - Advanced settings can be ear-specific in streaming program
 - Window-facing ear slide towards remote mic further than inward facing ear



GN Making Life Sound Better

ReSound GN

Programming Considerations for the Car: Continued

- **Seating arrangement**
 - Depending on seating, ear-specific programming could be beneficial
 - No Wind reduction in streaming program
- **App-based favorites**
 - DON'T HEARING AID PROGRAM AND DRIVE!
 - Make specific favorites that are easy to select prior to getting on the road
- **Can make very specific driving or seat-specific settings**
 - Mute left hearing aid mic for driver
 - Increase Multi Mic input
 - Save as labeled Favorites



GN Making Life Sound Better

ReSound GN



Lectures and the Classroom

GN Making Life Sound Better

Listening in the Classroom/Lecture: Challenges

Biggest Challenges

- Long distance between speaker and HA user
- Speaker may move around, increasing distance or turning away from HA user
- Signal to Noise Ratio (SNR)
 - High levels of background noise
 - Levels of reverberation

Why They Matter


- Increased distance reduces speaker volume at listener
 - Limits visual cues
- Moving speaker may negate conventional approach to noise reduction – may introduce distraction
- Noise and reverberation increase listening difficulty and decrease speech understanding
 - Directionality alone may not be enough

GN Making Life Sound Better

ReSound GN

Remote Microphone Benefits at a Distance

- **Directional Mics have a near-field effect**
 - Directional benefit diminishes quickly over longer distances
- **Remote Microphone technology allows for SNR benefits in noise to be preserved across distance**
- **Research by Jespersen & Laureyns (2012)⁵ measured SRT in noise through ReSound's Mini Mic**
 - Distances: 1.5m, 3m, & 6m
 - 3 Conditions: Directional Microphones alone, HA Mic + Mini Mic, Mini Mic with HA Mic Muted




GN Making Life Sound Better

ReSound GN

Remote Mic Benefits at a Distance

- **Results show significantly better SRT results using Mini Mic compared to directional microphones alone⁵**
 - Multi Mic improved range (90 feet)
- **Benefit shown with Mini Mic does not significantly decline with increased distance**
 - Directional microphone performance drops with increased distance




GN Making Life Sound Better

ReSound GN

Remote Mic Benefits at a Distance

- **HINT results show similar performance improvements with remote mic technology** (Rodemerk and Galster (2015)⁶)
- **Testing completed with different remote mic transmission methods**
 - ReSound Mini Mic Highlighted in red boxes
- **Testing performed at 6' and 12'**
 - Remote Microphone Improves performance in noise over HA alone even at 12 feet from target speaker




Rodemerk K, Galster J. The benefit of remote microphones using wireless protocols. *J Acad Audiol*. 2015;26:724-731.

GN Making Life Sound Better

ReSound GN

Programming Considerations for Worship Settings

- **Mic Balance**
 - Speaker often not sole sound source (ex: Music) – predictable turn taking
 - Slide Mic balance to center
- **Seating Arrangement**
 - Advice on seating arrangement to ensure hearing aid microphones can pick up music if desired
- **Noise reduction off**
- **BassBoost as needed**
- **Preach self-advocacy!**



GN Making Life Sound Better

ReSound GN

Programming for Lecture Settings

- **Mic Balance**
 - Slide mic balance towards Multi/Micro Mic
 - One target speaker
- **Expansion**
 - Increase expansion to reduce room noise
- **Seating arrangement**
 - Seating closer preferable if possible for visual cues as needed
- **Preach Self-Advocacy!**



GN Making Life Sound Better

ReSound GN

ReSound Pediatrics: Special Recommendations for Kids

- **Classroom environments will have similar challenges to church and lecture environments**
 - Reverberation
 - Background Noise
 - Good SNR critical for speech development
- **But also have unique considerations**
 - Benefits of incidental learning⁹
 - Moving target speaker
 - Transparency measures for remote mics




GN Making Life Sound Better


ReSound GN

ReSound Pediatrics: Special Considerations for Kids

- **Mic Balance**
 - For pediatrics, defaults to mic balance setting of 3 not 6 (adult)
 - May slide to middle (0 setting) to facilitate incidental learning
- **FM input option**
- **Multi Mic Volume Independent of HA gains¹⁰**
 - Allows wearer to adjust volume without impacting hearing aid gains
- **Volume Lock**
 - Volume settings on Multi Mic can be locked for mic-specific or FM-specific input



GN Making Life Sound Better



Listening in Meetings and Large Groups:
Multiple Speakers, in Noise, at a Distance – the perfect storm!

GN Making Life Sound Better

Challenges in Groups

Biggest Challenges

- **Multiple Noise Sources**
 - From different directions
- **Multiple target speakers**
 - May limit visual cues
- **Some distance from speakers**
 - Larger groups, bigger tables

Why They Matter?

- Noise sources are people – similar frequency content in noise compared to speakers
- Divided attention, overlapping speakers
- Decreased target speaker volume due to distance

Programming Considerations for Large Group Environments

- **Remote Microphone Placement**
 - Horizontal placement to enable Multi Mic table mode
 - Should be as close to equidistant from all speakers as possible from small tables
 - Longer table may require mic to be placed further from hearing aid user
- **Room and table dimensions**
 - For longer tables: Slide mic balance to 0 and place multi mic further away from listener
- **Presence of noise**
 - Slide balance more towards Multi Mic
 - Increase noise reduction






Other Usage Scenarios: Flexibility Through Functionality

GN Making Life Sound Better

Remote Microphone Benefits for Bimodal Patients

- **Bimodal patients struggle in difficult listening situations**
- **Research shows bimodal patients benefit from bilateral streaming**
 - 1st studies used FM
 - Proprietary connections proven effective
- **Benefits from Remote Microphones in noise¹¹**
 - Significant improvements in sentence recognition in noise
 - Speech: 65 dBA at test patients



Wolfe, J. Cochlear Wireless Accessories for Bimodal Users. Presentation, San Antonio, TX, 2015

GN Making Life Sound Better

ReSound GN

Programming Considerations for Bimodal Patients

- **Mic Balance**
 - Independent mic balance controls for HA & CI
 - Research proving effectiveness mic balance settings¹⁰
 - CI at 2:1 Accessory Mixing Ratio
 - HA Mic balance at 6 towards Remote Mic
- **Loudness balancing¹²**
 - Research shows best results for bimodal patients when ears are loudness balanced
- **Residual hearing in HA ear**



GN Making Life Sound Better

ReSound GN

Other Usage Scenarios

- **Tcoil**
 - For looped environments
 - Multi Mic should be vertically aligned (right)
- **Line In**
 - Conventional 9mm jack
 - Plug into headphone jack for easily listening to music and videos
- **FM**
 - Conventional 3-pin DAI port
 - Only requires one FM receiver to stream audio to both ears



GN Making Life Sound Better

ReSound GN

Other Usage Scenarios: Tips

- **Tcoil**
 - Multi Mic needs to be placed within the loop, even if listener is not
- **Line In**
 - Turn off Noise Reduction
 - No expansion for music
- **FM**
 - Multi Mic can remain with HA user, act as volume control for FM input



GN Making Life Sound Better

ReSound GN

Using the ReSound Smart App & Control App

- **Ear-specific volume control**
- **Hearing aid mic & Multi Mic independent volume controls**
 - Mute Controls
- **Discrete way to adjust devices on the fly**
- **Requires MFi connection to use Smart App**
 - MFi compatible devices required (iPhone 5 & beyond)
 - Now available for Samsung Galaxy devices (S5 & beyond)
- **Similar controls through ReSound Control App for non-MFi & Samsung smartphones**
 - Requires Phone Clip +



GN Making Life Sound Better

ReSound GN

Conclusions

- **Remote Microphone technology can help patients access target speech in more difficult listening situations**
 - Improved SNR
 - Improved listening at a distance
 - Improved listening in groups
- **Different environments require different programming considerations**
- **Utilizing clinical programming flexibility for best results**



Thanks for Listening!

Questions?

Email: nwright@gnresound.com



References

1. Abrams HB, Kihm J. An Introduction to MarkeTrak IX: A New Baseline for the Hearing Aid Market. *Hearing Review*. 2015;22(6):16.
2. Kochkin S. MarkeTrak VIII: Consumer satisfaction with hearing aids is slowly increasing. *Hear Jour*. 2010;63(1):19-27.
3. Bentler R. Effectiveness of directional microphones and noise reduction schemes in hearing aids: A systematic review of the evidence. *J Am Acad Audiol*. 2005;16(7):473-484.
4. Boothroyd A. Hearing aid accessories for adults: The remote FM microphone. *Ear Hear*. 2004;25:22-33.
5. Jespersen CT. A review of wireless hearing aid advantages. *Hearing Review*. 2012;19(2):48-55.
6. Rodemerk K, Galster J. The benefit of remote microphones using wireless protocols. *J Acad Audiol*. 2015;26:724-731.
7. Iwao, K., & Yamazaki, I. A study on the mechanism of tire/road noise. *JSAE review*, 1996;17(2), 139-144.
8. Summerfield Q. Lipreading and audio-visual speech perception. *Transactions Royal Soc London B: Biolog Sci*. 1992; 335(1273):71-78.
9. Saffran, J. R., Newport, E. L., Aslin, R. N., Tunick, R. A., & Barrueco, S. Incidental language learning: Listening (and learning) out of the corner of your ear. *Psychological science*, 1997;8(2), 101-105.
10. Kossek P, Jespersen CT. How to verify wireless technology. *Hearing Review*. 2013;20(11):16-20.
11. Wolfe, J. Cochlear Wireless Accessories for Bimodal Users. Presentation, San Antonio, TX, 2015.
12. Ching, T. Y., Incerti, P., & Hill, M. Binaural benefits for adults who use hearing aids and cochlear implants in opposite ears. *Ear and hearing*. 2004;25(1), 9-21.