



TruCore Introduction

For internal use only / © Resound 2016. All rights reserved.

TruCore: Truly Natural




The all TruCore platform provides unparalleled sound quality coupled with a new level of comfort for a more clear, natural, and relaxed hearing experience.

With QuadCore products we provided unprecedented speech understanding. Now with TruCore we complete the listening experience with outstanding comfort features.

For internal use only / © Resound 2016. All rights reserved.

Main Features in TruCore



For internal use only / © Resonance 2016. All rights reserved.

TruCore: Connect to Every Soundscape Naturally



For internal use only / © Resonance 2016. All rights reserved.

TruCore Binaural Processing



- **Benefit:**
 - **Creates an extremely accurate and natural binaural listening experience.**
- Combines directional microphone input of both hearing aids
 - Creates narrow directional beam
 - Engage directionality to front, back, left or right
 - Improves signal to noise ratio in more challenging situations
- Active whenever and wherever required
- Wireless platform

For internal use only / © Resound 2016. All rights reserved.

Comfort Features Voice Ranger



NEW!

- **Scenario:** distraction by competing conversations
- **Solution:** Voice Ranger
- **Benefit:** Reduces distractions and maintains the natural balance between the target speaker and competing speakers providing a more relaxed and enjoyable listening experience.
- Ideal for noisy environments where there are other voices behind the target speaker.
- Separates target speech from speech noise

For internal use only / © Resound 2016. All rights reserved.

Comfort Features Reverb Reducer



- **Scenario:** bothersome sound quality in environments with hard surfaces
- **Solution:** Reverb Reducer
- **Benefit:** Improves the wearer's listening comfort in reverberant environments
- Recommended use:
 - Places of worship
 - Museums
 - Airports
- A unique program option
- Minimizes distraction from reverberation
- Reverb Reducer cuts down echo effect
- Easily added to hearing instrument programming in Program Handling

For internal use only / © Resonance 2016. All rights reserved.

Comfort Features Music Enhancer



Live



Performing



Recorded

- **Scenario:** music lacks natural quality
- **Solution:** Music Enhancer
- **Benefit:** Preserves the emotion and richness of music for greater listening enjoyment
- Recommended use: patients with more discriminating music needs
 - Musicians
 - Music lovers
- **3 unique program options** easily added in Program Handling
 - Live music performances
 - Performing music
 - Recorded music

For internal use only / © Resonance 2016. All rights reserved.

Comfort Features XPhone




NEW!

- **Scenario:** trouble hearing the phone clearly
- **Solution:** XPhone
- **Benefit:** Binaural amplification of the phone for a louder, clearer signal.
- Recommended use: patients that use a hand-held phone.
- Unique phone program easily added in Program Handling
- Activates and de-activates automatically
- Simultaneously sends the phone signal to both ears

For internal use only / © Resonance 2016. All rights reserved.

Comfort Features Using the XPhone Magnet



There is a specially-designed XPhone magnet.

To apply the magnet to enhance the phone's magnetic field:

1. Remove the backing from the XPhone magnet,
2. Place it on the top of your phone, near the speaker,

For internal use only / © Resonance 2016. All rights reserved.

Comfort Features T-Coil Connectivity



- A built in Telecoil provides connectivity to:
 - Telephones
 - Landlines
 - Telecoil compatible cell phones
 - Hearing loops (Induction systems)
- The Autophone feature makes this function automatic when used with telephones.
 - Telecoil is not automatically activated in memory 1. It requires a dedicated memory.

For internal use only / © Rexton 2016. All rights reserved.

Speech Centered Features Auto Directional iLock

PROVEN!

Speech Understanding



PROVEN!

- **Scenario:** trouble hearing in restaurant, family gatherings, social events
- **Solution:** Auto Directional iLock
- **Benefit:**
 - **Better speech understanding in extremely challenging and noisy environments**
- Advanced adaptive directional system capable of creating directional beams as narrow as 10 degrees
- Default microphone in the first program
- Narrow directionality is activated automatically based on acoustic environment. App or remote are not required.

For internal use only / © Rexton 2016. All rights reserved.

Speech Centered Features Auto iFocus 360 & iFocus 360

Speech
Understanding



- **Scenario:** trouble hearing speech behind you
- **Solution:** Auto/iFocus 360
- **Benefit:**
 - Enhances speech understanding in environments where turning to see the speaker is difficult.
- Recommended use: to hear speech from front, back, left or right.
- Automatic in a car in the first program
- May also be added as a unique program for:
 - Lawyers
 - Medical professionals
 - People in wheelchairs
 - Teachers
- True directionality from any direction
 - Directionality is engaged automatically or with the use of apps

For internal use only / © Rexton 2016. All rights reserved.

Speech Centered Features Automatic Noise Reduction

Speech
Understanding



PROVEN!

- **Scenario:** speech drowned out by background noise
- **Solution:** Automatic Noise Reduction
- **Benefit:**
 - Provides improved comfort and reduces listening fatigue
- Ideal for anyone bothered by background noise.
- Automatically provides the appropriate amount of noise reduction
- Automatically reduces 3 types of noise:
 - Steady state noise (motors, fans, etc.)
 - Impact noise (dishes, papers, keys)
 - Speech noise (group noise, restaurants)

For internal use only / © Rexton 2016. All rights reserved.

Speech Centered Features Wind Noise Cancellation


Speech
Understanding




- **Scenario:** annoying wind noise
- **Solution:** Wind Noise Cancellation
- **Benefit:**
 - Increased user comfort and speech understanding in wind noise
- Ideal for active patients that need to be able to listen to conversations when outside.
- Automatic in the first program
- Effectively cancels wind noise without affecting speech signal
- Engages in less than 1 second of onset of wind noise

For internal use only / © Resound 2016. All rights reserved.


Automatic Optimization



Smart Automatic Equalizer



Bandwidth Compression




Smart Automatic Acclimatization




Feedback Preventer

For internal use only / © Resound 2016. All rights reserved.

Automatic Optimization Smart Automatic Equalizer




- Customize the response of the aids in different daily environments
- Flexible and customizable for an individual's preferences.
- Fast, easy, and convenient adjustments to personalize the Automatic Program.



For internal use only / © Rexton 2016. All rights reserved.


Automatic Optimization Smart Automatic Acclimatization



- Ensures for comfortable and positive first fit acceptance
- Reduce follow-up fittings
- Intuitive and easy to set

For internal use only / © Rexton 2016. All rights reserved.


Automatic Optimization Bandwidth Compression



- Bandwidth Compression takes un-aidable high frequencies and moves those frequencies to areas where amplification may be beneficial.
- Critical Speech elements are preserved.

For internal use only / © Resonance 2016. All rights reserved.

Automatic Optimization Feedback Preventer




- Incorporates new technologies under the umbrella of TruCore
 - Continuous Critical Gain Measurements using the sounds in the environment
 - Accounts for changes in posture or objects around the hearing aid
- Measurements and cancellation occur in parallel

For internal use only / © Resonance 2016. All rights reserved.

Automatic Optimization


Adjusting Feedback Preventer in Fine Tuning



For internal use only / © Rexton 2016. All rights reserved.

Automatic Optimization

Performing Critical Gain Measurements



For internal use only / © Rexton 2016. All rights reserved.

SecureTec Protection Package IP67 Rated



IP rating is a certification issued by independent laboratories.

IP 6X

- Dust tight: no ingress of dust

IP X7


- Protected against the effect of immersion between 15 cm and 1 m

Proven receiver technology



For internal use only / © Ranton 2016. All rights reserved.

Widest Fitting Range



General guidelines:

- Use open dome if low frequencies are 20 dB or better
- Use semiopen dome if low frequencies are 25-40 dB
- Use closed dome if low frequencies are 45-65 dB
- Use double dome if low frequencies are 65 dB or worse

For internal use only / © Ranton 2016. All rights reserved.

Rexton 6^c– Portfolio Fall 2016





XS = 10
S = 312
M = 13
P = 13
HP = 675

PL 80	Sterling 80 6i	inoX CIC 80 6i	Emerald XS 80 6i	Emerald S 80 6i	Emerald M* 80 4i	Finesse 18.2i S	Mosaic M 80 6i		
PL 60	Sterling 60 6i	inoX CIC 60 6i	Emerald XS 60 6i	Emerald S 60 6i	Emerald M* 60 4i	Accord 16.2i S	Mosaic M 60 6i		
PL 40	Sterling 40 6i	inoX CIC 40 6i	Emerald XS 40 6i	Emerald S 40 6i	Emerald M 40 4i	Joy S 40 4i	Mosaic M 40 6i		
PL 30	Sterling 30 4i		Emerald XS 30 4i	Emerald S 30 4i		Joy S 30 4i	Joy M 30 4i	Joy P 30 4i	
PL 20	Sterling 20					Joy S 20	Joy M 20	Joy P 20	Joy HP 20
PL 10	Sterling 10	inoX CIC 10 inoX ITC 10				Joy S 10	Joy M 10	Joy P 10	Joy HP 10
PL 5		inoX CIC 5 inoX ITC 5				Targa S 5		Targa P 5	Targa HP 5
PL 5A								Targa P 5A	Targa HP 5A
PL 3								Arena P 3	Arena HP 3
PL 1								Arena P 1	
Accessories	Connex Smart Key		Connex Smart Transmitter		Connex Fit2Go App		Connex CROS RIC		
Software / Launch	Connex 7.5 / Rextit 7.5 Previous		Connex 8.2 / Rextit 8.2 Nov 2016						

Tinnitus
*No Tinnitus
 HP Power
Customs

For internal use only / © Rexton 2016. All rights reserved.

Emerald XS 6^c – The Ultra-discreet RIC Hearing System



For internal use only / © Rexton 2016. All rights reserved.

Mosaic M 6^c— The Reliable Discreet Companion



Sterling 6^c – Optimum Coverage For All Patients Needs

- From tiny CIC to full shell
- Directionality for CIC PL 80
- Power matrices available up to 65dB gain



inoX CIC 6^c – Simplifying ITEs




- Click sleeve solution: natural sound, without feedback
- Wireless Synch
- Single mic directionality (PL 80 only)
- Audio streaming (with Smart Connect accessory)
- Also available in ITC size for PL 10 & 5



For internal use only / © Resonance 2016. All rights reserved.

inoX CIC 6^c – Click Sleeves



Click Sleeve Closed



Click Sleeve Vented



For internal use only / © Resonance 2016. All rights reserved.

Accessories



**CONNEXX
CROS RIC**

For unilateral hearing loss
Permanent **wireless** transmission
Compatible with
Emerald S 6° and Mosaic
M 6°

SMART KEY

Small and handy
Users without
smartphone
Acoustic signal



SMART TRANSMITTER



Easier connection to TV

LED indicators
On/off and volume
buttons
Compatible with
Smart Connect and mini Blu
RCU

For internal use only / © Resdon 2016. All rights reserved.

Connexx 8 Using Connexx 8



For internal use only / © Resdon 2016. All rights reserved.

Connexx 8 What's New


- **Mimic Fit**
 - Easily transfer demo settings into TruCore hearing aids.
- **Program Link**
 - Link programs together for efficient programming
- **Decouple Volume Control**
 - Personalize volume control settings
- **Copy to Program 1**
 - Flexible programming
- **Interactive Curve Display**
 - Quick adjustments
- **QR Code for Smart Remote app**
 - Simplifies setup of Smart Remote app



For internal use only / © Rexton 2016. All rights reserved.

Connexx 8 Mimic Fit

- Mimic Fit: transfer fitting session from any TwinCore and QuadCore device into the TruCore hearing aids.
 - Transfers gain settings of memory 1, program labels, and configuration settings. Other fine tuning changes are not transferred. Programs 2-6 are realigned with memory 1.
- Mimic Fit also allows receivers to be changed without requiring a new First Fit.



For internal use only / © Rexton 2016. All rights reserved.



Connexx 8 Program Link



- Multiple programs can now be linked. Setting changes will be applied to all programs simultaneously
- Can be applied in Basic Tuning or Fine Tuning/Amplification

For internal use only / © Raston 2016. All rights reserved.


Connexx 8 Decouple Volume Control



- Allows you to adjust the volume control of each instrument independently
- This feature may be used for those suffering with Meniere's

For internal use only / © Raston 2016. All rights reserved.


Using Connexx 8 Copy to Program 1




Copy Program allows you to move one of the programs, from 2 to 6, into program 1.

For internal use only / © Rexton 2016. All rights reserved.

Connexx 8 Interactive Curve Display




- Gain and output curves can be directly adjusted from the curve window making fine tuning quick and easy.




For internal use only / © Rexton 2016. All rights reserved.


Connexx 8 Fitting Formulas



SmartFit New



SmartFit Experienced





SmartFit Power

SmartFit is based on NAL-NL2 and is optimized for the new features found in TruCore hearing aids

For internal use only / © Rexton 2016. All rights reserved.

Personalizing the Hearing Aids Configuring the Rocker Switch



For internal use only / © Rexton 2016. All rights reserved.

Compatible Accesories



Apps



Smart Connect



Smart Power Charger



Smart Remote



BT Transmitter

For internal use only / © Rexton 2016. All rights reserved.

Compatible Accesories Connexx Smart Remote App



- Change programs
- Adjust volume
- Mute
- Treble control
- External audio source selection
- Battery life of instruments when set-up with QR Code
- Microphone Pattern Adjustment
- No intermediary device needed
- Bluetooth not required
- Wireless not required
- Uses Ultra High Frequency Coded Sound



For internal use only / © Rexton 2016. All rights reserved.

Compatible Accesories

QR Code in Documentation Tab for Smart Remote App

Connexx Smart Remote app

Did you know you can control your hearing aids with the Smart Remote app?

First download the app, which is available for Android and IOS. After installation open the Smart Remote App and scan the QR code. Now you are ready to use the app!



- New QR Code in Documentation Tab
- QR Code is read by the Smart Remote App
 - Transfer program configuration to the app for a faster and accurate setup
 - Enable battery status signal

For internal use only / © Rexton 2016. All rights reserved.

Compatible Accesories

Smart Connect



- Uses Bluetooth technology
 - Connects to most Bluetooth enabled cell phones
 - Encoded signal to prevent interference
- **Minimal effect on hearing aid batteries**
 - Connects to audio devices via Bluetooth or additional transmitter
- Simple pairing to hearing instrument and Bluetooth devices
- Connects to all wireless Rexton hearing aids including **CIC**

For internal use only / © Rexton 2016. All rights reserved.

Compatible Accesories Connexx Smart Connect App



- Change programs
- Adjust volume
- Mute
- Treble control
- External audio source selection
- Battery life of instruments and Smart Connect
- Quickly and easily pair Connexx Transmitter and Speech Connect to the Smart Connect
- Microphone Pattern Adjustment
- Requires the Smart Connect
- **App is NOT required for Bluetooth phone and audio streaming functions**
- Will act as Smart Remote App if Smart Connect is not present



For internal use only / © Roston 2016. All rights reserved.

Thank you!

For internal use only / © Roston 2016. All rights reserved.