

## A Deep Dive into Target for Pediatric Hearing Aid Fittings

Stacy Michels, M.A.



A Sonova brand



**PHONAK**  
life is on

### Disclosure in accordance with ASHA CEU Requirements

- I have a financial relationship to disclose:
  - Employee of Phonak who receives a salary
- I have no nonfinancial relationship to disclose



**Stacy Michels, AuD**  
Phonak Pediatrics  
IA, IL, KS, MN, ND, NE, SD, WI

## Agenda



1. Target Set Up
2. Fitting Session
3. Resources

3

Target Software for Phonak Pediatric Fittings


**PHONAK**  
life is on

# 1

## Target 5.1 Set up

## Set up for success

- From home page, go to Setup at bottom of page
- Allows you to customize your fitting session
- Show and allow changes to defaults
- Shows specific defaults for Junior mode
- Allows you to customize reports




5

Target Software for Phonak Pediatric Fittings

**PHONAK**  
life is on

## SoundRecover defaults



6

Target Software for Phonak Pediatric Fittings

**PHONAK**  
life is on

## Fitting session tab

- Allows you to set up fitting session windows
- Most important tab is presentation
- Can check show recommendation for Junior Mode
- Will create pop up window to show you that Junior Mode is on when you enter age into Noah or Client Info

The screenshot shows the 'Fitting session' tab selected in the 'Setup' menu. The 'Presentation' sub-tab is active, displaying a list of checkboxes for various presentation options. The 'Fitting session' menu item is highlighted in the top navigation bar.

Setup
General
My profile
<b>Fitting session</b>
Junior mode
Reports
Fitting device
Sound

Presentation
<input checked="" type="checkbox"/> Show recommendation for junior mode
<input checked="" type="checkbox"/> Show absolute tuning values
<input checked="" type="checkbox"/> Show compression values in tuning tools
<input checked="" type="checkbox"/> Show preparation information in AudiogramDirect assistant
<input type="checkbox"/> Show message on startup of training mode
<input checked="" type="checkbox"/> Show warning information when entering SoundRecover extended mode
<input checked="" type="checkbox"/> Show warning information when leaving SoundRecover extended mode

7

Target Software for Phonak Pediatric Fittings

**PHONAK**  
 life is on

## Create clinic defaults in software

The screenshot shows the 'Create clinic defaults' window. The 'Junior mode' tab is selected in the top navigation bar. The 'Defaults' section is expanded, showing a table of settings for different age groups: 0 - 3 years, 4 - 8 years, 9 - 12 years, and 13 - 18 years. The 'Fitting formula' section is also visible, showing a list of programs and their startup status.

Defaults	0 - 3 years	4 - 8 years	9 - 12 years	13 - 18 years
▶ Fitting formula				
▼ Program manager				
Automatic programs	Disabled	Disabled	Enabled	Enabled
1 First priority	Roger/DAI ± mic	Roger/DAI ± mic	Roger/DAI ± mic	Roger/DAI ± mic
2 Second priority	Acoustic phone	Speech in noise	None	StereoZoom
3 Third priority	None	None	None	None
RogerReady	Roger/DAI ± mic	Roger/DAI ± mic	Roger/DAI ± mic	Roger/DAI ± mic
EasyPhone	Disabled	Acoustic phone	Acoustic phone	Acoustic phone
Mute program	Disabled	Disabled	Disabled	Disabled

8

Target Software for Phonak Pediatric Fittings

**PHONAK**  
 life is on

## New Phonak Target Junior mode defaults - DSL

Setting	0-3 years	4-8 years	9-12 years	13-18 years
Startup	Roger/DAI+M	Roger/DAI+M	Roger/DAI+M	Roger/DAI+M
AutoSense Sky OS	Disabled	Disabled	Enabled	Enabled
Roger and directional	Omni	Real Ear Sound	Fixed directional	Fixed directional
Indicator light	Enabled	Enabled	Disabled	Disabled
Push button	Disabled	Disabled	Enabled	Enabled
Volume control	Disabled	Disabled	Enabled	Enabled
RogerReady	Enabled	Enabled	Enabled	Enabled

9

Target Software for Phonak Pediatric Fittings

**PHONAK**  
life is on

## Customize your reports

- Make your reports look professional
- Add your own logo
- Use preprinted letter head
- Choose your curve type
  - gain or output

A great way to communicate with school or other care providers.

Setup

General My profile Fitting session Junior mode Reports

Layout

Curve type

Layout

☐ Standard layout  
☒ Customized layout  
☒ Use own logo  
 Select file...  
☐ Use own text  
☐ Preprinted letterhead

10

Target Software for Phonak Pediatric Fittings

**PHONAK**  
life is on

## Sending client information

Setup		
General	My profile	Fitting session
Language	Import fitting	Export fitting
Presentation	The current Phonak Target setup settings will be overwritten.	Export the current setup settings to a file.
Import / Export fitting	Import...	Export...
Support	ⓘ Excluded settings: Language, my profile, protected fitting and sound system.	
Support tools		

# 2

## Fitting Session

## Fitting session at a glance

### Client

- Client details
  - Age
  - Will activate Junior Mode

### Audiogram

- Keeps history of audios
- Familiar sounds for counseling
- SPL view


### RECD

- Saves history of RECD

Will auto populate from NOAH

13


Target Software for Phonak Pediatric Fittings

**PHONAK**  
life is on


14

Target Software for Phonak Pediatric Fittings

**PHONAK**  
life is on



## Fitting session at a glance, *cont.*

🕒 Sky V70-P

Gain 66 (2cc)  
MPO 131 (2cc)

→ Acoustic parameters

Sky V70-P

POWER LEVELS

Standard 66/131 (2cc) SlimTube 58/126 (2cc)

DESCRIPTION

Advanced WaterResistant BTE (13) for mild to severe hearing loss. Wireless.

FEATURES

- AutoSense Sky OS advanced
- Speech in loud noise
- Speech in 360°
- UltraZoom advanced
- DuoPhone
- FlexControl
- SoundRecover2
- SoundRelax
- Roger and directional setting
- RogerReady
- Indicator light

▶ More information


### Instruments

- Hearing Instruments
  - Power levels
  - Features
  - Can help you decide between models
- Acoustic Parameters
  - Choose coupling
  - Save color choices
- Accessories
  - Identify accessories
  - Reprogram ComPilot for another hearing aid.

Hearing Instrument will auto populate with "Connect"



## Fitting session at a glance, *cont.*



### Fitting - Where the action happens!

- Gain and MPO changes
- Verification Assistant
- Sound Recover 2
- Program Manager
- Device options

17


Target Software for Phonak Pediatric Fittings

**PHONAK**  
life is on

## Gain and MPO

Most adjustments made on this screen

Access to Verification Assistant



18

Target Software for Phonak Pediatric Fittings

**PHONAK**  
life is on

## Verification assistant

Use to ensure all adaptive features are disabled including SoundRecover2



19

Target Software for Phonak Pediatric Fittings

**PHONAK**  
life is on

## Software update for audioscan

New /s/ and /sh/ signals can be downloaded

<http://software.audioscan.com>

### CURRENT SOFTWARE

VeriFit2 4.6.2  
Axiom 1.12.2  
RM500SL 3.14.4  
VeriFit 3.14.4  
NOAH Module 2.4.1  
Remote App 1.2.2

[Update now](#)  
[Register to stay connected!](#)

©2016, all rights reserved. Audioscan, VeriFit, RM500SL, Axiom, and Speechmap are all trademarks of Audioscan, a Division of Etymonic Design Incorporated.



20


Target Software for Phonak Pediatric Fittings

**PHONAK**  
life is on

## Verification of gain and MPO

Gain and MPO:

- Match prescribed targets
  - Average
  - Soft
  - Loud
  - MPO




21

Target Software for Phonak Pediatric Fittings

**PHONAK**  
life is on

## SoundRecover2 - intelligent adaptive frequency lowering algorithm




22

Target Software for Phonak Pediatric Fittings

**PHONAK**  
life is on

## SoundRecover vs. SoundRecover2




23

Target Software for Phonak Pediatric Fittings

**PHONAK**  
life is on

## Best match to prescribed targets without SoundRecover2




24

Target Software for Phonak Pediatric Fittings

**PHONAK**  
life is on

## Determine candidacy using /s/ stimulus




25

Target Software for Phonak Pediatric Fittings

**PHONAK**  
life is on

## SoundRecover2 fitting tool

- Adjustments made from
  - Verification Assistant
  - [Fine tuning], [SoundRecover2]**
- Two interdependent sliders are available:
  - Audibility ----- Distinction
  - Clarity ----- Comfort
- A change to top slider automatically resets second slider to optimal clarity




26

Target Software for Phonak Pediatric Fittings

**PHONAK**  
life is on

## Audibility – Distinction slider adjustments




27

Target Software for Phonak Pediatric Fittings

**PHONAK**  
life is on

## /s/ and /sh/ separation




28

Target Software for Phonak Pediatric Fittings

**PHONAK**  
life is on

## Clarity – Comfort slider adjustments



29

Target Software for Phonak Pediatric Fittings

**PHONAK**  
life is on

## SoundRecover2 fitting protocol summary is available in print or online

- Match gain and MPO for prescribed targets
- Assess candidacy for SoundRecover2 using /s/ stimulus
- Use Audibility --- Distinction slider to achieve audibility of /s/ stimulus
- Adjust for comfort if required using Clarity --- Comfort slider
- Check for 1/3 Octave separation between /s/ and /sh/

30

Target Software for Phonak Pediatric Fittings

**PHONAK**  
life is on

## Program manager

- Enable RogerReady
  - Automatic activation of Roger
  - In AutoSense Sky OS or startup program
- See programs within AutoSense Sky OS
  - The first automatic program just for kids.
- Can change Startup program
- Add additional programs

<< Program manager

AUTOMATIC PROGRAMS	Position	Startup	Easy access	Label for client
AutoSense Sky OS	A ▾	<input type="radio"/>		Automatic
Calm situation	– ▾	<input type="radio"/>		Calm situation
Speech in noise	– ▾	<input type="radio"/>		Speech in noise
Comfort in noise	– ▾	<input type="radio"/>		Comfort in noise
Music	– ▾	<input type="radio"/>		Music
+				
1 Roger/DAI + mic	1 ▾	<input checked="" type="radio"/>	<input checked="" type="checkbox"/> RogerReady	Roger/DAI + mic
2 Speech in loud noise	2 ▾	<input type="radio"/>		Speech in loud noise
Acoustic phone	– ▾	<input type="radio"/>	<input checked="" type="checkbox"/> EasyPhone	Acoustic phone
+				
STREAMING PROGRAMS				
Bluetooth audio + mic				Bluetooth audio + ..
Bluetooth phone / DECT + mic				Bluetooth phone / D...

☐ Enable mute program in last position

31

Target Software for Phonak Pediatric Fittings

**PHONAK**  
life is on

## RogerReady



100%  
of the  
time!

The Roger/DAI + Mic program  
is always in the program structure and  
cannot be deleted by accident

Automatic activation of the Roger/DAI + Mic program

No additional programming required

Every child's Sky V hearing aids are  
ready for Roger as soon as they are.

32

Target Software for Phonak Pediatric Fittings

**PHONAK**  
life is on



## RogerReady - Automatic activation of Roger

Transmission from the Roger microphone activates the Roger program

Automatic switching into the Roger program occurs in <400ms

Switching out of Roger program occurs after ~50 seconds of no transmission

Available in AutoSense Sky OS and also in startup program

RogerReady now enabled by default for all Junior mode fittings in Phonak Target

100%  
of the  
time!

33



Target Software for Phonak Pediatric Fittings

**PHONAK**  
life is on

## Roger and directional setting

Found in program options

- Three microphone modes available at programming
  - **Omni** – always omnidirectional microphone
  - **Real Ear Sound** – always Real Ear Sound
  - **Fixed directional** – adaptive activation of fixed beamformer



34

Target Software for Phonak Pediatric Fittings

**PHONAK**  
life is on

## Roger and Directional - **Omnidirectional**

- ✓ 0-3 Age group
- ✓ Hear from all around for incidental listening
- ✓ No changes from Sky Quest
- ✓ Recommended for infants and those with poor head control



35

Target Software for Phonak Pediatric Fittings

**PHONAK**  
life is on

## Roger and Directional - **Real Ear Sound**

- ✓ 4-8 yr old age group
- ✓ Hear from all around for incidental listening
- ✓ Essentially omnidirectional with small amount of directionality
- ✓ Helps with localization
- ✓ Helps with attenuation of high frequency sounds



36

Target Software for Phonak Pediatric Fittings

**PHONAK**  
life is on

## Roger and Directional - Fixed Directional

- ✓ 9-18 yr old age group
- ✓ Fixed null point allows her to hear her friends from the side
- Real Ear Sound is the default in quiet
- ✓ Fixed directional is automatically activated based on two parameters:
  - ☐ ~50% Speech in Noise classification is detected
  - ☐ Ambient noise level exceeds a moderate level



37


Target Software for Phonak Pediatric Fittings

**PHONAK**  
life is on

## Roger and Directional



### Three microphone modes:




38

Target Software for Phonak Pediatric Fittings

**PHONAK**  
life is on

## Push button

- Junior Mode
  - 0-3 and 4-8 yrs - disabled
  - 9-12 and 13-18 yrs - enabled
- Short and long press options
- QuickSync automatically activated



SHORT PRESS  
 Program change ☒ QuickSync  
 LONG PRESS  
 No function


39

Target Software for Phonak Pediatric Fittings

**PHONAK**  
 life is on

## Volume control

- Allows for increase and decrease in volume
- Junior Mode
  - 0-3 yrs and 4-8 yrs disabled
  - 9-12 yrs and 13-18 yrs enabled
- QuickSync automatically activated



FUNCTION  
 Volume ☒ QuickSync  
 No function

40

Target Software for Phonak Pediatric Fittings

**PHONAK**  
 life is on

## Volume behavior


Allows for setting of maximum volume both upper and lower limits

Individual options for both louder and softer:

- None
- 6 dB (3 steps)
- 10 dB (5 steps)

### Mic attenuation

Push button needs to be activated for long press, mic attenuation




41

Target Software for Phonak Pediatric Fittings

**PHONAK**  
life is on

## Signals and alerts

- Allows for manipulation of whether or not to enable/disable the signals
- Intensity and frequency are defaulted to medium
- Choose those alerts that apply to the hearing aid settings
- The low battery beep alert is audible to the child




42

Target Software for Phonak Pediatric Fittings

**PHONAK**  
life is on

## Indicator light

- 5 indication patterns
  - can be selected by the fitter
- Easily demonstrated on the hearing instrument when connected to Target 5.0
- Pattern explanation available in Target report
- Starting up in Roger please select top two to indicate that the hearing aid is working.




43

Target Software for Phonak Pediatric Fittings

**PHONAK**  
life is on

## Start up

- Start up program change option
- Start up melody on/off
- Additional start up delay option




44

Target Software for Phonak Pediatric Fittings

**PHONAK**  
life is on

## Datalogging

- Allows for datalogging on/off option
- Datalogging information can then be found under Datalogging tab



45

Target Software for Phonak Pediatric Fittings

# 3

## Resources

## Reports

Target audience	Documents for clients	Selected documents
Professional	INSTRUCTIONS	CLIENT
Client	Hearing aid instructions ✓	Hearing aid instructions
	Program description ✓	Program description
	INFORMATION	Understanding hearing loss
	SoundRecover	Speech, language & listening (7-12 months)
	Communication tips	
	Putting on hearing aids	
	Diary of observations	
	Understanding hearing loss ✓	
	Care and maintenance	
	SPEECH, LANGUAGE & LISTENING	
	0-6 months	
	7-12 months ✓	
	1-2 years	

Documents setup

español (Spanish) ▼

☐ Use large fonts

Preview Print Email Save

- For parent or professional
- 26 different languages
- Good reference for use and care
- Offers counselling opportunities
- Good way to communicate with school



47

Target Software for Phonak Pediatric Fittings

**PHONAK**  
life is on

## Did you notice?

Reports can be accessed in two different places



48

Target Software for Phonak Pediatric Fittings

**PHONAK**  
life is on



## Target 5.1 Demonstrator

Introduction

Conversation in restaurant

Phone call

LISTENING PERSON

Client  
Hearing loss: Mild

Edit

Significant other

Select



0:00 0:23

Hearing instruments None Basic Premium

Clients & sessions

Demonstrator

Media


News

Tools

Updates

Setup

## Pediatric fitting resources



### Desktop Fitting Guide

This guide provides you with a detailed instruction to fitting a Phonak CROS II system in Phonak Target 4.1 and later. For the fitting of CROS II it is recommended to use Phonak Target 4.1 and later. For the fitting of CROS II it is recommended to use Phonak Target 4.1 and later. For the fitting of CROS II it is recommended to use Phonak Target 4.1 and later.

### Best Practice Protocols

#### Electroacoustic verification of SoundRecover using the Audioscan VeriFit®

1. Create a SoundRecover in the hearing aid using the fitting software.
2. Increase the gain and gain of the hearing aid fitting using the fitting software.
3. Connect the hearing aid to the fitting software.
4. Connect the hearing aid to the fitting software.
5. Connect the hearing aid to the fitting software.
6. Connect the hearing aid to the fitting software.
7. Connect the hearing aid to the fitting software.
8. Connect the hearing aid to the fitting software.
9. Connect the hearing aid to the fitting software.
10. Connect the hearing aid to the fitting software.



### Pediatric Fitting Guide

1. Connect the hearing aid to the fitting software.
2. Connect the hearing aid to the fitting software.
3. Connect the hearing aid to the fitting software.
4. Connect the hearing aid to the fitting software.
5. Connect the hearing aid to the fitting software.
6. Connect the hearing aid to the fitting software.
7. Connect the hearing aid to the fitting software.
8. Connect the hearing aid to the fitting software.
9. Connect the hearing aid to the fitting software.
10. Connect the hearing aid to the fitting software.

**Together,  
we change lives**