



Sept. 28th, 2017
Audiology Online

BrainHearing and bone anchored devices: Why more output matters to all users

Marcus Holmberg, PhD
Senior Director, Clinical Affairs BAHS

Because sound matters

oticon MEDICAL

INIUM Sense Oticon BrainHearing™ Technology

Learning Objectives

- After the course, participants will be able to explain why listening is more effortful for those with hearing loss.
- After the course, participants will be able to describe a few different ways of measuring the effort used to understand speech, and why this is relevant for real life situations.
- After the course, participants will be able to discuss the limitations of maximum output of bone anchored devices, and how this affects different users of these devices

Because
sound matters

oticon
MEDICAL

This course is presented in partnership with

continued

BrainHearing

The reason we develop the products we do



Because
sound matters

oticon
MEDICAL

Outline

○ Part I

- Hearing and cognition
- Real world listening situations
- How we can rethink our measurement methods to include cognition
- Results using bone anchored devices and novel outcome measures including cognition

○ Part II

- Why maximum output of bone anchored devices affects all users
- Should this influence how we counsel patients?

○ Questions & comments

Because
sound matters

oticon
MEDICAL

This course is presented in partnership with

Part I: Overview

- Hearing and cognition
 - Why is listening more effortful for those with hearing loss?
 - Does a hearing loss give rise to long-term issues?
- Real-world daily life (ecological) listening situations
 - In which situations do hearing impaired typically communicate?
- Cognitive measures beyond speech performance testing
 - Effect of technology on cognition
 - Advanced hearing aid technology

Because
sound matters

oticon
MEDICAL

An example of processing

- You will be presented with two tasks to perform, then asked to compare your performance on each.
- Practice:

Red Purple ...

Because
sound matters

oticon
MEDICAL

This course is presented in partnership with

continued

Task 1: Read the word that is written

Red Green Black Green
Green Red Purple Blue

Because
sound matters

oticon
MEDICAL



Task 2: Read the color of the text

Red Green Black Green
Green Red Purple Blue

Because
sound matters


oticon
MEDICAL

This course is presented in partnership with



This course is presented in partnership with

The Listening Effort Framework




The Ease of Language Understanding -ELU- model. Adapted from Rönnberg et al., 2013

Because
sound matters

oticon
MEDICAL

The Listening Effort Framework



The Ease of Language Understanding -ELU- model. Adapted from Rönnberg et al., 2013


Because
sound matters

oticon
MEDICAL

This course is presented in partnership with

continued™

The Listening Effort Framework



Part I: Overview

- Hearing and cognition
 - Why is listening more effortful for those with hearing loss?
 - Does a hearing loss give rise to long-term issues?
- Real-world daily life (ecological) listening situations
 - In which situations do hearing impaired typically communicate?
- Cognitive measures beyond speech performance testing
 - Effect of technology on cognition
 - Advanced hearing aid technology

Because
sound matters

oticon
MEDICAL

This course is presented in partnership with


continued

Hearing Loss & Cognitive Load

Poorer hearing is associated with:


- A. Reduced language-driven activity in primary auditory pathways
- B. Increased compensatory language-driven activity in pre-frontal cortical areas

A Decreased language-driven speech activity in poorer hearers



Peele et al., J Neurosci, 2011

B




Grossman et al., Brain Lang, 2002

oticon
MEDICAL

Because
sound matters

Rhyme task (visual stimuli)
Normal hearing participants

Rhyme task (visual stimuli)
Deaf participants



Less Brain Activity

More Brain Activity


Adapted from Lazard & Giraud. Nature Communications, 2017.

This course is presented in partnership with

Rönnberg et al.: *Hearing Loss and LTM*

Journal of Speech, Language, and Hearing Research • Vol. 54 • 705-726 • April 2011


Figure 1. Visualized structural equation modeling (SEM) for the hearing losses of the better ear (BE; see upper panel) and worse ear (WE; see lower panel) and their relation to episodic long-term memory (ELTM) performance when taking age into account. SPTs = subject-performed tasks; Rec = recall.



Because
sound matters

oticon
MEDICAL

Conductive losses affect innervation already in the cochlea



Because
sound matters

Liberman et al., PlosOne, 2015

oticon
MEDICAL

This course is presented in partnership with

Summing up so far...

- Why is listening more effortful for those with hearing loss?
 - Less automatization compared to normal-hearing persons
 - More activation of working memory.
- Does the hearing loss give problems with memory and other long-term effects?
 - Yes, at least for untreated hearing loss
 - Work ongoing on the effect of different treatments

Because
sound matters

oticon
MEDICAL

Part I: Overview

- Hearing and cognition
 - Why is listening more effortful for those with hearing loss?
 - Does a hearing loss give raise to long-term issues?
- Real-world daily life (ecological) listening situations
 - In which situations do hearing impaired typically communicate?
- Cognitive measures beyond speech performance testing
 - Effect of technology on cognition
 - Advanced hearing aid technology

Because
sound matters

oticon
MEDICAL

This course is presented in partnership with

continued

Ecological communication situations

- In which situations are the hard of hearing people communicating in their daily life?
 - At which SNRs?
 - How much of the information is understood?


Work by Widex ORCA and Jade University of Applied Sciences, Oldenburg



Be
SOI

oticon
MEDICAL

People who are hard of hearingcommunication situations in daily life




Because
sound matters
Smeds et al. (2015)

oticon
MEDICAL


This course is presented in partnership with

People who are hard of hearingsound level variation in daily life



Subjective assessment of speech intelligibility under real life conditions


- Subjective intelligibility does not vary, mostly close to full intelligibility



This course is presented in partnership with

Subjective assessment of listening effort under real life conditions


- Subjective listening effort varies substantially



Summing up so far...

Ecological conditions (here):

- **Positive SNRs** (SNR + 5 dB or better)
- **Subjective effort varies**
- **Overall sound levels varies a lot**
- **Successful perception** (close to 100% correct performance).
- ... effort may vary although having high speech intelligibility.



McMaster University


Because sound matters

oticon MEDICAL


This course is presented in partnership with

continued

Speech in noise testing...





Speech in noise testing...



This course is presented in partnership with

...a new axis with ~100% speech intelligibility



Because
sound matters

Inspired by Pichora-Fuller (2006)

oticon
MEDICAL

Part I: Overview

- Hearing and cognition
 - Why is listening more effortful for those with hearing loss?
 - Does a hearing loss give rise to long-term issues?
- Real-world daily life (ecological) listening situations
 - In which situations do hearing impaired typically communicate?
- Cognitive measures beyond speech performance testing
 - Effect of technology on cognition
 - Advanced hearing aid technology

Because
sound matters

oticon
MEDICAL


This course is presented in partnership with

continued

Some examples of cognitive measures

Brain Activity

- Neuroimaging: fMRI/PET
- EEG/ MEG
- Optical Imaging fNIRS



Behavioral Response

- **Memory**
- Response Time
- Dual tasks

Physiological Correlates


- **Pupil dilation**
- Skin conductance
- Stress hormones

Because
sound matters

Adapted from: Peelle, 2017, manuscript in preparation

oticon
MEDICAL

Example: Pupillometry Algorithms for reduced listening effort



Because
sound matters

oticon
MEDICAL

This course is presented in partnership with

Example: SWIR test

Sentence – final Word Identification and Recall

... at an SNR which allows 95-100% intelligibility
(typically +5 ... +15 dB)

Task: **listen to speech in noise**

1. Her går alle med solbriller

Ecological ✓

Because
sound matters

oticon
MEDICAL

Method, SWIRL test

Task: **repeat** last word

1. Her går alle med solbriller

Because
sound matters

oticon
MEDICAL

This course is presented in partnership with

Method, SWIRL test

Task: **repeat** last word (of 7 sentences)

1. Her går alle med
2. Han ligger stadig i
3. Eleven skriver en lang
4. Hele byen kom til
5. Hans datter vil på
6. I går havde filmen
7. Fabrikkens port var ikke

?

Task: **recall** words

Because
sound matters



100% intelligibility
but << 100% recall

"Erm ...
lukket ...
premiere ...
sengen ...
sol..briller ..
... umm ..
hm!"

= Speech recall

oticon
MEDICAL

Method, SWIRL test

Task: **repeat** last word (of 7 sentences)

1. Her går alle med solbriller
2. Han ligger stadig i sengen
3. Eleven skriver en lang rapport
4. Hele byen kom til brylluppet
5. Hans datter vil på højskole
6. I går havde filmen premiere
7. Fabrikkens port var ikke lukket

Task: **recall** words

Because
sound matters



oticon
MEDICAL

This course is presented in partnership with

Experiment:

Sound quality of bone conduction devices

○ Softband needs higher drive force above 1 kHz


- => earlier into non-linear movement of transducer => sound quality/audibility issues
- => effect on memory?



Abutment



Softband



Because
sound matters

oticon
MEDICAL

Experiment:

Sound quality of bone conduction devices

- 16 Ponto Plus Power patients (10 women)
- Conductive and mixed hearing losses

Procedure:


○ Prepare Experiment

- Adjust SNR to ~95 % correct HINT sentences for Softband in ISTS background noise
- => average SNR 10 dB (SD = 4.8 dB)

○ SWIRL test

Contrast:

Abutment versus softband, 10 repetitions, randomized order



Because
sound matters


oticon
MEDICAL

This course is presented in partnership with

Experiment:

Sound quality of bone conduction devices

Recall is significantly better for abutment compared to softband fitting




Because sound matters

Lunner et al. (2016) Using speech recall in hearing aid fitting and outcome evaluation under ecological test conditions. Ear & Hearing, vol 37, supplement 1, 145S-154S

oticon
MEDICAL

Direct drive & bandwidth: Learning effects in children

Children learn new words
3 times quicker with extended bandwidth^{4*}




Because sound matters

^{4*}The study included both normal hearing and hearing impaired children who were exposed to limited and extended bandwidth. No conclusion can be drawn from this study with regards to Ponto 3.

oticon
MEDICAL

This course is presented in partnership with

Direct drive & bandwidth: Learning effects in children



Because
sound matters

**The study included both normal hearing and hearing impaired children who were exposed to limited and extended bandwidth. No conclusion can be drawn from this study with regards to Ponto 3.*

oticon
MEDICAL

This guides what products we bring to the market




Because
sound matters

oticon
MEDICAL

This course is presented in partnership with

continued

BrainHearing guides what products we bring to the market



Because
sound matters

Reinfeldt et al., 2015

oticon
MEDICAL

Ponto 3 – The world's most powerful family of abutment-level processors



Because
sound matters

oticon
MEDICAL

This course is presented in partnership with

continued

Part II: Outline



- Why all users benefit from higher maximum output
 - Technical limitations
 - Maximum output, gain, dynamic range, distortions – how does it all relate?
 - What will patients likely report back?
- Should this influence how we counsel patients?

Because
sound matters

oticon
MEDICAL

Part II: Outline




- Why all users benefit from higher maximum output
 - Technical limitations
 - Maximum output, gain, dynamic range, distortions – how does it all relate?
 - What will patients likely report back?
- Should this influence how we counsel patients?

Because
sound matters

oticon
MEDICAL

This course is presented in partnership with


The users are our starting point



Because
sound matters

oticon
MEDICAL

Users' challenges



Because
sound matters


oticon
MEDICAL

This course is presented in partnership with

continued

Challenge in BAHS

Maximum output is low



Device	Measured MPO	Reference	Manufacturer
Sophono Alpha 1	56 dB HL	Hol et al., 2013	Sophono, Boulder, US
Bonebridge	65 dB HL	Mertens et al, 2014	Med-El, Innsbruck, Austria
Standard Baha Divino/BP100	67-69 dB HL	Carlsson & Hakansson, 1997 Zwartenkot et al. 2014	Cochlear BAS, Goteborg, Sweden
Standard Ponto	67-69 dB HL	Zwartenkot et al. 2014	Oticon Medical, Askim, Sweden
Baha Cordelle, Baha 5 SuperPower*, Ponto 3 SuperPower*	80-85 dB HL	Idem	Cochlear BAS, Goteborg, Sweden Oticon Medical, Askim, Sweden

Objective measurement of the MPO of several hearing devices¹
(*updated spring 2017)

Because
sound matters

1) <http://www.snikimplants.nl>

oticon
MEDICAL

Part II: Outline


- Why all users benefit from higher maximum output
 - Technical limitations
 - Maximum output, gain, dynamic range, distortions – how does it all relate?
 - What will patients likely report back?
- Should this influence how we counsel patients?

Because
sound matters

oticon
MEDICAL

This course is presented in partnership with

Sound environments are dynamic




Because
sound matters

Wagener et al. (2008)

oticon
MEDICAL

Bone anchored users benefit from higher max. output




Because
sound matters


oticon
MEDICAL

This course is presented in partnership with

Bone anchored users benefit from higher max. output




Bone anchored users benefit from higher max. output




This course is presented in partnership with

Bone anchored users benefit from higher max. output




Bone anchored users benefit from higher max. output



This course is presented in partnership with


Max. output in devices with 45 dB HL fitting range and Ponto 3 SuperPower



Because sound matters

oticon
MEDICAL

Bone anchored users benefit from higher max. output




Because sound matters

oticon
MEDICAL


This course is presented in partnership with

continued

Bone anchored users benefit from higher max. output




Bone anchored users benefit from higher max. output




This course is presented in partnership with

The higher gain, the lower saturation level




The higher gain, the lower saturation level




This course is presented in partnership with

The higher gain, the lower saturation level



The higher gain needed, the lower saturation level



This course is presented in partnership with

Part II: Outline



- Why all users benefit from higher maximum output
 - Technical limitations
 - Maximum output, gain, dynamic range, distortions – how does it all relate?
 - What will patients likely report back?
- Should this influence how we counsel patients?

Because
sound matters

oticon
MEDICAL

Everyone deserves the best Sound Quality

Bone anchored users benefit from higher maximum output



Patients with

- Conductive hearing losses
- Mixed hearing losses
- SSD
- Softband fitting

Because
sound matters

oticon
MEDICAL

This course is presented in partnership with

continued

Ponto 3 SuperPower

Optimal audiological solution for all patients

- Higher MFO -> Increased dynamic range / headroom in the device
- **Patients with conductive or mixed hearing losses**
 - The patient's auditive dynamic range are better utilized
 - More sounds are reproduced naturally without being limited
- **Patients using softband and head band**
 - Higher gain is provided to adress skin attenuation -> smaller dynamic range. SP enlarges the dynamic range.
- **SSD patients**
 - Better ability to loudness match sounds from the device to the normal hearing ear



Ponto 3 SuperPower

Because
sound matters

oticon
MEDICAL

Will I be overfitting
the patients?

Associations to 'SuperPower' terminology




Because
sound matters

oticon
MEDICAL


This course is presented in partnership with

Using a SuperPower device, will it sound too loud/powerful?

- 1) BAHS's **MFO** is below patients' LDL
 - So no risk of "over-fitting" by using a BAHS SuperPower device



Warning when fitting AC power HA



Does not apply to BAHS

- 2) **Gain** is prescribed to HL and can be prescribed low enough
 - Patient can maximum turn up the volume control by 10 dB
 - So no risk of "over-amplification"

Because
sound matters

oticon
MEDICAL

Using a SuperPower device, what about annoyance to sounds

- Patients report same degree of aversiveness to sounds when fitted with Ponto 3 SuperPower compared to Ponto Pro Power (AV parameter in APHAB-questionnaire).¹
- Loudness discomfort levels were estimated (extrapolated) to be 82.8, 92.1 and 89.3 dB HL at 1, 2, and 4 kHz, respectively, for bone conduction thresholds ≤ 40 dB(HL).¹



Because
sound matters

1) Bosman et al (2017) On Maximum Force Output in Bone-Anchored Devices. Poster at Osseo2017 conference.

oticon
MEDICAL

This course is presented in partnership with

Will I be overfitting
the patients?

Associations to 'SuperPower' terminology

Patients' reactions:
do patients tell that it sounds
louder when fitted with
SuperPower?

Patients' reactions
are related
to their previous experience

Because
sound matters


oticon
MEDICAL

Patients' who are used to devices with lower MFO

Example: MFO as in Ponto 3



- Device is in saturation
- ↓
- Distortion
- ↓
- Sounds loud



Ponto 3



Because
sound matters

34989-00

oticon
MEDICAL


This course is presented in partnership with

Patients' who are used to devices with lower MFO

Example: now fitted with P3SP



- Less often in saturation
- Less distortion
- Sounds quieter, more calm



Ponto 3 SuperPower



Because
sound matters

34989-00


oticon
MEDICAL

Immediate reactions to the sound in Ponto 3 SP

Nijmegen study – Ponto 3 SP

- subjects
 - BAHs experiences: 5-year replacements
 - Ponto Pro Power (N=19), Baha Intenso (N=1)
 - unilateral (N=14), bilateral (N=6)
- inclusion criteria
 - compos mentis
 - age < 80 years (22-80 yrs; avg: 63.0 yrs)
 - substantial air-bone gap
 - BC-thresholds > 30 dB HL
- follow-up (N=17)
 - tests too complex/tiring (N=1), brain infarct (N=1)
 - appointment p

Pure-tone audiometry (N=20)



Some comments ...

- geluid is rustiger* calm, even sound
- my hearing is better* less startled by sudden loud sounds
- Ik krijg complimenten van mensen* gets compliments on his hearing
- Ik krijg complimenten van mensen* more natural sound
- geluid is natuurlijker* understands better in noise
- Versta beter in rumoerige omgeving* hears rear sounds better
- hoor beter wat en achter mij gebeurt* no reverberance in car
- more calm, even sound* more natural sound
- a. Geen galmend geluid in auto.....bij vorige wel.
b. Geeft over het algemeen een rustigere hoervervaring.
c. Alles klinkt natuurlijker...dit in vergelijking met het blikkende geluid van vorige toestel
d. Opvallend de andere draairichting van het volume-wiel.
e. Bij volume verhogen geen galm.....
f. Batterijen hebben opvallend korte levensduur.....????

Because
sound matters

34989-00

Bosman et al. (2017) Ponto 3 SuperPower test results. Radboud University, Nijmegen. Presented at Osseo2017

oticon
MEDICAL

This course is presented in partnership with

Part II: Outline




- Why all users benefit from higher maximum output
 - Technical limitations
 - Maximum output, gain, dynamic range, distortions – how does it all relate?
 - What will patients likely report back?
- Should this influence how we counsel patients?

Because
sound matters

oticon
MEDICAL

Choosing a sound processor



Patients don't ask for a SuperPower device



It is our joint task to make sure
the patients make an informed choice

Because
sound matters

34989-00

oticon
MEDICAL

This course is presented in partnership with

continued

Choosing a sound processor

Test BAHS in dynamic listening environments

Tangible

- Cosmetic
- Handling
- Wireless
- ...



Not tangible

- Sound quality
- Speech understanding
- Reliability
- ...

Let the patient experience the sound quality

Because
sound matters

34989-00

oticon
MEDICAL

Test BAHS in dynamic listening environments

A guide to help collecting experiences

- A test guide & note for the patient
 - to fill in during test at the hospital or at home.
 - it guides him/her to go to dynamic listening situations (e.g hospital café, outdoor)

Because
sound matters

oticon
MEDICAL

This course is presented in partnership with

Summary

- Fitting a product with higher max output ->
 - less distortion of sounds
 - the dynamic range of listening environments being naturally reproduced is enlarged
- Ponto 3 SuperPower has the highest MFO in any abutment level device, though its' MFO is below patient's Loudness-Discomfort-Level, therefore

All bone anchored users benefit from higher maximum output!
- Higher max output leads to better sound quality which has a positive effect on e.g. listening to speech in noisy environments and music.
- Let patients test Ponto 3 SuperPower in dynamic environments.
- Our aim is to take the effort out of listening supporting BrainHearing.

Because
sound matters

oticon
MEDICAL



Brain Hearing™ for Oticon Medical

Technology for Hearing Implants to reduce Listening Effort



Improved speech recognition

In all listening situations.
Technology that makes sense.



Reduced Listening Effort

More understanding with less effort.
Technology for intense living.



Increased recall abilities

Empowered cognition for enriched
Social Interactions and faster learning.

Because
sound matters

oticon
MEDICAL



This course is presented in partnership with