



2017

# BAHS Introduction

Barinder Samra BSc, MSc, CS (Aud)

Because sound matters

oticon  
MEDICAL

## Agenda

- BAHS Principles
- BAHS Candidates
  - *Advantages of BAHS device*
  - *Contraindications for surgery*
- Bilateral Conductive and Mixed Hearing Losses
  - *Which Side to Fit?*
- Ponto 3 Super Power
  - *The Importance of Increased MFO*
- Methods of Sound Transmission
  - *Skin Drive and Direct Drive*

Because  
sound matters

oticon  
MEDICAL

## Learning Goals

- To be able to calculate BAHS candidacy using the audiogram
- Identify the benefits of treating conductive/mixed and unilateral profound hearing losses using bone conduction devices.
- Define contraindications for BAHS surgery
- Ascertain whether a patient is a candidate for a bilateral or unilateral fitting using the audiogram.
- State the impact of low MFO on sound quality and explain why increasing the MFO of BAHS devices benefits all users.
- Indicate the affect different methods of sound transmission can have on output and outcome.

Because  
sound matters


**oticon**  
MEDICAL

## BAHS Principles




Because  
sound matters

**oticon**  
MEDICAL



## Osseointegration and Direct Bone Conduction



Osseointegration




Abutment (6, 9, 12, or 14 mm long)  
Implant (3 or 4 mm long)

Because  
sound matters

**oticon**  
MEDICAL


## Osseointegration



Because  
sound matters

**oticon**  
MEDICAL

## Direct Bone Conduction




Sounds are converted to vibrations, which are transmitted directly to both inner ears.

Because  
sound matters

**oticon**  
MEDICAL

## Ponto- The Bone Anchored System



Ponto (Latin)  
= Bridge

Cochlear owns the right to the name BAHA®;  
therefore we never write or call our system a BAHA  
(Bone-Anchored-Hearing-Aid).

Because  
sound matters

**oticon**  
MEDICAL

## Agenda

- BAHS Principles
- BAHS Candidates
  - *Advantages of BAHS device*
  - *Contraindications for surgery*
- Bilateral Conductive and Mixed Hearing Losses
  - *Which Side to Fit?*
- Ponto 3 Super Power
  - *The Importance of Increased MFO*
- Methods of Sound Transmission
  - *Skin Drive and Direct Drive*

Because  
sound matters

**oticon**  
MEDICAL

## BAHS Candidates



Because  
sound matters

**oticon**  
MEDICAL

## BAHS Candidacy Categories

- Conductive & mixed hearing losses
- Single-sided deafness (SSD)
- Other




Because  
sound matters

oticon  
MEDICAL

## Conductive and Mixed HL Candidacy

- Can the device compensate for hearing loss in the inner ear?
  - Average BC thresholds now up to 65 dB HL
- Is a bone anchored device a better solution than a traditional hearing aid?
  - If ABG > 30 dB, speech recognition with BAHS will likely be better than with conventional hearing aids.<sup>1, 2</sup>

Rule of  
Thumb



Because  
sound matters

1) Mylanus et al. Intraindividual comparison of the bone-anchored hearing aid and air-conduction hearing aids. *Otolaryngology-Head & Neck surgery*. 1998; 124(3): 271-6  
2) De Wolf et al. Better performance with bone-anchored hearing aid than acoustic devices in patients with severe air-bone gap. *The Laryngoscope* 121:613-616, 2011.


oticon  
MEDICAL

## Calculating Candidacy



### ○ Average BC threshold

- Calculate average BC response at .5 kHz, 1 kHz, 2 kHz & 3 kHz
- Average should be up to 65 dB HL



Average BC is less than 65 dB HL

Because  
sound matters


oticon  
MEDICAL

## Calculating Candidacy



### ○ Average air-to-bone gap

- Calculate average between AC and BC responses at .5 kHz, 1 kHz, 2 kHz & 4 kHz
- Average should be up to >30 dB HL



Av. ABG is larger than 30 dB

Because  
sound matters

oticon  
MEDICAL




## Conductive and Mixed Hearing Losses

Advantages of a BAHS device

Because  
sound matters

oticon  
MEDICAL

## Bone Anchored Sound Processor




The Ponto doesn't need to compensate for the air-to-bone gap (the conductive component of the hearing loss).

Because  
sound matters


oticon  
MEDICAL

## Traditional HA Versus Bone Anchored HA



### Traditional HA

- High gain
- Likely problems with feedback
- Tight earmould



### Bone anchored HA

- Low gain
- Better sound quality
- Ear canal is open

Because  
sound matters

**oticon**  
MEDICAL

## Example Candidates

- Chronic otitis media

Occluding earmoulds are contraindicated

Chronic ear discharge can lead to sensorineural loss<sup>3</sup>

BAHS helps to reduce discharge<sup>4</sup>

Reduction in infection = Fewer outpatient visits<sup>5</sup>



3) Papp, Z., et al. (2003). "Sensorineural Hearing Loss in Chronic Otitis Media." *Otology & Neurotology* 24(2): 341-344.

4) McDermott, A. L., et al. (2002). "An intra-individual comparison of the previous conventional hearing aid with the bone-anchored hearing aid: The Nijmegen group questionnaire." *J Laryngol. Otol Suppl* 28: 15-19.

5) Hol, M. S., et al. (2004). "The bone-anchored hearing aid: Quality-of-life assessment." *Archives of Otolaryngology-Head & Neck Surgery* 130(4): 394-399.

Because  
sound matters

**oticon**  
MEDICAL

## Example Candidates

- Otitis media w/ effusion
- Chronis otitis media
- Aural atresia/Microtia



Air-bone gap of 30 dB or less

Fewer follow up visits and fewer complications with BAHs when compared to surgery<sup>6</sup>

57% fulfilled this criteria after 1 year post op<sup>6</sup>

patients can still need amplification post surgery<sup>7</sup>



Because  
sound matters

6) Nadaraja, G. S., et al. (2013). "Hearing Outcomes of Atresia Surgery Versus Osseointegrated Bone Conduction Device in Patients With Congenital Aural Atresia: A Systematic Review." *Otology & Neurotology* 34(8): 1394-1399.

7) Evans, A. K. and K. Kazahaya (2007). "Canal atresia: "surgery or implantable hearing devices? The expert's question is revisited." *Int J Pediatr Otorhinolaryngol* 71(3): 367-374.

**oticon**  
MEDICAL

## Example Candidates

- Otitis media w/ effusion
- Chronis otitis media
- Aural atresia/Microtia
- Otosclerosis
- Traumatic injury to middle ear structures
- Ossicular disease
- Cholesteotoma


Because  
sound matters

**oticon**  
MEDICAL

## Case Study Patient BV

Long history of  
multiple ear  
infections

BAHS fitted in (L)  
2004 (R) 2005



Since BAHS fitted ear  
infections reduced, less  
visits to clinic

Because  
sound matters

**oticon**  
MEDICAL

## Candidacy Categories


- Conductive & mixed hearing losses
- Single-sided deafness (SSD)
- Other



Because  
sound matters

**oticon**  
MEDICAL

## Single-Sided Deafness (SSD) (Profound Unilateral Sensorineural HL)




Because  
sound matters

oticon  
MEDICAL

## Single Sided Deafness (SSD) Candidacy

- Does the better ear have normal or near-normal hearing?
  - Average air conduction threshold < 20 dB HL



- Or the patient chooses not to use an air conduction CROS system


Because  
sound matters

oticon  
MEDICAL

## Candidates for SSD

### ○ Single-sided deafness caused by:


- Acoustic neuroma
- Sudden deafness
- Congenital causes
- Cholesteatoma
- Ototoxic drugs
- Meniere's disease



Because  
sound matters


oticon  
MEDICAL

## CROS Versus Bone Anchored HA



### CROS

- Ear canal in good ear partially occluded by device
- Must wear two devices
- Poor acceptance<sup>8</sup>
- Can have high battery drain



### Bone anchored HA

- Ear canal is open
- Only one device
- Normal battery drain

Because  
sound matters


8.

Niparko, J. K., et al. (2003). "Comparison of the bone anchored hearing aid implantable hearing device with contralateral routing of offside signal amplification in the rehabilitation of unilateral deafness." *Otol Neurotol* 24(1): 73-78.

oticon  
MEDICAL

## Advantages of Treating SSD

High scores in subjective patient tests<sup>9</sup>



Better or equal speech understanding in noise compared to CROS aid<sup>8,9</sup>

8) Niparko, J. K., et al. (2003). "Comparison of the bone anchored hearing aid implantable hearing device with contralateral routing of offside signal amplification in the rehabilitation of unilateral deafness." *Otol Neurotol* **24**(1): 73-78.  
9) Kim, G., et al. (2017). "Efficacy of Bone-Anchored Hearing Aids in Single-Sided Deafness: A Systematic Review." *Otol Neurotol* **38**(4): 473-483.

Because  
sound matters

**oticon**  
MEDICAL

## Candidacy Categories

- Conductive & mixed hearing losses
- Single-sided deafness (SSD)
- Other




Because  
sound matters

**oticon**  
MEDICAL

## Other Medical Indications


- Maybe not a direct candidate, audiotically
- Physical constraints for fitting
  - Underlying skin condition
  - Inflammation
  - Reaction to earmoulds
  - Stenosis



Because  
sound matters

**oticon**  
MEDICAL

## Case Study




Because  
sound matters

**oticon**  
MEDICAL




## Candidacy - In Summary

Conductive/Mixed




Because  
sound matters

Single Sided Deafness



Medical Indications



oticon  
MEDICAL

## A Patient May Not be a Candidate for an BAHS Implant if:

- They have **poor hygiene** and are unable to care for their abutment site
- Their **bone quality** will not support the implant, for example, it is too thin, or it is brittle due to radiation.
- They have a **skin condition** that will reduce their ability to heal following surgery.
- They are under **5 years of age** (US, Canada and Singapore).

Because  
sound matters

oticon  
MEDICAL

## Materials

- Candidacy Guide
- Candidacy Quick Guide



Because  
sound matters

**oticon**  
MEDICAL

## Agenda

- BAHS Principles
- BAHS Candidates
  - *Advantages of BAHS device*
  - *Contraindications for surgery*
- Bilateral Conductive and Mixed Hearing Losses
  - *Bilateral or Unilateral Fitting- Making a Decision*
- Ponto 3 Super Power
  - *The Importance of Increased MFO*
- Methods of Sound Transmission
  - *Skin Drive and Direct Drive*

Because  
sound matters

**oticon**  
MEDICAL

## Bilateral Conductive and Mixed Hearing Losses

### *Bilateral or Unilateral Fitting - Making a Decision*






Because sound matters

**oticon**  
MEDICAL

## Bilateral Fitting- Making a Decision



**Bilateral  
Conductive/Mixed HL**

↓

**Look at  
Masked BC**

↓

**Symmetrical BC**

Less than 10  
dB HL on  
average

And less  
than 15 dB  
HL at  
individual  
frequencies


↓

**Fit Bilaterally= True  
Binaural Effect**

Because sound matters

**oticon**  
MEDICAL

## Bilateral Fitting



## Bilateral Fitting- Benefits

Skull provides enough attenuation to acoustically isolate BAHs output

Improvement in Speech Recognition in Noise

Binaural Loudness Summation

Utilise Binaural Squelch

Improved sound localisation

Because sound matters

oticon MEDICAL

## Unilateral Fitting-Making a Decision




Because  
sound matters

oticon  
MEDICAL

## Unilateral Fitting-Making a Decision?

Bilateral  
Conductive/Mixed HL



Asymmetrical BC


Fit BAHS to the side  
with best BC

If BAHS is fitted on  
poor BC side, then it  
will work like a CROS  
device

Because  
sound matters

oticon  
MEDICAL


## Unilateral Fitting-Making a Decision?



Because  
sound matters

oticon  
MEDICAL

## Unilateral Fitting-Making a Decision



Because  
sound matters

oticon  
MEDICAL

## Agenda

- BAHS Principles
- BAHS Candidates
  - *Advantages of BAHS device*
  - *Contraindications for surgery*
- Bilateral Conductive and Mixed Hearing Losses
  - *Which Side to Fit?*
- **Ponto 3 Super Power**
  - *The Importance of Increased MFO*
- Methods of Sound Transmission
  - *Skin Drive and Direct Drive*

Because  
sound matters

**oticon**  
MEDICAL



**Ponto 3 SuperPower**  
**-The Importance of Increased MFO**

Because  
sound matters

**oticon**  
MEDICAL

INTUM  
Sense  
Oticon BrainHearing™  
Technology


## Ponto 3 – The Definition of Power



Because  
sound matters

**oticon**  
MEDICAL

## Ponto 3 – BrainHearing at the Core



Because  
sound matters

**oticon**  
MEDICAL



## The Recipe for Ponto Sound Quality

Output –  
more loud  
and soft  
sounds



Bandwidth -  
greater range of  
high and low  
tones



Clarity -  
without  
distortion




These are all at core of our development and the fundamentals  
for delivering Ponto Sound Quality.

Because  
sound matters

<https://www.youtube.com/watch?v=kTMEydoxdpw>

**oticon**  
MEDICAL

## Challenge in BAHS Maximum Output is Low



Device	Measured MPO	Reference	Manufacturer
Sophono Alpha 1	56 dB HL	Hol et al., 2013	Sophono, Boulder, US
Bonebridge	65 dB HL	Mertens et al, 2014	Med-El, Innsbruck, Austria
Standard Baha Divino/BP100	67-69 dB HL	Carlsson & Hakansson, 1997 Zwartenkot et al. 2014	Cochlear BAS, Goteborg, Sweden
Standard Ponto	67-69 dB HL	Zwartenkot et al. 2014	Oticon Medical, Askim, Sweden
Baha Cordelle, Baha 5 SuperPower*, Ponto 3 SuperPower*.	80-85 dB HL	Idem	Cochlear BAS, Goteborg, Sweden Oticon Medical, Askim, Sweden

Objective measurement of the MPO of several hearing devices<sup>1</sup>  
(\*updated spring 2017)

Because  
sound matters

10) <http://www.snikimplants.nl>

**oticon**  
MEDICAL

## Challenge in BAHS

### Maximum Output is Low

- Transmitting sound via vibrations through the skull bone is not as efficient as transmitting sound via the ear canal.
- As such the maximum force output (MFO) of BAHS devices is much lower than the maximum power output (MPO) of hearing aids.
- The MFO of a BAHS device is much lower than uncomfortable loudness levels.

Because  
sound matters

**oticon**  
MEDICAL

## Ponto 3 – SuperPower made beautiful

### The world's first single-unit SuperPower

- Ponto 3 SuperPower combines a discreet-, small design and a battery that lasts.
- Abutment-level processor - no strings or need for any bulky ear- or body-worn devices



- Combine discretion with usability
- Follows the contours of the ear
- Left and right versions
- Available in six colours


Because  
sound matters

**oticon**  
MEDICAL

## Ponto 3 SuperPower

*Committed to Providing a Higher MFO*

Maximum Force Output level at 90 dB SPL input



Because  
sound matters

**oticon**  
MEDICAL

## Sound Environments are Dynamic




35190-00

Because  
sound matters


Wagener et al. (2008)

**oticon**  
MEDICAL


## Bone anchored users benefit from higher max. output




## Bone anchored users benefit from higher max. output




## Bone anchored users benefit from higher max. output




## Bone anchored users benefit from higher max. output




## Bone anchored users benefit from higher max. output



## Bone anchored users benefit from higher max. output



## Bone anchored users benefit from higher max. output



## Ponto 3 SuperPower optimal solution for all patient groups

- **Patients with conductive loss**
  - Higher MFO means better utilization of patients dynamic range
  - More natural sound in louder listening environments
- **Patients with mixed hearing losses**
  - Higher MFO gives larger dynamic range / headroom in the device, so more sounds are reproduced naturally without being limited
  - Higher gain needs excellent feedback management
- **Softband and head band users**
  - Higher MFO to address skin attenuation
- **SSD patients**
  - Better ability to loudness match sounds from the device to the normal hearing ear



Ponto 3 SuperPower

Because sound matters

**oticon**  
MEDICAL

## Agenda

- BAHS Principles
- BAHS Candidates
  - *Advantages of BAHS device*
  - *Contraindications for surgery*
- Bilateral Conductive and Mixed Hearing Losses
  - *Which Side to Fit?*
- *Ponto 3 Super Power*
  - *The Importance of Increased MFO*
- **Methods of Sound Transmission**
  - *Skin Drive and Direct Drive*

Because  
sound matters

**oticon**  
MEDICAL

## Methods of Sound Transmission

*-Skin Drive and Direct Drive*




Because  
sound matters

**oticon**  
MEDICAL



## Methods of Sound Transmission





Because  
sound matters

11) Briggs (2015) Clinical Performance of a New Magnetic Bone Conduction Hearing Implant System: Results from a Prospective, Multicenter, Clinical Investigation. [Otol Neurotol](#), 2015 Jun;36(5):834-41.

**oticon**  
MEDICAL

## Sounds are Attenuated in Skin Drive Solutions

- Physical facts with Skin Drive solutions:
  - There is 10-20 dB sound attenuation in mid to high frequency region<sup>2</sup>
  - These devices have lower perceived output in the mid to high frequencies



Ponto



Magnetic/Softband




Because  
sound matters

12) Verstraeten et al (2008) Comparison of the audiologic results obtained with the bone anchored hearing aid attached to the headband, the testband and to the 'snap' abutment. *Otolary & Neurology* 30: 70-75

**oticon**  
MEDICAL

## Skin Drive vs Direct Drive

**13%** improvement  
of words  
correctly recalled



Because  
sound matters

13) Lunner (2016) Using speech recall in hearing aid fitting and outcome evaluation under ecological test conditions. Ear & Hearing, vol 37, supplement 1, 145S-154S

**oticon**  
MEDICAL

Any Questions?

Because  
sound matters

**oticon**  
MEDICAL