

FEATURE AND PRODUCT QUICK GUIDE





HEARING **HAPPENS IN REAL** LIFE

ACHIEVING A PERFECT FIT IN A CLINIC CAN BE **DIFFICULT BECAUSE...**

- ... the wearers' environment is constantly changing
- ... wearers seek different levels of amplification depending on the circumstances
- ... wearers have different intentions. For example, they might intend to hear a conversation in one situation and tune out in others



CONSEQUENTLY. SOME WEARERS...

- ... return the hearing devices or seldom use them
- ... become frustrated with inadequate hearing
- ... come back to the clinic and ask for help, not always articulating their needs clearly.

WIDEX **EVOKETM EVOLVES** IN REAL LIFE

The new WIDEX EVOKE™, based on an advanced, dual-core chip, delivers intelligent automation as well as intuitive user control. It leverages industry-first machine learning technology to create the best real-life listening experience.

WIDEX EVOKE DISTRIBUTED COMPUTING SYSTEM









EVOKE PLATFORM

Widex Dual-Core Processor

Two core chips to provide the ideal balance of processing efficiency and flexibility.

CONNECTIVITY

Pure-link 2.4 GHz

Smartphone integration via Bluetooth expands processing capabilities.

MACHINE LEARNING

Smartphone app


EVOKE™ uses the computational power of the smartphone to operate machine learning features such as SoundSense Learn. WIDEX CLOUD

Cloud

The Widex platform spreads across many processing cores, with the potential to link even more cores into a network of extraordinary processing power to help all wearers.

MEETING THE WEARER WHERE THEY ARE

The WIDEX EVOKE™ offers additional and smarter sound class technology, new specialized programs, new fitting profiles for a more precise fit, and 30% more processing power to meet the needs of wearers more naturally than ever before.



KEY FEATURES	440	330	220	110	BENEFIT
Processing channels	15	12	10	6	For accurate processing and a more precise sound experience and speech intelligibility.
Fine-tuning channels	15	12	10	6	More options for personalized sound to cater to wearer preferences and hearing needs.
WidexLink wireless connectivity	•	•	•	•	Fast, echo-free transmission of sound between hearing devices for binaural coordination and to DEX accessories for streaming a variety of signals.
SoundSense Learn	•	•	•		Machine learning algorithm to effortlessly guide the wearer's to a better, more personalized hearing experience.
SoundSense Adapt	•	•	•		Adaptive learning to gradually remember the wearers preferred settings across multiple parameters for each sound class environment, resulting in better sound.
Fluid Sound Controller	•	•	•	•	Adjustment of more parameters in every sound class and in real time for more natural sound.
Fluid Sound Analyzer (Sound Classes)	11 (IE)	7 (IE)	4	3	Intuitive analysis of the sound environment in order to prioritize and select the most important sound sources.
Programs	5	4	3	3	Wide range of listening programs for specific and challenging situations.
Smartwind Manager	•				Keeps speech and other important sounds audible while significantly reducing the wind annoyance.
High-Frequency Boost	•				Clear and crisp sound in high frequencies, ideal for music.
Speech Enhancer RT	RT/IE	IE			Reduces noise and enhances speech in real time in response to individual sound situations and hearing loss.
Digital Pinna	•	•			Helps locate and focus naturally on the voice in front or behind by mimicking the normal pinna effect.
HD Locator	•	•	•		Better speech intelligibility in noise with channel-specific directionality for speech emphasis.
TruSound Softener	•	•	•		Preserves comfort for sudden sounds at any input level by taking the edge off impact/impulse sounds.
Soft-Level Noise Reduction	•	•	•	•	Unwanted soft background noises are reduced for greater listening comfort.
Noise Reduction	•	•	•	•	Listening comfort in noise while maintaining speech.
ZEN IE	•	•	•	•	A tinnitus management tool for reducing the effects of tinnitus using soothing tones and/or noise; additional programs with Zen+.
Audibility Extender	•	•	•	•	Improves perception and production of important high-frequency speech sounds, ideal for more significant high frequency losses.
Variable Speed Compression	•	•	•	•	Automatic gearing for smooth sound and stable audibility in each individual sound class.
EVOKE App	•	•	•		Easy connection and control from most smartphones.
TONELINK App	•	•	•	•	Easy connection and control from most smartphones.

THE WIDEX EVOKETM FAMILY

10 RIC







S-receiver

Yes

312

170 hrs

No

Optional/Via ToneLink or DEX

Standard/Via ToneLink or DEX

Yes

No

Yes

ΑII

63 dB

118 dB SPL

100 Hz - 10000 Hz

IP-68



M-receiver

Yes

312

170 hrs

No

Optional/Via ToneLink or DEX

Standard/Via ToneLink or DEX

Yes

No

Yes

ΑII

68 dB

124 dB SPL

100 Hz - 9200 Hz

IP-68

FUSION

312 RIC/RITE



P-receiver

No

312

160 hrs

No

Optional/Via ToneLink or DEX

Standard/Via ToneLink or DEX

Yes

No

Yes

ΑII

70 dB

131 dB SPL

100 Hz - 6700 Hz

IP-68



HP-receiver No

312

155 hrs

No

Optional/Via ToneLink or DEX

Standard/Via ToneLink or DEX

Yes

No

Yes

ΑII

84 dB

137 dB SPL

100 Hz - 6000 Hz


IP-68


PASSION STYLE


Receiver	S-receiver	M-receiver	
Open fit option	Yes	Yes	
Battery size	10	10	
Battery life	100 hrs	100 hrs	
Zpower rechargeable solution	No	No	
Preference control	Via ToneLink or DEX	Via ToneLink or DEX	
Program shift	Via ToneLink or DEX	Via ToneLink or DEX	
Telecoil	Via DEX	Via DEX	
2.4GHz mfi streaming	No	No	
Widex CROS compatible	Yes	Yes	
DEX devices*	All	All	
Max gain	63 dB	69 dB	
MPO	118 dB SPL	124 db SPL	
Bandwidth	100 Hz – 10000 Hz	100 Hz – 10000 Hz	
IP-rating	IP-68	IP-68	

^{*} COM-DEX, COM-DEX Remote Mic, CALL-DEX, UNI-DEX, RC-DEX. TV-DEX, FM+DEX, PHONE-DEX


IEC Standard
























312 BTE MINI





FASHION POWER

-	
9	
	- The second



FUSION2	
312 RIC/RITE 2	4


S-receiver	M-receiver	P-receiver	HP-receiver	
Yes	Yes	No	No	
312	312	312	312	
165 hrs	165 hrs	160 hrs	155 hrs	
Yes	Yes	Yes	Yes	
Optional/Via App or DEX				
Standard/Via App or DEX				
Yes	Yes	Yes	Yes	
Yes	Yes	Yes	Yes	
No	No	No	No	
All	All	All	All	
62 dB	69 dB	70 dB	84 dB	
118 dB SPL	124 dB SPL	132 dB SPL	138 dB SPL	
100 Hz - 10000 Hz	100 Hz - 10000 Hz	100 Hz - 6750 Hz	100 Hz - 7300 Hz	
IP-68	IP-68	IP-68	IP-68	

-
Yes
312
155 hrs
No
Optional/Via ToneLink or DEX
Standard/Via ToneLink or DEX
Via DEX
No
Yes
All
64 dB
135 dB SPL
100 Hz - 8450 Hz
IP-68


312 BTE
-
Yes
312
160 hrs
No
Standard/Via ToneLink or DE
Standard/Via ToneLink or DEX
Yes
No
Yes
All
71 dB
137 dB SPL
100 Hz - 7600 Hz
IP-68
·


13 BTE	
-	
No	
13	
285 hrs	
No	
Standard/Via ToneLink or DE	ΞX
Standard/Via ToneLink or DE	ΞX
Yes	
No	
Yes	
All	
77 dB	
138 dB SPL	
100 Hz - 6250 Hz	
IP-68	



















**Additional colors available for F2, FS, PA & FA: Metallic Blue (090), Silver White (112), Copper Brown (070), Shocking Pink (077), Lime Green (075), Mediterranean Turquoise (078)



























Tan















Blue 136



Sporty Red 076







Pearl White 068

Grey 119

Grey 121

Silk 081

Beige 123

Summer Gold 069

Capp. Brown 071

Tech Black 118

THE WIDEX EVOKETM FAMILY











	CIC-M
LE	10 CIC MIC

STYLE	10 CIC MICRO		
Receiver			
Open fit option	No		
Battery size	10		
Battery life	105 hrs		
ZPower rechargeable solution	No		
Preference control	-		
Program shift	-		
Telecoil	-		
Widex CROS compatible	-		
DEX devices*	-		
Max gain	60 dB		
MPO	118 dB SPL		
Bandwidth	100 Hz - 10000 Hz		
IP-rating	IP-68		

CIC 10 CIC

No

10

100 hrs

No Via ToneLink or DEX

Via ToneLink or DEX

Via DEX

Yes

ΑII

63 dB 124 dB SPL

100 Hz - 9700 Hz

IP-68

312 ITC 312 Half Shell

	312 Full Shell, Swing door		
	-		
	No		
	312		
	165 hrs		
	No		
-	Optional/Via ToneLink or DEX		
-	Optional/Via ToneLink or DEX		
	Via DEX		
	Yes		
	All		
	63 dB		
	123 dB SPL		


312 ITC, 312 Half Shell, 312 Full Shell, Swing door	312 ITC, 312 Half Shell, 312 Full Shell, Flip door	
-	-	
No	No	
312	312	
125 hrs	140 hrs	
No	No	
Optional/Via ToneLink or DEX	Via ToneLink or DEX	
Optional/Via ToneLink or DEX	Via ToneLink or DEX	
Via DEX	Yes	
Yes	Yes	
All	All	
67 dB	68 dB	
127 dB SPL	128 dB SPL	
100 Hz - 10000 Hz	100 Hz - 9800 Hz	





No 312 140 hrs neLink or DEX neLink or DEX Yes Yes ΑII 68 dB

* COM-DEX, COM-DEX Remote Mic, CALL-DEX, UNI-DEX, RC-DEX. TV-DEX, FM+DEX, PHONE-DEX

Custom models available in Summer 2018









100 Hz - 10000 Hz

IP-68



IP-68



IP-68

CIC/ITE models













Dark brown

COMPLETE DEX CONNECTIVITY

With WIDEX EVOKE™, connectivity doesn't stop with the apps. Our range of DEX communication solutions keeps wearers connected to a world of sounds.

TV-DFX

FEATURES

- · High-quality audio in stereo
- Ultra-low delay Echo-Free™
- Inputs for both TV and hi-fi systems
- Room Off feature
- 10 hours of non-stop capacity
- Rechargeable using the TV Base



RC-DEX

FEATURES

- Program toggle
- Volume up/down
- Keyring attachment
- Approximately 12 months of battery life
- Lock switch



COM-DEX

FEATURES

- True hands-free communication
- Compatible with most mobile phones
- Streams audio from most Bluetooth devices
- Available in three stylish colors
- 8 hours of streaming time
- 8 days of standby/control time*
- · Comfortable fabric neck loop
- Controls your hearing devices via COM-DEX app



UNI-DFX

FEATURES

- Connects to all devices with
 3.5 mm jack output
- Streams automatically
- Built-in microphone for handsfree talking
- Battery life 40 hours when streaming
- · Battery recharges in 1 hour
- · Room Off feature



FM+DEX

FEATURES

- A flexible device for FM, Audio and telecoil streaming
- Small and easy to handle
- Antenna range up to 30 meters
- 10 hours of battery capacity



COM-DEX REMOTE MIC

FEATURES

- Connects wirelessly via COM-DEX to all Widex wireless hearing devices
- · Convenient steel clip for easy attaching
- Mute/unmute microphone directly on the COM-DEX Remote Mic button
- Microphone can also be controlled using the COM-DEX app
- Up to 8 hours of streaming time



PHONE-DEX

FEATURES

- Range to base up to 900 ft.
- Phonebook for 50 numbers
- 100 hours of standby time and
 10 hours call time on battery
- Can also be used as a regular phone
 phone for the whole family
- a phone for the whole family




CALL-DEX

FEATURES

- Connects to most mobile phones with 3.5 mm jack output
- · Streams automatically
- No buttons
- Constant streaming up to 80 hours on size 10 battery
- Standby up to 3 weeks
- · Small and easy to handle



^{*12} hours per day



EVOKETM APP

Provides the ultimate flexibility for patients to control their hearing experience based on the environment and listening intent.

Wearers can benefit from binaural iPhone streaming through both speech and music templates for superior streaming sound quality.

Available Styles:

EVOKE FUSION2 312 RIC 2.4 Style



Features:

- SmartSense Learn Real-Time Machine Learning A/B profile comparisons
- Ability to create Personal Programs and select program name and picture
- · Program Control
- Right/Left Volume (Preference) Control
- Easy Mute Access
- Directional Focus to change from front, right, left, and back focal points.
- Equalizer for low, mid, and high modifications
- Sound Mixer to maximize streaming to environment preferences
- Find my Hearing Aids
- Geotagging
- · Works with most smartphones

TONELINK™ APP

Turns a smartphone into an EVOKE remote control to discreetly change programs, adjust volume, and change directional focus.

Available Styles:

All EVOKE styles except CIC-M



Features:

- Program Control
- · Volume (Preference) Control
- Easy Mute Access
- Directional Focus to change from front, right, left, and back focal points.
- · Works with most smartphones

Note: ToneLink app communicates through a high frequency signal from the phone speaker to matched EVOKE hearing devices. The phone cannot be paired to Bluetooth audio devices while using the TONELINK app.