



Closing the Gap: Mixed and Conductive Hearing Loss in Adults

Amy Isaacs Donaldson, AuD, CCC/A
Professional Education Manager


Hear now. And always.



Agenda



- Conductive and mixed hearing loss in children
- Treating conductive and mixed hearing loss in children with Baha® Technology
- Exploring Baha solutions for children



Conductive and Mixed Hearing Loss (CMHL) in Children



- Prevalence of conductive hearing loss in children may be up to 3.5% of the general population¹
- 29% of children with Down Syndrome have permanent mixed or conductive hearing loss; another 32% have fluctuating losses²
- 25% of children with cleft palate have conductive hearing loss of 25 dBHL or more³



1 - Feder KP, Michaud D, McNamee J, Fitzpatrick E, Ramage-Morris P, Beauregard Y. (2017) Prevalence of hearing loss among a representative sample of Canadian children and adolescents, 3 to 19 years of age. Ear & Hearing, 38(1):7-26.
2 - Nightingale E, Yoon P, Wolter-Warmuth K, Daniels D, Holey F (2017) Understanding hearing and hearing loss in children with Down Syndrome. American Journal of Audiology, 26:201-208.
3 - Gaddy B, Lott O, Canady J, Smith RJ (2008) Conductive hearing loss and otitis media in cleft palate patients. Otolaryngol Head Neck Surg. 134(6):946-8.

Congenital Aural Atresia & Microtia



- Prevalence of Microtia is 1.76 per 10,000 live births in North America¹
 - higher incidences in Hispanic populations (2.89) and Asian or Pacific Islander populations (2.08)¹
- Degree of hearing restoration by surgical correction is dependent on severity²
- Majority of pediatric patients who have reconstruction will still require some form of amplification after surgery³
- About 75% of patients who do not receive good hearing benefit from surgery also do not benefit from a revision²



1 – Lucetti DV, Leonelli E, Mastroiacovo P. (2011) Microtia-Anotia: A global review of prevalence rates. Birth Defects Res A Clin Mol Teratol. 91(8):813-822.
2 – Chang SO, Choi BY, Hur DG. (2008) Analysis of the long-term hearing results after the surgical repair of aural atresia. Laryngoscope. 118(10):1855-1861.
3 – Evans AC, Kozlowski K. (2007) Canal atresia: Surgery or implantable hearing device? The expert's question is evasive. International Journal of Pediatric Otorhinolaryngology. 71(3):367-374.

Treating CMHL Hearing Loss in Children



Intermittent/Fluctuating Hearing Loss



Permanent Hearing Loss



Why Consider Baha® Technology?



- ✓ Nothing worn in the ear canal, which is helpful in cases where ear drainage is present¹
- ✓ Baha does not need to be re-adjusted if the air conduction thresholds fluctuate because gain is only required for the bone conduction thresholds
- ✓ Pre-operative testing can be used to predict post-operative benefit²
- ✓ Use of Baha is associated with high user satisfaction and good long-term benefit^{3,4}



1 – Buchholz S, Arnold P, Saldaña L. (2010) Congenital aural atresia: bone-anchored hearing aid vs. external auditory canal reconstruction. Int J Pediatr Otorhinolaryngol. 74(2):270-7.
2 – Morris S, Filipp C, Almon P, Blagov M, Lucchini M, Barlow M. (2013) Individualized headband simulation test for predicting outcome after percutaneous bone conduction implantation. Acta Otolaryngol. 133:258-64.
3 – Giedd S, Andersen K, Faber CE, Wiersma J. (2011) Bone-anchored hearing aids are effective and associated with a high degree of satisfaction. Danish Medical Journal. 118(7):1010.
4 – Rasmussen J, Olsen S, Sørensen L. (2012) Evaluation of long-term patient satisfaction and experience with the Baha bone conduction implant. Int J Audiol. 51(3):164-8.

3

Baha Systems for Children: Surgical



Baha Attract

Baha Connect

- The Baha Attract provides clinically proven performance with nothing through the skin¹
- The Baha Connect's direct connection provides the maximum possible gain²
- Two systems with one strong foundation³

1. Briggs R, Van Hasselt A, Lutz M, Goyosolea M, Wignen S, Weber P, Smeets H, Flynn M, Cowan R. Clinical performance of a new magnetic bone conduction hearing implant system results from a prospective, multicenter, clinical investigation. *Choi Neurotol*. 2015;36(5):654-61.
2. Clinicaltrials.gov ID NCT01786238. Clinical and Health Economic Evaluation With a New Baha® Abutment Combined With a Minimally Invasive Surgical Technique.
3. Nalasek PC, Sathyan J, de Wolff MJ, Flynn MC, Wignen S, Day-Gordon M, Green K, Richman MP, Mylchreest SA, & Ho NW. Long-term stability, survival, and tolerability of a novel osseointegrated implant for bone conduction hearing: 5-year data from a multicenter, randomized, controlled, clinical investigation. *Choi Neurotol*. 2014; 35(8): 1486-91.

Baha 5 Sound Processors: Small. Smart. Powerful.



Baha 5 sound processors share the same unique technology building blocks to deliver a smart and seamless hearing experience to all patients.

- The industry's smallest sound processor.¹
- The only smart processors with direct-to-device wireless streaming and control.
- Two power sound processors including the industry's most powerful solution.²




1. Flynn MC. Smart and Small – Innovative technologies behind the Cochlear Baha 5 Sound Processor. Cochlear Bone Anchored Solutions AB. 629781, 2015.
2. Janssen, J. Review of Baha 5 series. Cochlear Bone Anchored Solutions AB. D773528, 2015.



Baha 5 Sound Processors for Children



- Made for iPhone compatibility to allow connection to iOS devices
- Baha 5 Smart App for iOS and Android™ allows parents to control and monitor their child's device
- Full FM compatibility with the amazing Mini Microphone 2+
- Various retention options and colors allow recipients to personalize their device for their needs
- Datalogging in all processors



Support and Resources




- Request free information packs for yourself or your candidates at our website
- Candidates may contact our innovative Concierge group for answers to their questions and to connect with a local volunteer if desired
- Parents and educators can find resources for rehabilitation online at our Communication Corner

www.cochlear.com/us

1-877-897-4474 or
concierge@cochlear.com

www.cochlear.com/us/Communication-Corner

Summary



- Mixed and conductive hearing loss in childhood can be common in patients with Down Syndrome, cleft palate, chronic otitis media and other congenital abnormalities
- Baha technology provides several options to treat both permanent and fluctuating hearing loss in children
- Baha 5 Sound Processors offer True Wireless and Made for iPhone technology to help children connect to the world around them
- Cochlear offers unparalleled support for candidates and recipients



Hear now. And always

©Cochlear Limited 2018. All rights reserved. Hear now, And always and other trademarks and registered trademarks are the property of Cochlear Limited or Cochlear Group. All other trademarks and registered trademarks are the property of their respective owners.

Android, Google Play and the Google Play logo are trademarks of Google Inc.

Apple, the Apple logo, Made for iPhone logo, iPhone, iPad Pro, iPad Air, iPad mini, iPad and iPad touch are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries.

Made for iPhone Apple Inc. Use Made for iPhone Hearing aids (Hearnow). Apple support. 2017 (last 24 February 2017). Available from: <https://support.apple.com/en-us/HT201446>. (Subject to change)

The Bluetooth® word mark and logo are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Cochlear is under license.
