

GSI AMTAS: Pro vs AMTAS Flex Operation and Integration

August 28, 2018

Laura Prigge, AuD



1



Learning Outcomes

- Describe the differences between tests available in AMTAS Flex and AMTAS Pro.
- Explain how to transfer to both the AMTAS Pro and Flex into GSI Suite.
- Define 2 scenarios for AMTAS Pro and 2 scenarios for AMTAS Flex.



2



Agenda

- Review of AMTAS
- Introduction of AMTAS Pro and AMTAS Flex
- Compare and Contrast AMTAS Pro and AMTAS Flex Operation
- Compare and Contrast AMTAS Pro and AMTAS Flex Integration
- Summary and Wrap up
- Discussion



3



Significant Audiology Changes Predicted

The Crisis Coming in Audiology. Barry Freeman, Audiology Today
Nov/Dec 2009

“A shortfall in the number of audiologists is predicted to occur in the not-so-distant future. If we do not address this impending crisis, any gaps in care will be filled by alternatives outside the profession.”



4



Incidence of Hearing Loss

- World Health Organization, March 2015
 - 360 million people
 - 328 million Adults and 32 million children
 - 1/3 of people over 65 have hearing loss
- US Department of Health and Human Services
 - About 2 to 3 out of every 1,000 children
 - 30 million aged 12 years or older has hearing loss



5



Flat Profession vs Growing Need

“Demand for Audiology Services: 30 Year Projection and Impact on Academic Programs” I. Windmill & B. Freeman, JAA 2013; 24, 407-413

- About 16,000 licensed audiologists
 - 20 % do not impact/participate in patient care.
 - Of the remaining 12,800, only 11,200 practice full time
 - About 600 Audiologist graduate each year (estimated that 40% never enter the field)
 - 400 retire each year
- Solution for demand
 - Number of audiologists must be increased
 - Increase productivity and efficiency.



6



Other Factors influencing change in Audiology

- Telehealth
- PSAPs
- Over the Counter Hearing Aids
- Demands for efficiency
- Hearing aid delivery models changing
- Big box stores selling hearing aids



7



Addressing the changing climate in Audiology.



8

amtas™



9



GSI AMTAS - Overview

- Automated **M**ethod for Testing Auditory Sensitivity
- Software program – works with GSI audiometers or Microsoft Tablets
 - **Software** is loaded on a PC or tablet
 - Connected to a GSI Audiometer
 - AudioStar Pro
 - Pello
 - Connected to a tablet
- Self administered automated test for obtaining a diagnostic or screening audiogram
- Patented algorithms ensure quality and reliability of evaluation.



10



AMTAS Pro and AMTAS Flex



11

amtas^{pro}

- Clinical Use:
 - Masked air and bone conduction thresholds
 - Masked SRT and WRS speech audiometry
- Required:
 - GSI Audiometer (AudioStar Pro or Pello with AMTAS License)
 - PC with AMTAS Software
- Optional:
 - Touch screen computer/monitor



amtas^{pro}


- Telehealth:
 - Masked air and bone conduction thresholds
 - Masked SRT and WRS speech audiometry
- Required:
 - GSI Audiometer (AudioStar Pro or Pello with AMTAS License)
 - PC with AMTAS Software and telehealth software
 - Telehealth camera with otoscope attachment
 - Quiet Testing Area



13

amtas^{flex}

- Screening Mode:
 - Pure tone screening – 4 frequencies
 - Adjustable level (limited)
- Required:
 - Windows Tablet with AMTAS Software
 - Calibrated Headphones
 - Quiet room



amtasflex™

- Threshold Mode:
 - Pure tone air conduction thresholds
 - Masking
- Required:
 - Windows Tablet with AMTAS Software
 - Calibrated Headphones
 - Quiet room



Key Differences between AMTAS Pro and Flex

AMTAS PRO

- Must be used with an audiometer
 - AudioStar Pro
 - Pello
- Diagnoses Hearing Loss
 - Site of lesion
 - Configuration
- Air/Bone/Speech

AMTAS Flex

- Used with a tablet
- Portable
- Air Conduction ONLY
 - Screening
 - Threshold
- Does not Diagnose Hearing Loss





Amtas Pro and Flex

Operation




17



Amtas General information


- Release version – 1.3
- Security options
 - **Windows Domain Integration**—No changes
 - **Single Password**—One admin account. Not connected to Windows
 - **No Security**—No passwords required; no restrictions
- Config App



18




AMTAS Pro Config App



19



AMTAS Flex Config App



20



Menu Bar: Flex Export



Data can be imported into GSI Suite 2.5.


Selecting Export from Menu bar performs the Batch export.



21




Reports: Flex and Pro



22




General Options: Flex and Pro



23




General Options-Flex and Pro



24




Workflow: Flex and Pro



25



Workflow: Before Test




26



Workflow: Headphone Check

Apply Headphones



To ensure an accurate evaluation, it is important that headphones are correctly placed on the left and right ears.

The left side has a blue indicator on it and the right side has a red indicator.

Quit Continue


27




Workflow: Headphone Check

Headphone Check


Select the ear where you hear the sound.



Left Ear



Right Ear



Quit Continue

28



Workflow: Text Instructions

Hearing Test Instructions

You will hear tones presented at different pitches. Some tones will be loud and others will be very quiet.

Press Yes, if you hear the tone, even if it is very quiet. Press No, if you did not hear the tone.

If you hear static or noise, ignore it, and only press Yes if you hear a tone.

Quit

Continue

29



Workflow: Test Options

Test Options

Test Type

Diagnostic Screening

Selected Tests



Pure Tone Air



Pure Tone Bone



Speech

Test Options

Test Type

Diagnostic Screening

Selected Tests



Pure Tone Air

30



Workflow: After Test

After Test

☒ Generate Report Setup

☐ Display Report

☒ Export Data Setup

☒ Return to Welcome Screen



31



Workflow: Generate Report

Allows user to name Report

Report File Name Setup

File Location: C:\Users\mms_kma\Desktop\Audiologist Testing\Bugs Browse

Available	Selected
Patient ID	Patient Last Name
DOB	Patient First Name
Gender	Session Date
Phone Number	
Email	
Session Time	

Delimiter: Underscore '_' Space ' ' Dash '-'

Example: Doe_Jane_20001231.pdf

Save Close



32



Workflow: After Test Flex Only

After Test

☒ Generate Report Setup

☐ Display Report

☐ Export Data Setup

☒ Return to Welcome Screen



33



Workflow: After Test Flex Export (EMR)

Data File Name Setup

File Location: C:\ProgramData\Grason-Stadler\GSI-AMTAS\ Browse

File Name Fields:

Available	Selected
Patient First Name	Session Date
Patient Last Name	Session Time
Patient ID	
DOB	
Gender	
Phone Number	
Email	

Delimiter: Underscore '_' Space ' ' Dash '-'

Example: 20001231_150102.xml

Save Close



34



Demographics: Flex and Pro



35




Screening: Flex and Pro



36




Air Diagnostic Flex



37



Air/Bone Diagnostic: Pro



38



Apply 4000 Hz bone RETSPL adjustment

Adds 15 dB to the 4000 Hz bone conduction threshold


This adjustment is provided to eliminate the air-bone gap at 4 kHz that has been reported in the literature.



39




Speech Diagnostic: AMTAS Pro



40



Config App vs Settings



41



Patient Set Up




42



Patient Set Up – Pro Only

- Bone Oscillator Placement
 - Forehead placement
 - AMBAND



43



Patient Set Up – Flex and Pro

- Headphone Placement
 - DD450 -Circumaural headphones




44

gsi
Grason-Stadler

Welcome to
GSI AMTAS™
Self Administered Automated Audiogram

≡




Get Started ▶

45


gsi
Grason-Stadler


Pure Tone – Flex and Pro

Listen for a TONE.




Did you hear the TONE?



YES



NO


Listen for a TONE.




Did you hear the TONE?


YES


NO



46



47




Speech - Pro Only



gsi



gsi




gsi



gsi




48




Thank you for completing the test.
Press the button below to conclude the test

Finish

* Contact your administrator




49



Amtas Pro and Flex

Integration



50




Automated Report


- Patient Information
- Audiogram
- Audiogram Symbol Legend
- Masking Level Table
- Quality Assessment Table
- Audiogram Classification Table
- Speech Recognition Table
- Comments




51



52





GSI AMTAS™ Hearing Screening Report

Results:


Frequency	Right	Left
1000 Hz	Yes	Yes
2000 Hz	Yes	Yes
4000 Hz	Yes	Yes
500 Hz	No	No
	Pass	Pass

Level: 25
Pass Criterion: At least 2 Yes responses
Predicted Accuracy: Poor

Interpretation:


Your Right Ear is within the normal hearing range.
Your Left Ear is within the normal hearing range.

This is a screening test.
If you have questions or concerns
contact an audiologist or a physician.



GSI AMTAS™ Your Hearing

AUDIOGRAM



The audiogram is a graph that shows what sounds you are able to hear and what sounds you cannot hear. Each mark on the audiogram is the softest sound you can hear - your threshold - for a particular frequency or pitch. Normal hearing people hear sounds that are 20 decibels (dB) or less at all the frequencies.

Hearing loss can be described by the degree of loss - how loud sounds have to be for you to hear them - and the pattern of the thresholds shown on the audiogram. The degree of loss can be mild, moderate, severe, or profound.

The speech area represents the sounds that make up everyday conversational speech. If your thresholds are all above the speech area, you hear all the sounds that are important for understanding speech. If your thresholds are all below the speech area, you hear none of the sounds in normal speech. If some of your thresholds are in or below the speech area, you are missing some of the sounds that are important for understanding everyday conversation.

Hearing loss is a symptom of a problem somewhere in the ear. It can be in the outer ear, the middle ear, or the inner ear. Some of these conditions can be treated with medication or surgery. Many people with hearing loss are helped by hearing aids. It is important to find out the cause of the hearing loss so that the appropriate treatment can be provided.

Your Hearing Report:

Your audiogram shows that the hearing in your RIGHT ear is A MILD, FLAT HEARING LOSS.


A flat hearing loss is one where the thresholds are about the same for all frequencies. When there is a mild, flat hearing loss some of the sounds of everyday speech can be heard and some can't be heard. This causes difficulty understanding speech especially when there is background noise, a soft speaker, or a reverberant room. People with mild, flat hearing loss usually can understand in quiet rooms with a clear speaker. Some people with mild, flat hearing loss are helped by hearing aids. Others get by without hearing aids.

Your audiogram shows that the hearing in your LEFT ear is NORMAL.

People with normal hearing understand everyday conversation without difficulty. When there is background noise, a soft speaker, or a reverberant room (one with a lot of echoes), some normal hearing people have trouble understanding what people are saying even though other people seem to understand.

It is possible to have ear diseases and still have normal hearing. If there are other symptoms of ear trouble, such as pain, ringing in the ear (tinnitus), liquid in the ear, or a plugged feeling, you should discuss this with your doctor.


53



Qualind: Interpretation

- Patented method for determining the accuracy of a test results
- Data was collected at three sites from a wide range of settings, patient demographics, and hearing loss characteristics.
- Large subject sample (n = 120), a strong relationship was found between predicted and measured accuracy.
- Method may be useful for automated test procedures to provide quality assurance.

Qualind™: A Method for Assessing the Accuracy of Automated Tests. Margolis, et al, 2007, JAAA, 18: 78-89.



54



Quality Assessment Table (Qualind)

- **Predicted Accuracy** (Overall Quality -Good, Fair, Poor)
- **Predicted Average Absolute Difference** –Difference between automated and manual thresholds
- **Masker Alerts**-thresholds where masking may have been too high or low
- **Time per Trial** - average time from stimulus to patient response
- **False Alarm Rate** – number of times patient responded with no stimulus presented divided by number of trials when no stim present
- **Average Test-Retest Difference** – average difference between 1 KHz test and retest in right and left hear
- **Quality Check Fail Rate** – number of times patient did not respond to stimulus above threshold divided by the number of measured thresholds
- **Number of Air / Bone Gap > 35 dB** – number of air bone gaps that exceed 35 dB
- **Number of Air / Bone Gap < -10 dB** – number of air bone gaps that are less than 10db



55



AMTAS Pro ONLY



56



Audiogram Classification Table (AMCLASS)

- The air- and bone-conduction thresholds are analyzed and the audiogram is given three descriptive categories based on configuration, severity, and site of lesion.
- In addition, the audiogram is analyzed for interaural asymmetry.
- 376 Billion unique audiograms



57



Configuration

- Normal Hearing
- Flat Hearing Loss
- Sloping Hearing Loss
- Rising Hearing Loss
- Trough-shaped Hearing
- Peaked Hearing Loss
- Other

Severity

- Normal
- Mild
- Moderate
- Severe
- Profound

Site of Lesion

- Conductive
- Sensorineural
- Mixed

Bilateral Symmetry

- Symmetric
- Asymmetric.



58



Additional Data Management

GSI Suite/Noah



59



Screening Audiometry – Flex and Pro

Screening Audiometry – Flex and Pro

Right Ear

Frequency [Hz]	Response
500	Yes
1000	Yes
2000	Yes
4000	Yes
	Pass

Left Ear

Frequency [Hz]	Response
500	Yes
1000	Yes
2000	Yes
4000	No
	Pass

Test Details

Screening Intensity Level	25
Screening Passing Criteria	2

Audiogram Classification

Ear	Hearing Range	Predicted Accuracy
Right	Within Normal	FAIR
Left	Within Normal	

60



Flex Data Transfer Options

1. **Local** → Available if Flex is on the same PC as GSI Suite
2. **Import** → Used to import data from a USB drive
3. **Wireless** → Used to transfer data from the Flex tablet to Suite PC via Bluetooth



61



Batch Transfer – Flex

Assign Tests

Patients				Unassigned Tests						
Last Name	First Name	Birth Date	Patient ID	Test Date	Last Name	First Name	Birth Date	Patient ID	Gender	Delete
GSI Suite	Test Patient	1/1/0001	0000140	1/29/2018 7:11:59 AM	R	D	1/28/2018		Unknown	Delete
				2/20/2018 8:13:00 PM	Morris	Karen	2/19/2018		Unknown	Delete
				2/20/2018 8:14:02 PM	Pawl	Karen	2/19/1969		Female	Delete
				2/20/2018 8:17:56 PM	Smith	Jonny	2/19/2018		Unknown	Delete
				3/5/2018 2:48:11 PM			1/1/0001		Unknown	Delete
				3/9/2018 10:22:31 AM	Morris	Karen	4/1/1989		Unknown	Delete
				3/9/2018 10:30:54 AM	Morris	Karen	1/1/0001		Unknown	Delete

Assigned Tests						
Test Date	Last Name	First Name	Birth Date	Patient ID	Gender	Unassign

OK Cancel

62



Import

1. Export Data from AMTAS Flex to USB drive



2. Launch Suite and select **Import** from the Batch dropdown
3. Locate the USB drive and import
4. Locate or Create GSI Suite patients
5. Assign tests to the appropriate patient



63



Wireless Transfer

1. Test patients with AMTAS Flex
2. Launch Suite and select **Wireless** from the Batch drop down menu
3. Imported data appears in the Assign Test dialog
4. Locate or Create GSI Suite Patient
5. Assign tests to the appropriate patient



64



After Flex Export

- Data that has been exported will be deleted from the Flex PC at the next export.
- Submenu: will restore the most recently exported test data.

This will restore 18 file(s) previously exported and marked for deletion.

Do you want to continue?

Yes

No

65





Transferring AMTAS *Pro* Data

1. Locate or create a patient in GSI Suite
2. Select the Patient button on the audiometer.
3. Select the Transfer Session menu option.
4. Select the AMTAS Pro data session for transfer. (Pro data is listed in the Session List with an "A" _Name_Date Time")



66





How does AMTAS fit with current practice?



69



Does AMTAS replace Audiologist? NO!!

AMTAS is only a hearing test.

Audiologist counsels, fits hearing aids, performs other necessary tests like Tympanometry and OAE, re-tests the patient that is inconsistent, masking dilemmas, makes recommendations and more.



70



Amtas Implementation

AMTAS Pro

- Basic Diagnostics at ENT
- Basic Diagnostics at satellite office
- Annual Evaluation
- Other?

AMTAS Flex

- Screen – hearing aid retail to get referrals
- Threshold – before patient moves to the booth for further testing
- Other?



71



Advantages of Automation

- Optimize use of audiologists' time
- Standardization
- Quantitative quality assessment
- Decrease errors
- Decrease cost
- Increase access
- Telemedicine



72



What patients would Benefit?

- Originally designed for Veterans Administration (VA)
- Busy audiology departments
- Busy ENT practices with limited staff
- School systems
- Satellite clinics
- Community outreach



73

AMTAS Publications



- Margolis, R.H. Automated Audiometry – Progress or Pariah? Audiology Online (www.audiologyonline.com). January 17, 2005.
- Margolis, R.H. Automated Audiometry: Progress or Pariah. Audiology Today 17:21, 2005.
- Margolis, R.H., Saly, G., Le, C., Laurence, J. Qualind™: A Method for Assessing the Accuracy of Automated Tests. J. Amer. Acad. Audiol., 18, 78-89, 2007.
- Margolis, R.H., Saly, G.S. Toward a standard description of hearing loss. Int. J. Audiology 46, 746-758, 2007.
- Margolis, R.H., Morgan, D.E. The Value of Automated Audiometry. Insights in Practice for Clinical Audiology, January 2008.
- Margolis, R.H., Saly, G.L. Prevalence of hearing loss types in a clinical population. Ear and Hearing, 29, 524-532, 2008.
- Margolis, R.H., Saly, G.L. Toward a standard description of hearing loss. Int. J. Audiology, 46, 746-758, 2007.
- Margolis, R.H., Saly, G.L. Asymmetrical Hearing Loss: Definition, Validation, Prevalence. Otology & Neurotology, 29, 422-431, 2008.
- Margolis, R.H., Morgan D.E. Automated Pure-Tone Audiometry – An Analysis of Capacity, Need, and Benefit. Amer. J. of Audiology, 17, 109-113, 2008.
- Margolis, R.H., Glasberg, B.R., Creeke, S., Moore, B.C.J. AMTAS – Automated Method for Testing Auditory Sensitivity: Validation Studies. Int. J. Audiology, 49, 185-194, 2010.
- Margolis, R.H., Frisina, R., Walton, J.P. (2011). Automated method for testing auditory sensitivity: II. Air Conduction Audiograms in Children and Adults. Int. J. Audiology, 50, 434-439, 2011.
- Margolis, R.H., Moore, B.C.J. Automated method for testing auditory sensitivity: III. Sensorineural hearing loss and air-bone gaps. Int. J. Audiology, 50, 440-447, 2011.
- Margolis, R.H., Eikelboom, R.H., Johnson, C., Ginter, S.W., Swanepoel, D.W., Moore, B.C.J. False Air-Bone Gaps at 4 kHz in Listeners with Normal Hearing and Sensorineural Hearing Loss. Int. J. Audiology 52:526-532, 2013.



74



Questions?



75