

- If you are viewing this course as a recorded course after the live webinar, you can use the scroll bar at the bottom of the player window to pause and navigate the course.
- This handout is for reference only. Non-essential images have been removed for your convenience. Any links included in the handout are current at the time of the live webinar, but are subject to change and may not be current at a later date.

© 2018 continued® No part of the materials available through the continued.com site may be copied, photocopied, reproduced, translated or reduced to any electronic medium or machine-readable form, in whole or in part, without prior written consent of continued.com, LLC. Any other reproduction in any form without such written permission is prohibited. All materials contained on this site are protected by United States copyright law and may not be reproduced, distributed, transmitted, displayed, published or broadcast without the prior written permission of continued.com, LLC. Users must not access or use for any commercial purposes any part of the site or any services or materials available through the site.

## Technical issues with the Recording?

- Clear browser cache using [these instructions](#)
- Switch to another browser
- Use a hardwired Internet connection
- Restart your computer/device

## Still having issues?

- Call 800-753-2160 (M-F, 8 AM-8 PM ET)
- Email [customerservice@AudiologyOnline.com](mailto:customerservice@AudiologyOnline.com)



## T Coils and Hearing Loops: Important to Client Satisfaction

Juliëtte Sterkens, AuD  
&  
Karen MacLennan, AuD/TSHH

### Learning Outcomes

After this course, participants will be able to:


- Describe telecoil orientation in custom and behind the ear hearing instruments.
- List three telecoil programs, and when and where to use them.
- Explain one in-office and one out-of-office telecoil validation measure.






CONTINUED<sup>®</sup>


## Hearing Loops: Perimeter Loop

CONTINUED<sup>®</sup>


## Hearing Loops: Phased Array

CONTINUED<sup>®</sup>

## Specifying T Coils



## T Coil Orientation in a Magnetic Field



CONTINUED

## The Importance of T-Coil Position in a Magnetic Field




Image © C. Compton-Conley, PhD

9

CONTINUED

## Programming T Coils For Use In a Hearing Loop



10


CONTINUED

## T-Coil Programming Tips

1. **Open/RIC/RITE:** consider manual T-coil *only* for loop program
2. **Closed/Occluded Fittings:** Consider manual M+T but take listening situation into account:
3. Considerations for HI/CI Users with FM & IR Receivers (Neckloops)
  - Amplified neckloops
  - Widening the neck loop
  - When a neckloop signal is just too low.....
  - explain clients have a right to a "reasonable accommodations" under the ADA




## Programming T Coils For Use In a Hearing Loop




## Comparing Mic & T Coil Gain/Frequency Using Spec Sheets

**Output - 2cc coupler**  
IEC 60118-7 / ANSI S3.22-2003



## Comparing Mic & T Coil Gain/Frequency Using Spec Sheets




Acoustic Input	Equivalent Magnetic Input
65 dB SPL	56.2 mA/m
60 dB SPL	31.6 mA/m
50 dB SPL	10 mA/m


CONTINUED<sup>®</sup>

## What is SPLITS?

### Coupler SPL for an Inductive Telephone Simulator




The ANSI SPLITS test allows the MFR to orient the instrument for optimum output – this may mean holding the instrument sideways to obtain best results



CONTINUED<sup>®</sup>

## What is SPLIV?

### SPLIV – Coupler SPL in Vertical Magnetic Field



The ANSI SPLIV test requires the instrument is oriented vertically (or as worn on the ear)

Technical Specifications	LT166H (H4P)		LT166H (L2P)	
	IEC 60318-4 EC T11 Ear simulator	IEC 60318-4 EC T11 Ear simulator	IEC 60318-4 EC T11 Ear simulator	IEC 60318-4 EC T11 Ear simulator
Reference test gain (50 dB SPL input)	1600 Hz/HFA	47	43	59
Full-on gain (50 dB SPL input)	1600 Hz/HFA	60	60	79
Maximum output (90 dB SPL input)	1600 Hz/HFA	130	121	137
Total harmonic distortion	500 Hz	0.6	0.4	0.5
	800 Hz	1.1	0.7	1.4
	1600 Hz	0.8	0.5	0.4
Equivalent input noise	1600 Hz/HFA	98	103	105
Frequency range (DIN 45605/ANSI)		22	20	24
Current drain		1.2	1.3	1.1


CONTINUED<sup>®</sup>

continued

## Determine if the HAs are capable of matching the Telecoil gain to the Mic gain

Look the ANSI Specs & compare to **Reference Test Gain (RTG)**


1. Do the SPLITS & RTG gain match?
2. **More important:** Do the SPLIV and the RTG match?  
(Ask MFR to provide in a format that makes comparing easy)
3. Note: If the T-coil is mounted vertically, or nearly so – (as is the case w/ most BTEs):  $SPLITS = SPLIV = RTG$



continued


However, ANSI Specs don't tell the whole story  
Verification is needed

(Some hearing aids do *not* perform in T-coil mode from what you see on your fitting screen)



Images used with permission  
Ingrid McBride, AuD

continued



**CONTINUED**

Oshkosh Convention Senior EXPO – Out of the Loop  
 – using the microphone on the iPhone



**CONTINUED**



CONTINUED

Oshkosh Convention Senior EXPO – in the Loop  
– using the telecoil equipped LoopBuds



CONTINUED


Verify Telecoil program  
in a HIT box  
once the fitting is *finalized*



audioScan  
VERIFIT  
The gold standard in digital instrument verification



CONTINUED



**continued**

## Telecoil Validation – is the T-coil setting right for your client?


- Call client at home and verify he/she is able to use telecoil function & hear adequately loud and clear enough
- Check telecoil setting in your calibrated waiting room loop
- Loop your receptionist's front-desk
- Install a floormat in reception area
- Refer clients to a local looped venue (or an ALS with neckloops)
- Refer clients to a local HLAA chapter Nearly all chapters utilize a hearing loop and offer CART where hearing aid users can try out their telecoil setting

Hearing Loop floor mat  
Note: A microphone on your receptionist's desk is required

Waiting room TV Loop

microphone  
Loop sign creates the magnetic field

## Hearing Loop vs. FM/IR system Logos



## Roadmap to a Looped Community

Sterkens J, Audiology Today 26 (3): 26-31

- Familiarize yourself the IEC60118-4 Hearing Loop Standard
- Fit clients with “transparent” telecoils & verify the telecoil response in a HIT box – query HA MFRs
- Loop your waiting room to validate T-coil setting & demo
- Educate clients about rights under the ADA, the different Assistive Listening Systems & refer to local (looped) venues – see [www.ALDLocator.com](http://www.ALDLocator.com) or [www.LoopFinder.com](http://www.LoopFinder.com)
- Find/foster local installers – [www.Hearingloop/vendors.htm](http://www.Hearingloop/vendors.htm)
- Identify places clients would benefit from loops. Make one hearing loop happen – and build from there. Engage clients.
- Build loop awareness via community lectures (e.g. Rotary clubs) social media and post loop info on your website.
- Empower clients by providing handouts that promote self-advocacy (next slide)



**CONTINUED**

*Share the gift of HEARING*

Hello, my name is \_\_\_\_\_

I use hearing devices to hear and to communicate. Like most people with hearing loss who use hearing devices, I often experience difficulty hearing in large gathering spaces even if a microphone and PA system are used.

On \_\_\_\_\_ (date) I experienced significant difficulty hearing what was said when I was at your facility.

**Lafayette Hearing Center**

Call us today at | 3122 South St. | Lafayette  
www.lafayettehearingcenter.com | 765-448-6226

**2** I'm one of 48 million Americans with some degree of hearing loss and had considerable difficulty hearing in your facility today.

May I recommend you consider installing a HEARING LOOP that would make your venue more accessible to those with hearing loss?

GITHInfo@hearingloss.org

**HLAA**  
GET IN THE HEARING LOOP

**Hearing Loops**

- transmit sound from public address systems wirelessly, much like Wi-Fi, to any hearing aid or cochlear implant with a telecoil
- deliver sound customized for one's own hearing loss pattern
- do not require those with hearing devices to seek, wear and return special equipment and are, for this reason, user-preferred
- are cost-effective and inconspicuous to use
- accommodate one or hundreds of users at the same time
- can also work with portable receivers for those who do not use a hearing device or Loopflute with a smartphone.

For more information:  
hearingloss.org/hearing/help/technology/what-hearing-loop-technology/

**Induction Hearing Loops**

May I recommend that your facility install an Induction Hearing Loop?

**Hearing Loops:**

- Broadcast sound from the PA system wirelessly to any hearing aid with a telecoil
- Deliver sound customized for one's own hearing
- Do not require those with hearing devices to seek, wear and return special equipment
- Are cost-effective and inconspicuous to use
- Reduce your need to purchase, maintain, sanitize and administer individual receivers
- Accommodate one or 100+ persons at the same time
- Can also work with portable receivers for non-hearing device users

For more information, encourage you to call:  
**Lafayette Hearing Center: 765-448-6226**

Ask for an informational packet on Induction Hearing Loops and learn how you can give persons with hearing loss equal access to sound in your facility.

**Hearing Loops**

Get more from hearing aids in a hearing loop

**A hearing loop provides significantly enhanced functionality for hearing aids and cochlear implants.** (Cons with the most up-to-date technology, these devices cannot completely replicate important sounds from background noise, nor do they pick up all sounds from a distant sound as found in a performance hall, a place of worship or even a home TV viewed from across the room. In such difficult listening settings, induction hearing loops are often a solution.)

A hearing loop is a wire connected to an electronic sound source that transmits that sound to the telecoil in the hearing aids. A loop can directly surround a room, a theatre, a hall or even be woven around the body. It can be connected to a public address system, a living room TV or computer.

**A hearing aid equipped with a manual telecoil is needed to hear in a hearing loop.** The telecoil also called T-coil, receives the signal from the loop and sends it back into sound to the hearing aid. This creates an alternative sound to the background noise. The telecoil means only some sound that is received whether it is speech from a speaker, a television, a computer.

**Hearing loops double hearing aid functionality.** Using the telecoil in conjunction with a hearing loop is the most cost-effective way to improve the quality of audio in various situations.

**Looped-receivers display this symbol:**

**How to "Get in the Loop":**

**Where are hearing loops used?**

Hearing loops are helpful in a variety of places. Some are used to improve access, and others are optional where a person is in need.

Examples of locations that use include:

- Theatres and performing arts centers
- Places of worship
- Board rooms and large meeting rooms
- High schools and college lecture halls
- Court rooms or City Council rooms
- Restaurants or sports facilities
- Following halls

Examples of manual telecoil include:

- Travel agencies and information booths
- Desktop phones and answering machines
- Once this and pick up wireless
- Electronic news and boards
- Museum exhibits

**A few common home devices that can be used with a neck or small room loop:**

- Television or computer
- Telephone (landline or cell phone)
- MP3 or speaker book CD player

For more information:  
www.hearingloss.org

**CONTINUED**

## Add'l Hearing Loops Resources

**The Telecoil**  
Connecting Directly to Sound

A telecoil can make a dramatic difference in your ability to hear speech and understand language.

**Get more from hearing aids or cochlear implants with hearing loops**

A hearing loop provides significantly enhanced functionality for hearing aids and cochlear implants.

**Hearing Loops**  
The Preferred Assistive Listening Technology

From Technologies For Hearing, Aug. 2014, with permission

**LET'S LOOP**  
America's WORSHIP CENTERS

By making assistive listening hearing aid compatible, churches are leading the way to doubled hearing aid functionality for people with hearing loss

by David G. Myers

**In the loop**  
(hearing)

Bestselling psychology textbook author David Myers is teaching his largest audience yet, through a national campaign to provide

**A Guide to Understanding Hearing Loops**

**HLAA**  
GET IN THE HEARING LOOP

"For the first time in my life, I was able to hear every single word."

Get in the Hearing Loop is a communication access program of the Hearing Loss Association of America

**America Is Getting in the Hearing Loop**

AN OFFICIAL HLAA CAMPAIGN

HLAA's public campaign, America Is Getting in the Hearing Loop, is a national effort to raise awareness of hearing loss and the benefits of hearing loops. The campaign is a part of the HLAA's ongoing efforts to improve the lives of people with hearing loss.

HLAA is a national organization that works to improve the lives of people with hearing loss. We provide a variety of services, including advocacy, education, and support. We also work to raise awareness of hearing loss and the benefits of hearing loops.

For more information, visit our website at [www.hearingloss.org](http://www.hearingloss.org).

continued

## Additional Resource Materials

**Test Loop on-Site Visit**  
Hearing loop systems are venue specific and almost always require an on-site visit ahead of time, to provide an accurate estimate of your installation cost. Most thorough site visits take two hours; more involved installations might require more time. Although some designs can be modeled on a computer, computer simulation cannot determine if magnetic background noise can be reduced on a computer. Your particular building has its own...

**Best Practices in Hearing Loop Procurement**  
Hearing loops are quickly becoming the hearing assistance hearing system of choice for bringing clear sound to people with hearing loss. From schools and houses of worship to concert venues, assisted living facilities and municipal buildings, hearing loops are a consistent preferred solution and the only system that is directly hearing-aid-compatible and will make your facility hearing friendly. This remarkable system will not only bring many of the 48 million Americans with hearing loss back to theaters, houses of worship and community activities, but it also is a proven vehicle to build businesses.

Hearing loop systems are used worldwide, and most installers follow the international standard IEC 60118-4 as developed under the auspices of the IEC (International Electrotechnical Commission).

This standard defines the strength of the magnetic field, frequency response and methods of measuring these requirements. It also specifies the maximum levels for electromagnetic background noise.

Compliance with the IEC standard means a hearing aid user can walk into Westminster Abbey in London, the Great Ford Airport in Grand Rapids, Michigan or the Robert Hoodson Theatre in New York City and hear the sound directly, and all at a comfortable level. If their hearing aid is equipped with a remote called a telecoil or telecoil, loop hearing devices are available for those who don't have telecoils or don't use hearing aids. People don't have telecoils or don't use hearing aids at all.

The checklist is meant to give you guidance in the due diligence process as you procure a hearing loop for your facility by choosing the right installer. In some geographic areas of the country, it might not be possible to find highly experienced installers. It is therefore recommended you choose an installer who has been trained in IEC standard verification, has technical support from the supplier and is legally allowed to carry out the installation in your geographic area.

**Hire a Knowledgeable and Committed Hearing Loop Installer**

- Who trained or certified the installer and is the installer available to provide references?
- What design, installation and audio experience does the installer or the supplier's technical support department have with the type of building that needs looping?
- Will the provider a Certificate of Conformity, and make actual hearing aid users in the first testing?
- Does the installer offer information about hearing loops and the IEC standard on their website?
- Does the installer list loop installations on their website or on national websites? If not, why not?

Two companies offer hearing loop certification: Certevo, Inc. and Loop Technologies.

It is advisable that the purchaser install on an IEC 60118-4 standard hearing loop installation in writing. This will not add to the cost of a loop installation but effectively guarantees a working system.

Buildings present many variables with regard to design and installation due to metal in the floors and ceilings. Electromagnetic interference due to radio, poorly installed wiring that might not meet the current electrical codes could cause interference, buzzing noises, that a hearing aid or loop device users can hear when they turn their head on. This magnetic noise must have been previously present in the facility, but might not have been of concern until now.

Note: If your facility is required to provide an Assistive Listening System (ALS) under the Americans with Disabilities Act (ADA) and magnetic background noise is determined to be of excessive levels during a hearing loop site visit, your facility, be advised that the ADA requires 25 percent of the ALS receivers be hearing and telecoil compatible or personal telecoils. Therefore magnetic background noise should always be investigated by a licensed electrician, even for the use of FM or infrared systems.

**Frequently Asked Questions About Hearing Loops**

- How many Americans live with hearing loss? In 610 million Americans adults report at least "a little hearing loss." According to a 2011 report based on nationwide testing of Americans 12 and older, the National Health and Nutrition Examination Survey (NHANES), 36 million Americans have at least a mild hearing loss, and 27 million have a moderate to severe hearing loss. About 18 million of those 36 million live in both ears and 48 million in one or both ears. About 18 million of those 36 million live in both ears and 48 million in one or both ears. About 18 million of those 36 million live in both ears and 48 million in one or both ears.
- Why are hearing loops needed? One hearing aid enables hearing? Today's digital hearing aids effectively enhance hearing in conventional settings. Yet the same people with today's digital hearing aids effectively enhance hearing in conventional settings. Yet the same people with today's digital hearing aids effectively enhance hearing in conventional settings. Yet the same people with today's digital hearing aids effectively enhance hearing in conventional settings.
- How many hearing aids have the national recall request for receiving hearing loop input? From its survey of hearing professionals, the Hearing Journal (April 2009) reported that 85% of hearing aid brands included a recall, as measured from 17% in 2007. In the 2009-2010 survey of hearing aid models, the Hearing Journal Product reported that 126 of 141 (89%) of 141 hearing aid models, including 10 of 10 of the hearing aid models included in the recall. In the 2009-2010 survey of hearing aid models, the Hearing Journal Product reported that 126 of 141 (89%) of 141 hearing aid models, including 10 of 10 of the hearing aid models included in the recall.
- Can hearing loops solve those without telecoils or without hearing aids? Yes, all forms of hearing loss, including hearing loss with hearing aids and hearing aids. Yes, all forms of hearing loss, including hearing loss with hearing aids and hearing aids. Yes, all forms of hearing loss, including hearing loss with hearing aids and hearing aids.
- How long does a hearing loop last? Cost ranges from \$250 to \$1000 for a standard hearing loop system. Most churches can install a professional installation in an average-sized auditorium or worship space. Most churches can install a professional installation in an average-sized auditorium or worship space. Most churches can install a professional installation in an average-sized auditorium or worship space.
- Is this a decade-old technology? Like electronic computers, magnetic induction loop technology began more than 70 years ago, and now is in newly developed forms (both low amplitude and telecoil technology), and now computer-aided designs for complex installations and with increasing applications.

continued

Thank You for Your attention  
Questions?

## References

- Caccavo M and Lopez S (2014). How to Win Patients and Influence Community Communication: A Step by Step Guide for Looping Your Community. *Audiology Practices* 6 (4): 26-27
- Kochkin S, Sterkens J, et al (2014) Consumer Perceptions of the Impact of Inductively Looped Venues on the Utility of Their Hearing Devices. *Hearing Review* 21(10): 16-26
- Kaufmann T, Sterkens J and Woodgate J (2014) Hearing Loops: The Preferred Assistive Listening Technology. *J Audio Eng. Soc.* 63 (4): 298-302
- Ledda K, et al (2018) Difference in Speech Recognition between a Default and Programmed Telecoil Program. *J Am Academy Audiology* – in press
- Remensnyder L. (2011) Invigorate Your Practice with Hearing Loops. *Audiology Practices* 3 (4): 32-35, 61
- Sterkens J (2011) Get your Practice and Patients in the Loop: Telecoil and Hearing Loop Essentials for Practitioners. *Audiology Practices* 3 (4): 36-39
- Sterkens J (2014) Roadmap to a Looped Community. *Audiology Today* 26 (3): 26-31.
- Valente M (2013) The Telecoil: The Lonely Transducer that Can Be a Big Producer. Retrieved at [www.audiologyonline.com/articles/telecoil-lonely-transducer-that-big-11911](http://www.audiologyonline.com/articles/telecoil-lonely-transducer-that-big-11911)

## Useful websites:

[www.HearingLoop.org](http://www.HearingLoop.org)  
[www.LoopWisconsin.com](http://www.LoopWisconsin.com)  
[www.hearingloss.org/programs-events/get-hearing-loop/](http://www.hearingloss.org/programs-events/get-hearing-loop/)

For further information: Contact Juliette Sterkens, AuD  
[jsterkens@hearingloss.org](mailto:jsterkens@hearingloss.org)

Karen MacLennan, AuD  
[drmaclellan@gmail.com](mailto:drmaclellan@gmail.com)



continued