



Evoked Potentials – Part 2

Electrocochleography, MLR, P300 and Cortical Responses

Kathleen Hill, Au.D.
Director of Commercial & Education Development
Otometrics, a division of Natus

natus.

1

otometrics
a division of natus.

Evoked Potentials

- Electrical signals generated by the nervous system in response to a stimulus
- Event related (evoked by onset of stimulus)
- Useful in diagnosing a variety of neurological disorders

natus.

2


otometrics
a division of natus.

Presented in Cooperation with

continued™

Auditory Evoked Potentials

OVERVIEW of AEP in the Time Domain



natus.

3

otometrics
a division of natus.

Recording Evoked Potentials

- Far Field recording-electrodes placed on the scalp
- Electrodes assigned as Active (+), Reference (-), and Ground
- Record from a pair of electrodes (i.e. Cz to A1)
- Recording from an non-inverting (+) and inverting (-) electrode
- Amplify the difference between the signals
- Proper care of the electrodes & impedance are important because this is the medium for collecting the EP data

natus.

4

otometrics
a division of natus.

Presented in Cooperation with

continued™

Recording EP Data

- Common mode rejection
 - The voltage is different between the active and reference electrodes
 - The voltage related to noise is similar at both electrodes and the response voltage has the greatest difference
 - The response at the reference electrode (inverting) is added to the response at the active (non-inverting) electrode
 - The components which are “common” to both electrodes are cancelled (i.e. noise - mains, biological or environmental)
- Inter-electrode impedance is the most important (common mode rejection does not work well if inter-electrode impedance varies)

Preparation


Presented in Cooperation with



10-20 System of Electrode Placement

True Cz = larger
amplitude
response

~15% reduction
in amplitude
using high
forehead




natus.

7

otometrics
a division of natus.

Impedance



natus.


8

otometrics
a division of natus.

Presented in Cooperation with

continued™

Checking Electrode Impedance



Goal: Low electrode readings of **<5K** Ohms each balanced with no more than **<2K** Ohms difference between electrodes


Inter-electrode impedance is the most important (common mode rejection does not work well if inter-electrode impedance varies)

natus.

9

otometrics
a division of natus.

ICS Chartr EP 200 Preamp



natus.

10

otometrics
a division of natus.

Presented in Cooperation with


continued™

Selecting a Protocol

Procedure
Selected
Protocol

The CHARTR EP default procedures are:

- **ABR (Auditory Brainstem Response)**
- **ASSR (Auditory Steady State Response) (optional)**
- **ECochG (Electrocochleography)**
- **VEMP (Vestibular Evoked Myogenic Response)**
 - **EMG monitor (optional)**
- **ALR (Auditory Late Response)**
- **AMLR (Auditory Middle Latency Response)**
- **P300 (optional)**



Listening Check

- Use subjects with hearing thresholds within normal limits
- Listen to the stimulus
- Should be able to barely hear it at 0-5 dB
- If the softest level anyone can hear the stimulus is 20 dB – there is a 20 dB correction factor
- Subtract the correction factor from the response obtained on the dial (i.e. patient's threshold is 60, 60-20=40; 40 is their true threshold)
- Call local representation or manufacturer to have system calibrated


Presented in Cooperation with

Electrocochleography

natus.

13

otometrics
a division of natus.




Electrocochleography (ECoChG)

Backed-up fluid in the sac and inner ear leads to swelling and pressure.

Swelling in the balance canals distorts balance information.

Swelling in the hearing canal distorts or blocks sound information.



Distorted information travels from the inner ear to the brain.

natus.


14

otometrics
a division of natus.

Presented in Cooperation with

continued™

Types of Electrodes - Tiptrode



natus.

15


otometrics
 a division of natus.

Application of Tiptrodes

- The ear canal must be scrubbed in order to obtain lower impedances. Be careful is using conductive paste on gold foil (may remove the gold).
- This electrode is further away from the generator site and results in less pronounced SP/AP ratio.
- Attach the metal clip making good contact with the gold foil near the base of the foam tip.
- Be careful when compressing the foam tip as not to remove the gold foil layer.
- Make sure to remove the red and blue tubing of the insert phones. Do not double the length of the tubing. Attach the tubing from the tiptrodes to the insert transducer boxes.
- Attach leadwire to the patient cable.

natus.

16



otometrics
 a division of natus.

Presented in Cooperation with



ECochG with Tiptrode

1-CHANNEL ECochG with Tiptrode (Horizontal Montage)




natus.

17

otometrics
a division of natus.

Extra tympanic electrode – TM ECochGtrode



Tymptrode and insert earphone in earcanal.
Closer to the neural generator gives a more
robust response.

natus.

18

otometrics
a division of natus.

Presented in Cooperation with


continued™

Application of TM ECochGtrode

- Perform Otoscopic exam first. Remove any cerumen obstructing the ear canal before placing the TM ECochGtrode.
- Do not place if there is a TM perforation.
- Place 2 cc of saline into the ear canal and then drain the ear. Saline will help reduce impedance. Make sure saline is body temperature (37 degrees Celsius or 98 degrees Fahrenheit).
- Place a small amount of conductive gel on the cover of the hydrogel end of the electrode. Too much gel can occlude the insert earphone.
- Secure the electrode with tape to prevent the cable from moving.
- Electrode impedance may be high (>20K ohms). Therefore, you may have a large interelectrode impedance difference.
- Might need to disable artifact rejection to collect the response.
- After testing, carefully remove insert eartip and electrode. TM may be red but inflammation will go away within a few hours.

ECochG with Tymptrode

1-CHANNEL ECochG with TM ECochGtrode (Horizontal Montage)




Presented in Cooperation with

continued™

Clinical Applications

- Assist with Meniere's diagnosis (endolymphatic hydrops) but can't diagnose Meniere's with this test alone.
 - Characteristics of Meniere's disease
 - Vertigo
 - Hearing loss (typically fluctuating)
 - Tinnitus
 - Aural fullness
- Improve identification of Wave I
 - Not always visible with high frequency hearing loss
 - Used for determining interwave latencies I-III and I-V
- Identification of cochlear microphonic
 - Auditory neuropathy
 - Must use rarefaction and condensation polarity

ECochG Test Parameters




Why use 90 dB? Ecog can be detected as low as 60dB but increases amplitude with increase in intensity. Easier to discern SP and AP.

Presented in Cooperation with

continued™

ECochG Test Parameters



Why use a slow rate? Response amplitude decreases with faster rates. Easier to discern SP and AP. Why alternating? Easier to discern SP and AP when CM is canceled out.

natus.

23

otometrics
a division of natus.

ECochG Components

CM (cochlear microphonic)

a stimulus-dependent cochlear response, which changes direction with changing polarity. Hence, it is not detected when averaging is performed to alternating polarity (**running a rare and con run will cancel out the CM and help locate baseline**)

SP (summating potential)

direct current response from the Organ of Corti hair cells. SP is often seen as a leading hump on the AP or wave I, although sometimes it can appear as a separate hump

AP (action potential)

alternating current response generated by the cochlear end of the 8th nerve (wave I). The AP represents the summed response of thousands of firing auditory nerve fibers

natus.

24


otometrics
a division of natus.

Presented in Cooperation with

continued™

ECochG Results and Interpretation

- **BL** – baseline
- **SP** – summing potential
- **AP** – action potential




natus.

25

otometrics
a division of natus.

Normal Range

- < 50% - Tiptrode
- < 25% - Tymptrode



Normal SP/AP ratio

Abnormal SP/AP ratio

natus.


26

otometrics
a division of natus.

Presented in Cooperation with

continued™

Calculation of SP/AP ratio



MLR (AMLR)

Presented in Cooperation with

continued™

MLR

- Derived from the medial geniculate, inferior colliculus and primary auditory cortex
 - Produces waveforms that are identifiable near threshold and may be useful in identifying low-frequency auditory sensitivity.
 - ABR latencies are more consistent
 - MLR has larger amplitude than ABR and less dependent on neural synchrony.
 - Amplitude is smaller if patient is asleep. Latency is unchanged based on patient state.
 - Anesthesia and light sedatives can eliminate MLR response.
-
- Montage – Cz/A1/A2/ground – Cz (vertex) gives best response.
 - C3 and C4 are also recommended instead of Cz.

Clinical Usage

- Determine auditory cortical function (presence/absence)
- If dysfunction is present determine side affected

Presented in Cooperation with



MLR Stimulus Factors

- Change in intensity does not effect latency
- Po, Pa and Pb are observed at higher intensities
- Pa will remain for lower intensities. With reduced intensity may see a slight increase in latency and decrease in amplitude.
- Not highly frequency dependent – increase stimulus frequency will result in slight reduction of latency and lower amplitude.
- Stimulus Rate should be between 6 -10 stimuli/sec


MLR Response

- MLR occurs at 10 to 80 msec
- Na – 12-18 msec
- Pa – 25-30 msec
- Pb - ~50 msec

Presented in Cooperation with



MLR Test Parameters




Why use a slow rate? Response amplitude decreases with faster rates.
Easier to the MLR.

natus.

33

otometrics
a division of natus.

MLR Response



natus.

34

otometrics
a division of natus.

Presented in Cooperation with

continued™

P300

P300

- Dependent on internal thought process – patient must be awake and alert to the stimuli
- Hippocampus generates the response (this is where short term memory functions)
- Montage – Cz/A1/A2/ground

Presented in Cooperation with



P300 Stimulus

- 2 stimuli (frequent and infrequent) (typically 1000 and 2000 Hz toneburst)
- Patient listens and counts the infrequent (rare) stimuli
- Ratio 5 frequent to 1 infrequent stimuli (~20 rare stimuli presented)
- Only the infrequent stimuli generate the P300
- More difficult the discrimination task the smaller the P300 response
- Can present stimuli binaurally


P300 Response

- Slow, Broad response with a peak at 250 to 300 msec
- Decrease in amplitude and increase in latency is indicative of dementia, Alzheimers, and various neurological and psychiatric diseases
- Diagnoses disorders of cognitive function

Presented in Cooperation with



P300 Test Parameters




natus.

39

 otometrics
a division of natus.

P300 Response



natus.


40

 otometrics
a division of natus.

Presented in Cooperation with

continued™

P300 Response



natus.

41

otometrics
a division of natus.

Cortical Response

natus.

42

otometrics
a division of natus.

Presented in Cooperation with

continued™


Cortical Response

- Derived from a slow vertex response
- Produces waveforms that are identifiable near threshold and may be useful in identifying hearing status using toneburst stimuli
- Used for medical-legal cases or patients difficult to test with audiometry
- Response is not fully mature until late teens
- Montage – Cz/A1/A2/ground – Cz (vertex) gives best response.
- Patient should be awake (reading a book)

<http://corticalera.com/> (for more info)

Cortical Response

- Trial Setting



Presented in Cooperation with

continued

Cortical Response

- P1, N1, P2 are marked
- Decrease in Intensity
 - Amplitude decreases
 - Latency does not change much
- Response habituates so very few sweeps (10 to 20) are collected


Cortical Response

- Cortical occurs at 50 to 300 msec
- N1 – ~100 msec
- P2 – ~200 msec
- Sum responses from 3 waveforms
- If response is $>5\mu\text{V}$ = threshold is 5 dB lower
- If response is $<5\mu\text{V}$ = threshold is same level

Presented in Cooperation with



Cortical Response



natus.

47

otometrics
a division of natus.

Thank You!

<https://otometrics.natus.com/en-us>

natus.

48

otometrics
a division of natus.

Presented in Cooperation with

continued™