




Learner Outcomes

- After this course learners will be able to describe the four new pieces of equipment in the Bio-logic range: the AudX, the AudX PRO, the AudX PRO FLEX and the NavPro ONE.
- After this course learners will be able to discuss the technology behind frequency modulated DPOAEs, as well as their benefits and clinical utility.
- After this course learners will be able to describe the technology behind binaural OAE, ABR and ASSR testing, their benefits and clinical utility.
- After this course learners will be able to discuss the technology behind multi-rate ASSR testing, its benefit and clinical utility.
- After this course learners will be able to describe the benefit of using spread spectrum technology for noise reduction in evoked potential testing.


THE NEXT GENERATION BIO-LOGIC BY OTOMETRICS (NATUS)



otometrics
a division of natus

3

The Bio-logic legacy lives on



Bio-logic[®]
otometrics

Introduction

Presented in Cooperation with

continued

INTERNAL COMMUNICATION

Meet Otometrics' newest family members



Screening

Diagnostic

otometrics
a division of natus

5

NEW BIO-LOGIC SOLUTIONS

Bio-logic® AuDX®



Screening OAE (TE or DP)

Optional modules:

- Audiometry Screening
- MAGIC

otometrics
a division of natus

6

Presented in Cooperation with

continued

NEW BIO-LOGIC SOLUTIONS

Bio-logic® AuDX® PRO

Diagnostic and screening DPOAE
Diagnostic and screening DPOAE & TEOAE

Optional Modules:

- Diagnostic and Screening Audiometry (Air/bone or Air/Bone/Speech)

otometrics
a division of natus

7

NEW BIO-LOGIC SOLUTIONS

Bio-logic® AuDX® PRO FLEX

Screening Tympanometry
Diagnostic Tympanometry

Optional modules:

- Diagnostic TE/DPOAE
- Screening TE/DPOAE
- Screening Audiometry
- Diagnostic Audiometry (Air/Bone/Speech)

otometrics
a division of natus

8

Presented in Cooperation with

continued

NEW BIO-LOGIC SOLUTIONS

Bio-logic® NavPRO ONE

Auditory Brainstem Response (click and chirp)

Optional modules:

- ASSR – Auditory Steady-state Response
- EABR – Electrical ABR for cochlear implant
- Electrocochleography (EcochG) (included in US base)
- Frequency specific stimulus package (Included in US base)
- Diagnostic and Screening TE/DPOAE
- Diagnostic Audiometry (Air/Bone/Speech) including Screening Audiometry
- Diagnostic Immittance



9

BIO-LOGIC

Full Modularity- Bio-logic family overview

Module	AuDX	AuDX PRO	AuDX PRO FLEX	AuDX NavPRO ONE
Scr DP	✓	Incl. With DX	+	Incl. With DX
Dx DP	or	✓	+	+
Scr TE	✓	and/or	+	Incl. With DX
Dx TE		✓	+	+
Scr Aud Air	+	Incl. With DX	+	Incl. With DX
Dx Ext HF Aud		+	+	+
Dx Aud Air		+	+	+
Dx Aud Air/Bone		+	+	+
Dx Aud Speech		+	+	+
MAGIC	+	+	+	+
Scr Tymp			✓	Incl. with DX
DX Tymp			+	+
Eustachian Tube Func.			+	+
ABR				✓
ASSR				+
EABR				+
ECoChG				+

✓ Included in base device

+ add-on as module or license



10

Presented in Cooperation with






Frequency modulated DPOAE

 **otometrics**
a division of natus


Presented in Cooperation with

continued




Frequency Modulated Distortion Product Otoacoustic Emission (FM-DPOAE™)

- When the ear is stimulated by two tones simultaneously, there are **two DPOAE sources**, one at the overlap of f1 and f2 close to the f2 place and the other at the 2f1-f2 place.



Frequency Modulated Distortion Product Otoacoustic Emission (FM-DPOAE™)


- There are several techniques to overcome this problem.
- One is to use a third tone to suppress the second source at the $2f_1-f_2$ place. However, for doing this, there is need for a third loudspeaker within the sound-probe.
- Windowing method and onset-decomposition technique
- Or a patented technique using frequency-modulated primary tones. To do this, primary tone frequencies are varied over time



otometrics
a division of natus

15

DPOAE fine structure



otometrics
a division of natus

http://www.otoemissions.org/guest_editorials/2009/dhar_2009.htm

Presented in Cooperation with

continued


Frequency Modulated Distortion Product Otoacoustic Emission (FM-DPOAE™)

- Primary tone frequencies are shifted between $\pm 100\text{Hz}$ with a modulation rate of about 1.5 Hz. Because of the associated phase shift the impact of the second source is reduced. As a consequence, DPOAE detection and hence hearing threshold estimation is significantly improved.
- Due to frequency modulation the number of stimulated OHCs is increased resulting in a higher DPOAE level. Thus, FM-DPOAE does not need any additional stimuli and does not extend test time.



Result examples DPOAE vs. FMDPOAE


$L_2 = 60, 45, 30 \text{ dB SPL}$



Presented in Cooperation with



Example data without (left figures for right and left ears) and with FM (right figures for right and left ears)



otometrics
a division of natus

Summary 1

- FMDPOAE can eliminate or reduce unwanted DPOAE fine structure
- Can reduce fine structure related refers---provides speed and efficiency by reducing the need for re-tests
- Threshold estimation is probably more accurate
- FMDPOAE could be used instead of standard DPOAE in virtually any protocol

otometrics
a division of natus

Presented in Cooperation with

continued

Binaural Stimulation and Multifrequency Stimulus

- Ability to test both ears at the same time
- Ability to test multiple frequencies at the same time

- **=TIME SAVINGS!**

- the binaural presentation of single-pair stimuli had subtle effects on DPOAE levels
- the effects of contralateral inhibition created by binaural stimulations are negligible

21



Binaural DPOAE test

Test both ears!
two frequencies per ear = four
frequencies **simultaneously!**



Presented in Cooperation with





THRESHOLD ESTIMATION USING OAES

23

 **otometrics**
a division of natus

DPOAE -Threshold

DPOAE Threshold (patented method by PATH)

→ Estimating hearing loss by means of extrapolated DPOAE I/O-functions

DPOAE Quick test

→ Indicating the presence of a valid DPOAE at one or more selected primary tone levels and frequencies.

Two different protocols are available:

Screening (at a selected screening level) and Diagnostic (multiple selections of stimulus parameters).

 **otometrics**
a division of natus

Presented in Cooperation with

continued[™]

DPOAE -Threshold

- DPOAE I/O functions are measured at 1000 (optional), 1500, 2000, 3000, 4000, 6000 Hz and 8000 Hz (optional).
- Before the measurement of each frequency starts, a fine adjustment is performed, in order to check at which frequency the emission can best be recorded (fine structure, jitter +/-100 Hz).



DPOAE THRESHOLD ESTIMATION

DPOAE thresholds = DPOAE audiograms


Quantitative and frequency-specific assessment of hearing loss by means of extrapolated DPOAE I/O-functions



Presented in Cooperation with




Estimated threshold from DPOAE vs. FMDPOAE



otometrics
a division of natus

DPOAE level as a function of L_1 and L_2 stimulus levels

65/55
55/45
45/35
35/25
25/15




$L_1 = 0.4 L_2 + 39$ „scissors paradigm“
Kummer et al., Hearing Research, 2000

otometrics
a division of natus

Presented in Cooperation with

continued™

Optimized stimulus levels for testing the DPOAE I/O function



$$L_1 = 0.4 L_2 + 39 \text{ „scissors paradigm“}$$

Scissor paradigm

Scissors paradigm	L1	L2	$f_2/f_1 = 1.2$
	65	65	
	63	60	
	61	55	
	59	50	
	57	45	
	55	40	
	53	35	
	51	30	
	49	25	
	47	20 dB SPL	

$$L_1 = 0.4 L_2 + 39 \text{ dB SPL}$$

Presented in Cooperation with


Correlation between behavioral pure tone threshold L_T and estimated DPOAE threshold L_{EDPT}

118 adults
with cochlear hearing loss

Boege and Janssen
JASA 2002

otometrics
a division of natus

Supporting DATA



Presented in Cooperation with

Study...in UK

- 20 brick workers assessed for NIHL
- Full PTA both ears
- DPOAE Threshold testing at 2 and 4 KHz
- 80 results to compare to PTA
- Average difference 8.7dB



DPOAE Threshold

Select the frequencies that
you want to test.

**Options for FMDPOAE and
Multi channel are available**



Presented in Cooperation with




DP Threshold in progress

A quick sweep of frequencies performed to find best start frequency

DPOAE I/O functions are measured from L2 = 65 dB SPL in 10-dB-steps decreasing to L2 = 15 dB SPL.



DP Threshold in progress – cartoon mode



An cartoon video can be shown to maintain the children's interest and avoid artifacts



Presented in Cooperation with




DP Threshold results

An arrow pointing down from the 50 dB line indicates that no sufficient emissions could be measured in order to fulfil the criterion. A hearing loss greater than 50dB could be the case.



Pressurized OAEs AuDX PRO FLEX

- OAEs conducted at peak admittance pressure to maximize amplitude of the response



Presented in Cooperation with



NEW TECHNOLOGY IN EP



39



Bio-logic®
NavPRO ONE

ABR

 **otometrics**
a division of natus

40

Presented in Cooperation with

continued

BIO-LOGIC

Binaural Recordings

- AEP in general can be recorded from both ears simultaneously.
- A traditional recording scheme is to just apply different but constant stimulus rates to both ears, e.g. 37 Hz and 41 Hz.
 - A more preferable choices would be stimulus rates that do not have a common period (1Hz in the example), such as 37.3 Hz and 41.1 Hz etc.
- If averaging is done in synch to each ear's stimulus rate, responses can be recorded independently

41



BIO-LOGIC

Spread Spectrum

ABR recording with spread spectrum technology

- Traditionally, ABR recording (time domains) is done with a constant stimulus rate, although it does not necessarily require a constant stimulus rate
- Modifying the rate during testing can drastically improve the robustness against artifacts caused by electric noise sources
- Moreover, binaural recording can be performed with equal average stimulus rates, which can be preferable over fixed different rates for both ears (such as 37 and 41 Hz)

42



Presented in Cooperation with



BIO-LOGIC

Weighted averaging

- Natural artifacts, such as myogenic activity, interfere with ASSR and ABR measurements
- Weighted averaging is the preferred method to implement artifact management

43

 **otometrics**
a division of natus

BIO-LOGIC

Statistical Analysis

Template matching

- ABR responses have a typical waveform e.g. I, III, V.
- This knowledge can be used to improve automatic detection of ABR responses.
- The recorded signal can be cross-correlated with a template that represents a typical waveform, and the correlation signal can then be statistically analyzed instead of the unprocessed signal.
- This technique can also be used to estimate the latency of the response, which can get a better estimation in noisy recordings than a plain peak search would.

44

 **otometrics**
a division of natus

Presented in Cooperation with

 **continued**

BIO-LOGIC

ABR Technology

- Using the techniques described above, the Bio-logic ABR module was designed for efficiency. It makes use of spread spectrum, template matching, weighted averaging and chirp stimuli.
- Tracings are shown as the original recording, while all statistics, if enabled, are performed on the template convolution signal. This combines detection performance with real trace view.

45



BIO-LOGIC

ASSR

Stimulus Rates

Variable rates

46




Presented in Cooperation with



BIO-LOGIC

Simple Cochlear Model




Cochlea response to a
500 Hz modulated
tone

 **otometrics**
a division of natus

47

Cochlea response to 500 Hz octave-wide chirp at 80 Hz vs. 40 Hz
stimulus rate



 **otometrics**
a division of natus


48

Presented in Cooperation with

continued

Type of Stimulus Affects Cochlear Response

2kHz narrow band chirp at a rate of 80Hz, compared to a modulated sine stimulus




49

otometrics
a division of natus

BIO-LOGIC

Binaural ASSR

Selecting different rates for the two ears also allows binaural ASSR recordings. This figure contains responses from eight different stimulus rates, four were presented to each ear.



50

Presented in Cooperation with

continued

BIO-LOGIC

Spread Spectrum ASSR

The rate variation during testing is moderate: a ± 1 Hz variation of the rate (centered at 37 to 160 Hz) does not impact AEP recording but improves robustness against interference significantly

51



BIO-LOGIC

ASSR Technology

- Since recording ASSR is a fully automated procedure, not too many parameters need to be configured.
- The main decision is to select a stimulus bandwidth, which is a trade over between test performance and frequency specificity.
 - Wider band stimuli excite larger portions of the cochlea and therefore generate stronger evoked responses.
- Spread spectrum is always enabled in the Bio-logic ASSR module.

52



Presented in Cooperation with



BIO-LOGIC

Thank you!



<https://otometrics.natus.com/bio-logic>

53

otometrics
a division of natus

Presented in Cooperation with

continued