

TympStar Pro 1.2:

Implementing New Features for Efficient Testing



1



Laura Prigge, AuD

- Part of a team of 3 three audiologists
- Doctorate of Audiology from AT Still University.
- 15 years experience include providing manufacturing support as well as managing educational audiology training for an international audiologic equipment company.
- Clinician in both ENT offices and a retail hearing center conducting audiologic evaluations and hearing aid fittings on all patient populations including infants, adults, geriatric, and special needs patients.
- Time at GSI focused on training and education. In addition to webinars and phone training, and speaks at audiology congresses and conferences around the world.



2



Course Outcomes

- After this course, the participant will be able to navigate the configuration application to customize settings in the TympStar Pro and implement new features into testing routine.
- After this course, the participant will be able to perform the Multi Hz test using the TympStar Pro.
- After this course, the participant will be able to perform the Patulous ET test using the TympStar Pro.



3



Agenda

- Overview
- Review Updates
 - *New Multi-Hz*
 - *New Patulous ET*
 - *New Seek Below Start*
 - *New Baseline Configuration Options*
- Config App Updates
- Summary



4



TympStar Pro™ – Overview

- Precise, Reliable, Efficient
- Large Color Touch Screen Display
- Ultra-light Probe
- Pre-Programmed Testing Protocols
- Multi-Frequency, Multi-Component Tympanometry
- Integrated Special Tests
 - Multi Hz
 - ARLT
 - Patulous ET



5

Front Panel



6

gsi Grason-Stadler

Touch Screen Display



7

Probe or Shoulder Box

- Change ear-side
- Launch test



8



Most Frequently Asked Questions


- Why is there a screening option on a clinical middle ear analyzer?
 - Fast method to evaluate middle ear status
 - Determine if more testing is appropriate
 - Infants, children and difficult to test patients
- What's the difference between screening and diagnostic?
 - Screening: As fast as possible (600/200 daPa per second), auto start, and recording will stop after a peak is detected and the tracing returns to baseline
 - Diagnostic: Full control of parameters such as pump speed, baseline, probe tone, start pressure, sweep direction, and admittance.



9

Most Frequently Asked Questions

- What is the fastest way to do reflex thresholds?
 - Fully Automatic: Auto Sequence
 - Semi-Automatic: Threshold Seek ON
 - Manual: Threshold Seek OFF




10





Most Frequently Asked Questions

- **How Do I Manage Direct Printing:**
 - **Screen Print** – prints the data displayed on the screen
 - **Summary Print** – prints the marked data to a one page report that includes Tymps and Reflex/Reflex Decay test results in tabular format (no graphs)
 - **Session Print** – prints all of the marked data from all test tabs and includes reflex and reflex decay graphs.
 - **Dialog Print** – displays the print dialog where you may choose what to print. “Ask me what to do”




11



Most Frequently Asked Questions

- **Can I View Stored Sessions**
 - Patient Screen → Session List
 - Load any previously saved test session
 - Current test session is labeled with an “*”

NOTE: data is active/not locked so it may be altered by the user.



12



Most Frequently Asked Questions

How do I set the time without a keyboard?



13



Most Frequently Asked Questions

Can I load a custom configuration without a computer?

- Now possible to load a custom configuration from a USB



14



TympStar Pro v1.2



15



What's New with 1.2?

- Multiple Frequency Tympanometry Test (Multi-Hz)
- Patulous Eustachian Tube Test
- Seek Below Start Level-Reflexes
- Normative data range displays for both compensated and uncompensated tymps for 1000 Hz probe tone.
- Set baseline configuration to Auto, Positive Tail, or Negative Tail.



16

Multi-Hz Tympanometry



17



Multi-Hz Test: Resonant Frequency of the ME

- The frequency at which the middle ear system is the most efficient for energy transfer
- The frequency at which the springy and mass components have equal contribution



18



Resonance Frequency

- Colletti (1975) showed that as the probe frequency increased the tympanometric shape changed in a predictable and systematic way
 - Moved from the familiar single peak shape to slightly notched, and ending up a deep V shape at higher frequencies
- The frequency where notching appears indicates the Resonant frequency
- Normal Ranges for Resonance Frequency
 - 600Hz to 1340 Hz (Colletti 1977)
 - 800Hz to 1200 Hz (Shanks 1984)



19



Resonance Frequency

- Low Frequency Probe Tone (226Hz)
 - Tympanogram will be \wedge
- Mid Frequency Probe Tone (678Hz & 1000Hz)
 - All tympanograms will have multiple peaks - $\wedge\wedge$, **inverted W**
- High Frequency Probe Tone (above resonance frequency)
 - Tympanograms will continue to evolve to a V shape – opposite of the 226Hz



20



Multi Hz Testing – TympStar V2 (Classic)

- First – Probe tone sweep at +200 daPa
- 250Hz to 2000Hz in 50Hz steps
- Then a Y, 226Hz tymp is taken to identify TPP
- A second probe tone sweep is performed at TPP
- Component (Y, B, G) and phase values are calculated for +200 and TPP sweeps and subtracted from each other to obtain the 'Delta B' curve

Test Time: 60-90 seconds plus; multiple button presses



21




Multi Hz Testing TympStar v2 Classic

Delta B curve (top) displays the total susceptance as it changes With probe tone frequency 250-2000Hz.

Resonance occurs when curve crosses 0 mmhos.

As frequency increases, the total susceptance progresses from positive values (stiffness controlled) toward zero (resonance) to negative value (mass controlled).



22



Multi-Hz Testing TymStar Pro

- Measurements and calculations are exactly the same as the Classic.
- A single sweep using a wideband 'probe tone' for the pressure range selected—default +200 to -400 daPa
- All calculations are performed from the single sweep results
- Automatically displays 226 Hz Y tymph, Delta B graphs, resonant frequency tymph.
- Can store only 1 Multi-Hz test per ear




23



Multi Hz TymStar Pro


Press start

1. Will run a tympanogram with the wideband stimulus
 2. Will analyze data
 - Test in Progress
 - Processing Data at Pressure
 - Processing Frequency
 - Completing Test
- Total time: ~25 seconds







24

Multi-Hz Results



Other Tabs





27



Multi Hz with TympStar Pro

- 5x faster – complete test in 13 seconds.
- Automatically displays 226 Hz tymp, Delta B, tymp at resonance frequency
- Allows you to scroll through the tymps from 250 Hz to 2000 Hz to verify resonance frequency.



28



Diagnostic Use of MultiHz

- Normal resonance frequency range: 800 to 1200 Hz
 - Resonance is point where starting admittance value equals the “tympanic peak” admittance value
- Disarticulation: decreases resonance Hz below 600 Hz
- Otosclerosis: increases the resonance Hz above 1300 Hz
- Documentation pre- and post- surgery—resonant freq should decrease post-stapedectomy
- Confirm success of surgery for physician and patient



29



Multi-Hz

- Multi Hz test type will be included with all TymStar Pros
- Upgrades Require:
 - Firmware Update
 - Probe Calibration*



30

gsi
Grason-Stadler

MultiHz

SYSTEM

Left Right Ear

mmho

Time 226 Hz

1.5

1.0

0.5

0.0

-0.5

-400 -300 -200 -100 0 100 200

daPa

Ready

0.00 mmho Admittance

-0.5 0.0 0.5 1.0 1.5

Pressure

-100 0 100 200

ip - 226 Hz

Peak

daPa mmho

Mode Not Calibrated

The system has not been calibrated for MultiHz testing. Contact your GSI service representative to calibrate.

Auto Start

Immittance Y

Baseline

Sweep Rate 50

Age Adult

7:57 AM
10/30/2018




31

gsi
Grason-Stadler

Test Results

- Print to PDF
- Print to MPT III thermal printer
- Print to LPQ80 printer
- Print to any compatible HP printer

NOTE: this test will be compatible with GSI Suite in the next release



32


Patulous ETF Test



33

What is a Patulous ET?

- Persistently open Eustachian Tube
- Symptoms
 - Aural fullness or plugged up feeling
 - Own voice is loud
 - Crackling sound when chewing
 - Loud breathing
 - Tinnitus
 - Hyperacusis



34



What causes a Patulous ET?

- Significant weight loss
- Hormonal changes--pregnancy
- Medications
- Scarring after surgery (adenoidectomy)
- Radiation
- Neurologic disorders causing muscle atrophy
- Stress



35

Patulous ETF on the TymStar Pro


- Measure admittance at peak pressure over time
- No stimulus
- Observe tracing for rhythmic movement
- Press Start to pressurize the ear
- Press Present to start the recording





Patulous ETF on the TympStar Pro

- Adjustable time base
- Customize labels for test conditions using Config App




37

Test Results-Normal

- Can save multiple tracings for the same test condition.
- Can save up to 8 tracings.




Zoom



39

Test Results Abnormal

- Rhythmic fluctuation of compliance (mmho)
- Fluctuation of greater than .07 mmho in any test condition is considered abnormal



40

References

J Am Acad Audiol 22:201–207 (2011)

Use of Middle Ear Immittance Testing in the Evaluation of Patulous Eustachian Tube

DOI: 10.3766/jaaa.22.4.2

Andrew P. McGrath*†
Elias M. Michaelides‡

J Am Acad Audiol 4: 53–57 (1993)

Patulous Eustachian Tube Identification Using Tympanometry

David F. Henry*
Joseph R. DiBartolomeo†



Test Results

- Print to PDF
- Print to MPT III thermal printer
- Print to LPQ80 printer
- Print to any compatible HP printer

NOTE: this test will be compatible with GSI Suite in the next release





Seek Below Start Level



43




Reflex Threshold Setting

- Applies when using Threshold Seek for reflexes.
- Seek Below Start level = On:
 - If a reflex is detected at the initial Start Level of Threshold Seek, the intensity level will decrease until no reflex is detected and then increase to find threshold.



44



Examples

Start Level: 100
Threshold level: .02
Verification: Repeat



Reflex

Right Ear

Reflex Results

500 Hz
Ipsi

1k Hz
Ipsi

2k Hz
Ipsi

Time
1.6 Sec

Admittance
-0.035 mmho

Mark Threshold

Mark No Response

Set Peak

Find Peak

Immittance ☒

Erase Selected

Close

Display Outcor

Auto Zero Probe Stim. Ear Threshold Scale

9:38 AM
12/21/2018



 **gsi**
Grason-Stadler

The screenshot displays the Reflex software interface for the left ear. The main graph shows a red line representing the reflex response over time. The x-axis is labeled 'Sec' and ranges from 0 to 20. The y-axis is labeled 'mmho' and ranges from -0.16 to 0.04. A yellow box highlights a region of the graph between approximately 14 and 18 seconds, with a peak value of 0.063 mmho. The interface includes a menu on the left with options for 500 Hz, 1k Hz, 2k Hz, and 500 Hz Contra. The control panel on the right includes buttons for Mark Threshold, Mark No Response, Set Peak, Find Peak, Immitance, Erase, and Close. The status bar at the bottom shows the current time as 1:33 PM on 12/1/2018.

Time	Amplitude (mmho)
0.021	-0.010
0.028	-0.015
0.010	-0.010
0.003	-0.010
0.012	-0.010
0.063	0.063
0.061	-0.010



continued®




Baseline Configuration Options



49

What is it?





3 Options

1. **Positive Tail (+200):** V_{ea} /Admittance for all probe tones is calculated at the positive tail.
2. **Negative Tail (-400):** V_{ea} /Admittance for all probe tones is calculated at the negative tail.
3. **Auto (default setting/recommended):** V_{ea} /Admittance for a 226Hz probe tone tympanogram is calculated at the positive tail and for a 1000Hz probe tone tympanogram is calculated from the negative tail.




51



Baseline Configuration

- Negative-tail compensated static admittance was selected because of the larger mean value on the negative tail side.
- The larger value makes it easier to distinguish between a normal tympanogram and one that is flat or nearly flat.




52



Baseline Configuration

- Negative Tail Calculation for 1000Hz
- Implemented so 1000 Hz norms from Margolis paper display correctly.



53



Config App Changes



54



Customizing TympStar Pro Configuration Application

- Set preferences
- Define Diagnostic and screening protocols
- Set Auto Sequences




55



Custom Protocols

- Adjust defaults
- Create protocols for specific patient populations
- Available for Screener, Diagnostic and Auto Sequence



56

gsi
Grason-Stadler

57

gsi
Grason-Stadler

Patulous ET Configuration

58



Seek Below Start

Stimulus Threshold Scaling

Threshold Seek: Minimum Change: 0.02 mmho Verification: **Allow Seek Below Start Intensity:**

dB Step: Start Level: 75 dB Min: 35 Max: 110 Stop Level: 100 dB Min: 35 Max: 110

Non-Threshold Seek Settings

dB Step: Start Level: 75 dB Min: 35 Max: 110



59



Baseline Configuration

Tympanometry

Probe Tone: 226 Hz 678 Hz 1000 Hz Immittance Component: Baseline: Tym Measure: Autostart: Peak Indicator: Sweep Quantity: 1 Min: 1 Max: 3 **Baseline Configuration:**


Sweep Rate: 600/200 daPa/s 200 daPa/s 50 daPa/s 12.5 daPa/s Pressure Range: Start Pressure: +400 daPa +300 daPa **+200 daPa** 0 daPa -200 daPa -300 daPa -400 daPa -500 daPa -600 daPa Normative Data:



60



Customizing TympStar Pro – Auto Sequence



61



Summary

- New TympStar Pro firmware v1.2 is released
 - New Multi-Hz
 - New Patulous ET
 - New Seek Below Start
 - New Baseline Configuration Options
 - Multi-Hz and Patulous ETF are included with all TSPs
 - Multi-Hz and Patulous ETF not yet implemented in GSI Suite
 - Probe must be calibrated for Multi-Hz
- Config App Updates – must use compatible config app for changes
- Advance Quick Videos



62



Questions?



63

Thank You!!!

lrpr@grason-stadler.com



64