

Real Ear Measurements Made Easy with TargetMatch and Phonak Marvel

A Sonova brand

PHONAK
life is on



PHONAK
life is on

Kailen Berry, Au.D.

Kailen joined Phonak in 2012 as a field-based representative. Since that time she worked as an Audiology Manager with Phonak headquarters in Illinois helping to create and deliver educational content for product launches. She is now working as a Clinical Trainer in the Pacific Northwest and fulfilling her true passion, providing meaningful education to providers and patients.

Financial – Phonak employee who receives a salary for employment

Non-financial – No relevant non-financial relationship exists

Presented in Partnership with

continued

Learner Outcomes




As a result of this activity, participants will be able to:

- Explain how to successfully use TargetMatch in a clinical environment.
- Discuss the importance of real ear measurement in clinical outcomes and patient satisfaction.
- Explain how TargetMatch works and is integrated into Target.

Why is verification so pivotal to the process?



- Individual ear canal variations are accounted for – patient ears are not 2cc couplers
- Actual instrument performance is ensured as compared to pre-calculated approximation
- The goals based on the prescriptive formula or algorithm are realized – which should blend both audibility and patient acceptance
- Provides objective guidance for further ‘fine-tuning’ rather than relying on the challenging task of interpreting perception into action



Presented in Partnership with



Measures of unaided threshold sensitivity

PHONAK
life is on


Do audiologic measures of pure tones, word recognition in quiet, etc. adequately predict success?

Research in this area (Gatehouse, Humes, Cox, Walden, Abrams, etc.) says: **No**


Why is this?

Variability in patient performance

PHONAK
life is on



Speech in quiet (% correct) vs pure tone average



Speech in noise (SNR) vs pure tone average

- Mild to Moderate
- Moderately Severe to Severe


Souza, P. (2009). Severe Hearing Loss Recommendations for Fitting Amplification. *Audiology Online*, January 19

Presented in Partnership with

continued

Variability in product performance

PHONAK
life is on




Difference between NAL-NL2 (0 dB line) and First-Fit REAR for 5 premium hearing aids @ 55, 65, and 75 dB SPL

Sanders, Stoodly, Weber and Mueller (2015)

Optimizing performance in noise

PHONAK
life is on

Mean SNR HL improvement of 6.6 dB following programming



Leavitt and Flexer (2012)

Presented in Partnership with

continued

Average Patient Visits (during first year) Influenced by Methods Utilized in Clinic

- 1 more visit per patient without verification or validation
- One practitioner sees 15 new patients a month
- That's 180 additional visits/year for that one practitioner!





FIGURE 2. Average patient visits to fit hearing aids based on use of verification (REM) & validation procedures (objective or subjective).

Source: Kochkin S. "MarkeTrak VIII: Reducing patient visits through verification and validation." *Hearing Review*. 2011;18(6):10-12.


How do REMs impact perceived value?

PHONAK
life is on


Quality Value




Emotional Value



Behavioral Intent



Willingness to Pay




Source: Amlani AM, Pumford J, Gessling E. Improving Patient Perception of Clinical Services Through Real-ear Measurements. *Hearing Review*. 2016;23(12):12.

Presented in Partnership with

continued

The Many 'REM Effects'

PHONAK
life is on



Why the gap between the evidence and routine practice?

PHONAK
life is on

- There are similarities between the reasons for patient resistance to treatment and clinicians resistance to REM
 - Cost, time, and complexity
 - Anecdotal evidence from listening to the experience of others
 - Past experience with products that failed to live up to expectations
- How do we work to gain trust and give them the experience of realizing what they are missing?
 - Blend evidence and patient stories
 - Work to ensure smooth transition
 - Align the efficiency of the process with the value that we know is there

Presented in Partnership with

continued

Closing the gap

- Efficiency and accuracy need to compliment value
- Inter Module Communication allows fitting and verification tools to exchange data as well as device control
- Experts in amplification and instrumentation have blended knowledge to bring even more value to providers and patients



TargetMatch

Presented in Partnership with

continued

TargetMatch – What is it?



Fully integrated verification system

- Seamless, step by step workflow to guide the HCP through verification
- Real ear data (REUG, REOG, MLE)
- Automatic match to targets
- Available for the first time in Phonak Target 4.3



TargetMatch – Why should I use it?



(Aazh & Moore, 2007; Amanli, 2016; ASHA, 2015; Bell, 2009; Cox, 2005; Jorgensen, 2014; Kochkin, 2011; Mueller, 2001 & 2005)

Presented in Partnership with



TargetMatch – Requirements

PHONAK
life is on

Requires:


- Fitting computer
- Phonak Target 4.3.1 or higher
- Otometrics Otosuite
 - Most current version
- AURICAL FreeFit Hardware
 - TargetMatch tab will appear once hardware is plugged in and turned on



TargetMatch - Benefits

PHONAK
life is on

- Guided Probe tube placement
- Fully controls AURICAL FreeFit
- Reproducible measurement results
- Automatic use of acoustic transformations
- Automatic or manual target matching
- Available for NAL, DSL, and Adaptive Phonak Digital fitting formulas




Presented in Partnership with

continued


TargetMatch – Where is it?

- Access through [Fitting] > [Global tuning] > [TargetMatch]



TargetMatch – Preparation

1. [REM]
2. [Measure new RECD]
3. Add [REUG measurement] (optional)




Presented in Partnership with

continued

TargetMatch supports verification of Adaptive Phonak Digital

- APD takes into account MLE, RECD, and REUG
 - RECD estimation within Phonak Target
- Gain and compression settings based on hearing loss configuration
- TargetMatch will also support verification of all prescriptive formulae found in Target software




RESULT
An individualized
fitting precalculation

TargetMatch: Guided probe tube placement

PHONAK
life is on

- On-screen instructions for each step to position the probe tube
- Use On/Off button on FreeFit device to proceed






Presented in Partnership with

continued

TargetMatch: Guided probe tube placement

PHONAK
life is on


- Move probe tube towards the eardrum as indicated
 - Values are in millimeters
 - Fit to the scale on probe tubes
- No reliable detection 
- Correct position 



TargetMatch: Real ear measurements

PHONAK
life is on

- Automatic measurements
 - RECD
 - REOG
 - Microphone check/MLE



Presented in Partnership with

continued

Why are these individual ear transformations important?

PHONAK
life is on



Acoustic transformations


- REUG
- RECD
- REOG
- MLE



TargetMatch: Measure & match targets

PHONAK
life is on

- Select **[Active program]**
- Measure & match targets
 - Automatic match
 - Manual match
 - Measure responses only




Presented in Partnership with

continued

TargetMatch: Automatic fitting

PHONAK
life is on

- Automatic application of RECD in the test box
- Response measurements
 - (50 dB, 65 dB, 80 dB)




TargetMatch: Automatic fitting

PHONAK
life is on

- Measured responses are automatically matched to the prescriptive target

— Prescribed target
- - - Individual target
— Measured response




Presented in Partnership with

continued

TargetMatch: Save results

- Apply changes to the fitting session
- Data storage in Noah database
 - All measurements
 - Notes



TargetMatch for 2cc test box

- Head control/body positioning
- Cooperation
- Environmental noise
- Lack of client feedback
- Fine-tuning without interaction or availability of the client

Presented in Partnership with

continued

TargetMatch

Real ear measures

Fully integrated workflow

- Based upon best practice guidelines
- Automatic match to targets with accuracy and precision *

2cc/test box


Enhanced workflow to include 2cc/test box measures

- Addresses a wider range of verification needs
- RECD for all couplings-including RIC, open, and vented fittings
- Automatic match to targets

Denys, S et al. (Sept. 2017) Real-Ear Measurements Integrated in the Fitting Software: Test-Retest Reliability, Matching Precision and Perceptual Outcomes. Hearing Review.

Evidence behind TargetMatch

- N = 10 hearing impaired participants (mean 44 years)
- Moderate/severe hearing loss bilaterally
- All experienced hearing aid wearers fit with Bolero V-70P



Denys, S et al. (Sept. 2017) Real-Ear Measurements Integrated in the Fitting Software: Test-Retest Reliability, Matching Precision and Perceptual Outcomes. Hearing Review.


Presented in Partnership with

continued

Evidence behind TargetMatch: Manual fit vs TargetMatch

PHONAK
life is on

- Targets are well matched within 3 dB across frequencies up to 4000 Hz
- *Statistically significant ($p < 0.05$) better target matching across most frequencies between 800-5000 Hz




Denys, S. et al. (Sept. 2017) Real-Ear Measurements Integrated in the Fitting Software: Test-Retest Reliability, Matching Precision and Perceptual Outcomes. Hearing Review.

Evidence behind TargetMatch: test/re-test reliability

PHONAK
life is on

- Test, re-test reliability is less than or equal to 2 dB for frequencies up to 4000 Hz across different acoustic couplings



Denys, S. et al. (Sept. 2017) Real-Ear Measurements Integrated in the Fitting Software: Test-Retest Reliability, Matching Precision and Perceptual Outcomes. Hearing Review.

Presented in Partnership with

continued

“The TargetMatch integrated/guided workflow **outperforms** traditional non-integrated REM fittings in terms of target matching and test, re-test reliability”

Denys, S (2016)-University of Leuven

Together,
we change lives

Presented in Partnership with

continued