

# ADHEAR Device Orientation Checklist



## Patient Details

Name:

Age:

Caregiver: Yes

No

## Fitting Details (enter details for side(s) fit with ADHEAR System)

	LEFT EAR	RIGHT EAR
Bilateral	Unilateral Left	Unilateral Right
Etiology of hearing loss	Conductive      SSD	Conductive      SSD
Audio Processor	ADHEAR Color	ADHEAR Color
	Serial #	Serial #
	Adapter Color	Adapter Color
	Accessory Cables	
	unilateral, long, audio unilateral, short, audio bilateral, long, audio bilateral, long, audio FM adapter cable angled adapter	
Date of fitting (mm/dd/yyyy)		
Audio Processor settings	VC Enabled? Yes      No	VC Enabled? Yes      No
	No. of Active Programs	No. of Active Programs
Aided testing comments		
Questionnaires administered		

## Patient Details

Name:

Age:

Caregiver: Yes

No

## Patient Orientation

Topic	Specifics
Battery	<p>Open and close door</p> <p>Change battery</p> <p>Understands batteries are toxic; keep out of reach of children, do not swallow, safe disposal</p>
Adhesive Adapter	<p>Correct placement</p> <p>Apply on clean, dry, hairless skin</p> <p>Apply without touching adhesive portion</p> <p>Apply light pressure for 30-60 seconds after placement</p> <p>Wait up to 15 minutes before attaching ADHEAR to adhesive</p>
Audio Processor	<p>Demonstrates correct attachment of audio processor to adhesive adapter</p> <p>Demonstrates correct detachment of audio processor from adhesive adapter</p> <p>Changes volume control (if enabled)</p> <p>Changes program (if enabled)</p> <p>Familiar with programs on device</p>
Security	Demonstrates use of retention clip
Sound Quality	Demonstrates ability to perform listening check (caregiver/teacher)
Care & Maintenance	Able to clean audio processor
Audio Sources	<p>Able to attach cable</p> <p>Able to select appropriate listening program for accessory cable</p>