

Learning Objectives

- After today's course, you will be able to:
 - Describe how SmartCompress technology can control compression in speech-in-noise environments and control gain in non-speech environments.
 - Manage acoustic feedback with the Adaptive Feedback Canceller Pro in Sonic Trek hearing aids.
 - List and describe all the available features and technology in Sonic's Trek Super Power and Ultra Power hearing aids.

www.sonici.com



The Hook of 4S Foundation

- Sonic is dedicated to improving lives through enhanced hearing by constantly focusing on the 4S Foundation
 - **Sound** that's natural
 - **Speech** understanding in noise
 - **Simplicity** in all we do
 - **Style** that stands out




www.sonici.com



Key Insight

- What are the market characteristics for these users?
 - Only an estimated 5-10% of all hearing impaired users are in the profound category*
 - Only an estimated 700,000 to 1.4 million power instruments are sold per year, based on 14 million units/year sold globally




(*Taylor, 2009)

www.sonici.com



Feature Benefit Value, Yet Again



www.sonici.com

- Higher market penetration than other segments
- Super Power | Ultra Power BTE users have a fundamental *need* for instrumentation

Image source: Niels Granholm-Leth, Carnegie Bank; WHO, Sonova via Hearing Review: HIA Comments on FTC "Now Hear This" Workshop and OTC Hearing Aids Published; May 24, 2017.



trekTM
Experience more



www.sonici.com





Trek. At-a-Glance

www.sonicl.com



See What's Inside

- Two performance levels
 - Trek⁸⁰ | Trek⁴⁰
- Two styles
 - Super Power (SP) | Ultra Power (UP)
- 2.4 GHz Bluetooth® Low Energy
- Near Field Magnetic Induction (NFMi)
- Telecoil
- LED Indicator
- Connectivity Solutions
- Pediatric friendly




Trek SP Trek UP

Made for
 iPhone | iPad | iPod


www.sonici.com

SONIC
 Everyday Sounds Better

TK 80 | 40 BTE SP Fitting Range




TK 80 | 40 BTE UP Fitting Range




	TK 80 40 BTE SP	TK 80 40 BTE UP
Battery Size	13	675
Peak Gain (dB/2cc)	79	83
Peak MPO (dB SPL/2cc)	139*	142*

*Special care should be taken when fitting and using a hearing instrument with maximum sound pressure capability in excess of 132 dB SPL (IEC 60318-4) since there may be a risk of impairing the remaining hearing of the hearing instrument user.



Sound That's Natural

www.sonici.com

 SONIC
Everyday Sounds Better

More Power, Control and Access to Everyday Sounds

- With Trek, Sonic delivers a natural sound with several technologies delivering power users need
 - Speech Variable Processing (SVP) with SmartCompress
 - Adaptive Feedback Canceller Pro (AFC Pro)
 - Frequency Transfer
 - Extended Dynamic Range
 - SmartMusic

www.sonici.com

 SONIC
Everyday Sounds Better

More Power for Optimized Performance



SoundDNA Platform

www.sonici.com

- **Phoneme Focus**
 - Offers maximum audibility of all the fine details in the speech signal rapidly adjusting the gain to apply the precise amount of amplification to each phoneme
- **Envelope Focus**
 - Designed to support patients who rely on information from the speech envelope applying less alteration to the amplification of the speech signal, preserving the envelope for greater contrast between phonemes



Intelligent Amplification in Noise



www.sonici.com

Configurable adaptive compression system for intelligent amplification in noise

Natural and detailed amplification for your patients – all with SmartCompress

Available in Trek⁸⁰



SmartCompress Delivers the Solution


- It's effective
 - Specifically addresses the problem of applying compression in noise
- It's adaptive
 - An adaptive compression system overcomes limitations of fixed Environment Classification system
- Works well with others 😊
 - Complements SVP, directionality and noise reduction, and any fitting rationale of your choice

www.sonici.com



SmartCompress - How does it do it?

- SmartCompress uses two SNR level estimators, compared to one in traditional systems
 - Green – fast phonemic estimator
 - Blue – slow, long-term estimator
- Provides accurate measurement of environmental changes in real time



www.sonici.com



SmartCompress

- Controls compression in speech-in-noise environments
 - Adaptively varies the compression in speech-in-noise environments
 - Determines how much the compression has to be decreased, or made more linear, for less noise in the output
- Controls gain in non-speech environments
 - Adaptively limits gain when speech is not present
 - Determines the occurrence of quiet and noise-only situations and applies less gain to the input

www.sonici.com




SONIC
Everyday Sounds Better

High SNR Speech in quiet



- SmartCompress does not activate!
- SVP applies the prescribed gain and compression as designated by the selected fitting rationale along with any fine tuning actions made by the hearing care professional

www.sonici.com





SONIC
Everyday Sounds Better

Very low SNR

Quiet and Noise-only environments (Without speech)

- SmartCompress via the **Gain Control** applies less gain to the incoming signal
- Offers greater listening comfort in non-speech listening environments






www.sonici.com

SONIC
Everyday Sounds Better

Low SNR

Speech-in-noise environments

- SmartCompress via the **Compression Control** reduces compression
- Minimizes the amplification of noise following short pauses in speech or even between phonemes
- Improves the output SNR
- Supports speech and listening comfort in complex environments



www.sonici.com

SONIC
Everyday Sounds Better

Better Target Match with Less Feedback Risk

Algorithm exclusive to 2019
Products on the SoundDNA
Platform


1 Feedback Management
system with 2 technologies



www.sonici.com



Robust Feedback Cancellation



- Standard system continually operates in background to quickly cancel in static conditions
- Second system quickly reacts to feedback-causing situations
- Speed allows more stable gain in the fitting offering higher feedback margins

www.sonici.com



Making Speech Cues Audible

Helps identify consonants for patients with severe high-frequency hearing losses

Access via Frequency Transfer control in EXPRESSfit® Pro



www.sonici.com



Frequency Transfer



www.sonici.com



- **10 Destination Regions**
 - Adjusts HF source and LF destination regions
- **7 Intensity Settings**
 - Adjusts the strength of the HF sounds (source) when mixed with the LF sounds (destination)
- **High Frequency Attenuation**
 - Turns gain in the HF bands above the destination on or off

Sounds, Loud and Clear



www.sonici.com



Advanced feature to keep loud sounds clear in challenging listening environment

Specifically designed for signals with high intensity

Extended Dynamic Range

- Listening environments with speech can get LOUD
 - Can quickly change in intensity
 - Presence of loud speech signals can have peaks that exceed > 95 dB A
- Digital hearing instruments limit inputs at 95 dB SPL
 - Creates distortion for higher input levels due to compression



www.sonici.com



Extended Dynamic Range - Summary

- Adaptively extends dynamic input range
- Expands to growing sound levels only as needed
 - Upper limit depends on input level of speech signal (from 95-113 dB SPL)
- Allows loud peaks to be amplified with less risk of distortion
- For everyday sounds
- Available in Trek⁸⁰ only

www.sonici.com



Live Music at its Best

Greater enjoyment of live music

Enthusiasts and musicians alike can enjoy clear sounding music again



www.sonici.com


SONIC
Everyday Sounds Better

4 Components to SmartMusic

- **Fixed Extended Dynamic Range**
 - Increases limit at the input to a fixed 113 dB SPL
- **Speech Variable Processing**
 - Maintains balance between harmonics
- **Wideband Frequency Response**
 - Important for perceived naturalness of music
- **Controls**
 - Fixed Hypercardioid polar setting
 - Position of null does not change
 - Eliminates noise from behind
 - Concentrate on what's on stage
- **Available in Trek⁸⁰**


www.sonici.com

SONIC
Everyday Sounds Better



Speech Understanding in Noise

www.sonici.com

 **SONIC**
Everyday Sounds Better

More Speech, Less Noise in More Places, and Attention to Detail

- With Trek, Sonic delivers what life sounds like in noise
 - SPiN Management
 - Binaural Noise Management
 - Impulse Noise Reduction
 - Soft Noise Reduction
 - Wind Noise Reduction
 - Tinnitus SoundSupport
 - VC Step Size

www.sonici.com

 **SONIC**
Everyday Sounds Better

More Speech, Less Noise

Speech in Noise (SPiN) Management introduces a trio of technologies

- SPiN Directionality
- SPiN Noise Reduction
- SPiN Engage



www.sonici.com



SPiN Management

- SPiN Directionality
 - Uses null steering to select the microphones' polar patterns that will produce the optimal SNR in 16 independent frequency bands
- SPiN Noise Reduction
 - Designed to reduce the gain of noise compared to speech operating in 16 independent frequency bands
- SPiN Engage
 - **Trek⁸⁰ only!**
 - Determines at which SNR the directional and noise reduction systems will adjust to the environment.

www.sonici.com



Less Noise in More Places



Support your patients who encounter situations with distracting, uneven noise from side to side with Binaural Noise Management

Available in Trek⁸⁰

www.sonici.com



Nobody Wants Noise

- Impulse Noise Reduction
 - Lessens the discomfort of sudden, loud sounds
- Soft Noise Reduction
 - Reduces those low-level noises without affecting speech
- Wind Noise Reduction
 - Adjusts to outside activities by monitoring the environment and preventing the rush of wind from being amplified

www.sonici.com



The Sounds of Relief



Aims to reduce your patient's perception of tinnitus by providing amplification and generating sound relief

Available in all Trek offerings

www.sonici.com



Tinnitus SoundSupport

- Offers the following tinnitus relief sound options:

- 4 broadband sounds

White Noise

Pink Noise

Red Noise

Shaped to
audiogram

- 3 nature-like sounds

Ocean 1

Ocean 2


Ocean 3

www.sonici.com



Control Tinnitus Sounds - App

- Sonic SoundLink 2 App
 - Adjust tinnitus sound volume or mute entirely
 - Modify tinnitus sound by changing its modulation rate or frequency response



www.sonici.com




Attention to Detail – VC Step Size

- Aims to help hearing aid wearers who have a reduced dynamic range
- Designed as a control for loudness growth issues, VC Step size makes finer volume adjustments in 1dB instead of 2.5 dB increments



www.sonici.com





Welcome to Wireless, Simple Essentials and Tools to Accomplish it All

- Trek delivers solutions that just make sense
 - Dual-Radio System
 - SoundLink 2 App
 - IFTTT
 - Low Frequency Enhancement
 - Binaural Coordination
 - IP Rating
 - EXPRESSfit® Pro

www.sonici.com

SONIC
Everyday Sounds Better

Welcome to Wireless



www.sonici.com




Wireless connectivity is a way of life for today's hearing aid users

Trek offers advanced wireless systems for their world, their way

Dual-Radio Wireless Technology

Two wireless systems integrated into one hearing aid

a) NFMI short distance ear-to-ear




b) 2.4 GHz long distance direct-to-hearing aid



www.sonici.com



SoundDNA Connectivity Portfolio



SoundLink 2 App

- Controls wireless operations like volume and program change and more
- Compatible with an iPhone® or Android™ smartphone¹
- Many conveniences such as “Find My Hearing Aid,” advanced Tinnitus SoundSupport controls, battery life percentage indicators, etc.



* Trek is a Made for iPhone® hearing aid that allows you to talk on the phone, stream stereo music, and stream video sound directly to your hearing aid. In addition, download the free SoundLink 2 app to control your hearing aid over your iPhone, iPad®, iPod touch®, and Android™ smartphones and tablets. For information on compatibility, please visit www.sonici.com/connectivity.

www.sonici.com

SONIC
Everyday Sounds Better

Make the Most of Smart Technology




www.sonici.com



If This Then That brings the Internet of Things to Trek

More than 600 internet-connected devices and services

How Trek is connected with IFTTT



Deep Bass Dynamics

Add low frequencies into direct audio streaming for a high fidelity signal




www.sonici.com

SONIC
Everyday Sounds Better

Low Frequency Enhancement

Available in Trek 80|40



www.sonici.us

SONIC
Everyday Sounds Better

Simple Essentials

- **Binaural Coordination**
 - Shares information between left and right ears
- **Binaural Synchronization**
 - Volume or program change made on one device is simultaneously made on the other
 - Push button mute with one touch
- **Non-Telephone Ear Control**
 - Reduce gain or mute the input on the opposite device
 - Reduce distractions caused by non-telephone sounds and increase perception of volume for the telephone

www.sonici.com



Personalized Sight-and-Sound

- **Multicolor LED and Audible Indicators**
 - On
 - Program
 - Volume
 - Mute
 - Accessories
 - Flight Mode
 - Battery Warning
 - Microphone Service Check



www.sonici.com



Impressive Protection



Trek stands up to everyday exercise, spontaneous splashes, and the imperfections of day-to-day life

www.sonici.com



More Protection than Ever

- Compared to previous devices in this segment
- Ingress Protection (IP) rating
- IP68 = Dust-tight and protection against continuous immersion in three feet of water or more.*

* Although highly protected, Trek is not intended for swimming and showering. Encourage patients to use commonsense measures to preserve their electronic devices.




www.sonici.com



Tools of the Trade

- Fitting Rationales
 - Default DSL 5
- Adaptation Manager
 - Gradually increase gain
- Data Logging
 - Records information about patient's listening behaviors
- Transfer Fitting
 - Transfer settings from one instrument to another
- Real Ear Fit
 - Uses Inter Module Communication 2 - IMC2 to communicate to REM system
- Noahlink Wireless
 - Best in class wireless programmer
- SoundStudio
 - Experience the sound of Trek right in your office

www.sonici.com



Style that
Stands Out

www.sonici.com



Trek. Stylish by Design.



- 1. Single Push Button**
- » Activates Program change and Flight Mode
 - » Provides easy access for up to four listening environments
 - » Configurable alerting tones provide audible cue as program is changed

- 2. Double Push Button**
- » Activates Volume Control and Mute
 - » More sensitive 1 dB volume step size
 - » Configurable alerting tones provide audible cue as volume is changed

- 3. Multicolor LED Indicator**
- » Allows visual confirmation of Program and Volume change, Battery status, etc.
 - » Useful for caregivers, parents, or teachers to check/confirm status of operation

- 4. Dual Covered Microphones**
- » Allow for sophisticated directional features to enhance performance in noise
 - » Protection from debris and moisture

- 5. Battery Compartment**
- » Easy-open design
 - » Integrated on/off switch
 - » Size 13 (SP) or 675 (UP) battery
 - » Tamper-resistant lock keeps the battery secure from children

- 6. DAI/FM Compatible**
- » Provides connectivity for Direct Audio Input or FM audio sources
 - » DAI adapter 1000
 - » FM adapter 10

- 7. Acoustic Options**
- » Ships with undamped earhook
 - » Damping element available
 - » Pediatric earhook available
 - » Custom molds available in a variety of styles and configurations

- 8. Telecoil**
- » Provides clear sound quality for use with telephones or in looped environments

- 9. Non-Telephone Ear Control**
- » Reduces gain or mutes the input from the non-telephone ear

- 10. Wireless Connectivity**
- » Near Field Magnetic Induction enables binaural communication between devices
 - » 2.4 GHz provides connectivity to external audio sources (mobile phone, television, computer, etc.)

- 11. If This Then That (IFTTT)**
- » Connects Trek instruments to other internet-connected smart devices via Sonic Go service

- 12. Ingress Protection Rating**
- » Hydrophobic coating and overall design protects against dust and continuous immersion in three feet of water or more
 - » IP68

Model Overview



	SP	UP
Battery Size	13	675
Power Model	■	■
Directional Microphones	■	■
Program Button	■	■
Volume Control	■	■
Telecoil	■	■
Wireless Connectivity	■	■
Earhook	■	■
IP Rating	IP68	IP68
Tinnitus SoundSupport	■	■
IFTTT Internet Connection*	■	■
DAI/FM	■	■
LED	■	■

■ Standard

* Requires SoundLink 2 App.



Pediatric Care Kit

Updated Pediatric Kit includes a variety of accessories for the caregiver and child:

- Bag
- Drying and cleaning set
- Cleaning cloth
- Listening tube
- Hand puppet Lion
- Safety cord
- Care guide booklet with instructions for caregivers



www.sonici.com

SONIC
Everyday Sounds Better

EduMic

- Trek supports all wireless accessories available for the SoundDNA product families **as well as** a new remote microphone – called EduMic
- EduMic is designed for use in modern classrooms
- EduMic offers a transmission range for classroom coverage, one-to-many communication and has a long-lasting battery life to last an entire day in school.



www.sonici.com



Reasons to Believe

www.sonici.com



Trek's Take Home Message

- Increased power with less audible feedback helps patients hear sounds they were previously unable to hear
- So patients can engage in more conversations, connections, and experiences
- And with advanced connectivity, reliability, and style, Trek gives users the confidence they need to experience more of everyday life

www.sonici.com




Thank you!

- As always, thank you for joining me for today's course,
Trek: Experience More
- Any Questions?
 - 888.423.7834
 - support@sonici.com
 - www.sonici.com




www.sonici.com



The logo for the SoundDNA Platform features a central circular emblem with a stylized DNA double helix. This emblem is surrounded by a semi-circular arc of 18 small, white icons representing various audio and technology concepts, such as a speaker, a microphone, a musical note, a bar chart, a signal wave, a gear, a magnifying glass, a document, a person, a house, a car, a plane, a ship, a train, a bus, a truck, a boat, and a cell tower. Below the emblem, the text "SoundDNA Platform" is written in a white, sans-serif font.

SoundDNA Platform

www.sonici.com

 **SONIC**
Everyday Sounds Better

The SONIC logo consists of a stylized graphic of three curved lines on the left, resembling sound waves or a speaker's output. To the right of this graphic, the word "SONIC" is written in a large, bold, white, sans-serif font. Below "SONIC", the tagline "Everyday Sounds Better" is written in a smaller, white, sans-serif font.

 **SONIC**
Everyday Sounds Better