

## Dual Radio: Stay Connected

www.sonici.com



1


## Learning Outcomes

- After this course, learners will be able to:
  - Describe how the Dual-Radio System uses 2 separate wireless technologies and the purpose of each
  - Pair compatible Sonic hearing aids containing Dual Radio technology to a variety of wireless products
  - Explain how Sonic hearing aids with Dual Radio technology can access the Internet of Things via a smartphone

www.sonici.com



2




# Advanced Wireless Capabilities with Dual-Radio System

[www.sonicl.com](http://www.sonicl.com)

**SONIC**  
Everyday Sounds Better

3




Fast and direct wireless transmissions with optimized battery life

**Dual-Radio System**

4

## Benefits of New Chip in SoundDNA



www.sonici.com

**SONIC**  
 Everyday Sounds Better


5

## Dual-Radio Wireless Technology

Two wireless systems integrated into one hearing aid

NFMI short distance ear-to-ear

2.4 GHz long distance direct-to-hearing aid




www.sonici.com

**SONIC**  
 Everyday Sounds Better

6

## Ear-to-Ear Wireless Transmission

- Near Field Magnetic Induction signal (NFMI)
- Transmit only over a short distance
- Higher rate of information exchange




www.sonici.com

**SONIC**  
Everyday Sounds Better

7

## Direct-to-Hearing Aid Wireless Transmission

- 2.4 GHz without an intermediary device
- Transmit over a long distance
- “Line-of-sight” connection to avoid interference



www.sonici.com

**SONIC**  
Everyday Sounds Better

8



# Connectivity Portfolio

The diagram illustrates the Connectivity Portfolio, centered around a smartphone displaying the SoundLink 2 app. A jagged starburst shape radiates from the phone, connecting to various products and services:

- Phone Adapter 2**: A small black rectangular adapter.
- TV-A**: A black rectangular adapter for TV connectivity.
- Remote RC-A**: A black remote control with three buttons.
- Wireless Programming**: A black wireless programming device and a USB dongle.
- Sound Clip-A**: A small black clip-on device.
- Internet Connection IFTTT**: A logo for IFTTT (If This Then That) for internet connectivity.
- iPhone® direct streaming**: Text indicating direct streaming capability for iPhones.
- SoundLink 2 app for Apple® and Android™**: Text describing the app available on both platforms.
- Download on the App Store** and **GET IT ON Google Play**: Badges for app availability.
- Made for iPhone | iPad | iPod**: A logo indicating compatibility with Apple devices.

9

## Technology Levels

### Five Wireless Technology Levels

100 Premium	80 Advanced	60 Mid-Level	40 Basic	20 Entry
----------------	----------------	-----------------	-------------	-------------

**captivate<sup>100</sup>**  
**enchant<sup>100</sup>**

**captivate<sup>80</sup>**  
**enchant<sup>80</sup>**  
**trek<sup>80</sup>**

**captivate<sup>60</sup>**  
**enchant<sup>60</sup>**

**enchant<sup>40</sup>**  
**trek<sup>40</sup>**

**enchant<sup>20</sup>**

[www.sonicl.com](http://www.sonicl.com)

**SONIC**  
Everyday Sounds Better

10

## Captivate Styles

Captivate offers an array of wireless  
behind-the-ear styles



2.4 GHz & NFMI Standard

www.sonici.com

**SONIC**  
Everyday Sounds Better

11

## Enchant Styles

Enchant offers a full line of wireless styles\*



\*NOTE: 2.4 GHz not available in CIC (NFMI only)

www.sonici.com

**SONIC**  
Everyday Sounds Better

12

# Trek Styles

Trek offers robust power with wireless capabilities



2.4 GHz & NFMI Standard

[www.sonici.com](http://www.sonici.com)

**SONIC**  
Everyday Sounds Better

13

## Connectivity with SoundDNA

Streaming Phone Calls

[www.sonici.com](http://www.sonici.com)

**SONIC**  
Everyday Sounds Better

14

## Made For iPhone Hearing Aid

- Streaming of phone calls, music and videos from iPhone®, iPad®, and iPod®.



32 ft | 10 m transmission range

[www.sonici.com](http://www.sonici.com)

**SONIC**  
Everyday Sounds Better

15

## Made For iPhone Hearing Aid


- 3 Transmission Protocols
  - Standard Bluetooth® Low Energy (BLE) protocol
  - 2.4 GHz proprietary long-distance protocol
  - Apple Bluetooth® Low Energy protocol (Made for iPhone®)

[www.sonici.com](http://www.sonici.com)

**SONIC**  
Everyday Sounds Better

16

# Three Bluetooth® Protocols



www.sonicl.com

SONIC  
Everyday Sounds Better

17

## Standard Bluetooth® Low Energy (BLE) Protocol




18

## 2.4 GHz Proprietary Long Distance Protocol

Data and Sound  
→

Data  
→



TV-A

50 ft | 15 m transmission range

[www.sonici.com](http://www.sonici.com)

**SONIC**  
Everyday Sounds Better

19

## Apple Bluetooth® Low Energy Protocol

Data and Sound  
→

Data  
→



Made for  
**iPhone | iPad | iPod**



32 ft | 10 m transmission range

[www.sonici.com](http://www.sonici.com)

**SONIC**  
Everyday Sounds Better

20



## Android™ Smartphones and Tablets

- Volume and program changes via the SoundLink 2 app
- No streaming of sound or music
- Enter the SoundClip-A



www.sonicl.com

SONIC  
Everyday Sounds Better

21

## SoundClip-A

Connect to Bluetooth devices

www.sonicl.com

SONIC  
Everyday Sounds Better

22

# SoundClip-A


- SoundClip-A is a wireless device for audio streaming from any Bluetooth 2.1 device to compatible hearing aids
- Remote mic functionality
- Basic remote control functions
- No lanyard required

www.sonici.com



23


## SoundClip-A Overview



www.sonici.com



24



25

## TV-A

- 2.4 GHz direct streaming to hearing aids
- Unlimited number of pairings
- Streaming to multiple instruments at the same time
- Digital stereo sound supported


www.sonici.com

SONIC  
Everyday Sounds Better

26

## TV-A Direct Streaming

- Not compatible with products supporting SoundGate
- 50 ft | 15 m transmission range
- Regulate the volume of the streamed signal with
  - RC-A or SoundClip-A
  - SoundLink 2
  - Ear Level Volume Control




www.sonici.com

**SONIC**  
Everyday Sounds Better

27

## TV-A Streams in Stereo



www.sonici.com

**SONIC**  
Everyday Sounds Better

28

## TV Program

- Available once device is paired to the TV-A
- Access the TV Program to hear streamed signal
  - Local push button
  - RC-A
  - SoundLink 2 app
- Program cycle depends on the programs already available in your hearing aids
  - Dedicated TV Program indicator
- Mute hearing aid microphones while TV Program is active

[www.sonici.com](http://www.sonici.com)



29

## RC-A

### Remote Control


[www.sonici.com](http://www.sonici.com)



30

## RC-A

- Operates using 2.4 GHz BLE
- Not compatible with products supporting SoundGate



[www.sonici.com](http://www.sonici.com)



31

## RC-A

- RC-A will be printed on the back of the device
- 2x AAAA batteries
  - RC-A battery life is up to one year with normal use
- 6 ft / 1.8 m transmission range
- Keylock functionality available

[www.sonici.com](http://www.sonici.com)




32





33



34

## Non-Rechargeable

- Ensure batteries are fresh and hearing aids have been programmed together in fitting software
- Open and close battery doors and wait for start up indicator
- Hearing aids will be in pairing mode for approximately 3 minutes



www.sonicl.com

**SONIC**  
Everyday Sounds Better

35


## Rechargeable (Li Ion)

- Ensure batteries are charged and hearing aids have been programmed together in fitting software
- Two options for rebooting hearing aids
  - Place hearing aids in and out of charging unit
  - Press and hold push button (top or bottom) for 5 seconds to shut aids down, to reboot press and hold bottom push button for 5 seconds
- Hearing aids will be in pairing mode for approximately 3 minutes


www.sonicl.com

**SONIC**  
Everyday Sounds Better

36




37



38

## Unpairing - iPhone



Scroll to bottom of screen – ‘Forget this device’ – removed from both MFI settings and SoundLink 2 app


[www.sonici.com](http://www.sonici.com)



39

## Pairing- Android

**\*\*Devices are paired directly through the app\*\***



Once found device name will be listed- Touch on ‘Pair’ – Devices will then connect


[www.sonici.com](http://www.sonici.com)



40

# Unpairing- Android

**\*\*Devices are removed in Bluetooth Settings\*\***



[www.sonicl.com](http://www.sonicl.com)

**SONIC**  
Everyday Sounds Better

41


## SoundClip-A

[www.sonicl.com](http://www.sonicl.com)

**SONIC**  
Everyday Sounds Better

42

## Pairing Hearing Aids to SoundClip-A




www.sonici.com

**SONIC**  
Everyday Sounds Better

43

## Pairing SoundClip-A to Mobile Phone



www.sonici.com

**SONIC**  
Everyday Sounds Better

44





45

## Pairing with the TV-A

- Close proximity pairing (initial pairing)
- Fast and easy
- Make sure TV-A is powered on
- Place hearing aids into pairing mode and place on top of TV-A
  - Both power and TV indicators will be fading Blue
- Pairing is complete when power and TV indicators are blue for about 30 seconds
  - Power and TV indicators will be solid Blue for 30 seconds and then Power indicator will turn solid Red



www.sonici.com

46

## How to pair the TV-A



www.sonici.com


 **SONIC**  
Everyday Sounds Better

47

## Removal of Paired Instrument

- Must perform a system reset to remove all pairings
- Press and hold the Setup button for more than 5 seconds
- Indicators fade to blue to confirm that all pairings have been removed

TV-A	Description
 	No pairings

 Fading


www.sonici.com

 **SONIC**  
Everyday Sounds Better

48



49



50

## Pairing with the RC-A

- Close proximity pairing
- No need to pair in EXPRESSfit® Pro
- Place hearing aids into pairing mode and position 6-8 inches from RC-A
- Press and release any button to complete pairing
  - Status light will turn solid Green for 10 seconds



www.sonici.com

**SONIC**  
Everyday Sounds Better

51

## How to pair the RC-A



www.sonici.com

**SONIC**  
Everyday Sounds Better

52

## Removal of Paired Instrument

- RC-A can only be paired to one set of devices at a time
- Must clear RC-A if necessary to pair to different devices
- To clear pairing press and hold all three buttons for 5 seconds
  - Status light will turn solid Green for 10 seconds when cleared



www.sonici.com

**SONIC**  
Everyday Sounds Better

53


## How to unpair the RC-A



www.sonici.com

**SONIC**  
Everyday Sounds Better


54



# SoundLink 2 App

Android and iOS App

[www.sonici.com](http://www.sonici.com)




55

## SoundLink 2 Overview

Modify Tinnitus Sound Support program

Available for iOS and Android™ devices



Remote control functionality




Activate streaming from TV-A  
Change volume of the TV-A

Activate streaming from SC-A  
Change volume of the SC-A

Access the Internet of Things (IFTTT)

[www.sonici.com](http://www.sonici.com)




56




## Select Programs


Touch 'Remote Control' from Home Screen



Current program is highlighted



Select the desired program



www.sonici.com

**SONIC**  
Everyday Sounds Better

57

## Volume Options

- App will show the current program and current volume status
- Volume Increase
- Volume Decrease
- Mute
- Adjust together or independently




www.sonici.com

**SONIC**  
Everyday Sounds Better

58

## SoundClip-A



- Place near desired sound source
- Will transmit the sound directly to the hearing instruments
- Adjust volume of both SC-A and HA microphone


www.sonicl.com

**SONIC**  
Everyday Sounds Better

59

## Tinnitus Options

- Available in app once tinnitus program added within software
- Adjust volume of both tinnitus sound support and HA microphones
- Tinnitus options




www.sonicl.com

**SONIC**  
Everyday Sounds Better

60

## Tinnitus Options

- Advanced Settings
  - Noise variations
  - Noise equalizer
- Restore defaults
- Adjustments within the app are based on programming




www.sonici.com

**SONIC**  
Everyday Sounds Better

61

## Access IFTTT


- From Home screen select IF to access IFTTT
  - Connect to hundreds of internet based services
  - Create Applets
  - <https://ifttt.com/go>



www.sonici.com

**SONIC**  
Everyday Sounds Better

62




# Internet Connectivity

If This Then That (IFTTT)  
Internet Service

[www.sonicl.com](http://www.sonicl.com)

**SONIC**  
Everyday Sounds Better

63



Connect hearing devices to other  
devices via the Internet of Things

**Internet Connection**

64

## If This Then That (IFTTT)

- Internet-based service
- Hundreds of internet-connected devices and services
  - 'Internet of Things'
- Connect various services, apps and devices to create 'recipes'
- Every day a new company enters the IFTTT world, offering greater opportunities



www.sonici.com

**SONIC**  
Everyday Sounds Better

65

## Sonic devices and IFTTT


- All devices with 2.4 GHz
- Capable of connecting to various services and products
- Compatible iPhone®, iPad®, iPod® and Android™ devices
- IFTTT is **only** accessible via SoundLink 2 app on the smartphone or tablet
  - Without the app, device cannot be connected to the internet
  - Must set up an account within IFTTT
  - Smartphone/tablet needs to be connected to the internet

www.sonici.com

**SONIC**  
Everyday Sounds Better


66

## How are hearing aids connected with IFTTT?



67

## If This Then That



www.sonici.com

**SONIC**  
 Everyday Sounds Better

68

# Sonic Devices and GO Service

- GO service is the interface between Sonic devices and other services in IFTTT



69

## GO Service: Available Triggers and Actions

<https://ifttt.com/go>

GO Triggers	GO Actions
<ul style="list-style-type: none"> <li>Program has changed</li> <li>Battery is low</li> <li>Connection to Bluetooth® device is lost</li> <li>Connection to Bluetooth® device is established</li> <li>Mute state is changed</li> <li>Volume is changed</li> <li>Power is on</li> </ul>	<ul style="list-style-type: none"> <li>Set a program</li> <li>Increase the volume</li> <li>Decrease the volume</li> <li>Play the text to speech</li> <li>Set the volume</li> </ul>

www.sonici.com

**SONIC**  
Everyday Sounds Better

70

# Triggers and Actions

## Applet using a GO Trigger

- If the hearing aid batteries are low, then send me a text message
- If the hearing aid is turned on, then send me a text message

## Applet using a GO Action



- When I arrive at work, then change the program in my hearing aids
- If motion is detected on my video doorbell, then play text to speech

www.sonici.com



71

# Triggers and Actions




www.sonici.com



72



## Example: Applet available on GO Service



www.sonici.com



73

## Benefits of IFTTT

- Growing number of devices are internet connected
- Hearing aid is a part of this internet-based world
- Offers convenience and individualization

www.sonici.com



74

# Summary

www.sonicl.com



75

## Connectivity Portfolio



76

## Sonic and Connectivity

- Full line of both custom and behind-the-ear products boasting full Bluetooth connectivity options
- Standard WL in BTE products, optional WL in custom products
- 2.4 GHz Direct-to-ear streaming to iPhone with no intermediary device
- SoundClip-A available for Android streaming, hands free use of phone and remote mic functionality
- Easy pairing between accessories and hearing aids
- Hassle-free connectivity
- Convenience and individualization

[www.sonici.com](http://www.sonici.com)



77

## Thank You!



SoundDNA Platform

- Thank you so much for joining us today!
- If you have any questions...
  - 888.423.7834
  - [www.sonici.com](http://www.sonici.com)
  - [support@sonici.com](mailto:support@sonici.com)

[www.sonici.com](http://www.sonici.com)



78



Bluetooth is a registered trademark owned by Bluetooth SIG, Inc, USA.

Apple, the Apple logo, iPhone, iPad, iPod touch, and Apple Watch are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc. Android, Google Play, and the Google Play logo are trademarks of Google LLC.

© 2020 Sonic Innovations, Inc. All rights reserved.