

Do Hearing Aids Prevent Cognitive Decline?

A Sonova brand

PHONAK
life is on

Disclosure in accordance with CEU Requirements

PHONAK
life is on

Paul Kreimer, M.A., CCC-A

Manager, Clinical Training, Phonak, US

- I have a financial relationship to disclose:
 - Employee of Phonak who receives a salary
- I have no nonfinancial relationship to disclose

Presented in Partnership with

continued

Learning Outcomes



After this course, learners will be able to:

- describe the associations between cognitive decline and hearing loss.
- differentiate between different hypotheses about the origin of these associations.
- evaluate the current evidence related to this topic, cognitive decline and hearing loss.

Introduction



- Paul J Kreimer, M.A., CCC-A
 - Clinical Training Manager – Central Region
 - 11 years in a clinical setting
 - Hospital
 - ENT & Otology Practice
 - VA Hospital – Indianapolis
 - Private Practice
 - 11 years with Phonak as a Clinical Trainer and manager

Presented in Partnership with



PHONAK
life is on

Connectivity?

OTC?

Streaming?

eSolutions?

Hearables?

Framing?

Healthy Aging?



What does Healthy Aging Look Like to You?

PHONAK
life is on



Presented in Partnership with

continued

Well hearing = Well-being



Agenda


PHONAK
life is on

1. What we know about hearing loss and cognition
2. What we know about hearing aids and cognition
3. What opportunities do we have?

Presented in Partnership with

continued

What we know about hearing loss and cognition



Cognitive abilities

- Attention
- Processing speed
- Short-term memory
- Working memory
- Long-term memory
- Executive functions
- Semantic and language knowledge

Dementia


PHONAK
life is on

- Deterioration in cognitive function beyond normal aging.
- It affects:
 - Memory
 - Thinking
 - Orientation
 - Comprehension
 - Calculation
 - Learning capacity
 - Language and judgement

Presented in Partnership with

continued



PHONAK
life is on



Presented in Partnership with



6

Brain aging – normal and abnormal



- Degeneration in brain physiology
- Reduction in cognition functioning

Brain aging – normal and abnormal



- *MCI is estimated to effect 22.6% of adults over 60*
- *The number of individuals with dementia is expected to double every 20 years, resulting in over 65 million individuals affected by 2030*

Presented in Partnership with

CONTINUED


Hearing Loss: Precursor & potential risk factor for cognitive decline

PHONAK
life is on



Hearing Loss: Precursor & potential risk factor for cognitive decline

PHONAK
life is on




Presented in Partnership with

continued


Hearing Loss: Precursor & potential risk factor for cognitive decline

PHONAK
life is on



Hearing Loss: Precursor & potential risk factor for cognitive decline

PHONAK
life is on



Hearing loss was identified as the largest modifiable risk factor for dementia


Presented in Partnership with

continued

Correlation between hearing loss & various cognitive functions



- A small but significant association was found for age related hearing loss (ARHL) within all domains of cognitive function
- ARHL precedes the onset of clinical dementia by 5 to 10 years making it a possible non-invasive biomarker




Loughrey, D. G., Kelly, M. E., Kelley, G. A., Brennan, S., & Lawlor, B. A. (2018). Association of Age-Related Hearing Loss With Cognitive Function, Cognitive Impairment, and Dementia: A Systematic Review and Meta-analysis. *JAMA otolaryngology-- head & neck surgery*, 144(2), 115–126. doi:10.1001/jamaoto.2017.2513

Hearing Loss is Correlated with Incident Dementia



Health, Aging and Body Composition (Health ABC) study (2017):

- N = 1889
- age: 70 – 79 years
- no dementia at baseline
- adjusted for demographic and cardiac factors
- Normal/mild/moderate-severe HL



Conclusion:

- HL is associated with increased risk of developing dementia in older adults
- Randomized trials are needed to determine possible dementia onset postponement through treatment of HL

Presented in Partnership with



Summary




- Hearing loss and cognitive decline are associated with each other.
- This holds for specific cognitive functions and dementia
- People with HL also have accelerated cognitive decline compared to NH
(Amieva et al., 2015)
- HL is linked to accelerated brain atrophy (Gurgel et al., 2014; Lin et al., 2011, 2013, 2014)

Note:

- Study samples are usually large and effects are small
- Not all studies observe an association

**What we know about
the underlying
mechanisms**




Presented in Partnership with



Auditory Event-Related Potentials and Cognition

PHONAK
life is on

- Older adults with sloping HL and MCI have longer latencies and lower amplitude P1-N1-P1 complex, compared to without MCI ^{1,2,3}
- Hemispheric activation pattern more diffuse with MCI
- “Results suggest auditory signal traveling in the superiorly in the system appears to be degraded in individuals with MCI compared to those without”
Fulton, Seminars in Hearing, 2015
- ERPs may be a useful indicator of emergent MCI or AD




1. Fluegel B A, Matthews C, Harrison Bush A L, Lister J J, Edwards J D, Andel R. Neurophysiological indicators of early-stage cognitive decline. Paper presented at: American Academy of Audiology Conference; March 27–29 2014; Orlando, FL.
2. Golob E J, Irimajiri R, Starr A. Auditory cortical activity in amnesic mild cognitive impairment: relationship to subtype and conversion to dementia. *Brain*. 2007;130(Pt 3):740–752. [PubMed: 17293359]
3. Golob E J, Johnson J K, Starr A. Auditory event-related potentials during target detection are abnormal in mild cognitive impairment. *Clin Neurophysiol*. 2002;113(1):151–161. [PubMed: 11801437]

Potential Mechanisms behind the Association

PHONAK
life is on

Common Cause vs Cascade Effect



(Lin et al., AAIC 2016)

Presented in Partnership with

continued

What we know about hearing aids and cognition




Untreated Hearing Loss Drives the Association

PHONAK
life is on

“Cognitive decline associated with ARHL is probably preventable by early rehabilitation and increased opportunistic screening for the elderly”

Ray et al. *JAMA Otolaryngol Head Neck Surg.* 2018;144(10):876-882.

- English longitudinal study of aging
- Sample: N = 7385; 50+ years; no dementia, Alzheimer, Parkinsons, ear infections, CI
- Cross-sectional analyses



Presented in Partnership with

continued

Does hearing aid use alter cognitive trajectories?

PHONAK
life is on



Rate of cognitive decline
faster than controls



Rate of cognitive decline
similar to controls

- Amieva et al 2015: Over 25-yr follow-up, observational study,

Self-Reported Hearing Loss, Hearing Aids, and Cognitive Decline in Elderly Adults: A 25-Year Study

Hélène Amieva, PhD, Caroline Cherad, MSc, Caroline Girosch, MSc, Céline Méillon, MSc, Larissa Ballieu, PhD, and Jean-François Dartigues, MD, PhD

OBJECTIVES: To investigate the association between hearing loss, hearing aid use, and cognitive decline.

DESIGN: Prospective population-based study.

SETTING: Data gathered from the French Agir cohort study, a cohort study begun in 1989-90.

PARTICIPANTS: Individuals aged 65 and older (N = 1,670).

MEASUREMENTS AND MAIN RESULTS: At baseline, hearing loss was determined using a questionnaire assessing self-perceived hearing loss. 1,157 subjects reported no hearing loss, 1,159 reported moderate hearing loss, and 1,159 reported severe hearing loss. Cognitive decline was measured using the Mini-Mental State Examination (MMSE), administered at follow-up visits over 25 years.

RESULTS: Self-reported hearing loss was significantly associated with lower baseline MMSE scores ($\beta = -0.47$, $P < .001$) and greater decline during the 25-year follow-up period ($\beta = -0.04$, $P = .01$) independent of age, sex, and education. A difference in the rate of change in MMSE scores over the 25-year follow-up was observed between participants with hearing loss not using hearing aids and controls ($\beta = -0.06$, $P = .003$). In contrast, subjects with hearing loss using a hearing aid had no difference in cognitive decline ($\beta = 0.07$, $P = .08$) from controls.

CONCLUSIONS: Self-reported hearing loss is associated with accelerated cognitive decline in older adults; hearing aid use attenuates this decline. J Am Geriatr Soc 63:2009-2014, 2015.

Key words: hearing loss; hearing aids; cognitive decline; elderly.

From Epidemiology and Biostatistics, Institut National de la Santé et de la Recherche Médicale (INSERM) U1058, University of Bordeaux, Bordeaux, France. Address correspondence to Prof. Hélène Amieva, INSERM U1058, Epidemiology and Biostatistics, University of Bordeaux, 146 Rue Lyautey, 33076 Bordeaux Cedex, France. E-mail: helen.amieva@inserm.fr

DOI: 10.1111/gps.12407

© 2015 The Authors. Journal of Clinical Gerontology. © 2015 The American Geriatrics Society

0898-2643/15/0000-0000

Hearing loss is the third most common chronic health condition affecting older adults. Approximately 10% of individuals aged 65 and older have some degree of hearing loss, with estimates ranging from 70% to 90% of those aged 85 and older.¹⁻³ Individuals with hearing loss often experience depressive symptoms and social isolation. There is also evidence that older adults with hearing loss have poorer cognitive performance.⁴⁻⁷ Two longitudinal studies showed an association between hearing loss and cognitive decline over a range of follow-up.^{8,9} In particular, the Health, Aging and Body Composition (Health ABC) Study, conducted in a sample of 1,984 community-dwelling individuals aged 70 to 79, showed that hearing loss measured using audiometric testing was independently associated with accelerated cognitive decline and incident cognitive impairment during the 6 years of follow-up.⁸

Despite the high prevalence and consequences for health outcomes, hearing loss is largely underdiagnosed and thus undermanaged.¹⁰ Recent research of older adults with hearing impairment do not use hearing aids.¹¹ Little is known about the effect of hearing aids on health outcomes in older adults, in particular cognitive decline. In one study,¹² hearing aid use was associated with slightly lower rates of cognitive decline and risk of incident cognitive impairment in individuals with hearing loss, but the results were not statistically significant, possibly because of the short follow-up and lack of statistical power. In the other longitudinal study,¹³ because of small sample size, the association between hearing aid use and cognitive decline could not be properly examined.

With the ongoing follow-up of the French Agir cohort study, the 25-year follow-up of the cohort was used to assess the relationship between hearing loss and long-term age-related cognitive decline. Using self-reported hearing loss, the association between hearing loss and cognitive trajectories over 25 years was investigated in a community-based cohort of older adults, and the impact


Does hearing aid use alter cognitive trajectories?

PHONAK
life is on

“Providing hearing aids or other rehabilitative services for hearing impairment much earlier in the course of hearing impairment may stem the worldwide rise of dementia”

- Health and Retirement Study (HRS)
 - Repeated cognition measurements every 2nd year 1996 – 2014
- Sample: N = 2040; 50+ years
 - Started using hearing aids during evaluation period
- Longitudinal analyses
- Better memory scores in HA users
- Slower decline in memory scores after than before using HAs

(Maharani et al. 2018. JAGS.)




Presented in Partnership with


continued

Hearing Aids and Listening Effort


PHONAK
life is on



- Noise reduction reduces cognitive load.
- Listeners have free capacity to do other things in parallel with listening.



- DuoPhone reduces cognitive load as if listening at a 3 dB higher SNR



- StereoZoom increases performance and reduces subjective and objective listening effort, noise annoyance, and memory effort in adverse listening conditions

Lower risk with use of hearing instruments!

PHONAK
life is on

Journal of the American Geriatrics Society

AGS
American Geriatrics Society
Leading Through Improving Care for Older Adults

Clinical Investigation

Can Hearing Aids Delay Time to Diagnosis of Dementia, Depression, or Falls in Older Adults?

Elham Mahmoudi PhD, Tanima Basu MS, Kenneth Langa MD, PhD, Michael M. McKee MD, MPH, Philip Zazove MD, Neil Alexander MD, Neil Kamdar MA

First published: 04 September 2019 | <https://doi.org/10.1111/jgs.16109>

- Diagnosed with dementia – 18% lower
- Diagnosed with depression – 11% lower
- Treated for fall-related injuries – 13% lower

Listen Up: Hearing Aids May Do More Than Just Treat Hearing Loss

Kristen Geranchar Contributor @healthcare
I write about chronic diseases, treatments and patients' experiences.

Age-related hearing loss is rarely dramatic, happening slowly over time, robbing older people of snippets of conversation and leaving them unable to filter out background noise as well as they once could. Early signs can include routinely cranking up the TV volume and having trouble with phone dialogue. Left untreated, hearing loss can lead to social isolation and mental and physical decline.

Presented in Partnership with


CONTINUED

So can we conclude that hearing aids prevent cognitive decline?

Longitudinal randomly controlled trials (RCTs) are needed

Cognitive decline - Achieve Study

ACHIEVE Trial Conceptual Model



- ACHIEVE Study (started early 2018)
 - Multicenter RCT investigating whether hearing loss treatment could reduce cognitive decline; duration 5 years
 - N = 850; age 70–84 years; cognitively normal with hearing loss
 - Arm1 hearing intervention: hearing needs assessment, fitting of HAs, education/counseling
 - Arm2 successful aging intervention: individual sessions with a health educator on healthy aging

Presented in Partnership with




Evidence suggests that...



- Hearing loss and cognitive decline are associated
 - This holds for specific cognitive functions and dementia
- Untreated hearing loss drives the association
 - Better cognition with HA use
 - Slower cognitive decline with HA use
 - Reduced listening effort with HA use
- Hearing loss is a **modifiable** risk factor. Early intervention is strongly recommended
- **Results from longitudinal RCTs needed to confirm the causal pathway!**

What opportunities do we have to become patient advocates?



Presented in Partnership with





Create Awareness

Optimized Amplification Leads to Reduced Listening Effort

Access to Technology



Exceptional sound quality



Reduced listening effort in noise


RogerDirect™



Presented in Partnership with

continued

Why is auditory training needed?



- Auditory training can help speed up this process

It Can Start with YOU!

- Hearing Care Professionals are in an ideal position to observe the longitudinal cognitive trajectory of the older adult and recognize pathological cognitive function in the preclinical phase
 - Nature of the service delivery model
 - Counseling segment




Presented in Partnership with

continued



Test Yourself!

Rules:
(1 minute)



“Draw a clock. Put in all the numbers and set the time to 10 after 11”.

How did you do?



Three criteria:



1. Contour (1 pt.): the clock face must be a circle with only minor distortion acceptable (e.g., slight imperfection on closing the circle);
2. Numbers (1 pt.): all clock numbers must be present with no additional numbers; numbers must be in the correct order and placed in the approximate quadrants on the clock face; Roman numerals are acceptable; numbers can be placed outside the circle contour;
3. Hands (1 pt.): there must be two hands jointly indicating the correct time; the hour hand must be clearly shorter than the minute hand; hands must be centred within the clock face with their junction close to the clock centre.

Presented in Partnership with



It Can Start with YOU!

- Consider the nature of your hearing care service delivery model
- Counseling segment



- AAA and ASHA have both included cognitive screening in the defined scope of practice for audiologists for over a decade

Integration Tips and Resources



- Understand your scope of practice
- Be sensitive to the “sensitive”
- Be prepared to explain why



Suggested language for referral:

“Our clinic has been utilizing the MoCA to screen patients for possible cognitive impairment. _____ recently scored ____ out of 30 on his/her exam which indicates a need for referral. Please evaluate at your discretion.”

Presented in Partnership with



Family Centered Care



Shared **understanding**, **empathy**, and **responsibility** for managing hearing difficulties for the patient and the family

Family Centered Care



*FCC will help develop a **shared understanding** and **shared responsibilities** for treating the communications difficulties both parties have.*

Presented in Partnership with

continued

Summary

PHONAK
life is on

- Hearing loss and cognitive decline are associated
- Untreated hearing loss drives the association
- Further research needed
- HCPs are uniquely positioned to help with early detection
- Advocate for your patients' cognitive health and overall well-being
- Family unit involvement is integral to your patients' socio-emotional health



Together,
we change lives

Presented in Partnership with

