



<p>MED-EL Cochlear Implants & Electrical Stapedial Reflex Thresholds (ESRTs)</p>	
<p>Protocol and Case Studies from Chattering Children/The River School</p>	

<div data-bbox="1128 1136 1386 1182" data-label="Image"> </div> <div data-bbox="232 1199 589 1249" data-label="Section-Header"> <h2>Learning Outcomes</h2> </div> <div data-bbox="232 1295 1344 1535" data-label="List-Group"> <ul style="list-style-type: none"> – After this course, learners will be able to diagram the auditory mechanism that elicits ESRT measures. – After this course, learners will be able to list the equipment required to perform ESRTs on a cochlear implant (CI) recipient. – After this course learners will be able to describe some clinical tips for working obtaining ESRTs in the pediatric population. </div>	
© MED-EL	

Presented in partnership with



Fitting Fundamentals with MED-EL




Clinicians' priority (in general) should be to focus on measuring accurate Maximum Comfortable Loudness (MCL) levels.

- All other parameters can potentially be left at default
- For some patients, other parameters may be important

© MED-EL

Question:



“If MCL is set to Maximum does that mean the patient is hearing at Maximum loudness?”

© MED-EL

Presented in partnership with

continued



With MED-EL, patients hear at MCL level when...


- User volume setting is set to 100%
- AND
- Input level is ~105 dB SPL

At all other times, stimulation occurs below MCL

(Stoebich et. al, Ear & Hearing 1999; 20; 104 – 116)

© MED-EL

AGC/Adaptive Sound Window




© MED-EL

Presented in partnership with



Setting Threshold (THR)

- Not a critical parameter for speech perception performance
 - **Should be inaudible**
- Setting options
 - Measure behaviorally and set below audibility
 - Quick fit approaches
 - Leave at default of 8 - 10%



© MED-EL


General Fitting Considerations

Average MCL level:

- **8 to 15 qu** (first fitting)
- **15 to 25 qu** (stabilized map)

~Be careful with high stimulation

- Overstimulation
- Poor speech discrimination/performance
- Facial Nerve Stimulation (FNS)





© MED-EL

Presented in partnership with

continued

MED-EL Fitting Reminders

- MCL = most important parameter
- THR = must be inaudible
- Default volume at 100%
- Watch for outliers 
- Audiogram results at ~25 - 30 dB are normal



Fitting goal


Provide access to sound
with the appropriate
balance of:

- ✓ Comfort
- ✓ Satisfaction
- ✓ Performance



Presented in partnership with

continued



ESRT

- Not always present
- High correlation between ESRT and MCLs
- Verifies appropriately-fit maps

Kosaner, et. al (2009, 2018)

© MED-EL

Presented in partnership with

continued



ESRT in the Trenches: *A Practical Guide to Implementing ESRT in a Pediatric Setting*

Sydney Bednarz, AuD
Julia Reid, AuD

CHATTERING
CHILDREN



THE
RIVER
SCHOOL

Presented in partnership with

continued[®]

Table of Contents


- Our Clinic Protocol
- ESRT at Our Clinic
- Benefits & Challenges
- Billing & Coding
- Case Studies

Our Clinic Protocol

Presented in partnership with



The River School & Chattering Children



Our Typical CI Protocol

- During the First Year
 - Activation
 - 2 weeks
 - 4 weeks
 - 9 weeks
 - 3 months
 - 6 months
 - 9 months
 - 1 year
- After the First Year
 - If concerns, every 3 months
 - 1-8 years old: Biannually
 - >9-10 years old: Annually
- Booth testing
 - Warble tones starting at 4 weeks
 - Speech as soon as possible

Presented in partnership with




When We Use ESRT

- Begin at 2-week or 1-month follow-up
 - Fill in electrodes as needed at further mapping visits
- If it has never been measured on an existing patient
- Transfer patients
- Patient complains that things are too loud overall
- Patient's dynamic range is really narrow
- Annually *
 - ESRTs increase over the first year of implant use (Gordon et al, 2004; Kosaner et al, 2009)
 - ESRTs stabilize over time (Kosaner et al, 2009)

Verification of Mapping

Threshold Determination

- Warble Tones
- Ling 6 Sounds



www.JaneMadell.com

Source: Jane Madell, used with permission.

Presented in partnership with

continued

Verification of Mapping

Speech Perception

- Words
- Sentences – Quiet and Noise
- Closed Set → Open Set
- Picture Pointing → Speech Production
- Scoring
 - Word Score
 - Phoneme Score

Details

- 1 meter away from speaker
- 60 dB SPL
- Speech front, Noise front
- Recorded Materials (if possible)
- Based off of the Pediatric Minimum Speech Test Battery (PMSTB)

(Ulher et al, 2017)

Pediatric Minimum Speech Test Battery

- Parental Questionnaires
- VRISD
- ESP Patterns, Spondees, Monosyllables
- PSI Words & Sentences
- MLNT/LNT
- CNC
- BKB Quiet
- Baby Bio
- BKB SIN
- Baby Bio Noise

MLNT/LNT

If <25% STOP

If 25-79%, CNC,
repeat at follow-up

If >80% twice:
CEILING, CNC

(Ulher et al, 2017)


Presented in partnership with



ESRT at Our Clinic

Room Set-up


- Child-sized chair for patient
- TV Screen for distraction
- GSI Tymstar Pro
- Programming Laptop
 - on a child sized table
- Clinician's Chair
 - May or may not be an adult-sized chair
 - ... OR sitting on the ground




Presented in partnership with

CONTINUED

Parameter Adjustments





ESRT Step by Step




Presented in partnership with

CONTINUED

ESRT Step by Step



ESRT Step by Step



(Wolfe et al, 2017; Carranco et al, 2019)

Presented in partnership with

continued

ESRT Step by Step



ESRT Step by Step


- Connect the child's processor to the programming computer using the Max interface
- If the child is bilateral or bimodal, the opposite ear's device should be removed
- Be sure patient/child is still and avoids excessive movements such as swallowing



Presented in partnership with


continued

Using the Fitting Tab



ESRT Step by Step

- Look for a time-locked response on the tympanometer
- Response may be up-ward deflecting or down-ward



Presented in partnership with

CONTINUED

Benefits & Challenges

Benefits of ESRT

- Helps set M levels during programming
- Ensures not over or under-stimulating children who cannot other report loudness
- Correlates to target thresholds for Ling sounds for MED-EL
- Can be used with any manufacturer
 - MED-EL typically has a 1:1 ratio for where you measure ESRT and MCL levels
- Mostly stable over time
 - Annually tested to verify stability

Presented in partnership with



Common Challenges

- Keeping the child still or from talking
 - “Catch a bubble!”
 - Screen time
- Limited ESRT responses across the array
 - Fill in the gaps at subsequent programming sessions
- Admittance Issues
 - Maintaining a seal
 - Maintaining zero pressure
 - Adjust pressure dial
 - Middle Ear Issues, PE Tubes, perforation, etc.

Billing & Coding

Presented in partnership with



National Correct Coding Initiative (CCI)

Column 1 CPT Procedure (one)	Column 2 Paired With (one)	Can be performed on same date? Yes/No		If so, use what modifier? [2]
		Office setting	Hospital outpatient setting ^[1]	
92601	92507, 92508, 92521, 92522, 92523, 92524, 92550, 92567, 92568, 92570, 92585, 92586, 92602, 92626	Y	Y	-59
92601	92552, 92553, 92555, 92556, 92557, 92558, 92561, 92562, 92563, 92564, 92565, 92571, 92572, 92575, 92576, 92577, 92579, 92582, 92583, 92584, 92587, 92588, 92596, 92597, 92603, 92604, 97755, 0208T, 0209T, 0210T, 0211T, 0212T	N	N	N/A

ASHA (n.d.)

National Correct Coding Initiative (CCI)

Column 1 CPT Procedure (one)	Column 2 Paired With (one)	Can be performed on same date? Yes/No		If so, use what modifier? [2]
		Office setting	Hospital outpatient setting ^[1]	
92601	92507, 92508, 92521, 92522, 92523, 92524, 92550, 92567, 92568, 92570, 92585, 92586, 92602, 92626	Y	Y	-59
92601	92552, 92553, 92555, 92556, 92557, 92558, 92561, 92562, 92563, 92564, 92565, 92571, 92572, 92575, 92576, 92577, 92579, 92582, 92583, 92584, 92587, 92588, 92596, 92597, 92603, 92604, 97755, 0208T, 0209T, 0210T, 0211T, 0212T	N	N	N/A

ASHA (n.d.)

Presented in partnership with



Case Studies

Case Study # 1


- 4 year, 8 month old male
- Underwent cochlear implant surgery for the RIGHT ear on July 15, 2016 and received a MED-EL SYNCHRONY cochlear implant
 - Activated on September 16, 2016 with the SONNET BTE cochlear implant processor
- Underwent cochlear implant surgery for the LEFT ear on December 30, 2016 and received a MED-EL SYNCHRONY cochlear implant cochlear implant
 - Activated on February 14, 2017 with a SONNET BTE cochlear implant processor
- Left revision surgery on August 21, 2019
 - Activated on June 9, 2019
- Uses RONDO processors as back ups

Presented in partnership with




Case Study #1

- Pre-ESRT map (71)
 - Measured Ts
 - Flat map as per recommended



Case Study #1

- Pre-ESRT map (71)
- NU-CHIPS AU- 80%
 - DNT individual ear
- Ling thresholds
 - mm- 35 dB
 - oo- 65 dB
- Datalogging- 7 hours/day
- Prior to re-implantation




Presented in partnership with

CONTINUED


Case Study #1

- ESRT map (79)
- Ts locked at 10%



Case Study #1

- ESRT based map (80)
- Ms reduced 3%
- Ts locked at 10%




Presented in partnership with

continued



Case Study #1

- ESRT based map (81)
- Ms reduced 15%
- Ts locked at 10%



Case Study #1

- 4 months post re-activation
- ESRT based map (81)
- NU-CHIPS Left-20%
- Datalogging- 5.5 hours/day



Presented in partnership with


continued

Case Study #2

- 31-year-old female
- Enlarged vestibular aqueducts (EVA)
- Spoken English and American Sign Language
- Between 2-3 years: Diagnosed with bilateral hearing loss
- 3 years old: Fit with bilateral HAs
- 7 years old: Discontinued right HA due to limited benefit
- 25 years old: Received a left MED-EL CONCERT FLEX28 internal electrode

Case Study #2

- Pre-ESRT based map
- Consistent CI User
- CNC = 74%
- BKB-SIN = 9 dB SNR




Presented in partnership with

continued


Case #2

- Pre-ESRT map



Case #2

- ESRT map
- Did NOT like
- Facial stimulation noted on high frequency channels




Presented in partnership with

continued


Case #2

- ESRT-based map
 - Able to bring down MCL levels, but still higher than ESRT
 - Measured Ts and decreased/swept to ensure inaudibility
- Recommendation:
 - Sound Diet



Case Study #2

- In a live mode, patient thought it was too soft and felt like everyone was whispering but gradually adjusted.
- After the appointment, patient reported much more comfort with the sound of running water.
- CNC = 72%



Presented in partnership with



continued

Case Study #3

- 4-year-old female (currently)
- Enlarged Vestibular Aqueducts (EVA)
- Asymmetric, progressive hearing loss
- 7 mos: Right hearing aid
- 23 mos: Left hearing aid
- 27 mos: Received a right MED-EL FLEX28 internal electrode

Case #3

- Pre-ESRT




Presented in partnership with


continued

Case #3

- ESRT levels
- Note:
 - Higher stimulation levels in the high frequencies than a typical map. ESRT is “typically” flat across the array.



Case #3




Right CI	NUCHIPS MLNT	80% 79%
Bimodal	LNT	74%

Presented in partnership with




Case #3



Case Study #3

- Pediatric AzBio
- Quiet = 96%; Noise +5 SNR = 64%

Right CI	PBK	76%
Left HA	PBK	84%
Bimodal	PBK BKB-SIN	80% 3.5 SNR



Presented in partnership with



Case Study #3

CASL-2 Subtest	Standard Score	Percentile Rank	Test-age Equivalent
Receptive Vocabulary	119	90	4;9-4;11
Expressive Vocabulary	110	75	4;6-4;8
Sentence Expression	119	90	5;3-5;5
Grammatical Morphemes	114	82	5;0-5;2
Sentence Comprehension	124	95	6;0-6;2
Inference	108	70	4;3-4;5
Pragmatic Language	108	70	5;0-5;2

CASL-2 Indexes	Standard Score	Percentile Rank
General Language Ability	112	79
Receptive Language Index	119	90
Expressive Language Index	109	73
Lexical Semantic Index	113	81
Syntactic Index	119	90

Conclusion

- Great objective measure
 - To verify where you've set your MCL levels
 - When patients aren't reliable at loudness scaling
 - When patients are loudness junkies
- ESRT can be completed on children
- Better patient outcomes and comfort
- Validation and verification of programming is essential

Presented in partnership with



References

- American Speech Language Hearing Association. (n.d). National Correct Coding Initiative (CCI) and Outpatient Code Edits (OCE) Edit Tables. Retrieved from <https://www.asha.org/Practice/reimbursement/coding/CCI-Edit-Tables-Audiology/>
- Carrano Hernandez, L., Cristerna Sanchez, L., Camacho Olivares, M., Rodriguez, C., Finley, C.C., Saoji, A. A. (2019). Effect of Probe-Tone Frequency on Ipsilateral and Contralateral Electrical Stapedial Reflex Measurement in Children with Cochlear Implants. *Ear Hear*, 40(3): 732-740.
- Gordon, K.A., Papsin, B.C., & Harrison, R.V. (2004). Toward a battery of behavioral and objective measures to achieve optimal cochlear implant stimulation levels in children. *Ear Hear*, 25(5):447-63.
- Kosaner, J., Anderson, I., Turan, Z., & Deibl, M. (2009). The Use of ESRT in Fitting Children with Cochlear Implants. *Int. Adv. Otol.*, 5(1): 62-71.
- Kosaner, J., Spitzer, P., Bayguzina, S, Gultekin, M., & Behar, L.A. (2018). Comparing eSRT and eCAP measurements in pediatric MED-EL cochlear implant users. *Cochlear Implants Intl*, 19(3): 153-161.
- Madell, J.R. (n.d.). Are Children Hearing What They Need to Hear? [PDF file]. Retrieved from www.earfoundation.org.uk/files/download/1630
- Uhler, K., Warner-Czyz, A., Gifford, R., & PMSTB Working Group.(2017). Pediatric Minimum Speech Test Battery. *J Am Acad Audiology*, 28:232-247
- Wolfe, J., Gilbert, M., Schafer, E., Litvak, L. M., Spahr, A. J., Saoji, A., & Finley, C. (2017). Optimizations for the Electrically-Evoked Stapedial Reflex Threshold Measurement in Cochlear Implant Recipients. *Ear Hear*, 38(2) 255-261. ASHA

Thank you!

Presented in partnership with

