

1. This document was created to support maximum accessibility for all learners. If you would like to print a hard copy of this document, please follow the general instructions below to print multiple slides on a single page or in black and white.
2. If you are viewing this course as a recorded course after the live webinar, you can use the scroll bar at the bottom of the player window to pause and navigate the course.
3. This handout is for reference only. Non-essential images have been removed for your convenience. Any links included in the handout are current at the time of the live webinar, but are subject to change and may not be current at a later date.
4. Copyright: Images used in this course are used in compliance with copyright laws and where required, permission has been secured to use the images in this course. All use of these images outside of this course may be in violation of copyright laws and is strictly prohibited.

How to print Handouts

- On a PC
 - Open PDF
 - Click Print
 - Choose # of pages per sheet from dropdown menu
 - Choose Black and White from “Color” dropdown
- On a Mac
 - Open PDF in Preview
 - Click File
 - Click Print
 - Click dropdown menu on the right “preview”
 - Click layout
- Choose # of pages per sheet from dropdown menu
- Checkmark Black & White if wanted.
- If more details needed please visit our FAQ page: <https://www.audiologyonline.com/help>

No part of the materials available through the continued.com site may be copied, photocopied, reproduced, translated or reduced to any electronic medium or machine-readable form, in whole or in part, without prior written consent of continued.com, LLC. Any other reproduction in any form without such written permission is prohibited. All materials contained on this site are protected by United States copyright law and may not be reproduced, distributed, transmitted, displayed, published or broadcast without the prior written permission of continued.com, LLC. Users must not access or use for any commercial purposes any part of the site or any services or materials available through the site.

Technical issues with the Recording?

- Clear browser cache using [these instructions](#)
- Switch to another browser
- Use a hardwired Internet connection
- Restart your computer/device

Still having issues?

- Call 800-753-2160 (M-F, 8 AM-8 PM ET)
- Email customerservice@AudiologyOnline.com

Oticon CROS

A revolutionary sound experience



Maegan Mapes, Au.D
Education & Training Specialist

oticon
life-changing technology

Course presented in partnership with

continued

Learning Objectives

Participants will be able to:

- describe the acoustic limitations of individuals with single-sided deafness
- describe how an unprecedented signal processing technology can provide better conditions for speech understanding and localization
- identify the full portfolio of devices available to support Oticon CROS
- list the steps to setup and program Oticon CROS

oticon

Course presented in partnership with

continued



Oticon CROS

Outstanding sound quality for
those with single-sided deafness

SSD can result in poor speech
understanding especially in noise
and reduced localization ability


Oticon CROS offers outstanding
sound quality and **OpenSound
Navigator**

World's first **TwinLink dual-
streaming** solution for people with
SSD

Course presented in partnership with

continued

Background on SSD



How many of you...

**...have seen a
client with single-
sided deafness?**



**...in your mind,
what constitutes
single-sided
deafness?**


oticon

Course presented in partnership with

continued

What is single-sided deafness?

*" ...single-sided deafness (SSD) is defined as **unaidable hearing in one ear and normal hearing or an aidable hearing loss in the other ear** "*



Taylor 2010


oticon

Course presented in partnership with

continued

Do you know...

**...how many
suffer from
SSD?**



oticon

Course presented in partnership with

continued

Every year...

60,000
USA




≈ 1% of the general population

oticon

Course presented in partnership with

continued

**What has caused
the SSD cases the
you have seen?**



Voelker and Chole 2010

oticon

Course presented in partnership with

continued

Who are the clients?


SSD can occur in all population groups



oticon

Course presented in partnership with

continued


A photograph of three people sitting at a table in a restaurant. A woman with curly hair is smiling and looking towards a man in a dark blazer who has a serious expression. Another man is partially visible on the right, holding a piece of food. The background is a blurred restaurant interior.

What do they
struggle with?

oticon

Course presented in partnership with

continued



Reduced localization
and ability to judge

Poor Speech
Recognition in noise

oticon

Course presented in partnership with

continued

The importance of binaural hearing

2 ears are necessary in order to obtain the benefits of...

Binaural loudness summation

Binaural redundancy



Binaural squelch

Head shadow effect

Dillon 2012



oticon

Course presented in partnership with

continued

Treatment options for SSD



Benefits and Considerations

	CROS 	BAHS 
Pros	<ul style="list-style-type: none">• Non-surgical• Can accommodate for hearing loss on better ear	<ul style="list-style-type: none">• Wears one device, nothing in better ear
Cons	<ul style="list-style-type: none">• The better ear may be partially occluded• User must wear two devices• Does not restore binaural hearing	<ul style="list-style-type: none">• Minimally invasive surgery• Does not restore binaural hearing

oticon

CROS/BiCROS – contralateral routing of signal

Two methods for wireless transmission from a poorer ear to a better ear



Taylor 2010

oticon

Candidacy



Course presented in partnership with

continued

Candidacy

**...can the user
benefit from
conventional
amplification?**



**...is one ear normal
or not?**

oticon

CROS for kids?

- CROS is intended for children over 36 months
- No tamper resistant solution available
- Available in 7 standard colors only



oticon

Course presented in partnership with

continued

Oticon CROS



Course presented in partnership with

continued

Velox S™ – power for BrainHearing™

The most advanced and powerful platform ever

- Increased processing power
- 56,000 measurements per second (30 times faster)



oticon

Course presented in partnership with

continued



Oticon CROS

Outstanding sound
quality for those with SSD

World's first CROS featuring **OpenSound Navigator™**

World's first CROS with **TwinLink™ dual-streaming**


oticon

Course presented in partnership with

continued

Featuring OpenSound Navigator


Optimized for CROS



Analyze



Balance



Noise removal


oticon

Course presented in partnership with

continued

Industry first - TwinLink dual streaming

Receiver can simultaneously receive 2.4GHz and NFMI streamed signals




oticon

Course presented in partnership with

continued

What is Oticon CROS?

Contralateral Routing Of Signal



 **TwinLink**
NFM + 2.4 GHz

oticon

Course presented in partnership with

continued

Evidence

WHITEPAPER
2019

Oticon CROS

SUMMARY

Oticon now introduces Oticon CROS, the first ever wireless CROS/BiCROS solution using an open sound paradigm and with dual-streaming capability using TwinLink™.

Oticon CROS is a unique solution for persons with single-sided deafness and it now makes the advanced signal processing of the OpenSound Navigator™ feature available to even more people with hearing loss.

A 2019 internal study shows the benefit of having dual-streaming capability within a CROS solution. While streaming sound using 2.4 GHz Bluetooth Low Energy technology, the CROS transmission was switched on and off in order to determine whether or not transmission of sound from the test subject's poorer ear side affected a person's awareness of speech in the environment. Results showed a 50% improvement in speech awareness from 33% to 49% awareness when transmission was active during streaming from an external source.

This tells us that CROS/BiCROS users can benefit greatly from having access to speech information in their environment, also when streaming.



Susanna Love
Callaway, Au.D.
Director of Clinical Audiology,
Oticon A/S, Denmark



Pernille Aaby Gade
Clinical Research
Audiology Assistant
Oticon A/S, Denmark

oticon
PEOPLE FIRST

Study purpose



- Participants will perform better hearing speech from the poorer ear side when wearing transmitter vs not wearing transmitter.
- Simultaneous streaming using wireless blue-tooth low energy connection to the hearing aid is possible while CROS transmits signal from the impaired side to the hearing aid via NFMI.

oticon


Study methodology



- Internal lab test
- Participants = 8
 - 4 CROS
 - 4 BiCROS
- Dual-task paradigm
 - Task related to 2.4 GHz streaming content
 - Task related to NFMI transmission content
- Outcome measure
 - Speech awareness: number of words heard with transmission ON vs. OFF

oticon

Study setup




- Target speech - side
 - 208 common two syllable words & phrases
 - Recorded in-house
- Target speech – front
 - Directly streamed into hearing aid
 - 1 minute news stories
- Noise source – behind
 - 4 speaker array
 - ICRA speech-shaped noise

oticon

Study results


50% improvement in speech awareness while streaming



(Callaway & Aaby Gade, 2019)

oticon

Practicalities



The Oticon CROS solution

An overview

Transmitter



Receiver



oticon

Course presented in partnership with

continued

How to tell the difference



oticon

Course presented in partnership with

continued

The transmitter

Transmitter

**Speaker is real but
has no sound or
acoustic purpose**

**10 kHz broad
bandwidth signal
is transmitted**

**Sound processing
before sending
to receiver**



oticon

The transmitter

Transmitter



**Volume control of
transmitted sound**

**No Bluetooth® low
energy technology**

**IP 68
dust & water
resistant**

oticon

The transmitter

Sound processing

Transmitter



OpenSound Navigator settings:

- Simple environments: 0 dB NR
- Complex environments: -5 dB NR
- Transition: High


oticon

Course presented in partnership with

continued

The transmitter

Sound processing



High pass filter with 1500 Hz cut-off:

- Sounds above 1500 Hz are transmitted via NFMI
- Sounds below 1500 Hz passes the head naturally

oticon

Flexible controls

Transmitter

Mute/unmute
Press at least 4 sec
(Short press to unmute)



oticon

Transmitter



Receiver



oticon

Course presented in partnership with

continued

The receiver

**Receive NFMI signal
and mix with own
mic signal**

**Monaural
hearing aid
(normal hearing aid)**

Receiver




Oticon ON App

**Control of
overall volume**

oticon

Oticon CROS

Compatibility



Oticon Xceed / Xceed Play



Oticon Opn S™ / Opn Play™ / Oticon Ruby



Please note: Oticon CROS is not compatible with Oticon Opn S 3 and Oticon Opn Play 2

oticon

Course presented in partnership with

continued

Oticon CROS

Compatibility



- For new CROS/BiCROS fitting: ideal to recommend miniRITE T for receiver to match transmitter style
- For existing wearer of other style: counsel on physical differences between devices
- Receiver must have Firmware 8.0

oticon

What if the receiver is a different style?



MiniRITE
no volume control



oticon

Oticon CROS

Find the client's favorite look



CO90
Chroma Beige



CO94
Terracotta



CO93
Chestnut Brown



CO63
Diamond Black



CO92
Steel Grey



CO91
Silver Grey




CO44
Silver

oticon

Course presented in partnership with

continued


Programming



Course presented in partnership with

continued

Pairing transmitter and receiver for the first time



oticon

Course presented in partnership with

continued



Remember when programming...

Wireless fitting is recommended for best performance during fitting

CROS transmitter is not visible in Genie 2

Transmitter is active during fitting

It will not be possible to firmware update CROS Transmitter

oticon

Course presented in partnership with

continued

Objective verification

Verification with REMs are considered Best Practice


oticon

Course presented in partnership with

continued

Balancing transmitter and receiver

The transferred sound should not be louder than other sound entering the better ear



oticon

Course presented in partnership with

continued

A photograph of three people (two men and one woman) sitting at a table in a restaurant, smiling and looking at each other. The background is a blurred bar area.

In summary

SSD can result in poor speech understanding especially in noise and reduced localization ability

Oticon CROS offers outstanding sound quality and **OpenSound Navigator**

World's first **TwinLink dual-streaming** solution for people with SSD

oticon

Course presented in partnership with

continued



Thank you!

Course presented in partnership with

continued