

1. This document was created to support maximum accessibility for all learners. If you would like to print a hard copy of this document, please follow the general instructions below to print multiple slides on a single page or in black and white.
2. If you are viewing this course as a recorded course after the live webinar, you can use the scroll bar at the bottom of the player window to pause and navigate the course.
3. This handout is for reference only. Non-essential images have been removed for your convenience. Any links included in the handout are current at the time of the live webinar, but are subject to change and may not be current at a later date.
4. Copyright: Images used in this course are used in compliance with copyright laws and where required, permission has been secured to use the images in this course. All use of these images outside of this course may be in violation of copyright laws and is strictly prohibited.

## How to print Handouts

- On a PC
  - Open PDF
  - Click Print
  - Choose # of pages per sheet from dropdown menu
  - Choose Black and White from “Color” dropdown
- On a Mac
  - Open PDF in Preview
  - Click File
  - Click Print
  - Click dropdown menu on the right “preview”
  - Click layout
- Choose # of pages per sheet from dropdown menu
- Checkmark Black & White if wanted.
- If more details needed please visit our FAQ page: <https://www.audiologyonline.com/help>

No part of the materials available through the continued.com site may be copied, photocopied, reproduced, translated or reduced to any electronic medium or machine-readable form, in whole or in part, without prior written consent of continued.com, LLC. Any other reproduction in any form without such written permission is prohibited. All materials contained on this site are protected by United States copyright law and may not be reproduced, distributed, transmitted, displayed, published or broadcast without the prior written permission of continued.com, LLC. Users must not access or use for any commercial purposes any part of the site or any services or materials available through the site.

## Technical issues with the Recording?

- Clear browser cache using [these instructions](#)
- Switch to another browser
- Use a hardwired Internet connection
- Restart your computer/device

## Still having issues?

- Call 800-753-2160 (M-F, 8 AM-8 PM ET)
- Email [customerservice@AudiologyOnline.com](mailto:customerservice@AudiologyOnline.com)



# Paradise in Target 7.0

The fitting software of choice

Alyssa Ricevuto, AuD

## Learning Outcomes

- After this course learners will be able to describe the controls that are available to program and adjust Phonak Paradise.
- After this course learners will be able to explain how to adjust Phonak Paradise settings based upon patient needs.
- After this course learners will be able to identify situations when adjustments will be beneficial



# Key innovations of Audéo Paradise

## Unrivaled sound quality

- New sound processing chip - PRISM™
- New motion sensor
- New Adaptive Phonak Digital 2.0
- New AutoSense OS™ 4.0
  - Speech Enhancer
  - Dynamic Noise Cancellation
  - Motion Sensor Hearing



## Universal connectivity


- Connects to smartphones, TV & more
- Made For All
- Multiple Bluetooth® connections
- New Tap Control for easy access to Bluetooth functionalities

## Personalized digital solutions


- New Hearing Screener 2.0
- New myPhonak 4.0 app
  - Hearing Diary 2.0
  - AudiogramDirect (via Remote Support)
  - Remote Control
- New Client Dashboard 2.0

Paradise is available with groundbreaking lithium-ion rechargeable technology

# Unrivaled sound quality



## New hardware



### **PRISM™**

Processing Real-time Intelligent Sound Management

Double  
memory

More connectivity  
options



### **Motion sensor**

Usability  
options

Hearing  
performance feature

# A Paradise for every patient



P-312



P-13T



P-R



P-RT

Available in:



Sand Beige



Beige



Chestnut



Champagne



Silver Gray




Graphite Gray



Velvet Black


# Fitting range

**S**




Mild to moderately-severe hearing loss.

**M**




Mild to severe hearing loss.

**P**



Mild to severe hearing loss.

**UP**



Moderate to profound hearing loss.

# Phonak Paradise Feature Summary

## 90 Premium

Comfort in Echo  
Speech in Car  
Speech in Loud Noise  
Music  
Comfort in Noise  
Speech in Noise  
Calm Situation  
Streaming Speech  
Streaming Music  
Comfort in Echo (manual)  
Speech in Loud Noise (manual)  
Speech in 360° (manual)  
UltraZoom Premium  
Dynamic Noise Cancellation  
Speech Enhancer  
EchoBlock  
Tap Control\*  
StereoZoom  
WindBlock  
SoundRelax  
Motion Sensor Hearing\*  
DuoPhone  
Real Ear Sound  
SoundRecover2  
WhistleBlock  
RogerDirect™  
Tinnitus Balance  
20 Channels

## 70 Advanced


Music  
Comfort in Noise  
Speech in Noise  
Calm Situation  
Streaming Speech  
Streaming Music  
Speech in Loud Noise (manual)  
Speech in 360° (manual)  
UltraZoom Advanced  
EchoBlock  
Tap Control\*  
StereoZoom  
WindBlock  
SoundRelax  
Motion Sensor Hearing\*  
DuoPhone  
Real Ear Sound  
SoundRecover2  
WhistleBlock  
RogerDirect  
Tinnitus Balance  
20 Channels

## 50 Standard

Comfort in Noise  
Speech in Noise  
Calm Situation  
Streaming Speech  
Streaming Music  
UltraZoom Standard  
Motion Sensor Hearing\*  
DuoPhone  
Real Ear Sound  
SoundRecover2  
WhistleBlock  
RogerDirect  
Tinnitus Balance  
16 Channels

## 30 Essential

Speech in Noise  
Calm Situation  
Streaming Speech  
Streaming Music  
UltraZoom Essential  
SoundRecover2  
WhistleBlock  
RogerDirect  
Tinnitus Balance  
12 Channels

 AutoSense OS\* 4.0


\*Rechargeable models only

## Phonak Paradise charging options

Batteries engineered to last 6 years. Your clients can enjoy a full day of hearing.




Charger Case Combi



Power pack option



Charger BTE RIC



Mini Charger case

# Paradise has a space for all of your devices-






# Paradise has **two times** the Bluetooth!



## Target 7.0 workflow



# Landing screen, setup and navigation

The screenshot displays the Phonak Target 7.0.5 software interface. At the top, a dark grey header bar contains the text "Phonak Target 7.0.5" and menu options "File", "eService", and "Help". Below this is a light grey bar labeled "Client & session" which contains the text "Paradise Welcome", highlighted by a red rectangular box. The main content area is titled "Noah audiogram" and features two line graphs for the right (R) and left (L) ears, both dated 8/31/2020. To the right of these graphs, a message states "No fitting data available" with an orange button labeled "New fitting session". Below the audiogram section, there is a section for "eSolutions" with a warning icon, a "PhonakPro login..." button, and a "CLIENT DASHBOARD" section with an "Open client dashboard" button. At the bottom of the interface is a dark grey navigation bar with several icons and labels: "Client & session" (highlighted in green), "Demonstrator", "Media", "News", "Trial & tools", "Updates" (with a blue notification bubble), and "Setup". In the top right corner of the main content area, there are additional buttons for "New fitting session" and "Reports...".

# Landing screen, setup and navigation

The screenshot displays the Phonak Target 7.0.5 software interface. At the top, a dark header bar contains the text "Phonak Target 7.0.5" and menu items "File", "eService", and "Help". Below this is a light gray bar with the text "Client & session".

The main content area is divided into sections. The first section is titled "Paradise Welcome". The second section is titled "Noah audiogram" and contains two line graphs. The left graph is labeled "R" and the right graph is labeled "L", both dated "8/31/2020".

Below the audiogram section, there is a "CLIENT DASHBOARD" section with a button labeled "Open client dashboard". To the right of this section is a "PhonakPro login..." button.

A red rounded rectangle highlights a message box on the right side of the screen. The message box contains the text "No fitting data available" and a button labeled "New fitting session".


At the bottom of the interface is a dark navigation bar with several icons and labels: "Client & session" (highlighted in green), "Demonstrator", "Media", "News", "Trial & tools", "Updates", and "Setup".

In the top right corner of the main content area, there are two buttons: "New fitting session" and "Reports...".

# Landing screen, setup and navigation

The screenshot displays the Phonak Target 7.0.5 software interface. At the top, a dark header bar contains the text "Phonak Target 7.0.5" and menu items "File", "eService", and "Help". Below this, a light gray bar shows "Client & session" with window control icons on the right. The main content area is divided into sections: "Paradise Welcome" with a horizontal line, "Noah audiogram" with two line graphs for the right (R) and left (L) ears dated 8/31/2020, and "eSolutions" with a warning icon. A "PhonakPro login..." button is positioned to the right of the eSolutions section. Below the audiogram, a "CLIENT DASHBOARD" label is next to an "Open client dashboard" button. On the right side, a message states "No fitting data available" above a "New fitting session" button. At the bottom, a dark navigation bar features icons for "Client & session" (highlighted in green), "Demonstrator", "Media", "News", "Trial & tools", "Updates", and "Setup". A red rectangle highlights the "Client & session" icon and the "New fitting session" button. In the top right corner of the main content area, there are dropdown menus for "New fitting session" and "Reports...".

# Trial and tools device setup



# Trial and tools device setup

Phonak Target 6.1

Trial & tools

Trial hearing aids

Assistant to change the performance level of trial hearing aids.

Configure...

### Configure trial hearing aid

1 Preparation 2 Selection 3 Configure device 4 Finish

Select desired performance level.

	30 Essential	50 Standard	70 Advanced	90 Premium
AutoSense OS 3.0	✓	✓	✓	✓
RogerDirect	✓	✓	✓	✓
Speech in loud noise			✓	✓
Speech in car				✓
Speech in echo				✓
Speech in 360°			✓	✓
Tinnitus balance	✓	✓	✓	✓
SoundRecover2	✓	✓	✓	✓
UltraZoom	✓	✓	✓	✓
DuoPhone		✓	✓	✓
Real ear sound		✓	✓	✓

Select Select Select Select (Current level)

Start over Continue Cancel

17

Client & session Demonstrator Media News **Trial & tools** Updates Setup

# Landing screen, setup and navigation

The screenshot displays the Phonak Target 7.0.5 software interface. At the top, a dark header bar contains the text "Phonak Target 7.0.5" and menu items "File", "eService", and "Help". Below this is a light gray bar with the text "Client & session". The main content area is divided into sections: "Paradise Welcome" with a horizontal line, "Noah audiogram" with two line graphs for the right (R) and left (L) ears, both dated 8/31/2020. To the right of the audiogram is a message "No fitting data available" and an orange button labeled "New fitting session". Below the audiogram is an "eSolutions" section with a warning icon and a "PhonakPro login..." button. At the bottom of the main area is a "CLIENT DASHBOARD" section with an "Open client dashboard" button. A dark footer bar at the bottom contains several navigation icons: "Client & session" (highlighted in green), "Demonstrator", "Media", "News", "Trial & tools", "Updates", and "Setup" (highlighted with a red box). On the right side of the footer bar, there are buttons for "New fitting session" and "Reports...".



# Fitting session presets

Phonak Target 6.1 File eService Help

Setup

General My profile **Fitting session** Junior mode Reports Fitting device Sound system Internet

**Diagnostics**

Instruments

Fitting

Easy navigation

AudiogramDirect

Presentation

Diagnostics

Default measurement condition Foam tip

Bone conduction symbol

Left -->-- Right --<--<

Left --<--< Right -->-->

Use Auditory Brainstem Response (ABR)

DSL v5a - ABR nHL correction values

250	500	750	1k	1.5k	2k	3k	4k	6k	8k	Hz
30	20	17	15	12	10	7	5	5	5	dB

Reset

Client & session Demonstrator Media News Trial & tools Updates **Setup**



# Fitting workflow and navigation

Phonak Target 7.0.5 File eService Help

Client & session


Paradise Welcome

Noah audiogram

No fitting data available

[New fitting session](#)

eSolutions  [PhonakPro login...](#)

CLIENT DASHBOARD [Open client dashboard](#)


[New fitting session](#) [Reports...](#)

**Client & session** Demonstrator Media News Trial & tools Updates Setup

# Target software training option

The screenshot displays the Phonak Target 6.1 software interface. At the top, the title bar reads "Phonak Target 6.1" with menu options "File", "eService", and "Help". Below this, the "Client & session" section shows the name "Smith Jenny". The "Noah audiogram" section contains two small line graphs, one for the right ear (R) and one for the left ear (L), both dated "12/27/2019". In the "eSolutions" section, there is a "PhonakPro login..." button and a message "No fitting data available" with a "New fitting session" button. At the bottom, a navigation bar features icons for "Client & session", "Demonstrator", "Media", "News", "Trial & tools", "Updates", and "Setup". A dropdown menu is open over the "New fitting session" button, showing "New fitting session" and "New training session", with the latter option highlighted by a red rounded rectangle.

# Basic navigation



# Fitting Navigation

The screenshot shows the Phonak Target 7.0.5 software interface. The top menu bar includes 'File', 'Fitting', 'eService', and 'Help'. The main window is divided into several sections:

- Client Information:** Displays 'Client Paradise, Welcome'.
- Instruments:** A green header section with sub-tabs for 'Details' and 'Accessories'. Below this, there are two horizontal lines representing hearing aid data, with the text 'No hearing aid' displayed on the right.
- Styles List:** A list of hearing aid styles with their corresponding page numbers:

▶ Ambra	12
▶ Audéo P	16
▶ Audéo M	20
▶ Audéo B	22
▶ Audéo V	16
▶ Audéo Q	12
▶ Audéo S	10
▶ Bolero M	8
▶ Bolero B	15
▶ Bolero V	12
▶ Bolero Q	16
▶ Cassia	12
▶ CROS	13
- Preview:** A section labeled 'Preview' with the text 'No device to preview'.

A callout box on the left side of the interface contains the text: 'All client information, such as personal data and audiogram, are found on this tab.'

# Fitting Navigation

The screenshot shows the Phonak Target 7.0.5 software interface. The top menu bar includes 'File', 'Fitting', 'eService', and 'Help'. The main window is divided into several sections: 'Client' (Paradise, Welcome), 'Instruments' (highlighted in green), and 'Fitting'. A tooltip is displayed over the 'Instruments' section, stating: 'All hearing aids, acoustic parameters, remote controls and other accessories are found here. Note: Hover over hearing aid icon to get more information about battery state of charge (rechargeable only) and Roger license.' Below the tooltip, a list of hearing aid models is shown under the 'Family' dropdown menu. The list includes models like Ambra, Audéo P, Audéo M, Audéo B, Audéo V, Audéo Q, Audéo S, Bolero M, Bolero B, Bolero V, Bolero Q, Cassia, and CROS, each with a corresponding number. A 'Show only new hearing aids' checkbox is also visible.

Family	Count
▶ Ambra	12
▶ Audéo P	16
▶ Audéo M	20
▶ Audéo B	22
▶ Audéo V	16
▶ Audéo Q	12
▶ Audéo S	10
▶ Bolero M	8
▶ Bolero B	15
▶ Bolero V	12
▶ Bolero Q	16
▶ Cassia	12
▶ CROS	13

# Fitting Navigation

The screenshot displays the Phonak Target 7.0.5 software interface. The top menu bar includes 'File', 'Fitting', 'eService', and 'Help'. A secondary bar contains 'Client view', 'Select media...', and 'Save & close session'. The main interface is divided into several sections:

- Client:** Shows 'Paradise, Welcome'.
- Instruments:** The active section, containing sub-tabs for 'Hearing aids', 'Acoustic parameters', 'Details', and 'Accessories'. The 'Hearing aids' sub-tab is selected, displaying 'No hearing aid'.
- Preview:** Shows 'No device to preview'.
- Device List:** A table listing various hearing aid models and their corresponding page numbers.



A callout box on the right side of the 'Instruments' section contains the text: "All device adjustments are done here."

Family	Page
▶ Ambra	12
▶ Audéo P	16
▶ Audéo M	20
▶ Audéo B	22
▶ Audéo V	16
▶ Audéo Q	12
▶ Audéo S	10
▶ Bolero M	8
▶ Bolero B	15
▶ Bolero V	12
▶ Bolero Q	16
▶ Cassia	12
▶ CROS	13

Show only new hearing aids



- Feedback & real ear test
- AudiogramDirect
- Global tuning
- Fine tuning
- DataLogging
- Device options



No "Feedback & real ear test" measured.

Start both measurements

No "Feedback & real ear test" measured.

Overtuning provides no additional benefit.

OVERTUNING

Higher gain limit  
Higher feedback risk

Better sound quality  
Low battery consumption

Off

Overtuning provides no additional benefit.





Connect 75% 75%

Navigation icons: Home, Ear, Microphone, Speaker, Headphones, Hand, Headset

Client Another day, in Paradise Instruments Audéo P90-R

Fitting Calm situation

Feedback & real ear test AudiogramDirect Global tuning Fine tuning DataLogging Device options



No "Feedback & real ear test" measured.

Start both measurements

No "Feedback & real ear test" measured.

Overtuning provides no additional benefit.

OVERTUNING Off

Higher gain limit  
Higher feedback risk

Better sound quality  
Low battery consumption

OVERTUNING Off

Overtuning provides no additional benefit.



Connect



Client Paradise, Welcome

Instruments Audéo P90-R

Fitting Calm situation

Hearing aids Acoustic parameters Details Accessories

Audéo P90-R (P)



Gain 59 (2cc)  
MPO 122 (2cc)  
Receiver P

Acoustic parameters

Audéo P90-R (P)



Gain 59 (2cc)  
MPO 122 (2cc)  
Receiver P

Acoustic parameters

Family

Styles

Features

Search...

Ambra 12

Audéo P 16

Audéo P90-R Premium 71/132 (2cc)

Audéo P70-R Advanced 71/132 (2cc)

Audéo P50-R Standard 71/132 (2cc)


Audéo P30-R Essential 71/132 (2cc)

Audéo P90-RT Premium 71/132 (2cc)

Show only new hearing aids

Audéo P90-R

RECEIVERS




DESCRIPTION

Premium WaterResistant external receiver instrument (rechargeable) for mild to profound hearing loss. Direct connectivity (Bluetooth®).

FEATURES

- AutoSense OS 4.0 premium
- Speech enhancer
- Comfort in echo
- Rechargeable
- Own voice pick-up
- Tinnitus balance
- Dynamic noise cancellation
- Environmental balance
- SoundRecover2
- Tap control
- RogerDirect
- UltraZoom premium
- Motion sensor hearing
- StereoZoom
- DuoPhone





Client Paradise, Welcome

Instruments Audéo P90-R Fitting Calm situation

Hearing aids Acoustic parameters

**(R)** Audéo P90-R (P)



Gain 59 (2cc)  
MPO 122 (2cc)  
Receiver P

Family  
▶ Ambra  
▼ Audéo P

-   Audéo P90-R  
Premium 71/132 (2cc)
-   Audéo P70-R  
Advanced 71/132 (2cc)
-   Audéo P50-R  
Standard 71/132 (2cc)
-   Audéo P30-R  
Essential 71/132 (2cc)
-   Audéo P90-RT  
Premium 71/132 (2cc)





Show only new hearing aids

### Connect


**?** Select hearing aids

**(R)** No hearing aid **(L)** No hearing aid

**(C)** Detecting...

	Audéo P90-R SN: SN69390788		<b>(R)</b> Select	<b>(L)</b> Select
	Audéo P90-R SN: SN36310893		<b>(R)</b> Select	<b>(L)</b> Select

**(i)** If your hearing aid is not detected, turn it off and on again.

 Press the button on the hearing aid to highlight it in the list.

Continue Cancel

→ Acoustic parameters




UP  
32 (2cc)

rechargeable) for mild to profound hearing loss.

- Autosense OS 4.0 premium
- Speech enhancer
- Comfort in echo
- Rechargeable
- Own voice pick-up
- Tinnitus balance
- Dynamic noise cancellation
- Environmental balance
- SoundRecover2
- Tap control
- RogerDirect
- UltraZoom premium
- Motion sensor hearing
- StereoZoom
- DuoPhone





Client Paradise, Welcome

Instruments Audéo P90-R Fitting Calm situation

Hearing aids Acoustic parameters

**(R)** Audéo P90-R (P)

Gain 59 (2cc)  
MPO 122 (2cc)  
Receiver P

Family ▾  
▶ Ambra  
▼ Audéo P

- Audéo P90-R Premium 71/132 (2cc)
- Audéo P70-R Advanced 71/132 (2cc)
- Audéo P50-R Standard 71/132 (2cc)
- Audéo P30-R Essential 71/132 (2cc)
- Audéo P90-RT Premium 71/132 (2cc)

Show only new hearing aids

### Connect

**Select hearing aids**

<b>(R)</b> Audéo P90-R SN: SN69390788	<b>(L)</b> Audéo P90-R SN: SN36310893
<b>Detecting...</b>	
Audéo P90-R SN: SN69390788	<b>(R) Selected</b> <b>(L) Select</b>
Audéo P90-R SN: SN36310893	<b>(R) Select</b> <b>(L) Selected</b>

**(i)** If your hearing aid is not detected, turn it off and on again.

Press the button on the hearing aid to highlight it in the list.

**Continue** **Cancel**

- AutoSense OS 4.0 premium
- Speech enhancer
- Comfort in echo
- Rechargeable
- Own voice pick-up
- Limbus balance
- Dynamic noise cancellation
- Environmental balance
- SoundRecover2
- Tap control
- RogerDirect
- UltraZoom premium
- Motion sensor hearing
- StereoZoom
- DuoPhone



Client Paradise, Welcome

Instruments Audéo P90-R Fitting Calm situation (0)

Hearing aids Acoustic parameters

### Connect

✓ Connection successful ✓ Connection successful

Ⓡ  Audéo P90-R SN: SN69390788 Battery level: Good	Ⓛ  Audéo P90-R SN: SN36310893 Battery level: Good
---------------------------------------------------------	---------------------------------------------------------

Unmute both

Close

Ⓡ Audéo P90-R (P) SN: SN69390788

Gain 59 (2cc)  
MPO 122 (2cc)  
Receiver P

Family ▾

► Ambra

▼ Audéo P

Audéo P90-R  
Premium 71/132 (2cc)

- Audéo P70-R  
Advanced 71/132 (2cc)
- Audéo P50-R  
Standard 71/132 (2cc)
- Audéo P30-R  
Essential 71/132 (2cc)
- Audéo P90-RT  
Premium 71/132 (2cc)

Show only new hearing aids

S	M	P	UP
47/111 (2cc)	51/114 (2cc)	59/122 (2cc)	71/132 (2cc)


**DESCRIPTION**  
Premium WaterResistant external receiver instrument (rechargeable) for mild to profound hearing loss. Direct connectivity (Bluetooth®).

- FEATURES**
- |                            |                              |                         |
|----------------------------|------------------------------|-------------------------|
| • AutoSense OS 4.0 premium | • Tinnitus balance           | • RogerDirect           |
| • Speech enhancer          | • Dynamic noise cancellation | • UltraZoom premium     |
| • Comfort in echo          | • Environmental balance      | • Motion sensor hearing |
| • Rechargeable             | • SoundRecover2              | • StereoZoom            |
| • Own voice pick-up        | • Tap control                | • DuoPhone              |



Client Paradise, Welcome Instruments Fitting

Hearing aids Acoustic parameters Details Accessories



Family Styles Features

- Audéo P
- Audéo P90-R Premium 71/132 (2cc)
- Audéo P70-R Advanced 71/132 (2cc)
- Audéo P50-R Standard 71/132 (2cc)
- Audéo P30-R Essential 71/132 (2cc)
- Audéo P90-RT Premium 71/132 (2cc)
- Audéo P70-RT Advanced 71/132 (2cc)
- Audéo P50-RT Standard 71/132 (2cc)

Show only new hearing aids


Phonak Target

Create new fitting

Optimize gain for

First time user 0-3 months

OK



DESCRIPTION  
Premium WaterResistant external receiver instrument (rechargeable) for mild to profound hearing loss. Direct connectivity (Bluetooth®).

FEATURES

- AutoSense OS 4.0 premium
- Speech enhancer
- Comfort in echo
- Rechargeable
- Own voice pick-up
- Tinnitus balance
- Dynamic noise cancellation
- Environmental balance
- SoundRecover2
- Tap control
- RogerDirect
- UltraZoom premium
- Motion sensor hearing
- StereoZoom
- DuoPhone

More information

Connect 75% 75%

Client view Select media... Close training

Client Paradise, Welcome



Instruments Audéo P90-R

Fitting Calm situation

Hearing aids Acoustic parameters Details Accessories

R Audéo P90-R (P) SN: SN69390788

L Audéo P90-R (P) SN: SN36310893




Acoustic code	XXXXXX	
Earpiece	Power dome	
Dome size	Small	
Receiver	P	Check
Wire length	OP R	

Acoustic code	XXXXXX	
Earpiece	Power dome	
Dome size	Small	
Receiver	P	Check
Wire length	OP L	



# Target 7.0: Acoustic parameters update



Client fitting, transfer   Instruments Audéo P90-R   Fitting Calm situation

Hearing aids   **Acoustic parameters**   Details   Accessories

Audéo P90-R (M) SN: 2002H0612      Audéo P90-R (M) SN: 2002H0612

Acoustic code: XXXXXX

Earpiece: SlimTip

Vent: Cavity vent Ø 0.9 - 1.1 mm

Receiver: M      Check

Wire length: OM R

Occluded

Open

**NORMAL VENT**

Ø 0.9 - 1.1 mm

Ø 1.2 - 1.6 mm

Ø 1.7 - 2.2 mm

Ø 2.3 - 3.0 mm

Ø 3.1 - 4.5 mm

**CAVITY VENT**


Ø 0.9 - 1.1 mm

Ø 1.2 - 1.6 mm

Ø 1.7 - 2.2 mm


Ø 2.3 - 3.0 mm





No "Feedback & real ear test" measured.

**R** Start both measurements **L**



No "Feedback & real ear test" measured.

OVERTUNING  Higher gain limit Higher feedback risk  OVERTUNING

Off  Better sound quality Low battery consumption  Off




Client Paradise, Welcome

Instruments Audéo P90-R

Fitting Calm situation

Feedback & real ear test AudiogramDirect Global tuning Fine tuning DataLogging Device options



No "Feedback & real ear test" measured.

Start both measurements

No "Feedback & real ear test" measured.

OVERTUNING

Higher gain limit  
Higher feedback risk



Better sound quality  
Low battery consumption

Off

Connect 75% 75%

Client Paradise, Welcome Instruments Audéo P90-R Fitting Calm situation

Feedback & real ear test AudiogramDirect Global tuning Fine tuning DataLogging Device options



Use test result to predict vent  Start both measurements  Use test result to predict vent

OVERTUNING


Higher gain limit  
Higher feedback risk

Better sound quality  
Low battery consumption


OVERTUNING



Successfully measured



Successfully measured



Start both measurements

More than 5 dB overtuning can lead to distortion.

OVERTUNING

Higher gain limit  
Higher feedback risk

Better sound quality  
Low battery consumption

OVERTUNING

Off

+5 steps

Off


# Target 7.0: Interaction with WhistleBlock

The screenshot displays the Phonak fitting software interface. At the top, it shows 'Client: Target, Phonak', 'Instruments: Audéo P90-R', and 'Fitting: Calm situation'. The main navigation bar includes 'Feedback & real ear test', 'AudiogramDirect', 'Global tuning', 'Fine tuning', 'DataLogging', and 'Device options'. On the left, the 'Program manager' lists various programs, with 'WhistleBlock' highlighted in a red box. The central area shows the 'WhistleBlock' feature with a value of 20, labeled 'Strong'. Below it are sliders for 'Speech enhancer' (14, Moderate), 'SoundRelax' (8, Weak), 'NoiseBlock' (8, Weak), 'Soft noise reduction' (0, Off), and 'WindBlock' (16, Moderate). The right side shows sliders for 'Real ear sound' (4), 'Speech enhancer' (14, Moderate), 'NoiseBlock' (13, Moderate), 'SoundRelax' (8, Weak), 'NoiseBlock' (8, Weak), 'Soft noise reduction' (0, Off), and 'WindBlock' (16, Moderate). The bottom navigation bar includes 'Gain & MPO', 'Audibility fine tuning', 'Program options', 'SoundRecover2', 'Automatic fine tuning', and 'Tinnitus balance'.


# AudiogramDirect

The screenshot displays the Phonak Target 7.0.5 software interface. At the top, the title bar reads "Phonak Target 7.0.5" and includes menu items: File, Fitting, eService, and Help. On the right side of the title bar, there are buttons for "Client view", "Select media...", "Save & close session", and window control icons. Below the title bar is a toolbar with various icons, including a "Connect" button and two "100% TRIAL" indicators. The main interface area has a header with "Client Alyssa, Paradise" and "Instruments Audéo P90-RT". A green navigation bar contains tabs: "Feedback & real ear test", "AudiogramDirect" (which is selected), "Global tuning", "Fine tuning", "DataLogging", and "Device options". The "AudiogramDirect" tab is active, showing a "History" dropdown menu with options for "AC" (red circle) and "UCL" (red square), and a "Start" button. Below this, there are two large empty boxes, each containing the text "No AudiogramDirect measured" and a "Start" button. The "Fitting" section on the right shows "Calm situation" with a speaker icon.

# AudiogramDirect



# AudiogramDirect





# AudiogramDirect

Phonak Target 7.0.5 File Fitting eService Help Client view Select media Save & close

Client Alyssa Paradise

Connect 100% TRIAL 100% TRIAL

**AudiogramDirect**

1 Preparation 2 Acoustic parameters 3 Measure AudiogramDirect

Muted

0 20 40 60 80 100 120

125 250 500 1k 2k 4k 8k

Mandatory frequency

AC UCL

Play signal [Space]

Save point [Enter]

Talk to client [T]

Enable HL > 100 dB [E]

Pulsed tone (AC only)

Save and proceed with other side Save & close Cancel

# AudiogramDirect

The screenshot shows the Phonak Target 7.0.5 software interface. The main window is titled "AudiogramDirect" and is divided into three steps: 1. Preparation, 2. Acoustic parameters, and 3. Measure AudiogramDirect. The graph displays the following data:


Frequency (Hz)	AC (dB)	UCL (dB)
125	-	~50
250	~25	~20
500	~25	~20
1k	~25	~20
2k	~20	~20
4k	~25	~20
8k	~25	~20

The graph also shows a "Muted" area and a "No Audiogram" area. The control panel on the right includes the following options:

- AC
- UCL
- Play signal [Space]
- Save point [Enter]
- Talk to client [T]
- Enable HL > 100 dB [E]
- Pulsed tone (AC only)

At the bottom of the window, there are buttons for "Save and proceed with other side", "Save & close", and "Cancel".

# AudiogramDirect





Connect

Client Paradise, Welcome Instruments Audéo P90-R

Fitting Calm situation

Feedback & real ear test AudiogramDirect Global tuning Fine tuning DataLogging Device options

FITTING FORMULA  
Adaptive Phonak Digital 2.0



GAIN LEVEL  
80%

Auto acclimatization

OCCUSION COMPENSATION  
Off

COMPRESSION  
Prescribed compression

Initial fitting SoundRecover Tinnitus balance



Connect 75% 75%



Client Paradise, Welcome



Instruments Audéo P90-R

Fitting Calm situation

Feedback & real ear test AudiogramDirect Global tuning Fine tuning DataLogging Device options

FITTING FORMULA

- Adaptive Phonak Digital 2.0
- DSL [i/o]
- DSL v5a Adult
- DSL v5a Pediatric
- NAL-NL 1
- NAL-NL 2
- NAL-NL 2 Tonal
- Adaptive Phonak Digital 2.0
- Adaptive Phonak Digital Contrast 2.0



GAIN LEVEL  
80%

Auto acclimatization

OCCUSION COMPENSATION  
Off

COMPRESSION  
Prescribed compression

Initial fitting SoundRecover Tinnitus balance

## Adaptive Phonak Digital 2.0



**New**  
fitting formula



### Hearing inspired by nature

1

Adaptive  
compression  
speeds

2

Linearized  
high level gain

3


Pre-calculation  
for mild-moderate  
hearing losses



# 1. Adaptive compression speeds – Intelligent compression


## Advantages

- More dynamic range<sup>1</sup>
- Reduction of perceived reverberation<sup>1</sup>
- Comfort in noise<sup>2</sup>
- Comfort for loud impulse sounds<sup>3</sup>



1. Hassager, H. G., Winberg, A., and Dau, T. (2017a). "Effects of hearing-aid dynamic range compression on spatial perception in a reverberant environment," J. Acoust. Soc. Am. 141, 2556–2568.
2. Kowalewski, B., Zaar, J., Fereczkowski, M., MacDonald, E., Strelcyk, O., May, T., & Dau, T. (2017). Effects of slow-and fast-acting compression on hearing impaired listeners' consonant-vowel identification in interrupted noise. In *International Symposium on Auditory and Audiological Research*. The Danavox Jubilee Foundation.
3. Lopez-Poveda, E. A., Johannesen, P. T., Perez-González, P., Blanco, J. L., Kalluri, S., & Edwards, B. (2017). Predictors of hearing-aid outcomes. *Trends in hearing*, 21, 1-28.

## 2. Linearized high level gain



### Advantages

- Loudness comfort at high input levels<sup>1</sup>
- Speech intelligibility for loud speech in noise<sup>1</sup>
- More fine-tuning flexibility

— APD  
— APD 2.0

1. Lopez-Poveda, E. A., Johannesen, P. T., Perez-González, P., Blanco, J. L., Kalluri, S., & Edwards, B. (2017). Predictors of hearing-aid outcomes. Trends in hearing, 21, 1-28.



### 3. New pre-calculation

- Mild to moderate hearing losses
- G50 targets remains unchanged
- G65 and G80 targets are reduced for mid-loud inputs
- All experience levels

Senn. M., (2020). An analysis of – fine-tuning applied at first fit. Phonak Target Fitting Data. Unpublished data. Contact [claims@phonak.com](mailto:claims@phonak.com) if you are interested in further information.



Connect 75% 75%



📍 🧠 🗂️ 🛠️ 📄 🗑️ 🗨️

Client Paradise, Welcome Instruments Audéo P90-R

Fitting Calm situation

Feedback & real ear test AudiogramDirect Global tuning Fine tuning DataLogging Device options

FITTING FORMULA  
Adaptive Phonak Digital 2.0



GAIN LEVEL  
80%

Auto acclimatization

OCCUSION COMPENSATION  
Off

COMPRESSION  
Prescribed compression

Initial fitting SoundRecover Tinnitus balance



Connect

Client Paradise, Welcome Instruments Audéo P90-R

Fitting Calm situation

Feedback & real ear test AudiogramDirect Global tuning Fine tuning DataLogging Device options

FITTING FORMULA  
Adaptive Phonak Digital 2.0



GAIN LEVEL  
84%

Auto acclimatization  
Target 100 %  
Remaining 100 days

OCCUSION COMPENSATION  
Off

COMPRESSION  
Prescribed compression


Initial fitting SoundRecover Tinnitus balance




Program manager >>

ALL PROGRAMS

- AUTOMATIC PROGRAMS 9
  - A** AutoSense OS 4.0
    - Calm situation
    - Speech in noise
    - Speech in loud noise
    - Speech in car
    - Comfort in noise
    - Comfort in echo
    - Music
  - AutoSense OS 4.0 (streaming)
    - Media speech + mic
    - Media music + mic
- STREAMING PROGRAMS 3
  - PartnerMic + mic
  - Phone call + mic
  - RogerDirect + mic
- ADDITIONAL PROGRAMS +



MPO	92	105	109	113	118	124	123	120	114	106
All	250	500	750	1k	1.5k	2k	3k	4k	6k	8k
G80	4	15	19	23	29	39	41	39	39	35
G65	8	21	26	31	37	46	47	45	44	41
G50	11	26	32	37	44	50	49	46	45	41
CR	1.5	1.6	1.7	1.8	1.9	1.6	1.4	1.3	1.3	1.3



MPO	91	105	110	114	117	122	123	120	115	107
All	250	500	750	1k	1.5k	2k	3k	4k	6k	8k
G80	3	15	19	23	28	37	40	38	37	33
G65	8	21	26	30	36	45	46	44	40	36
G50	11	26	31	36	42	49	49	45	40	36
CR	1.5	1.6	1.7	1.8	1.9	1.7	1.4	1.3	1.1	1.1









Client Paradise, Welcome

Instruments Audéo P90-R

Fitting Calm situation

Feedback & real ear test AudiogramDirect Global tuning Fine tuning DataLogging Device options

- Program manager >>
- ALL PROGRAMS
    - AUTOMATIC PROGRAMS (9)
      - AutoSense OS 4.0
        - Calm situation
        - Speech in noise
        - Speech in loud noise
        - Speech in car
        - Comfort in noise
        - Comfort in echo
        - Music
      - AutoSense OS 4.0 (streaming)
        - Media speech + mic
        - Media music + mic
    - STREAMING PROGRAMS (3)
      - PartnerMic + mic
      - Phone call + mic
      - RogerDirect + mic
  - ADDITIONAL PROGRAMS (+)



MPO	92	105	109	113	118	124	123	120	114	106
All	250	500	750	1k	1.5k	2k	3k	4k	6k	8k
G80	4	15	19	23	29	39	41	39	39	35
G65	8	21	26	31	37	46	47	45	44	41
G50	11	26	32	37	44	50	49	46	45	41
CR	1.5	1.6	1.7	1.8	1.9	1.6	1.4	1.3	1.3	1.3

MPO	91	105	110	114	117	122	123	120	115	107
All	250	500	750	1k	1.5k	2k	3k	4k	6k	8k
G80	3	15	19	23	28	37	40	38	37	33
G65	8	21	26	30	36	45	46	44	40	36
G50	11	26	31	36	42	49	49	45	40	36
CR	1.5	1.6	1.7	1.8	1.9	1.7	1.4	1.3	1.1	1.1

Gain & MPO Audibility fine tuning Program options SoundRecover2 Automatic fine tuning Tinnitus balance





Client Paradise, Welcome

Instruments Audéo P90-R

Fitting Calm situation

Feedback & real ear test AudiogramDirect Global tuning Fine tuning DataLogging Device options

- Program manager >>
- ALL PROGRAMS
    - AUTOMATIC PROGRAMS (9)
      - AutoSense OS 4.0
        - Calm situation
        - Speech in noise
        - Speech in loud noise
        - Speech in car
        - Comfort in noise
        - Comfort in echo
        - Music
      - AutoSense OS 4.0 (streaming)
        - Media speech + mic
        - Media music + mic
    - STREAMING PROGRAMS (3)
      - PartnerMic + mic
      - Phone call + mic
      - RogerDirect + mic
    - ADDITIONAL PROGRAMS (+)



3 6 10 20

Verification assistant

MPO	92	105	109	113	118	124	123	120	114	106
All	250	500	750	1k	1.5k	2k	3k	4k	6k	8k
G80	4	15	19	23	29	39	41	39	39	35
G65	8	21	26	31	37	46	47	45	44	41
G50	11	26	32	37	44	50	49	46	45	41
CR	1.5	1.6	1.7	1.8	1.9	1.6	1.4	1.3	1.3	1.3

MPO	91	105	110	114	117	122	123	120	115	107
All	250	500	750	1k	1.5k	2k	3k	4k	6k	8k
G80	3	15	19	23	28	37	40	38	37	33
G65	8	21	26	30	36	45	46	44	40	36
G50	11	26	31	36	42	49	49	45	40	36
CR	1.5	1.6	1.7	1.8	1.9	1.7	1.4	1.3	1.1	1.1

Gain & MPO Audibility fine tuning Program options SoundRecover2 Automatic fine tuning Tinnitus balance





Client Paradise, Welcome



Instruments Audéo P90-R

Fitting Calm situation

Feedback & real ear test AudiogramDirect Global tuning Fine tuning DataLogging Device options

Program manager >>

- ALL PROGRAMS
- ▼ AUTOMATIC PROGRAMS 9
  - AutoSense OS 4.0
    - Calm situation
    - Speech in noise
    - Speech in loud noise
    - Speech in car
    - Comfort in noise
    - Comfort in echo
    - Music
  - AutoSense OS 4.0 (streaming)
    - Media speech + mic
    - Media music + mic
- ▼ STREAMING PROGRAMS 3
  - PartnerMic + mic
  - Phone call + mic
  - RogerDirect + mic
- ADDITIONAL PROGRAMS +



3 6 10 20

Verification assistant

MPO	92	105	109	113	118	124	123	120	114	106
All	250	500	750	1k	1.5k	2k	3k	4k	6k	8k
G80	4	15	19	23	29	39	41	39	39	35
G65	8	21	26	31	37	46	47	45	44	41
G50	11	26	32	37	44	50	49	46	45	41
CR	1.5	1.6	1.7	1.8	1.9	1.6	1.4	1.3	1.3	1.3

MPO	91	105	110	114	117	122	123	120	115	107
All	250	500	750	1k	1.5k	2k	3k	4k	6k	8k
G80	3	15	19	23	28	37	40	38	37	33
G65	8	21	26	30	36	45	46	44	40	36
G50	11	26	31	36	42	49	49	45	40	36
CR	1.5	1.6	1.7	1.8	1.9	1.7	1.4	1.3	1.1	1.1

Gain & MPO Audibility fine tuning Program options SoundRecover2 Automatic fine tuning Tinnitus balance







Connect 75% 75%



Client Paradise, Welcome



Instruments Audéo P90-R

Fitting Calm situation

Feedback & real ear test AudiogramDirect Global tuning Fine tuning DataLogging Device options

Program manager

- ALL PROGRAMS
- ▼ AUTOMATIC PROGRAMS 9
  - A AutoSense OS 4.0
    - Calm situation
    - Speech in noise
    - Speech in loud noise
    - Speech in car
    - Comfort in noise
    - Comfort in echo
    - Music
- AutoSense OS 4.0 (streaming)
  - Media speech + mic
  - Media music + mic
- ▼ STREAMING PROGRAMS 3
  - PartnerMic + mic
  - Phone call + mic
  - RogerDirect + mic
- ADDITIONAL PROGRAMS +



Select type Own voice - Own voice

Select issue Overall too loud

NEXT CHANGE Reduce gain and increase compression for own voice at low frequencies.

NEXT CHANGE Reduce gain and increase compression for own voice at low frequencies.

Fix issue Undo

Gain & MPO Audibility fine tuning Program options SoundRecover2 Automatic fine tuning Tinnitus balance



Client Paradise, Welcome Instruments Audéo P90-R Fitting Calm situation

Feedback & real ear test AudiogramDirect Global tuning Fine tuning DataLogging Device options

- Program manager >>
- ALL PROGRAMS
- ▼ AUTOMATIC PROGRAMS 9
    - A** AutoSense OS 4.0
      - Calm situation
      - Speech in noise
      - Speech in loud noise
      - Speech in car
      - Comfort in noise
      - Comfort in echo
      - Music
    - AutoSense OS 4.0 (streaming)
      - Media speech + mic
      - Media music + mic
  - ▼ STREAMING PROGRAMS 3
    - PartnerMic + mic
    - Phone call + mic
    - RogerDirect + mic
  - ADDITIONAL PROGRAMS +

Transition speed 2  
Balanced

SPEECH IN NOISE & LOUD NOISE

Motion sensor hearing 3  
Wide

SPEECH IN LOUD NOISE

Program activation  Enable as automatic program

Activation level 10  
Loud noise



Connect 75% 75%

Client Paradise, Welcome Instruments Audéo P90-R Fitting Media speech + mic

Feedback & real ear test AudiogramDirect Global tuning Fine tuning DataLogging Device options

Program manager >>

ALL PROGRAMS

- ▼ AUTOMATIC PROGRAMS 9
  - AutoSense OS 4.0
    - Calm situation
    - Speech in noise
    - Speech in loud noise
    - Speech in car
    - Comfort in noise
    - Comfort in echo
    - Music
    - AutoSense OS 4.0 (streaming)
      - Media speech + mic
      - Media music + mic
- ▼ STREAMING PROGRAMS 3
  - PartnerMic + mic
  - Phone call + mic
  - RogerDirect + mic
- ADDITIONAL PROGRAMS +

"Media streaming" activation  Automatic

---

"TV Connector" activation
 

- Automatic
- Manual
- Manual (with beep)

Gain & MPO Audibility fine tuning Program options SoundRecover2 Automatic fine tuning



Connect 75% 75%



Client Paradise, Welcome

Instruments Audéo P90-R

Fitting PartnerMic + mic

Feedback & real ear test AudiogramDirect Global tuning Fine tuning DataLogging Device options

Program manager >>

- ALL PROGRAMS
- ▼ AUTOMATIC PROGRAMS 9
  - AutoSense OS 4.0
    - Calm situation
    - Speech in noise
    - Speech in loud noise
    - Speech in car
    - Comfort in noise
    - Comfort in echo
    - Music
  - AutoSense OS 4.0 (streaming)
    - Media speech + mic
    - Media music + mic
- ▼ STREAMING PROGRAMS 3
  - PartnerMic + mic
  - Phone call + mic
  - RogerDirect + mic
- ADDITIONAL PROGRAMS +

Program activation

- Automatic
- Manual
- Manual (with beep)

---

Hearing aid microphone

Enable microphone -6 dB -6 dB

---

WhistleBlock *i*

20 Strong

---

SoundRelax

8 Weak

---

NoiseBlock

8 Weak

---

Soft noise reduction

0 Off

---

WindBlock

0 Off

Gain & MPO Audibility fine tuning Program options SoundRecover2 Automatic fine tuning

Connect

Client Paradise, Welcome Instruments Audéo P90-R Fitting Phone call + mic

< >
Feedback & real ear test
AudiogramDirect
Global tuning
Fine tuning
DataLogging
Device options

- Program manager >>
- ALL PROGRAMS
- ▼ AUTOMATIC PROGRAMS 9
    - AutoSense OS 4.0
      - Calm situation
      - Speech in noise
      - Speech in loud noise
      - Speech in car
      - Comfort in noise
      - Comfort in echo
      - Music
    - AutoSense OS 4.0 (streaming)
      - Media speech + mic
      - Media music + mic
  - ▼ STREAMING PROGRAMS 3
    - PartnerMic + mic
    - Phone call + mic
    - RogerDirect + mic
  - ADDITIONAL PROGRAMS +

Hearing aid microphone -6 dB

---

SoundRelax 8 Weak

NoiseBlock 12 Moderate

Soft noise reduction 0 Off

WindBlock 16 Moderate

---

Bluetooth (i) Bandwidth can be changed here → Bluetooth

Gain & MPO
Audibility fine tuning
Program options
SoundRecover2
Automatic fine tuning



Connect 75% 75%



Client Paradise, Welcome

Instruments Audéo P90-R

Fitting RogerDirect + mic

Feedback & real ear test AudiogramDirect Global tuning Fine tuning DataLogging Device options

Program manager >>

- ALL PROGRAMS
- ▼ AUTOMATIC PROGRAMS 9
  - AutoSense OS 4.0
    - Calm situation
    - Speech in noise
    - Speech in loud noise
    - Speech in car
    - Comfort in noise
    - Comfort in echo
    - Music
  - AutoSense OS 4.0 (streaming)
    - Media speech + mic
    - Media music + mic
- ▼ STREAMING PROGRAMS 3
  - PartnerMic + mic
  - Phone call + mic
  - RogerDirect + mic
- ADDITIONAL PROGRAMS +

Program activation

- Automatic
- Manual
- Manual (with beep)

---

Adaptive microphone behavior

- Standard
- Dual

Enable microphone 0 dB 0 dB

Hearing aid microphone

Fixed directional 12

WhistleBlock 20 Strong

SoundRelax 8 Weak

NoiseBlock 0





Connect 75% 75%



Client Paradise, Welcome

Instruments Audéo P90-R

Fitting RogerDirect + mic

Feedback & real ear test AudiogramDirect Global tuning Fine tuning DataLogging Device options

Program manager >>

- ALL PROGRAMS
- ▼ AUTOMATIC PROGRAMS 9
  - AutoSense OS 4.0
    - Calm situation
    - Speech in noise
    - Speech in loud noise
    - Speech in car
    - Comfort in noise
    - Comfort in echo
    - Music
  - AutoSense OS 4.0 (streaming)
    - Media speech + mic
    - Media music + mic
- ▼ STREAMING PROGRAMS 3
  - PartnerMic + mic
  - Phone call + mic
  - RogerDirect + mic
- ADDITIONAL PROGRAMS +

Program activation

- Automatic
- Manual
- Manual (with beep)

Adaptive microphone behavior

- Standard
- Dual

- Calm situation
- Speech in noise
- Speech in loud noise
- Music
- Acoustic phone
- Comfort in noise
- Speech in 360°
- Comfort in echo
- Custom program
- 1:1 copy of...
- Add multiple programs...

Enable microphone 0 dB 0 dB



12 Fixed directional

20 Strong

8 Weak

NoiseBlock

Gain & MPO Audibility fine tuning Program options SoundRecover2 Automatic fine tuning



Connect 75% 75%



Client Paradise, Welcome

Instruments Audéo P90-R

Fitting RogerDirect + mic

Feedback & real ear test AudiogramDirect Global tuning Fine tuning DataLogging Device options

Program manager >>

- ALL PROGRAMS
- ▼ AUTOMATIC PROGRAMS 9
  - AutoSense OS 4.0
    - Calm situation
    - Speech in noise
    - Speech in loud noise
    - Speech in car
    - Comfort in noise
    - Comfort in echo
    - Music
  - AutoSense OS 4.0 (streaming)
    - Media speech + mic
    - Media music + mic
- ▼ STREAMING PROGRAMS 3
  - PartnerMic + mic
  - Phone call + mic
  - RogerDirect + mic
- ADDITIONAL PROGRAMS +

Program activation

- Automatic
- Manual
- Manual (with beep)

---

Adaptive microphone behavior

- Standard
- Dual

Enable microphone 0 dB 0 dB

Hearing aid microphone

Fixed directional 12

WhistleBlock 20 Strong

SoundRelax 8 Weak

NoiseBlock 0



Connect 75% 75%



Client Paradise, Welcome

Instruments Audéo P90-R

Fitting RogerDirect + mic

Feedback & real ear test AudiogramDirect Global tuning Fine tuning DataLogging Device options

<< Program manager

AUTOMATIC PROGRAMS		9	POSITION	STARTUP	EASY ACCESS	LABEL FOR REPORTS
<b>A</b>	AutoSense OS 4.0		A	<input checked="" type="radio"/>		Automatic
	Calm situation		-	<input type="radio"/>		Calm situation
	Speech in noise		-	<input type="radio"/>		Speech in noise
	Speech in loud noise	<input checked="" type="checkbox"/>	-	<input type="radio"/>		Speech in loud noise
	Speech in car		-	<input type="radio"/>		Speech in car
	Comfort in noise		-	<input type="radio"/>		Comfort in noise
	Comfort in echo		-	<input type="radio"/>		Comfort in echo
	Music		-	<input type="radio"/>		Music
<b>B</b>	AutoSense OS 4.0 (streaming)				<input checked="" type="checkbox"/>	
	Media speech + mic					Media speech
	Media music + mic					Media music
STREAMING PROGRAMS		3	POSITION	STARTUP	EASY ACCESS	LABEL FOR REPORTS
<b>C</b>	PartnerMic + mic				<input checked="" type="checkbox"/>	PartnerMic
<b>D</b>	Phone call + mic				<input checked="" type="checkbox"/>	Phone call
<b>E</b>	RogerDirect + mic				<input checked="" type="checkbox"/>	RogerDirect
ADDITIONAL PROGRAMS		+	POSITION	STARTUP	EASY ACCESS	LABEL FOR REPORTS

- Automatic
- Manual
- Manual (with beep)

- Standard
- Dual

Enable microphone 0 dB 0 dB



12

Fixed directional

20

Strong

8

Weak

SoundRecover2 Automatic fine tuning



Connect 75% 75%



Client Paradise, Welcome

Instruments Audéo P90-R


Fitting Calm situation

Feedback & real ear test AudiogramDirect Global tuning Fine tuning DataLogging Device options

Program manager >>

- ALL PROGRAMS
- ▼ AUTOMATIC PROGRAMS 9
  - A AutoSense OS 4.0
    - Calm situation
    - Speech in noise
    - Speech in loud noise
    - Speech in car
    - Comfort in noise
    - Comfort in echo
    - Music
  - AutoSense OS 4.0 (streaming)
    - Media speech + mic
    - Media music + mic
- ▼ STREAMING PROGRAMS 3
  - PartnerMic + mic
  - Phone call + mic
  - RogerDirect + mic
- ADDITIONAL PROGRAMS +

- Hearing aid microphone
- Speech enhancer
- WhistleBlock
- SoundRelax
- NoiseBlock
- Soft noise reduction
- WindBlock







Connect 75% 75%


Client Paradise, Welcome Instruments Audéo P90-R Fitting Calm situation

Feedback & real ear test AudiogramDirect Global tuning Fine tuning DataLogging Device options

Program manager >>

ALL PROGRAMS

- ▼ AUTOMATIC PROGRAMS 9
  - A AutoSense OS 4.0**
    - Calm situation
    - Speech in noise
    - Speech in loud noise
    - Speech in car
    - Comfort in noise
    - Comfort in echo
    - Music
- AutoSense OS 4.0 (streaming)
  - Media speech + mic
  - Media music + mic
- ▼ STREAMING PROGRAMS 3
  - PartnerMic + mic
  - Phone call + mic
  - RogerDirect + mic
- ADDITIONAL PROGRAMS +



Noise generator is not enabled  
→ Enable in "Global tuning"

Gain & MPO Audibility fine tuning Program options SoundRecover2 Automatic fine tuning **Tinnitus balance**



Connect 75% 75%



Client Paradise, Welcome

Instruments Audéo P90-R

Fitting Calm situation



Feedback & real ear test AudiogramDirect Global tuning Fine tuning DataLogging Device options

FITTING FORMULA

Adaptive Phonak Digital 2.0

TINNITUS BALANCE SHAPE

- Set to hearing loss
- Set to hearing loss
- Set to white noise
- Set to pink noise



Enable noise generator

Noise generator maximum: 77 dB(A)

▲ 0 ▼

▶ 0 ◀

Noise generator maximum: 77 dB(A)

Initial fitting

SoundRecover

Tinnitus balance






- STATISTICS
  - Overview
  - Program activity
- HISTORY
  - Automatic programs
  - Automatic programs (streaming)
  - Streaming programs
- ADAPTATION OF LOGGED DATA
  - Automatic programs
  - Automatic programs (streaming)
  - Streaming programs

**R** Last session

DataLogging was enabled → DataLogging options

Total days	30	Start	12/13/2020
Total hours	322	End	1/12/2021
Hours/day	10.7		

Average daily usage time (h/d) **i**




**L** Last session

DataLogging was enabled → DataLogging options

Total days	30	Start	12/13/2020
Total hours	369	End	1/12/2021
Hours/day	12.3		

Average daily usage time (h/d)



Client Paradise, Welcome Instruments Audéo P90-R Fitting Calm situation

Feedback & real ear test AudiogramDirect Global tuning Fine tuning DataLogging Device options

HEARING AIDS R: SN69390788 L: SN36310893

SERVICES Phone compatibility Available apps Check warranty Reports... Product support

Behavior Standard Streaming Phone call



SHORT PRESS Volume ±10 dB (±5 steps)

LONG PRESS Program change AutoSense OS 4.0 (Startup) Button down has same function

DOUBLE TAP R: Start/stop voice assistant L: No function in Standard

BUTTON FUNCTIONS Volume Program change Accept/end call Mute/unmute mic (not for phone call)

TAP CONTROL Accept/end call Pause/resume streaming Start/stop voice assistant Tap control training...

Manual controls Bluetooth Signals & alerts Startup DataLogging



# Behavior configurations

The screenshot shows the 'Behavior' configuration screen in the Phonak software. The interface includes a top navigation bar with 'Client Target, Phonak', 'Instruments Audéo P90-R', and 'Fitting Calm situation'. Below this is a menu with options like 'Feedback & real ear test', 'ECLD', 'AudiogramDirect', 'Global tuning', 'Fine tuning', 'DataLogging', and 'Device options'. The left sidebar lists 'HEARING AIDS' (R: 2002H06GJ, L: 2002H06H2), 'ROGERDIRECT' (Roger receiver status), and 'SERVICES' (Phone compatibility, Available apps, Check warranty, Reports..., Product support). The main area is titled 'Behavior' and has three tabs: 'Standard' (selected), 'Streaming', and 'Phone call'. An image of the hearing aid is shown. The 'SHORT PRESS' setting is 'Volume ±10 dB (±5 steps)'. The 'LONG PRESS' setting is 'Program change AutoSense OS 4.0 (Startup) Button down has same function'. The 'DOUBLE TAP' setting is 'R: Start/stop voice assistant L: No function in Standard'. At the bottom, there are sections for 'BUTTON FUNCTIONS' (Volume, Program change, Accept/end call, Mute/unmute mic) and 'TAP CONTROL' (Accept/end call, Pause/resume streaming, Start/stop voice assistant, Tap control training...). The bottom navigation bar includes 'Manual controls', 'Bluetooth', 'Signals & alerts', 'Startup', and 'DataLogging'.

# Behavior configurations

Client Target, Phonak Instruments Audéo P90-R Fitting Calm situation

Feedback & real ear test ECLD AudiogramDirect Global tuning Fine tuning DataLogging Device options

**HEARING AIDS**  
R: 2002H06GJ  
L: 2002H06H2

**ROGERDIRECT**  
Roger receiver  
R: (02) installed  
(03) installed  
L: (03) installed  
(02) installed

**SERVICES**  
Phone compatibility  
Available apps  
Check warranty  
Reports...  
Product support

**Behavior** Standard Streaming Phone call

**SHORT PRESS**  
Volume  
±10 dB (±5 steps)

**LONG PRESS**  
Program change  
AutoSense OS 4.0 (Startup)  
Button down has same function

**DOUBLE TAP**  
R: Start/stop voice assistant  
L: No function in Standard

**BUTTON FUNCTIONS**  
 Volume  
 Program change  
 Accept/end call  
 Mute/unmute mic (not for phone call)

**TAP CONTROL**  
 Accept/end call  
 Pause/resume streaming  
 Start/stop voice assistant  
 Tap control training...

Manual controls Bluetooth Signals & alerts Startup DataLogging

# Behavior configurations

Client Target, Phonak Instruments Audéo P90-R Fitting Calm situation

Feedback & real ear test ECLD AudiogramDirect Global tuning Fine tuning DataLogging Device options

HEARING AIDS  
R: 2002H06GJ  
L: 2002H06H2

ROGERDIRECT  
Roger receiver  
R: (02) installed  
(03) installed  
L: (03) installed  
(02) installed

SERVICES  
Phone compatibility  
Available apps  
Check warranty  
Reports...  
Product support

Behavior Standard Streaming Phone call

SHORT PRESS  
Volume  
±10 dB (±5 steps)

LONG PRESS  
Program change  
AutoSense OS 4.0 (Startup)  
Button down has same function

DOUBLE TAP  
R: Start/stop voice assistant  
L: No function in Standard


BUTTON FUNCTIONS  
 Volume  
 Program change  
 Accept/end call  
 Mute/unmute mic (not for phone call)

TAP CONTROL  
 Accept/end call  
 Pause/resume streaming  
 Start/stop voice assistant  
 Tap control training...

Manual controls Bluetooth Signals & alerts Startup DataLogging



# Tap Control



# Tap Control

The screenshot displays the Phonak software interface for configuring a hearing aid. The main window is titled 'Device options' and shows various settings for the device 'Audéo P90-R'. The 'TAP CONTROL' section is highlighted with a red box, indicating the focus of the presentation.

**Client:** Target, Phonak  
**Instruments:** Audéo P90-R  
**Fitting:** Calm situation

**Navigation:** Feedback & real ear test | AudiogramDirect | Global tuning | Fine tuning | DataLogging | **Device options**

**HEARING AIDS:**  
 R: 2002H06GJ  
 L: 2002H06H2

**ROGERDIRECT:**  
 Roger receiver  
 R: (02) installed  
 (03) installed  
 L: (03) installed  
 (02) installed

**SERVICES:**  
 Phone compatibility  
 Available apps  
 Check warranty  
 Reports...  
 Product support

**Behavior:** Standard | Streaming | Phone call

**SHORT PRESS:**  
 Volume  
 ±10 dB (±5 steps)

**LONG PRESS:**  
 Program change  
 AutoSense OS 4.0 (Startup)  
 Button down has same function


**DOUBLE TAP:**  
 R: Start/stop voice assistant  
 L: No function in Standard

**BUTTON FUNCTIONS:**  
 Volume  
 Program change  
 Accept/end call  
 Mute/unmute mic (not for phone call)

**TAP CONTROL (highlighted):**  
 Accept/end call  
 Pause/resume streaming  
 Start/stop voice assistant  
 Tap control training...

**Bottom Navigation:** Manual controls | Bluetooth | Signals & alerts | Startup | DataLogging


myPhonak app



**My Hearing Aids**

Updated with new Paradise features

# Tap Control



# Tap Control

The screenshot displays the Phonak software interface for a client named 'Target, Phonak'. The main window shows the 'Fitting' tab for the 'Audéo P90-R' instrument. A 'Tap control training' dialog box is open, providing instructions for the client to perform tap control training by doing two quick taps at the ear. The dialog indicates that successful double taps are indicated by a green check mark. The dialog shows a counter for 'Successful double taps = 5' with a green check mark, and a 'Stop' button at the bottom right.

**Tap control training for client**

The client can perform tap control training now by doing two quick taps at the ear. Successful double taps are indicated by a green check mark.

Successful double taps = 0

Successful double taps = 5

Stop

**Behavior** Standard Streaming Phone call

**SHORT PRESS**

Volume ±10 dB (±5 steps)

**LONG PRESS**

Program change AutoSense OS 4 Button down has s

**BUTTON FUNCTIONS**

- Volume
- Program change
- Accept/end call
- Mute/unmute mic (not for phone call)

**TAP CONTROL**

- Accept/end call
- Pause/resume streaming
- Start/stop voice assistant
- Tap control training...

Manual controls Bluetooth Signals & alerts Startup DataLogging





Connect



Client Alyssa, Paradise

Instruments Audéo P90-RT

Fitting Calm situation

Feedback & real ear test AudiogramDirect Global tuning Fine tuning DataLogging Device options

HEARING AIDS

R: 2031N2XU1  
L: 2031N2XU5

SERVICES

Phone compatibility

Available apps

Check warranty

Reports...

Product support

Behavior

Standard

Streaming

Phone call



SHORT PRESS

Volume  
±10 dB (±5 steps)

LONG PRESS

Program change  
AutoSense OS 4.0 (Startup)  
Button down has same function

DOUBLE TAP

R: Start/stop voice assistant  
L: No function in Standard

BLUETOOTH NAME

Alyssa Paradise RT

BLUETOOTH SIDE

Right

BLUETOOTH NAME

Alyssa Paradise RT

BLUETOOTH PHONE

Fixed bandwidth

Adaptive bandwidth

Delete pairings

Manual controls

Bluetooth

Signals & alerts

Startup

DataLogging



Connect



Client Alyssa, Paradise

Instruments Audéo P90-RT

Fitting Calm situation

Feedback & real ear test AudiogramDirect Global tuning Fine tuning DataLogging Device options

HEARING AIDS

R: 2031N2XU1  
L: 2031N2XU5

SERVICES

- Phone compatibility
- Available apps
- Check warranty
- Reports...
- Product support

Behavior

Standard

Streaming

Phone call



SHORT PRESS

Volume  
±10 dB (±5 steps)

LONG PRESS

Program change  
AutoSense OS 4.0 (Startup)  
Button down has same function

DOUBLE TAP

R: Start/stop voice assistant  
L: No function in Standard



INTENSITY

Medium

FREQUENCY

Medium

ENABLE

- General alerts
- Volume control
- Program changes
- Ring tones
- Connection alerts

DEMO

- Low battery
- Pairing successful
- Streaming program available
- End of trial period

Manual controls

Bluetooth

Signals & alerts

Startup

DataLogging



Connect



Client Alyssa, Paradise

Instruments Audéo P90-RT

Fitting Calm situation

Feedback & real ear test AudiogramDirect Global tuning Fine tuning DataLogging Device options

HEARING AIDS

R: 2031N2XU1  
L: 2031N2XU5

SERVICES

- Phone compatibility
- Available apps
- Check warranty
- Reports...
- Product support

Behavior

Standard

Streaming

Phone call



SHORT PRESS

Volume  
±10 dB (±5 steps)

LONG PRESS

Program change  
AutoSense OS 4.0 (Startup)  
Button down has same function

DOUBLE TAP

R: Start/stop voice assistant  
L: No function in Standard

STARTUP PROGRAM

AutoSense OS 4.0

STARTUP DELAY

12 seconds additional de. Startup melody

Automatically turn on when removed from charger

Manual controls

Bluetooth

Signals & alerts

Startup

DataLogging



Connect



Client Alyssa, Paradise

Instruments Audéo P90-RT

Fitting Calm situation

Feedback & real ear test AudiogramDirect Global tuning Fine tuning DataLogging Device options

HEARING AIDS

R: 2031N2XU1  
L: 2031N2XU5

SERVICES

- Phone compatibility
- Available apps
- Check warranty
- Reports...
- Product support

Behavior

Standard Streaming Phone call



SHORT PRESS

Volume  
±10 dB (±5 steps)

LONG PRESS

Program change  
AutoSense OS 4.0 (Startup)  
Button down has same function

DOUBLE TAP

R: Start/stop voice assistant  
L: No function in Standard

DATALOGGING

On  
Off

Manual controls

Bluetooth

Signals & alerts

Startup

DataLogging



Client Target, Phonak Instruments Audéo P90-R Fitting Calm situation


Feedback & real ear test ECLD AudiogramDirect Global tuning Fine tuning DataLogging Device options

**HEARING AIDS**  
 R: 2002H06GJ  
 L: 2002H06H2

**ROGERDIRECT**  
 Roger receiver  
 R: (02) installed (03) installed  
 L: (03) installed (02) installed

**SERVICES**  
 Phone compatibility  
 Available apps  
 Check warranty  
 Reports...  
 Product support

**Behavior** Standard Streaming Phone call



**SHORT PRESS**  
 Volume  
 ±10 dB (±5 steps)

**LONG PRESS**  
 Program change  
 AutoSense OS 4.0 (Startup)  
 Button down has same function

**DOUBLE TAP**  
 R: Start/stop voice assistant  
 L: No function in Standard

**BUTTON FUNCTIONS**  
 Volume  
 Program change  
 Accept/end call  
 Mute/unmute mic (not for phone call)

**TAP CONTROL**  
 Accept/end call  
 Pause/resume streaming  
 Start/stop voice assistant  
 Tap control training...

88

Manual controls Bluetooth Signals & alerts Startup DataLogging



# RogerDirect™ status

Client: Paradise, Phonak

Connect

**Audéo P90-R (M)**

SN: 2002H06GJ  
Battery level: High  
Roger receiver: (02) installed  
(03) installed

Connected

Feedback & real ear test

HEARING AIDS

R: 2002H06GJ  
L: 2002H06H2

ROGERDIRECT

Roger receiver  
R: (02) installed  
(03) installed  
L: (03) installed  
(02) installed

SERVICES

- Phone compatibility
- Available apps
- Check warranty
- Reports...
- Product support


Fitting (Calm situation)

Fine tuning DataLogging Device options

ROGER RECEIVER

R: (02) installed  
(03) installed  
L: (03) installed  
(02) installed

# Printing a customized report



The screenshot shows the Phonak Roger app interface. At the top, there are navigation icons and status indicators for 'Connect', battery levels (73% and 74%), and various settings. Below this, the 'Client' is identified as 'Target, Phonak' and 'Instruments' as 'Audéo P90-R'. The 'Fitting' status is 'Calm situation'. A navigation bar includes options like 'Feedback & real ear test', 'ECLD', 'AudiogramDirect', 'Global tuning', 'Fine tuning', 'DataLogging', and 'Device options'. The left sidebar contains sections for 'HEARING AIDS' (listing R: 2002H06GJ and L: 2002H06H2), 'ROGERDIRECT' (listing Roger receiver status for R and L), and 'SERVICES'. The 'SERVICES' menu has a red box around the 'Reports...' option. The main content area is titled 'Behavior' and shows settings for 'Standard', 'Streaming', and 'Phone call'. It includes a diagram of the Roger receiver and lists functions for 'SHORT PRESS' (Volume), 'LONG PRESS' (Program change), and 'DOUBLE TAP' (Start/stop voice assistant). At the bottom, there are sections for 'BUTTON FUNCTIONS' and 'TAP CONTROL' with checkboxes and dropdown menus.



Client Paradise, Welcome Instruments Audéo P90-R Fitting (Calm situation)

Feedback & real ea

HEARING AIDS

R: SN69390788  
L: SN36310893

SERVICES

- Phone compatibility
- Available apps
- Check warranty
- Reports...
- Product support

### Reports

Target audience	Documents for clients	Selected documents
Professional	INSTRUCTIONS	CLIENT
<b>Client</b>	Hearing aid <input checked="" type="checkbox"/>	Hearing aid instructions
	Program description <input checked="" type="checkbox"/>	Program description

Documents setup

English (English)

Use large fonts

Preview Print Email Save

Report setup... Close

Mute/unmute mic (not for phone call) Tap control training...

Hearing aid instructions

Program description



Zoom 55 %




1-2 of 5 pages

Paradise, Welcome  
Client number: 000015

Audio P90-R (SN69307768)  
Audio P90-R (SN69310860)

### Hearing aid instructions

Overview



- 1 Multi-function button
- 2 Microphone input with wind and weather protector
- 3 Battery indicator light
- 4 Right or left indicator (right=red, left=blue), placed in the name plate

---

#### Using the hearing aid

1 Multi-function button

<p>Standard</p> <p>Short press</p> <p><b>Volume</b> ±10 dB (±5 steps)</p> <p>Double tap</p> <p><b>R:</b> Start/stop voice a... <b>L:</b> No function in Sta</p>	<p>Long press</p> <p><b>Program change</b> AutoSens... (Startup) <small>Button down has same function</small></p>
<p>Streaming</p> <p>Short press</p> <p><b>Streaming volume*</b> ±10 dB (±5 steps) <small>*An increase will reduce environmental sounds</small></p> <p>Double tap</p> <p><b>R:</b> Start/stop voice a... <b>L:</b> Pause/resume str...</p>	<p>Long press</p> <p><b>Program change</b> AutoSens... (Startup) Streaming <small>Button down has same function</small></p>

1 / 5 9/14/2020 Read the user guide for further information PHONAK

Paradise, Welcome  
Client number: 000015

Audio P90-R (SN69307768)  
Audio P90-R (SN69310860)

### Hearing aid instructions

Phone call

<p>Short press</p> <p><b>Accept call</b> (during ringing) <small>Button down has same function</small></p> <p><b>Phone volume*</b> ±10 dB (±5 steps) <small>*An increase will reduce environmental sounds</small></p> <p>Double tap</p> <p><b>Accept/end call</b></p>	<p>Long press</p> <p><b>Reject/end call</b> <small>Button down has same function</small></p>
---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	--------------------------------------------------------------------------------------------------

2 / 5 9/14/2020 Read the user guide for further information PHONAK

CLIENT



Zoom 55 %

3-5 of 5 pages

Hearing aid instructions

Program description

Paradise, Welcome  
Client number: 0000015

Audio P90-R (SN6930768)  
Audio P90-L (SN6311083)

Hearing aid instructions

Switching on/off

The hearing aid automatically turns on when removed from the charger. This behavior can be changed with the remote app.

To manually switch the hearing aid on, press the down button for 3 seconds. The indicator light will blink green.

To switch the hearing aid off, press the down button for 4 seconds. The indicator light will turn red.

If the hearing aid does NOT turn on, place it in the charger.

Available programs

**Automatic (startup program) Melody**  
This is the program that automatically adapts to changing environments to provide optimum clarity and comfort in different listening situations.

AutoSense OS 4.0 (streaming)

Automatically identifies when speech or music is present. It seamlessly switches between media speech and media music to optimize your listening experience.

**"TV Connector" and streaming activation**  
When streaming is detected, it is activated automatically.

Streaming programs

**PartnerMic + mic**  
When the PartnerMic is detected, it is activated automatically.

**Phone call + mic**  
This program is optimized for Bluetooth phone calls. The microphone of the hearing aid is used to pick up your voice like a headset.  
A short press on the up or down button on the hearing aid will accept the call. A long press (> 2 seconds) will reject the call. Alternatively, you can accept and reject the call directly on your mobile phone.

**RogerDirect + mic**  
When the RogerDirect is detected, it is activated automatically.

**Pairing mode**  
After turning on, the hearing aid will be in pairing mode for 3 minutes.

**Bluetooth name:**  
L-Phonak hearing aid / R-Phonak hearing aid

**Flight mode**  
To enter flight mode turn the hearing aid off. Then press and hold the down button until the indicator light turns orange.  
To exit flight mode turn the hearing aid off and on again.

Paradise, Welcome  
Client number: 0000015

Audio P90-R (SN6930768)  
Audio P90-L (SN6311083)

Hearing aid instructions

Indicator light

Hearing aid status:	Indicator light:
Turning on	Green blinking
Turning off	Red solid
Charging, battery low	Red blinking
Charging, battery medium	Orange blinking
Charging, battery high	Green blinking
In charger, battery fully charged	Green solid

Streaming functions

**Phone call**  
In combination with your mobile phone, a phone call can be accepted, rejected or ended on the hearing aid.

**Media streaming**  
Streaming from a Bluetooth device to the hearing aid starts automatically.

**Tap control**  
In combination with your hearing aid, a double tap is performed by two quick taps at the ear. This allows you to accept/end a phone call, pause/resume streaming and start/stop voice assistant.

Phonak accessories

Streaming from a Phonak accessory to the hearing aid can be started automatically or manually.

For more information, visit the support page:

<https://bluetooth.phonak.com>



Available apps

For detailed information about available apps, please visit:

[www.phonak.com/availableapps](http://www.phonak.com/availableapps)





# Client notes

The screenshot shows the Phonak software interface with a 'Phonak Target' dialog box open. The dialog box contains the following text:

**Phonak Target**

Notes

+ New note

15.06.2020 09:32: First follow-up appt

First follow-up appt  
Changed power dome from small to medium.

Discuss manual activation behavior for TV Connector at the next appt.


Remind me again (in next session)

Close

The background interface shows two graphs: 'Calm situation' and 'Muted'. The 'Calm situation' graph shows a black line representing the target and several red lines representing the current fit. The 'Muted' graph shows a black line representing the target and several blue lines representing the current fit. The interface also includes a 'Program manager' on the left and a 'Verification assistant' on the right.



Connect



Client Alyssa Paradise

Instruments Audéo P90-RT

Fitting Calm situation

Feedback & real ear test AudiogramDirect Global tuning Fine tuning DataLogging Device options

HEARING AIDS

R: 2031N2XU1  
L: 2031N2XU5

SERVICES

- Phone compatibility
- Available apps
- Check warranty
- Reports...
- Product support

Behavior Standard Streaming Phone call

SHORT PRESS LONG PRESS DOUBLE TAP

Phonak Target

Notes + New note

New session note...

9/2/2020 1:45 PM: Go over wax filters next visit

Remind me again (in next session)

Close

- BUTTON FUNCTIONS
- Volume
  - Program change
  - Accept/end call
  - Mute/unmute mic (not for phone call)

- TAP CONTROL
- Accept/end call
  - Pause/resume streaming
  - Start/stop voice assistant
  - Tap control training...



Connect 100% TRIAL 99% TRIAL



Client Alyssa Paradise



Instruments Audéo P90-RT

Fitting Calm situation

Hearing aids Acoustic parameters Details Accessories

® Audéo P90-RT (M) SN: 2031N2XU1

Ⓛ Audéo P90-RT (M) SN: 2031N2XU5



Ⓛ

Acoustic code	<input type="text" value="XXXXXX"/>	<input type="text" value="XXXXXX"/>
Earpiece	<input type="text" value="Vented dome"/>	<input type="text" value="Vented dome"/>
Dome size	<input type="text" value="Small"/>	<input type="text" value="Small"/>
Receiver	<input type="text" value="M"/> <input type="button" value="Check"/>	<input type="text" value="M"/> <input type="button" value="Check"/>
Wire length	<input type="text" value="OM R"/>	<input type="text" value="OM L"/>



**Notes reminder**  
 You wished to be reminded of the notes from a previous session.  
 → Notes

Client Alyssa, Paradise



Hearing aids



Fitting Calm situation

Ⓡ Audéo P90-RT (M) SN: 2031N2XU1

Ⓛ Audéo P90-RT (M) SN: 2031N2XU5



Acoustic code:

Earpiece:

Dome size:

Receiver:

Wire length:

Acoustic code:



Earpiece:

Dome size:

Receiver:

Wire length:





Phonak Target

Notes + New note

9/2/2020 1:45 PM: Go over wax filters next visit

Remind me again (in next session)

Close

Acoustic code	XXXXXX	XXXXXX
Earpiece	Vented dome	Vented dome
Dome size	Small	Small
Receiver	M Check	M Check
Wire length	0M R	0M L





Connect 100% TRIAL 100% TRIAL

Save and return to Phonak Target start screen

Client Alyssa, Paradise

Instruments Audéo P90-RT

Fitting Calm situation

- Feedback & real ear test
- AudiogramDirect
- Global tuning
- Fine tuning
- DataLogging
- Device options

HEARING AIDS

R: 2031N2XU1  
L: 2031N2XU5

SERVICES

- Phone compatibility
- Available apps
- Check warranty
- Reports...
- Product support

Behavior

- Standard
- Streaming
- Phone call



SHORT PRESS

Volume  
±10 dB (±5 steps)

LONG PRESS

Program change  
AutoSense OS 4.0 (Startup)  
Button down has same function

DOUBLE TAP

R: Start/stop voice assistant  
L: No function in Standard

BUTTON FUNCTIONS

- Volume
- Program change
- Accept/end call
- Mute/unmute mic (not for phone call)

TAP CONTROL

- Accept/end call
- Pause/resume streaming
- Start/stop voice assistant
- Tap control training...

Manual controls

Bluetooth


Signals & alerts

Startup

DataLogging



Connect



Client Alyssa, Paradise

Instruments Audéo P90-RT

Fitting Calm situation

Feedback & real ear test

Close session

Select the items to save

Hearing aids

- L: Audéo P90-RT (SN: 2031N2XU5) Trial period: 6 weeks
- R: Audéo P90-RT (SN: 2031N2XU1) Trial period: 6 weeks

Accessories

No configurable accessory

Database

Optional session note...

Save Cancel

- Accept/end call
- Mute/unmute mic (not for phone call)
- Start/stop voice assistant
- Tap control training...

Manual controls

Bluetooth

Signals & alerts

Startup

DataLogging


HEARING AIDS

R: 2031N2XU1  
L: 2031N2XU5

SERVICES

- Phone compatibility
- Available apps
- Check warranty
- Reports...
- Product support

DOUBLE TAP  
R: Start/stop voice assistant  
L: No function in Standard



Client Alyssa, Paradise

Instruments Audéo P90-RT

Fitting Calm situation

Feedback & real ear test

HEARING AIDS

R: 2031N2XU1  
L: 2031N2XU5

SERVICES

- Phone compatibility
- Available apps
- Check warranty
- Reports...
- Product support

Close session

✓ Successfully saved

Hearing aids

- ✓ L: Audéo P90-RT (SN: 2031N2XU5) Battery level: High
- ✓ R: Audéo P90-RT (SN: 2031N2XU1) Battery level: High

Accessories

No configurable accessory

✓ Database

Optional session note...

Close session

- ✓ Accept/end call
- ✓ Start/stop voice assistant
- ☐ Mute/unmute mic (not for phone call)
- Tap control training...

Manual controls

Bluetooth


Signals & alerts

Startup

DataLogging

## Audeo Paradise and unrivaled sound quality

**PHONAK**  
life is on




## Speech Enhancer

More gain of soft speech inputs





# Target 7.0: Speech Enhancer



The screenshot displays the Phonak fitting software interface. At the top, the status bar shows 'Client: Target, Phonak', 'Instruments: Audéo P90-R', and 'Fitting: Calm situation (R)'. The 'Fine tuning' tab is selected and highlighted with a red box. Below the tab, the 'Program manager' on the left lists various programs, with 'Calm situation' selected and 'Speech enhancer' highlighted in a red box. The main area shows several sliders for different features: 'Real ear sound' (set to 4), 'Speech enhancer' (set to 14, Moderate), 'WhistleBlock' (set to 13, Moderate), 'SoundRelax' (set to 8, Weak), 'NoiseBlock' (set to 8, Weak), 'Soft noise reduction' (set to 0, Off), and 'WindBlock' (set to 16, Moderate). The bottom navigation bar includes 'Gain & MPO', 'Audibility fine tuning', 'Program options' (which is highlighted), 'SoundRecover2', 'Automatic fine tuning', and 'Tinnitus balance'.

## Dynamic Noise Cancellation

Up to 4 dB improvement in  
SNR in challenging environments

Client-driven personalization  
with  
myPhonak 4.0 app



# Target 7.0: Dynamic Noise Cancellation

The screenshot displays the Phonak fitting software interface. At the top, the 'Fine tuning' tab is highlighted and circled in red. Below it, the 'Program manager' on the left lists various programs, with 'Speech in noise' selected and highlighted in green. The main area shows the 'Dynamic noise cancellation' settings, which are also highlighted with a red box. The settings include a slider for 'UltraZoom' set to 20, and sliders for 'Weak' (8), 'Moderate' (13), and 'Weak' (8). Other settings like 'WhistleBlock', 'SoundRelax', 'NoiseBlock', 'Soft noise reduction', and 'WindBlock' are visible but not selected. The bottom navigation bar includes 'Gain & MPO', 'Audibility fine tuning', 'Program options', 'SoundRecover2', 'Automatic fine tuning', and 'Tinnitus balance'.

## The myPhonak app: What's new?




### Remote Control

Hearing aid personalization  
in real time

# myPhonak: Remote Control- Speech Focus Slider



## Paradise behavior in Speech in Noise





## Motion Sensor Hearing

Higher preference rating for:

Speech intelligibility

Localization

Overall listening experience

Sound quality



# Target 7.0: Motion Sensor Hearing

The screenshot displays the Phonak fitting software interface. At the top, the client is identified as 'Target, Phonak' and the instrument as 'Audéo P90-R'. The 'Fitting' tab is active, showing 'Calm situation -D'. The navigation bar includes 'Feedback & real ear test', 'AudiogramDirect', 'Global tuning', 'Fine tuning', 'DataLogging', and 'Device options'. The 'Fine tuning' tab is highlighted with a red box. On the left, the 'Program manager' shows a list of programs under 'AUTOMATIC PROGRAMS'. 'A AutoSense OS 4.0' is selected and highlighted with a red box. Under it, 'Motion sensor hearing' is also highlighted with a red box. The main area shows sliders for 'Transition speed' (set to 2, 'Balanced'), 'SPEECH IN NOISE & LOUD NOISE' (set to 3, 'Wide'), 'SPEECH IN LOUD NOISE' (checkbox checked, 'Enable as automatic program'), and 'Activation level' (set to 10, 'Loud noise'). The bottom navigation bar includes 'Gain & MPO', 'Audibility fine tuning', 'Program options', 'SoundRecover2', 'Automatic fine tuning', and 'Tinnitus balance'.

# Fitting Paradise in 7 steps





- 1) Connect
- 2) Select gain level
- 3) Check that acoustic parameters match in Target and in patients ears
- 4) Run feedback test
- 5) Increase MPO overall
- 6) Devices options tab-
  - Manual controls tab: confirm button is doing what you would like
  - Bluetooth tab: Bluetooth name- type in name for right and left ears
- 7) Save and close



Connect



Client Paradise, Welcome

Instruments

Fitting

Hearing aids Acoustic parameters Details Accessories

R

L

No hearing aid

No hearing aid

Family

Styles

Features

Search...


Preview

No device to preview.

▶ Audéo P	16
▶ Audéo M	4
▶ Bolero M	8
▶ CROS	9
▶ Lyric	1
▶ Naída Link M	1
▶ Sky Link M	1
▶ Naída M	4
▶ Sky M	9
▶ Virto M	14


Show only new hearing aids





Client Paradise, Welcome Instruments Fitting

Hearing aids Acoustic parameters Details Accessories



- Family ▾ Styles ▾ Features
- ▼ Audéo P
- Audéo P90-R**  
Premium 71/132 (2cc)
  - Audéo P70-R**  
Advanced 71/132 (2cc)
  - Audéo P50-R**  
Standard 71/132 (2cc)
  - Audéo P30-R**  
Essential 71/132 (2cc)
  - Audéo P90-RT**  
Premium 71/132 (2cc)
  - Audéo P70-RT**  
Advanced 71/132 (2cc)
  - Audéo P50-RT**  
Standard 71/132 (2cc)
- Show only new hearing aids


Phonak Target

**① Create new fitting**

Optimize gain for

First time user 0-3 months ▾

**OK**



**DESCRIPTION**  
Premium WaterResistant external receiver instrument (rechargeable) for mild to profound hearing loss. Direct connectivity (Bluetooth®).

**FEATURES**

- AutoSense OS 4.0 premium
- Speech enhancer
- Comfort in echo
- Rechargeable
- Own voice pick-up
- Tinnitus balance
- Dynamic noise cancellation
- Environmental balance
- Rechargeable
- Tap control
- RogerDirect
- UltraZoom premium
- Motion sensor hearing
- StereoZoom
- DuoPhone

► More information



Connect 75% 75%



Client Paradise, Welcome

Instruments Audéo P90-R

Fitting Calm situation <|)



Hearing aids

Acoustic parameters

Details Accessories

Ⓜ Audéo P90-R (M) SN: SN31751500

Ⓛ Audéo P90-R (M) SN: SN89365078



Acoustic code: XXXXXX

Earpiece: Open dome

Dome size: Small

Receiver: M

Wire length: OM R

Acoustic code: XXXXXX

Earpiece: Open dome

Dome size: Small

Receiver: M

Wire length: OM L



Connect 75% 75%




Client Paradise, Welcome

Instruments Audéo P90-R

Fitting Calm situation

Feedback & real ear test AudiogramDirect Global tuning Fine tuning DataLogging Device options


Estimated feedback threshold



No "Feedback & real ear test" measured.

Overtuning provides no additional benefit.

Estimated feedback threshold



No "Feedback & real ear test" measured.

Overtuning provides no additional benefit.

**R** Start both measurements **L**

OVERTUNING



Higher gain limit  
Higher feedback risk

Off

Better sound quality  
Low battery consumption

Off

- Program manager >>
- ALL PROGRAMS
    - AUTOMATIC PROGRAMS 9
      - A** AutoSense OS 4.0
        - Calm situation
        - Speech in noise
        - Speech in loud noise
        - Speech in car
        - Comfort in noise
        - Comfort in echo
        - Music
      - AutoSense OS 4.0 (streaming)
        - Media speech + mic
        - Media music + mic
    - STREAMING PROGRAMS 3
      - PartnerMic + mic
      - Phone call + mic
      - RogerDirect + mic
  - ADDITIONAL PROGRAMS +



3 6 10 20

MPO	65	80	88	94	101	105	107	104	97	87
All	250	500	750	1k	1.5k	2k	3k	4k	6k	8k
G80	1	1	2	3	10	17	20	16	9	2
G65	1	1	2	3	12	20	23	19	11	5
G50	1	1	2	4	12	21	25	20	13	6
CR	1	1	1	1	1.1	1.2	1.2	1.2	1.1	1.2

Verification assistant

MPO	65	80	88	94	101	105	107	104	97	87
All	250	500	750	1k	1.5k	2k	3k	4k	6k	8k
G80	1	1	2	3	10	17	20	16	9	2
G65	1	1	2	3	12	20	23	19	11	5
G50	1	1	2	4	12	21	25	20	13	6
CR	1	1	1	1	1.1	1.2	1.2	1.2	1.1	1.2



Client Paradise, Welcome Instruments Audeo P90-R Fitting Calm situation (L)

Connect 75% 75%

Feedback & real ear test AudiogramDirect Global tuning Fine tuning DataLogging **Device options**

HEARING AIDS

R: SN31751500  
L: SN89365078

SERVICES

- Phone compatibility
- Available apps
- Check warranty
- Reports...
- Product support

Behavior Standard Streaming Phone call



SHORT PRESS

Volume  
±10 dB (±5 steps)

LONG PRESS

Program change  
AutoSense OS 4.0 (Startup)  
Button down has same function

DOUBLE TAP

R: Start/stop voice assistant  
L: No function in Standard

BUTTON FUNCTIONS

- Volume
- Program change
- Accept/end call
- Mute/unmute mic (not for phone call)

TAP CONTROL

- Accept/end call
- Pause/resume streaming
- Start/stop voice assistant
- Tap control training...



Manual controls Bluetooth Signals & alerts Startup DataLogging





Connect 75% 75%



Client Paradise, Welcome

Instruments Audeo P90-R

Fitting Calm situation

Feedback & real ear test AudiogramDirect Global tuning Fine tuning DataLogging Device options

HEARING AIDS

R: SN31751500  
L: SN89365078

SERVICES

- Phone compatibility
- Available apps
- Check warranty
- Reports...
- Product support

Behavior Standard Streaming Phone call



SHORT PRESS  
Volume  
±10 dB (±5 steps)

LONG PRESS  
Program change  
AutoSense OS 4.0 (Startup)  
Button down has same function

DOUBLE TAP  
R: Start/stop voice assistant  
L: No function in Standard

BLUETOOTH NAME  
My Paradise

BLUETOOTH SIDE  
Right

BLUETOOTH NAME  
L-Phonak hearing aid

BLUETOOTH PHONE  
 Fixed bandwidth  
 Adaptive bandwidth  
 Delete pairings

Manual controls

Bluetooth

Signals & alerts

Startup

DataLogging



Connect



Client Paradise, Welcome

Instruments Audéo P90-R

Fitting Calm situation

Feedback & real ear test AudiogramDirect Global tuning Fine tuning DataLogging Device options

HEARING AIDS

R: SN31751500  
L: SN89365078

SERVICES

Phone compatibility

Available apps

Check warranty

Reports...

Product support

Behavior

Standard

Streaming

Phone call



SHORT PRESS

Volume  
±10 dB (±5 steps)

LONG PRESS

Program change  
AutoSense OS 4.0 (Startup)  
Button down has same function

DOUBLE TAP

R: Start/stop voice assistant  
L: No function in Standard

BLUETOOTH NAME

My Paradise

BLUETOOTH SIDE

Right

BLUETOOTH NAME

My Paradise

BLUETOOTH PHONE

Fixed bandwidth

Adaptive bandwidth

Delete pairings

Manual controls

Bluetooth

Signals & alerts

Startup

DataLogging

## Key innovations of Audéo Paradise

### Unrivaled sound quality

- **New sound processing chip - PRISM™**
- **New motion sensor**
- **New Adaptive Phonak Digital 2.0**
- **New AutoSense OS™ 4.0**
  - Speech Enhancer
  - Dynamic Noise Cancellation
  - Motion Sensor Hearing



### Universal connectivity

- Connects to smartphones, TV & more
- Made For All
- Multiple Bluetooth® connections
- New Tap Control for easy access to Bluetooth functionalities

### Personalized digital solutions

- **New Hearing Screener 2.0**
- **New myPhonak 4.0 app**
  - Hearing Diary 2.0
  - AudiogramDirect (via Remote Support)
  - Remote Control
- **New Client Dashboard 2.0**

Paradise is available with groundbreaking lithium-ion rechargeable technology

Together,  
we change lives