

1. This document was created to support maximum accessibility for all learners. If you would like to print a hard copy of this document, please follow the general instructions below to print multiple slides on a single page or in black and white.
2. If you are viewing this course as a recorded course after the live webinar, you can use the scroll bar at the bottom of the player window to pause and navigate the course.
3. This handout is for reference only. Non-essential images have been removed for your convenience. Any links included in the handout are current at the time of the live webinar, but are subject to change and may not be current at a later date.
4. Copyright: Images used in this course are used in compliance with copyright laws and where required, permission has been secured to use the images in this course. All use of these images outside of this course may be in violation of copyright laws and is strictly prohibited.

How to print Handouts

- On a PC
 - Open PDF
 - Click Print
 - Choose # of pages per sheet from dropdown menu
 - Choose Black and White from “Color” dropdown
- On a Mac
 - Open PDF in Preview
 - Click File
 - Click Print
 - Click dropdown menu on the right “preview”
 - Click layout
- Choose # of pages per sheet from dropdown menu
- Checkmark Black & White if wanted.
- If more details needed please visit our FAQ page: <https://www.audiologyonline.com/help>

No part of the materials available through the continued.com site may be copied, photocopied, reproduced, translated or reduced to any electronic medium or machine-readable form, in whole or in part, without prior written consent of continued.com, LLC. Any other reproduction in any form without such written permission is prohibited. All materials contained on this site are protected by United States copyright law and may not be reproduced, distributed, transmitted, displayed, published or broadcast without the prior written permission of continued.com, LLC. Users must not access or use for any commercial purposes any part of the site or any services or materials available through the site.

Technical issues with the Recording?

- Clear browser cache using [these instructions](#)
- Switch to another browser
- Use a hardwired Internet connection
- Restart your computer/device

Still having issues?

- Call 800-753-2160 (M-F, 8 AM-8 PM ET)
- Email customerservice@AudiologyOnline.com



VEMP Testing and Analysis with the GSI Audera Pro

Karen Morris, Product Manager GSI





MEET THE TEAM

Our Audiology team has over 65 years of experience presenting audiology and product knowledge to professionals around the world. Our ADVANCE team enjoys all things audiology and are ready to assist you with any questions. Please reach out to our education team at audiology@grason-stadler.com.



LAURA PRIGGE, AUD
GSI CLINICAL APPLICATION
SPECIALIST



KAREN MORRIS, MS, CCC-A
PRODUCT MANAGER



TONY LOMBARDO, MS
GSI CLINICAL APPLICATION
SPECIALIST





Agenda

- Overview of cVEMP testing
- Overview of oVEMP testing
- Utilizing the VEMP Analysis Module
- 2021 VEMP CPT Codes





Learning Objectives

- After this course, participants will be able to describe the patient task when collecting cVEMP and oVEMPs.
- After this course, participants will be able describe 2 differences between oVEMP and cVEMP recording.
- After this course, participants will be able to list 2 features of the GSI VEMP Analysis module that can assist in test interpretation.






What is a cVEMP?

- Part of the test battery for patient's with dizziness or vertigo.
- Generated from the inner ear (sacculle-in the cochlea) in response to a very *loud* sound
- Recorded from the *contracted* sternocleidomastoid (SCM) muscle.
- Amplitude of the VEMP scales in proportion to the tonic EMG activity
- Calculate the amplitude asymmetry between R/L sides





Electrode Placement*



Active/+ : Sternal notch
Ground: Forehead
Right/- : Belly of the Right SCM
Left/- : Belly of the Left SCM

SCM = Sternocleidomastoid muscle

*P1 is downward, N1 upward





Protocol Parameters: cVEMP

- Stimulus: 500Hz TB; 4 cycles, Blackman, Inserts
- Phase: Alternating
- Intensity: 90-95 dBnHL*
- Rate: 5.1
- Sweeps: 50-100
- Gain: 5000
- Filters: 10-1000Hz
- Window: 50-100 msec

cVEMP Threshold: 80-85 dBnHL





Patient Task: Optimal (SCM) Activation Technique

Stronger the contraction; larger the VEMP response



The Effects of Amplitude Normalization and EMG Targets on cVEMP Interaural Amplitude Asymmetry

Devin L. McCaslin, Gary P. Jacobson, Kelsey Hatton, Andrea P. Fowler, and Andrew P. DeLong

EAR & HEARING, 2013 , VOL. 34, NO. 4, 482–490



cVEMP Amplitude with Optimal Activation Technique

The Effects of Amplitude Normalization and EMG Targets on cVEMP Interaural Amplitude Asymmetry

Devin L. McCaslin, Gary P. Jacobson, Kelsey Hatton, Andrea P. Fowler, and Andrew P. DeLong

EAR & HEARING, 2013 , VOL. 34, NO. 4, 482-490

Pediatric Ages: 5-17
 Young Adult Ages: 19-40
 Older Adult Ages: 41-70





EMG Expected Ranges with Optimal Activation Technique

The Effects of Amplitude Normalization and EMG Targets on cVEMP Interaural Amplitude Asymmetry

Devin L. McCaslin, Gary P. Jacobson, Kelsey Hatton, Andrea P. Fowler, and Andrew P. DeLong




EAR & HEARING, 2013 , VOL. 34, NO. 4, 482-490

Pediatric Ages: 5-17
Young Adult Ages: 19-40
Older Adult Ages: 41-70





GSI VEMP Biofeedback EMG Monitoring





Step by Step

- Prepare the software
- Prepare the patient-scrub skin, apply electrodes, check impedance
- Position the patient and provide instructions: lift and hold head off the table slightly when they hear the stimulus. (consider fatigue)
- Start the collection
- Stop the collection when a consistent waveform appears or max sweeps reached
- Replicate



Marking Waveforms: cVEMP


Latency

P1: 13 - 20.5 ms

N1: 22 - 29.3 ms

cVEMP Threshold

~80-85dBnHL






What is an oVEMP?

- Part of the test battery for patient's with dizziness or vertigo.
- Generated from the inner ear (utricle/superior vestibular nerve) in response to *loud* sound.
- Recorded from inferior oblique muscle of the eye
- Recorded from the eye *contralateral* to the stimulated ear
- Calculate the amplitude asymmetry between R/L sides





Electrode Placement



Active/+ : Low Forehead
Ground: High Forehead
Right/- : Under the Left Eye
Left/- : Under the Right Eye

Recorded from the eye *contralateral* to the stimulated ear





Protocol Parameters: AC oVEMPs


- Stimulus: 500Hz TB; 4 cycles, Blackman, Inserts
- Phase: Alternating
- Intensity: 95 dBnHL*
- Rate: 5.1
- Sweeps: 50-100
- Gain: *100,000*
- Filters: 1-1000Hz
- Window: 50-100 msec

oVEMP Threshold: 90-94dBnHL






Patient Task: Look Upward



Place a spot on the wall for the patient to focus on
Relax forehead/look up with eyes only





Step by Step

- Prepare the software
- Prepare the patient-scrub skin, apply electrodes, check impedance
- Position the patient and provide instructions: Look upward about 35-40 degrees/gaze at spot on the wall when they hear the loud sound.
- Start the collection
- Replicate



Marking Waveforms: oVEMP

Latency

N1: 10 - 12 ms

P1: 15 - 17 ms

Amplitude* (decreases with age)


N1-P1: ~5.0 uV

(Range 1.0-25uV)

*Varies depending on study

Threshold

90-94 dBnHL





How are cVEMP/oVEMPs analyzed?

- Amplitude Asymmetry Ratio (between sides)

$$\frac{\text{Right (P1/N1 amplitude)} - \text{Left (P1/N1 Amplitude)}}{\text{Right (P1/N1 amplitude)} + \text{Left (P1/N1 Amplitude)}} \times 100$$

- Amplitude asymmetry greater than 35 - 45% is considered significant (varies depending on the study.)
- Obtain cVEMP threshold if SSCD is suspected





Causes of Absent VEMPs

1. Conductive component as minimal as 5dB can cause NR
2. Stimulus intensity is too low
3. oVEMP absence with AC stimulus; consider BC stimulation, age of patient
4. cVEMP muscle contraction is too weak (should be minimum of 30-50uV)





References: cVEMP

- Akin, F., Murnane, O., Profitt, T. The Effects of Click and Tone Burst Stimulus Parameters on the Vestibular Evoked Myogenic Potential (VEMP). Journal of the American Academy of Audiology, 2003, 14(9), 500-509.
- Janky, K. L., & Shepard, N. (2009) Vestibular-Evoked Myogenic Potential Testing: Normative Threshold Response Curves and Effects of Age. Journal of the American Academy of Audiology, (20)8, 514-522.
- McCaslin, D.L., Jacobson, G.P., Hatton, K, Fowler, A.P., DeLong, A.P. The Effects of Amplitude Normalization and EMG Targets on cVEMP Interaural Amplitude Asymmetry. Ear & Hearing, 2013, Vol. 34, No. 4, 482-490.
- Rosengren, S.M., Colebatch, J.G., Young, A.S., Govender, S., Welgampola, M.S. (2019). Vestibular Evoked Myogenic Potentials in practice: Methods, Pitfalls and Clinical Applications. Clin. Neurophysiol Practice (4) 47-68.





References: oVEMP

- Piker, Erin G, Jacobson, Gary P., McCaslin, Devin L., Hood, Linda. Normal Characteristics of the Ocular Vestibular Evoked Myogenic Potential. J Am Acad Audiol 2011. 22:222-230.
- Murnane, O. D., Akin, F. W., Kelly, K. J., & Byrd, S. (2011). Effects of stimulus and recording parameters on the air conduction ocular vestibular evoked myogenic potential. Journal of the American Academy of Audiology, 22(7), 469-480. doi: 10.3766/jaaa.22.7.7.
- McCaslin, D.L., Piker, E.G. A Quick Look at Ocular Vestibular-Evoked Myogenic Potentials, Audiology Today, Sept Oct 2011, p 28-35.
- Rosengren, S.M., Colebatch, J.G., Young, A.S., Govender, S., Welgampola, M.S. (2019). Vestibular Evoked Myogenic Potentials in practice: Methods, Pitfalls and Clinical Applications. Clin. Neurophysiol Practice (4) 47-68.





VEMP Analysis Module





VEMP *Analysis* Module Features

- Mark oVEMP or cVEMP Waveforms
- Automatically calculates the asymmetry ratio
- Saving automatically creates a report in the EP Module
- Rectification of cVEMP waveforms (a way to correct for asymmetrical muscle contraction using pre-stimulus EMG activity)
- Allows for viewing each sweep in a recording
- Allows for 'culling' or excluding sweeps from the grand average whose EMG activity is greater than or less than 1, 1.5 or 2 SD's





How does it work?

- Collect VEMPs in EP module
- Open VEMP Analysis module
- Select the waveforms
- Mark N1, P1
- Asymmetry Ratio is automatically calculated





2021 VEMP CPT Codes*

92517	Cervical VEMP only (cVEMP)	\$81.35
92518	Ocular VEMP only (oVEMP)	\$75.51
92519	Cervical VEMP AND Ocular VEMP (cVEMP and oVEMP)*	\$126.72

*2021 Medicare Fee Schedule for Audiologists; ASHA





Questions

krmr@grason-stadler.com

