

1. This document was created to support maximum accessibility for all learners. If you would like to print a hard copy of this document, please follow the general instructions below to print multiple slides on a single page or in black and white.
2. If you are viewing this course as a recorded course after the live webinar, you can use the scroll bar at the bottom of the player window to pause and navigate the course.
3. This handout is for reference only. Non-essential images have been removed for your convenience. Any links included in the handout are current at the time of the live webinar, but are subject to change and may not be current at a later date.
4. Copyright: Images used in this course are used in compliance with copyright laws and where required, permission has been secured to use the images in this course. All use of these images outside of this course may be in violation of copyright laws and is strictly prohibited.

How to print Handouts

- On a PC
 - Open PDF
 - Click Print
 - Choose # of pages per sheet from dropdown menu
 - Choose Black and White from “Color” dropdown
- On a Mac
 - Open PDF in Preview
 - Click File
 - Click Print
 - Click dropdown menu on the right “preview”
 - Click layout
- Choose # of pages per sheet from dropdown menu
- Checkmark Black & White if wanted.
- If more details needed please visit our FAQ page: <https://www.audiologyonline.com/help>

No part of the materials available through the continued.com site may be copied, photocopied, reproduced, translated or reduced to any electronic medium or machine-readable form, in whole or in part, without prior written consent of continued.com, LLC. Any other reproduction in any form without such written permission is prohibited. All materials contained on this site are protected by United States copyright law and may not be reproduced, distributed, transmitted, displayed, published or broadcast without the prior written permission of continued.com, LLC. Users must not access or use for any commercial purposes any part of the site or any services or materials available through the site.

Technical issues with the Recording?

- Clear browser cache using [these instructions](#)
- Switch to another browser
- Use a hardwired Internet connection
- Restart your computer/device


Still having issues?

- Call 800-753-2160 (M-F, 8 AM-8 PM ET)
- Email customerservice@AudiologyOnline.com

March 24, 2021

Counseling Patients through the Ponto Trial Journey: Benefits of Direct Sound Transmission

Carissa Moeggenberg, MA, CCC-A
Training Manager, Oticon Medical



Learner Outcomes

- Following completion of this course the participants will be able to describe the benefits of direct sound transmission.
- Following completion of this course the participants will be able to list the steps involved in the trial process of a Ponto System.
- Following completion of this course the participants will be able to illustrate effective counseling tips to use during the patient journey process.

Agenda

Bone Conduction
Hearing

The Ponto System

The Trial Process

Benefits of Direct
Sound
Transmission


Transitioning to the
Abutment and
Implant

Clinical Evidence
Supporting Direct
Sound
Transmission



Bone Conduction Hearing and The Ponto System

Components of a Bone Anchored Hearing System




Osseointegration and Direct Bone Conduction



Abutment (6, 9, 12, or 14 mm long)


Implant (3 or 4 mm long)



Osseointegration





Direct Bone Conduction



Sounds are converted to vibrations, which are transmitted directly to both inner ears.

The Ponto System Today

- Ponto 4 and Ponto 3 Superpower share the same basic foundation built on *BrainHearing*TM 
 - a. Offer the widest frequency bandwidth in the market 
 - b. Durable design and coupling
 - c. Incorporate Oticon, Inc. sound processing and feedback management



World's most powerful
abutment-level sound
processor



World's smallest &
advanced sound
processing



The Trial Process

The Trial Process

Different for Adults and Children

Adults

- Acute simulation completed in one appointment
- Patient fitting on test band
- Verification of benefit using a test protocol
- Subjective questionnaire administered to patient
- Counseling



Pediatrics

- Not an acute trial but a non-surgical process until child meets minimum age requirements
- Softband fitting
- Counseling
- Parent report measure
- Baseline testing
- Follow-up testing and monitoring

Clinical Tools and Resources

To have available for Patient Simulation

- Ponto Sound Processor (i.e. Ponto 3 SP)
- Verification/Test Protocol
- Subjective Questionnaire/Parent Report Measure
- Test band (preferred), Head band, or Soft band
- Ponto Counseling Kit
- Ponto Care APP – Trial Companion
- Ponto Candidate Materials



Why the Non-Surgical Process is Important



Opportunity to
hear with the
device



Verification of
benefit



Comparison of
different solutions



Assessment of
subjective
impressions



Assist with
counseling



Provides best
practices



Direct Sound Transmission

Counseling Patients

Making an informed decision



- More options than ever for bone anchored patients
- There are both benefits and limitations in the different solutions for the individual patient
 - Does the solution provide enough amplification / output?
 - Is it comfortable enough to wear for a full day?
 - Cosmetically, how does it compare?
 - **Are there any benefits available beyond better hearing?**
- As professionals we must ensure that the patients make informed decisions before any surgical intervention

Direct Sound Transmission vs. Skin Transmission Systems


There are two types of bone conduction devices:

1. Direct Sound Transmission



Devices send vibrations via direct route to bone

2. Skin Transmission



Devices send vibrations through the skin to bone

Bone Conduction Devices

Direct Sound Transmission
Vibrations directly to the bone

Skin Transmission
Vibrations through the skin

Active transcutaneous
Transducer in implant

Percutaneous
Transducer in external sound
processor

Non-surgical solutions

Passive transcutaneous
Implanted magnet



Bone Conduction Devices

Direct Sound Transmission
Vibrations directly to the bone



Skin Transmission
Vibrations through the skin

Active transcutaneous
Transducer in implant

Percutaneous
Transducer in external sound
processor

Non-surgical solutions

Passive transcutaneous
Implanted magnet




Benefits of Direct Sound Transmission

- Sounds are attenuated in Skin Transmission solutions – even with Softband correction.




Skin Transmission




Direct Sound Transmission



- Direct Sound Transmission systems can provide an extra 10-20 dB of output in the mid to high frequency range.*



- The mid to high frequency range contains the most important sound information for speech understanding.




* MFO of Ponto 3 corrected by Skin Transmission dampening as measured in ref 5.

Benefits of Direct Sound Transmission


- Consequence of skin attenuation of speech phonemes




Skin
Transmission



Direct Sound
Transmission



Another Way to Look at it !



Counseling Considerations

Benefits of Direct Sound Transmission

- Hearing thresholds obtained are 5-20 dB lower (better) between 600 – 6000 Hz
- Speech reception thresholds are 4-7 dB lower (better) than conventional devices

Benefit to your patients:

- Access to more sounds
- Access to softer sounds



Clearer, more natural sound with fewer distortion and less feedback

Other Important Considerations



Wearing Comfort

Ponto

Minimal physical sensation of the device

Device stays attached to abutment

Nothing prevents the user from wearing it all waking hours

90% of users wear 8 hours or more

Magnetic Solution

One-third of patients experienced pain

High risk of retention issues

Tissue reduction if skin over 6 mm

22% of users wear 8 hours or more



Cosmetic Considerations



Magnetic Solution

Ponto



Surgical Considerations


- A surgery should be as minimally invasive as possible for less numbness and pain
- Numbness is virtually gone as complication after MIPS surgery³
- With Attract, tissue reduction needed if skin thickness over 6mm
- Tissue preservation with MIPS is therefore reversible without side effects
- Average 45 minute Attract surgery⁴ vs. 7 min for MIPS⁴



Post-Surgical Complications

Percutaneous Devices

- Accepted classification system
 - Holgers Score




Passive Transcutaneous

- No accepted classification system
- Systematic review (Cooper et al.)
 - 26 articles covering 482 implantations
 - Major complications was reported in 25 cases
 - Minor complications mostly consisting of pain or erythema due to magnet strength

MRI Considerations

- MRI is widely used procedure with the likelihood that a patient will need on in their lifetime
- The implant may cause an artifact
 - Artifact shadowing will obstruct the scan
- Acoustic neuroma patients undergo several MRI's post surgery



<http://www.magnetic-resonance.org/ch/21-01.html>



Ponto

- 3 Tesla
- 10 mm artifact





Magnetic system

- 1.5 Tesla
- 115 mm artifact

*Oticon Medical MRI Security card
Surgery guide, Cochlear Baha 4 Attract System,
Surgical procedure*

Note! The illustrations are not exact. Just a schematic view of the artifact effect

Who are the Candidates if Sound Matters?



* Above the age of 5 years



Transitioning to the Abutment and Implant

Transitioning to the Implant

Minimally Invasive Ponto Surgery




80% of all Ponto operations in the US are done with MIPS!

- Minimizing post-op complications
 - Less tissue disturbance, numbness and irritation
- Suture free operation
 - No scarring = Better cosmetic outcomes
 - Less surgical time
- Tailor-made operational components
 - Allows for minimal and safe procedure
 - Long-term goal is to further reduce the need for post-surgical treatment



Transitioning to the Implant

Osseointegration and Direct Transmission



*Ponto is an auditory
osseointegrated bone
conduction implant
system.*

Osseointegration




Because
sound matters

oticon
MEDICAL


Transitioning to the Implant

Bone bonding – fast and strong osseointegration

Implant




Bone structure



First laser-ablated titanium surface

Site-specific modification at the root of the threads


Matching the natural bone structure at macro-, micro-, and nano-scale.



Stronger than bone

The bone bonding to the BHX implant is stronger than the bone itself.⁵

Ponto BHX implant



“


Nanometer roughness plays an important role in osseointegration. The improvement in biomechanical capacity is even greater than I imagined.”

R. Brånemark, Associate Professor MD Msc. PhD

Ponto Abutments



A unique design for tissue preservation



The proven surface

- Ponto abutments feature the proven surface for tissue preservation surgery.⁶

Perfect fit

- Even supportive surface and perfectly matched interface is ideal for MIPS/tissue preservation

OptiFit™ -
Perfect fit with
no pockets or
pathways



I would highly recommend it. It has had a huge impact on my life... I don't know it's there. Once it's healed, you don't really need to do anything different."

Louise

Transitioning to the Implant

The MIPS Procedure

1. Preparing the Site
2. Punching and Inserting the Cannula
3. Initial Drilling with Cannula Guid Drill
4. Drilling with the Cannula Widening Drill
5. Implant Installation
6. Attach the Healing Cap and Dressing






Safe Outcomes

Safe Outcomes

Safe Outcomes: Technology and Procedure go hand and hand



Safe Outcomes

Quality of Life

- Five studies evaluating 176 patients with conductive or mixed hearing loss or single-sided deafness.
- Conditions:
Change in quality of life measured after Ponto surgery.
- Task:
The Glasgow Benefit Inventory (GBI) health-related quality of life questionnaire.

98% of users reported improved quality of life after Ponto surgery¹




1. Lagerkvist, H, et al. (2020). Ten years' experience with the Ponto bone anchored hearing system – a systematic literature review. *Clinical Otolaryngology* (under review).

Safe Outcomes

Quality of Life

- Published Ponto data
 - 176 patients (5 studies)
 - 98% of patients reported improvement in quality of life
 - Average improvement 33 points
- Improvement after Ponto surgery larger than for middle ear implants, and almost as large as for cochlear implants¹




1. Hendry, Chin, Swan, Akeroyd and Browning. *The Glasgow Benefit Inventory: a systematic review of the use and value of an otorhinolaryngological generic patient-recorded outcome measure.* Clin Otolaryngol. 2016; 41(3): 259–275.

Safe Outcomes

Surgical

- **Participants:**
Twenty-seven studies evaluating more than 1,100 patients with conductive or mixed hearing loss and single-sided deafness.
- **Conditions:**
Adverse skin reactions according to *Holgers (across visits) up to five years follow-up (mean follow-up time 16 months).
- **Task:**
*Holgers score ≥ 2 indicates need for skin-related aftercare treatment.

In 95% of follow-up visits no skin-related aftercare treatment* is required¹




1. Lagerkvist, H, et al. (2020). Ten years' experience with the Ponto bone anchored hearing system – a systematic literature review. *Clinical Otolaryngology* (under review).

Safe Outcomes



Surgical



A systematic literature review reporting on 10 years' experience with the Ponto™ System demonstrated **an overall implant survival rate of 98%.¹**



1. Lagerkvist, H, et al. (2020). Ten years' experience with the Ponto bone anchored hearing system – a systematic literature review. *Clinical Otolaryngology* (under review).



Ponto Care App

Safe Outcomes

Ponto Care App

A tool for you and your patients
during social distancing times

- ✓ Developed with usability across all age groups in




- ✓ Adhering to the highest privacy standards- photos are saved on patient's phone, never on the app



- ✓ Free for everyone



- ✓ Does not require log in / account



Trial Section

1. Flow
2. Benefits & service evaluation
3. Best practices

Aftercare section

1. Flow
2. Benefits
3. Best practices

Ponto Care App

What are the benefits?

Patient



- Guidance and engagement during the trial process
- Aftercare information and ability to self-monitor their implant site and their hearing

Professional



- Provides the ability to better engage with patients
- Offers information that can assist with patient management and counseling

- Can have a structured approach to the post trial visit, assist with counseling and overall management
- Patients and HCPs have a better understanding of the evolution of the implant site and their hearing; they can refer to precise and concrete data



Clinical Evidence

Introducing ground-breaking evidence

The choice of system and sound processor will impact the user's everyday life in fundamental ways:

The choice of **System** will

- Affect how fast a child will **learn**
- Have an impact on the ability to **remember**

The choice of **sound processor** will

- Impact the **level of effort** needed for understanding speech

1. Pittman, A. L. Bone conduction amplification in children: Stimulation via a percutaneous abutment vs. a transcutaneous softband. *Ear Hear* (under review).
2. Lunner, T., Rudner, M., Rosenbom, T., Ågren, J., and Ng, E.H.N. (2016) Using Speech Recall in Hearing Aid Fitting and Outcome Evaluation Under Ecological Test Conditions. *Ear Hear* 37 Suppl 1: 145S-154S.
3. Bianchi, F., Wendt, D., Wassard, C., Maas, P., Lunner, T., Rosenbom, T., and Holmberg, M. Benefit of higher maximum force output on listening effort in bone-anchored hearing system users: a pupillometry study. *Ear Hear* (in press).

Clinical Evidence

- We use BrainHearing™ as our guiding star and strive to develop products that make it as easy as possible for the brain to make sense of sound.
- With the Ponto™ System more cognitive resources can be used for understanding, remembering and enjoying life.



LEARN FASTER



REMEMBER MORE



LESS EFFORT

Clinical Evidence

Learn Faster


ASU Arizona State University



- Principal investigator:
Prof. Andrea Pittman
- Participants:
17 children (mean age: 11 years old);
16 with conductive hearing loss,
1 SSD.




- Conditions:
Power version of Ponto
optimally fitted on
abutment and softband.



- Task:
Listen to and learn six
new words.

Clinical Evidence

Learn faster



Children learn new words
2.5 x faster
using a system with
Direct Sound Transmission
compared to a solution
with Skin Transmission.
(Pittman AL, 2019)

Clinical Evidence

Remember more



- Principal investigators:
Prof. Thomas Lunner & Oticon Medical
- Participants:
16 adults (mean age: 58 years old)
with conductive or mixed hearing loss




- Conditions:
Power version of Ponto optimally
fitted on abutment and softband

1. Everybody wears **sunglasses**.
2. He is still lying on the **sofa**.
3. The student will write a long **report**.
4. The whole town came to the **wedding**.
5. His daughter wants to go to **college**.
6. Yesterday was the film's **premiere**.
7. The factory port was not **closed**.


Translated from Danish.

- Task:
Recall words after listening to 7
sentences (SWIR test⁴)



Clinical Evidence

Remember More



13% relative improvement in recall performance with **Direct Sound Transmission**. Fewer resources are needed to process the signal, and more resources can be used for remembering.

Clinical Evidence


Less Effort



- Principal investigator:
Oticon Medical
- Participants:
21 adults (mean age: 59 years old)
with conductive or mixed hearing
loss.



- Conditions:
Sound processors with different
maximum output: Ponto Pro, Ponto
3, Ponto 3 SuperPower.



- Task:
Listen to and repeat sentences, while
an eye-tracking camera monitors
pupil dilation, as a measure of
listening effort.



Clinical Evidence

Less Effort

Listening
window



Retaining
window



A significant decrease in listening effort can be achieved with the Ponto 3 SuperPower, as indicated by reduced pupil dilation.³



Pediatric Case Study


The Importance of Providing Early Amplification

- Age at fitting of amplification is predictive for speech perception, speech production, and spoken language skills. *(Sininger et al 2010)*
- Auditory system development, and particularly development of speech perception, is guided by access to relevant acoustic and linguistic information early in life. *(Kuhl PK, 2000)*



Pediatric Case Study: Importance of Timely Transition

- Transition to Direct Sound Transmission was necessary to maintain normal language development
- "Language development testing showed an accelerated improvement in speech development after implantation"



Pediatric Summary

- Children need to wear a hearing solution:
 - With maximum audibility in the high frequency region in order to facilitate word learning
 - That can be worn all day to develop language and communication skills, both in school and during social activities
 - That is not easily lost
- Soft band is the only choice before age of surgery
- After minimum age for surgery, for normal development:
 - Best access to high frequencies must be provided in order to facilitate word learning
 - Full day of use must be secured



Ponto evidence – Impact beyond better hearing



LEARN FASTER

System level
Direct Sound Transmission vs.
Skin Transmission
(Children)

2.5 x faster learning



REMEMBER MORE

System level
Direct Sound Transmission vs.
Skin Transmission
(Adults)

13 % better recall



LESS EFFORT


Product level
Ponto Pro, Ponto 3 vs.
Ponto 3 SuperPower

Significantly less effort




Ponto - A direct way to Easier Listening

- Results show - Direct Sound Transmission and sound processor with higher max output provides best outcomes.



The Ponto System uses Direct Sound Transmission with an abutment.



Ponto 3 SuperPower is the most powerful abutment-level bone anchored sound processor in the market.

Take-aways

- Direct sound transmission delivers the optimum sound quality
- Increased sound quality leads to a reduction in listening effort
- This is especially important in the pediatric population as it can have a significant effect on other areas such as language development
- There are a wide array of considerations to take into account that the patients may be aware of, or misconceptions they may hold
- It is important the patient makes an informed decision whichever device they choose; they can only do this when armed with the correct information

References

- Briggs, R., et al. (2015). "Clinical Performance of a New Magnetic Bone Conduction Hearing Implant System: Results From a Prospective, Multicenter, Clinical Investigation." *Otology & Neurotology* 36(5): 834-841.
- Dimitriadis PA, Farr MR, Allam A, Ray J. Three year experience with the cochlear BAHA attract implant: a systematic review of the literature. *BMC Ear, Nose, and Throat Disorders*. 2016;16:12. doi:10.1186/s12901-016-0033-5.
- Steehler, Mark W et al. "A Comparison of the Operative Techniques and the Postoperative Complications for Bone-Anchored Hearing Aid Implantation" *International archives of otorhinolaryngology* vol. 22,4 (2018): 368-373.
- Calon et al. Minimally Invasive Ponto Surgery Versus the Linear Incision Technique With Soft Tissue Preservation for Bone Conduction Hearing Implants: A Multicenter Randomized Controlled Trial. *Otol Neurotol*. 2018;39(7):882-893.
- Cooper at., 2017. Passive Transcutaneous Bone Conduction Hearing Implants: A Systematic Review *Otology & Neurotology*: [October 2017 - Volume 38 - Issue 9 - p 1225–1232](#)

References

- Pittman, A. L. Bone conduction amplification in children: Stimulation via a percutaneous abutment vs. a transcutaneous softband. Ear Hear (under review).
- Lunner, T., Rudner, M., Rosenbom, T., Ågren, J., and Ng, E.H.N. (2016) Using Speech Recall in Hearing Aid Fitting and Outcome Evaluation Under Ecological Test Conditions. Ear Hear 37 Suppl 1: 145S-154S.
- Bianchi, F., Wendt, D., Wassard, C., Maas, P., Lunner, T., Rosenbom, T., and Holmberg, M. Benefit of higher maximum force output on listening effort in bone-anchored hearing system users: a pupillometry study. Ear Hear (in press).
- Verhagen et al (2008) The Baha Softband a new treatment for young children with bilateral congenital aural atresia. International Journal of Pediatric Otorhinolaryngology 72, 1455-1459

References

- Sininger et al (2010) Auditory development in early amplified children: Factors influencing Auditory-based communication outcomes in children with hearing loss. Ear and Hearing 2010, 31(2), pp166-185
- Kuhl, P. K. (2000). A new view of language acquisition. Proc Natl Acad Sci, 97, 11850–11857.

Thank you