

House keeping

I am speaking now....Can you hear me?

If you are having technical problems, please stay logged on and call Audiology Online at 1-800-753-2160



House keeping

This session is available for 1/.1 CEU.

Must stay logged on for the full session.

Must successfully complete a short quiz.



Learning objectives:

- 1. Able to identify the product lines utilizing the new Synergy platform
- 2. Able to identify the algorithms that make up the Acuity Lifescape Analyzer
- 3. Know what products are available in the CROS/BiCROS configuration

















Made for Life Video



Patient Feedback

"This is what I've been looking for for year strain to be a continue of the property of the pr



"The quality and performance is really beyond anything I've ever seen."

Patients and professionals agree – Muse hearing aids are made for life.











Muse Platform

What is a platform?

The underlying hardware and computer system on which application programs run

The base configuration that an application, program, or process is integrated on











Muse Operating System

What is an operating system?

The most important software that runs on a computer

An interface between the computer and the user and also between the computer and the processing





Made to deliver Superior Sound Acuity Speech Optimization Acuity Lifescape Analyzer Music Enhancement Multiflex Tinnitus Technology Limitless Wireless









Combination of speech audibility, overall comfort, and enjoyment to deliver enhanced Starkey Sound











MPO Boost




Additional maximum output for severe to

Input and audiogram dependent

Available on all technology tiers

MPO Boost





profound hearing loss

Input and audiogram depender

Available on all technology tier







Muse



Micro RIC 312t Custom

Wireless

i2400 i2000 i1600

Muse Mini BTE 312 Models Available: i2400, i2000, i1600 Wireless 900 MHz

Features and Options

Telecoil and Autophone
Programmable Push Button
Adult and Pediatric Earhook
Thin Tube Option
Universal Right/Left Indicators
Removable MIC Cover (no tool)
Flex Strip Wired Programming
IP 57 Quality Rating

Matrix Options: 130/60 Battery Size: 312

Standard Colors:



Muse Micro RIC 312t



Features and Options

Telecoil and Autophone Programmable Push Button SnapFit Receivers Hear Clear Universal Right/Left Indicators Removable MIC Cover (no tool) CS44 Wired Programming IP 57 Quality Rating

Stock Receiver: 115/50, 120/60 Custom AP Receiver: 115/50, 123/60, 130/70 Matrix Options:

Battery Size: 312



Muse Custom



Wireless 900 MHz

Same user control options as Z Series Hear Clear CIC – Now with Autophone Flex Strip Wired Programming Features and Options

Matrix Options: 110/40 (except ITE), 115/50, 120/60, 130/70 Battery Size: 13, 312, 10a

Standard Colors: . . .



SoundLens Synergy



Invisible In Canal

Non-Wireless Wireless 2400 i2400 2000 1600

SoundLens Synergy



Features and Options No user controls onboard device Hear Clear Autophone Flex Strip Wired Programming Matrix Options: 110/35, 110/40, 115/50

Battery Size: SurfLink Power 10a



Halo 2



RIC 13

i2400 Wireless i2000 i1600



Halo 2 RIC 13

Features and Options

Telecoil and Autophone Programmable Push Button SnapFit Receivers Hear Clear Universal Right/Left Indicators Removable MIC Cover (no tool) CS44 Wired Programming IP 57 Quality Rating

Stock Receiver: 115/50, 120/60 Custom AP Receiver: 115/50, 123/60, 130/70

Battery Size: 13 Standard Colors:



New SnapFit Smart Receivers More accurate acoustic option modeling for earbuds New Comfort Earbuds

Technology Tier Portfolio

Feature	Premium i2400	Advanced i2000	Select i1600
Sound Imaging: Channels Bands	24	20	16
Speech Compression: Acuity Speech Optimization	•	•	•
Music Optimization: Music Enhancement	•	•	•
Music Memory Music Adaptation EZE Music Adaptation	1234	1 2 5 4	1 2 3
Ear-to-Ear Technology: Acuity Binaural Imaging	•	•	0
E2E Wind Noise Management	1 2 3 4	1 2 3 4	1 2 3
E2E Machine Noise Adaptation	1 2 3 4	1 2 5 4	1 2 3
E2E Directionality	•	•	•
E2E Phone Streaming*	•	•	•











Acuity Lifescape Analyzer

Improved **classification** performance
Algorithms run in **parallel** more seamlessly
Includes an additional sound classification for **music**





























Ear-to-Ear Phone

Hearing aid automatically streams phone signal to other ear











TRULINK 3.0



















