


#### LEARNER OBJECTIVES


- Attendees will be able to list the models available with the Widex Unique
- Attendees will be able to describe the new features in the Widex Unique and how they benefit patients
- Attendees will be able to explain how to determine the best product recommendation based on improved signal processing

W

















# UNIQUE **NUMBERS**



	<b>U</b> -PLATFORM
LINEAR INPUT LEVEL	113 dB SPL
SAMPLING FREQUENCY	33.1 kHz
BIT DEPTH	18 bit
LINEAR DYNAMIC RANGE	108 dB
PROCESSING POWER	600 Million Operations/sec
MEMORY	1Mbit
INTEREAR	21 Hz
DELAY	1.3-5 ms
CURRENT CONSUMPTION	< 1.00mA





# INTRODUCING WIDEX UNIQUE













### **PURIFY ALL**


Bringing sound quality to the next level

















# SOFT LEVEL NOISE REDUCTION BENEFITS FOR HEARING AID USERS



- Audibility of the entire input range
- Not bothered by undesirable noise in quiet environments
- Experience enhanced sound quality
- Greatest benefit will be experienced by those with mild to moderate hearing losses






#### **PURIFY ALL**

Bringing sound quality to the next level




























# WIND NOISE MANAGER BENEFITS FOR HEARING AID USERS

- Improvement speech intelligibility (backed by internal trial results)
- No need to remove HAs when in windy environments
- Benefits monaural hearing aid users











W

**GOVERNMENT\*SERVICES** 







Various sound environments



Optimal sound for specific situations

# WIDEX SOUND CLASS SYSTEM




**GOVERNMENT\*SERVICES** 


# UNIQUE INTRODUCES: WIDEX **SOUND CLASS TECHNOLOGY**







**GOVERNMENT\*SERVICES** 









# **UNIQUE INTRODUCES:** WIDEX SOUND CLASS TECHNOLOGY









### THE CONTROLLER


- Different Sound Environments = Different Settings within the Hearing Aids
- Continuous Smooth Transitions
- Seamless Processing
- Providing Audibility and Comfort











GOVERNMENT\*SERVICES




















#### COM-DEX AND COM-DEX APP



- TRUE hands-free solution
- Stream from all Bluetooth devices
- I-Phone & Andriod App
- Compatible with:
  - UNIQUE
  - DREAM











## COM-DEX FEATURES AND BENEFITS

## **FEATURES**


- Plug and Play
- Stylis Compact Design
- Hands-Free
- Discreet Control via APP
- Compatible Android and iOS Smartphones
- Connect to other BT Devices
- Three Color Options

#### **BENEFITS**

- Seamless Connectivity
- Easy to Demonstrate
- Easy to Set-up and Use
- Intuitive
- Long Battery life
- Solution for all WidexLink customers
- Handsfree Safe Communication









# WIDEX UNIQUE"

# **BREAKING THE SOUND BARRIER**





## ADDITIONAL INFORMATION

- Widex VA Pro
  - www.vawidexpro.com
- Widex Learn
  - www.widex.pro/en/learn
- Audiology Online
  - aoinquiry@widex.com

