

If you are viewing this course as a recorded course after the live webinar, you can use the scroll bar at the bottom of the player window to pause and navigate the course.

This handout is for reference only. It may not include content identical to the powerpoint. Any links included in the handout are current at the time of the live webinar, but are subject to change and may not be current at a later date.

**School of Audiology and Speech Sciences**
Blending real-life learning with sound science.


Auditory brainstem responses (ABR) to brief-tone bone-conducted (BC) stimuli



*Susan Small, PhD
Associate Professor
Hamber Professor of Clinical Audiology
University of British Columbia
AudiologyOnline, June 15, 2016*





**School of Audiology and Speech Sciences**
Blending real-life learning with sound science.

TOPIC AREAS TO BE ADDRESSED

BC ABR to brief-tones -- preamble

Overview of methodology

Estimation of infant hearing thresholds

Isolation of test cochlea

Case Study



School of Audiology and Speech Sciences

Blending real-life learning with sound science.

LEARNING OUTCOMES

As a result of this Continuing Education Activity, participants will be able to:

- 1) Explain how to set-up clinical equipment to conduct brief-tone bone-conduction auditory brainstem testing in infants
- 2) Estimate bone-conduction hearing thresholds in infants with normal hearing and hearing loss using bone-conduction auditory brainstem testing
- 3) Explain how to isolate the test cochlea using brief-tone bone-conduction auditory brainstem testing



School of Audiology and Speech Sciences

Blending real-life learning with sound science.

BC ABR to brief tones-- preamble

Clinical goal for BC testing?

Accurate estimation of BC thresholds to determine type of hearing loss responsible for elevated air-conduction (AC) thresholds

Conductive? Sensorineural? Mixed?
- How much is conductive?

- Standard practice for pure-tone audiometry
- Should be standard practice for infant ABR testing



School of Audiology and Speech Sciences

Blending real-life learning with sound science.

Very Brief History of BC ABR testing:

- In the late 1970s and 1980s, BC ABR research emerged (brief tones and clicks) – some technical issues arose but research continued

Examples of early studies:

- ❖ Mauldin & Jerger (1979) found that adult wave V latencies to BC clicks were longer than AC clicks
- ❖ Boezeman et al. (1983) found the same for 2000-Hz brief tones
- ❖ Cornacchia et al. (1983) compared AC & BC ABR wave V latencies in infants & adults; found that infant wave V latencies to BC stimuli were prolonged relative to adults

-- Differences in AC vs BC ABR results and maturational effects emerging in ABR research



School of Audiology and Speech Sciences

Blending real-life learning with sound science.

Examples of early studies that focussed on clinical use of brief-tone BC ABR testing in infants:

Clicks:

- ❖ Yang et al. (1987), Stuart et al. (1993, 1994), Gorga et al. (1993)

Brief-tones:

- ❖ Gravel et al. (1989), Stapells (1989), Stapells & Ruben (1989), Cone-Wesson (1995), Cone-Wesson & Ramirez (1997), Foxe & Stapells (1993) & Noursak & Stapells (1992); Sininger et al. (1997)

Examples of more recent infant brief-tone BC ABR research:

- ❖ Ferm et al. (2014), Hatton et al. (2012), Vander Werff et al. (2009), & Valeriotte & Small (2015)

- (i) feasible to record brief-tone BC ABRs clinically
- (ii) frequency- and mode- (AC vs BC) dependent infant-adult differences to be accounted for in their interpretation



School of Audiology and Speech Sciences

Blending real-life learning with sound science.

Joint Committee on Infant Hearing (2007)

“The audiological assessment should include:


... A frequency-specific assessment of the ABR using air-conducted tone bursts and bone-conducted tone bursts when indicated. When permanent hearing loss is detected, frequency-specific ABR testing is needed to determine the degree and configuration of hearing loss in each ear for fitting of amplification devices.”



School of Audiology and Speech Sciences

Blending real-life learning with sound science.

Overview of methodology



UBC School of Audiology and Speech Sciences
Blending real-life learning with sound science.

Air- and bone-conduction brief tones
Stimulus parameters

"2-1-2" (cycles) linearly-gated tones; 5-cycle Blackman tones (no plateau)

500 Hz:	AC/BC	4 ms rise/fall, 2 ms plateau; 10 ms
1000 Hz:	AC/BC	2 ms rise/fall, 1 ms plateau; 5 ms
2000 Hz:	AC/BC	1 ms rise/fall, 0.5 ms plateau; 2.5 ms
4000 Hz:	AC	0.5 ms rise/fall, 0.25 ms plateau; 1.25 ms
	* BC	1 ms rise/fall, 0.25 ms plateau, 2.25 ms

** to reduce ringing*

(Small & Stapells, 2003)



School of Audiology and Speech Sciences

Blending real-life learning with sound science.

Calibration

Supra-aural TDH49/ER3-A insert earphones/B71 transducer

UNITS: dB peak (peak hold) minus 3 dB = dB ppe x

AC: dB ppe SPL

BC: dB ppe re: 1 μ N

Acoustic calibration for 0 dB nHL

500 Hz			1000 Hz			2000 Hz			4000 Hz		
AC	AC	BC	AC	AC	BC	AC	AC	BC	AC	AC	BC
TDH	ER3-		TDH	ER3-		TDH	ER3-		TDH	ER3-	
49	A		49	A		49	A		49	A	
25	22	67	23	25	54	26	20	49	29	26	46

(Small & Stapells, 2003)

ASSR: Bone oscillator coupling method in infants

- Small et al. (2007) compared infant ASSR & adult behavioural thresholds obtained with elastic headband and hand-held coupling method (for trained individuals)



hand-held



Elastic head band

- **No significant differences 500-4000 Hz**

Caveat: used trained individuals for both methods


-- We do not recommend that a parent couple the bone oscillator to their child's skull

- BC EHP: clinicians use hand-held method
– least likely to wake up infant

ASSRs: Bone oscillator placement

- Small et al. (2007) compared infant ASSR thresholds at different positions on the skull

No difference
for
T versus M
position



Significantly
poorer
for
F versus T or M
position

- BC EHP: clinicians use the T position method
 - greatest range of intensities available
 - easier to maintain firm consistent placement than M position


School of Audiology and Speech Sciences
Blending real-life learning with sound science.

EEG recording parameters


AC: 1-channel recording

➤ Consider two-channel AC recordings if large asymmetry between ears


Adapted from Stapells & Oates (1997)


School of Audiology and Speech Sciences
Blending real-life learning with sound science.


EEG recording parameters



BC: 2-channel recording


School of Audiology and Speech Sciences
Blending real-life learning with sound science.

EEG recording parameters




BC: 2-channel recording

➤ ***Always use two-channel BC recordings***

UBC School of Audiology and Speech Sciences
Blending real-life learning with sound science.

EEG recording parameters



EEG filter: High Pass: 20-30 Hz
Low Pass: 1500-3000 Hz (Slope: 6 or 12 dB/octave, analog)

Artifact reject: Trials exceeding $\pm 25 \mu\text{V}$ (equals Nicolet "50 μV ")
-- can reduce this to optimize recording (balance with rejection rate)

of trials: Typically 2000 per average

of sweeps: Two or more as needed to obtain good signal-to-noise ratio

Display Scale: Avoid too large a scale close to threshold

Adapted from Stapells & Oates (1997)

UBC School of Audiology and Speech Sciences
Blending real-life learning with sound science.

Air- and bone-conduction brief tones

Setting latency window for signal-to-noise ratio (SNR) & residual noise (RN)

	500 Hz	1000 Hz	2000 Hz	4000 Hz
	AC & BC	AC & BC	AC & BC	AC & BC
Start	10.5	7	6.5	5
End	20.5	17	16.5	15

*Note: should shift 10 ms window later for higher presentation levels to avoid stimulus artifact

500 Hz BC stimuli > 30 dB re: 1 μN → 14 – 24 ms (BCEHP, 2012)



School of Audiology and Speech Sciences

Blending real-life learning with sound science.

Interpretation of waveforms



Response present:

- ✓ SNR > 1
- ✓ Wave V visually replicates



Response absent:

- SNR << 1
- Visually flat
- RN < 0.05-0.08 μ V




Cannot evaluate:


- SNR << 1
- RN > 0.08 μ V
- No repeatable peaks
- Not visually flat

(BCEHP, 2012)

SNR latency windows

Infant: 12 months
BC left mastoid





School of Audiology and Speech Sciences

Blending real-life learning with sound science.

What does the absence or presence of a response mean re: the infant's hearing?

- need to relate these results to what is “normal” or “near normal” for AC & BC stimuli for infants
- need to know how these elevated responses predict the degree & type of hearing loss



School of Audiology and Speech Sciences

Blending real-life learning with sound science.

Definition of terms (BCEHP, 2012)

Normal behavioural threshold:

- 25 dB HL

Normal ABR maximum level:

- ABR presentation level at which the majority of normal-hearing infants have a response present

normal ? → response must be present at normal ABR

eHL correction:

max

- Correction factor used to estimate behavioural hearing threshold (dB HL) from the ABR threshold (dB nHL)

ABR threshold
(dB nHL)

—

eHL correction
(dB)

=

estimated behavioural
threshold
(dB HL)

Normal ABR maximum levels & eHL correction for infants *Air- and bone-conduction ABR*

	500 Hz	1000 Hz	2000 Hz	4000 Hz
	AC	AC	AC	AC
BC EHP				
Normal ABR Max (dB nHL)	35	35	30	25
<i>Range in literature</i>	30-35	30-35	20-30	20-25
BC EHP				
eHL correction (dB)	10	10	5	0
<i>Range in literature</i>	10-15	5-10	0-5	-5-0

(BC-EHP 2012, 2015; Small & Stapells, Ch. 21, 2017)

Normal ABR maximum levels & eHL correction for infants
Air- and bone-conduction ABR

	500 Hz		1000 Hz		2000 Hz		4000 Hz	
	AC	BC	AC	BC	AC	BC	AC	BC
BC EHP								
Normal ABR Max (dB nHL)	35	20	35	na	30	30	25	na
<i>Range in literature</i>	30-35	20	30-35	na	20-30	30	20-25	na
BC EHP								
eHL correction (dB)	10	5	10	na	5	5	0	na
<i>Range in literature</i>	10-15	-5	5-10	na	0-5	5	-5-0	na



(BC-EHP 2012, 2015; Small & Stapells, Ch. 21, 2017)




School of Audiology and Speech Sciences

Blending real-life learning with sound science.

Estimation of infant hearing thresholds





UBC School of Audiology and Speech Sciences
Blending real-life learning with sound science.

- If audiologist conducts only AC ABR testing and tympanometry & otoacoustic emissions (OAEs) to identify a conductive component
➤ *May lead to error*
- Tympanometry in very young infants:
 - may fail to identify middle-ear involvement
 - flat tympanogram does not assess amount of hearing loss attributed to the conductive component
- OAEs:
 - sensitive to middle-ear involvement but only helpful if present

Only BC thresholds can distinguish between sensorineural, conductive and mixed losses
AND
determine magnitude of conductive loss



School of Audiology and Speech Sciences

Blending real-life learning with sound science.

How well do BC ABR results predict the nature of the hearing loss (conductive versus sensorineural loss?)

Data collected from BC EHP diagnostic follow up:

	Nature of loss is certain	All data (includes cases where assumptions made)
500 Hz	91.9% (65 cases)	81.2% (126 cases)
2000 Hz	94.2% (37 cases)	93.7% (49 cases)

(Hatton, Janssen & Stapells, 2012)



School of Audiology and Speech Sciences

Blending real-life learning with sound science.

Where does BC ABR testing fit in the diagnostic protocol?



- ☐ After AC thresholds are established for both ears?
- ☐ After AC threshold search, tympanometry & OAEs?
- ☐ Before AC testing?
- ☒ As soon as AC thresholds are determined to be elevated?

➤ **Want BC ABR thresholds early in diagnostic testing to avoid delays in medical follow-up or intervention**

✓ BC testing occurs after AC thresholds are shown to be elevated in at least one ear at 2000 Hz (> normal max) -- *before AC threshold search*

BCEHP, 2012







Isolation of test cochlea



BC ABR: Isolation of test cochlea

INFANTS

- Clinical masking?
 - IA for AC stimuli are not known
 - IA for BC have been approximated with indirect measures (ABR & ASSR data)
 - effective masking levels for BC not known for ABR (BC ASSR data available)
 - are corrections for occlusion effect needed? (BC ASSR data available)



Fig. 2. Illustration of the effect of occlusion on the ear. (Source: Biele, G., Mäkelä, J., & Mäkelä, R. (1990). Effects of occlusion on the ear. In: J. Mäkelä, G. Biele, & R. Mäkelä (Eds.), *Acoustic and Auditory Physiology*. Berlin: Springer-Verlag, 1990, pp. 101-112.)

ADULTS

- Use masking to isolate test ear as needed
 - interaural attenuation (IA) & effective masking levels for AC & BC stimuli are well established
 - corrections for occlusion effect are known

Interaural attenuation of BC stimuli – indirect measures

Study	Method	Indirect measure	Age	Interaural Attenuation (dB)
Yang & Stuart 1987	ABR clicks	Wave V latency	Adult	0-10
			Neonate	25-35
			12 months	15-25

Interaural attenuation of BC stimuli – indirect measures

Study	Method	Indirect measure	Age	Interaural Attenuation (dB)
Yang & Stuart 1987	ABR clicks	Wave V latency	Adult	0-10
			Neonate	25-35
			12 months	15-25
Small & Stapells 2008	ASSR- AM/FM 500-1000 Hz Fc	Ipsi/contra asymmetries	Adult	0-10
			0-6 months	10-30

Interaural attenuation of BC stimuli – indirect measures


Study	Method	Indirect measure	Age	Interaural Attenuation (dB)
Yang et al 1987	ABR clicks	Wave V latency	Adult	0-10
			Neonate	25-35
			12 months	15-25
Small & Stapells 2008	ASSR- AM/FM 500-1000 Hz Fc	Ipsi/contra asymmetries	Adult	0-10
			0-6 months	10-30
Hansen 2010 (M.Sc. Thesis)	ASSR- AM/FM 1000 Hz	Effective masking levels (Binaural AC)	Adult	0
			0-7 months	10-15

➤ Interaural attenuation for BC stimuli in infants is a minimum of 10-35 dB depending on the age

Utilize ipsilateral/contralateral asymmetries?

- Expected pattern for normal cochleae up to 102 years of age -- normal hearing or conductive loss (e.g., aural atresia)
[e.g., Foxe & Stapells, 1993; Stapells & Ruben, 1989; Stapells & Mosseri, 1991]


BC left mastoid 2000 Hz @ 40 dB nHL



Amplitude: contra smaller than ipsi **Latency:** contra later than ipsi


Utilize ipsilateral/contralateral asymmetries?

- Expected pattern for normal cochleae up to 102 years of age -- normal hearing or conductive loss (e.g., aural atresia)
[e.g., Foxe & Stapells, 1993; Stapells & Ruben, 1989; Stapells & Mosseri, 1991]



Amplitude: contra smaller than ipsi **Latency:** contra later than ipsi



Bone-conduction ABR



Contra >> IPSI

SEVERE UNILATERAL (RIGHT EAR) SNHL

(Stapells, personal communication)





School of Audiology and Speech Sciences

Blending real-life learning with sound science.

Factors contributing to ipsi/contra asymmetries?

1. Greater IA compared to adults due to unfused cranial sutures

2. Infant-adult differences in positioning of neural generators



Evidence: infant AC ABR/ASSRs show consistent ipsi/contra asymmetries; adult AC ABR/ASSRs do not show these patterns

(Reviewed in Small & Stapells, 2017)

➤ **Two-channel recordings are routinely used by our provincial program (BCEHP) for BC brief-tone ABRs**

NOTE: Can also use ipsi and contra EEG channel for AC if a large difference in thresholds between ears exists (and contra masking not used)



School of Audiology and Speech Sciences

Blending real-life learning with sound science.

What if ipsi/contra asymmetries in BC ABRs are ambiguous?

➤ **Need clinical masking**

Main reason masking not routinely used clinically for infant BC ABRs:

-- *effective masking levels (EMLs) for BC ABR stimuli in young infants have not been measured directly*

➤ **What do we know about EMLs for BC auditory evoked potentials?**

-- EMLs for infant BC ASSR stimuli were estimated for 500-4000 Hz using binaural AC masking (Hansen & Small, 2012; Small, Smyth & Leon, 2014)



School of Audiology and Speech Sciences

Blending real-life learning with sound science.

Recommended EMLs (dB SPL) for BC ASSR stimuli presented at 35 dB HL

	Frequency (Hz)				
	500	1000	2000	4000	
Infant	81	68	59	45	-10
Adult	66	63	59	55	

* Significant infant minus adult EML difference (dB)

- Frequency-dependent infant-adult differences in EMLs except at 2000 Hz

(Hansen & Small, 2012; Small, Smyth & Leon, 2014)



School of Audiology and Speech Sciences

Blending real-life learning with sound science.

Is there an occlusion effect (OE) in infants?

- Adults with normal hearing or a sensorineural hearing loss: occluding the ear canal results in a significant improvement in pure-tone BC thresholds
- Do we need to correct for an OE in infants when we obtain BC thresholds with earphones in place?

We investigated this phenomenon in infants (2 studies):

- (i) Small, Hatton & Stapells, 2007
 - no occlusion effect for BC ASSR thresholds 500-4000 Hz
- (ii) Small & Hu, 2011
 - Sound pressure ↑ in ear canal when occluded: infants >> adults
 - % occurrence of OE:
 - Older infants: OE emerging at 500 & 1000 Hz
 - Young infants: OE absent at 1000 Hz (very small at 500 Hz)


**School of Audiology and Speech Sciences**
Blending real-life learning with sound science.





***Earphones in or out
during BC testing?***



➤ **Recommendations (conservative):**



- (i) **Young infants: leave earphones in place**
- (ii) **Older infants: remove earphones**

**School of Audiology and Speech Sciences**
Blending real-life learning with sound science.

Case Study







Stimulus	ABR (dB nHL)		Estimated Behavioural Threshold (dB eHL)	
	RIGHT	LEFT	RIGHT	LEFT
AC – 500 Hz	> 100	80	> 110	70
BC – 500 Hz	> 20	> 20	> 25	>25
AC – 1000 Hz	> 100	≤ 55	> 110	≤ 45
AC – 2000 Hz	> 100	40	> 105	35
BC – 2000 Hz	> 60	35-60	> 65	30-55
AC – 4000 Hz	> 90	25	> 90	25

+ ipsi/contra asymmetries (BC & AC) support left ear responding

- 1st appointment: BC ABR established nature & severity of loss L & R
- 2nd and 3rd appointment completed AC ABR testing:
 - L: thresholds at 500, 1000, 2000 & 4000 Hz
 - R: established profound loss
- MRI/CT: confirmed absence of cochlear nerve on the R
(click ABR– no clear signs of ANSD)



School of Audiology and Speech Sciences

Blending real-life learning with sound science.



Thoughts?

References

- Boezeman EHJF, Kapteyn TS, Visser SL, Snel AM (1983) Comparison of the latencies between bone and air conduction in the auditory brainstem evoked potential. *Electroencephalography and clinical Neurophysiology* 56:244-247.
- Cone-Wesson B (1995) Bone-conduction ABR tests. *American Journal of Audiology* 4:14-19.
- Cornacchia L, Martini A, Morra B (1983) Air and bone conduction brain stem responses in adults and infants. *Audiology* 22:430-437.
- Davis, H., Hirsh, S. K., Popelka, G. R., & Formby, C. (1984). Frequency selectivity and thresholds of brief stimuli suitable for electric response audiometry. *Audiology*, 23, 59-74.
- Ferm, I., Lightfoot, G., & Stevens, J. (2013). Comparison of ABR response amplitude, test time, and estimation of hearing threshold using frequency specific chirp and tone pip stimuli in newborns. *Int J Audiol*, 52(6), 419-423. doi:10.3109/14992027.2013.769280
- Foxe JJ, Stapells DR (1993) Normal infant and adult auditory brainstem responses to bone-conducted tones. *Audiology* 32:95-109.
- Gorga MP, Kaminski JR, Beauchaine KL, Bergman BM (1993) A comparison of auditory brain stem response thresholds and latencies elicited by air- and bone-conducted stimuli. *Ear and Hearing* 14:85-93.
- Gravel JS, Kurtzberg D, Stapells DR, Vaughan HGJ, Wallace IF (1989) Case studies. *Seminars in Hearing* 10:272-287.
- Hansen, E.E. (2010) Effective masking levels for bone-conduction auditory steady-state response thresholds in infants. University of British Columbia, M.Sc. Thesis.
- Hansen, E. E., & Small, S. A. (2012). Effective masking levels for bone conduction auditory steady state responses in infants and adults with normal hearing. *Ear Hear*, 33(2), 257-266. doi:10.1097/AUD.0b013e31822f67f6
- Hatton, J. L., Janssen, R. M., & Stapells, D. R. (2012). Auditory Brainstem Responses to Bone-Conducted Brief Tones in Young Children with Conductive or Sensorineural Hearing Loss. *International Journal of Otolaryngology*, 2012, 284864. <http://doi.org/10.1155/2012/284864>
- Joint Committee on Infant Hearing (JCIH). (2007). Year 2007 position statement: Principles and guidelines for early hearing detection and intervention programs. *Pediatrics*, 120(4), 898-921.
- Mauldin L, Jerger J (1979) Auditory brainstem evoked responses to bone-conducted signals. *Archives of Otolaryngology* 105:656-661.

References cont'd

- Nousak JK, Stapells DR (1992) Frequency specificity of the auditory brain stem response to bone-conducted tones in infants and adults. *Ear and Hearing* 13:87-95
- Sininger YS, Abdala C, Cone-Wesson B (1997) Auditory threshold sensitivity of the human neonate as measured by the auditory brainstem response. *Hearing Research* 104:27-38
- Small, S. A., & Hansen, E. E. (2012). Effective masking levels for bone-conducted amplitude- and frequency-modulated tones in adults with normal hearing: A behavioural study. *Int J Audiol*, 51(3), 216-219. doi:10.3109/14992027.2011.622303
- Small, S. A., Hatton, J., Stapells, D.R. (2007). Effects of bone oscillator coupling method, placement location, and occlusion on bone-conduction auditory steady-state responses in infants. *Ear & Hearing*, 28(1), 83-98.
- Small, S. A., & Hu, N. (2011). Maturation of the occlusion effect: a bone conduction auditory steady state response study in infants and adults with normal hearing. *Ear Hear*, 32(6), 708-719. doi:10.1097/AUD.0b013e31821de1b0
- Small, S. A., Smyth, A., & Leon, G. (2014). Effective masking levels for 500 and 2000 Hz bone conduction auditory steady state responses in infants and adults with normal hearing. *Ear Hear*, 35(1), 63-71. doi:10.1097/AUD.0b013e31829f2657
- Small, S. A., & Stapells, D. R. (2003). Normal brief-tone bone-conduction behavioral thresholds using the B-71 transducer: three occlusion conditions. *J Am Acad Audiol*, 14(10), 556-562.
- Stapells DR (1989) Auditory brainstem response assessment of infants and children. *Seminars in Hearing* 10:229-251

References cont'd

- Small, S. A., & Stapells, D. R. (2008). Normal ipsilateral/contralateral asymmetries in infant multiple auditory steady-state responses to air- and bone-conduction stimuli. *Ear Hear*, 29(2), 185-1
- Stapells, D.R., & Mosseri, M. (1991). Maturation of the contralaterally recorded auditory brainstem response. *Ear Hear*: 12, 167-173.
- Stapells, D. R., & Oates, P. (1997). Estimation of the pure-tone audiogram by the auditory brainstem response: a review. *Audiology & Neuro Otology*, 2, 257-280.
- Stapells DR, Ruben RJ (1989) Auditory brain stem responses to bone-conducted tones in infants. *Annals of Otology, Rhinology and Laryngology* 98:941-949
- Stuart A, Yang EY, Green WB (1994) Neonatal auditory brainstem response thresholds to air- and bone-conducted clicks: 0 to 96 hours postpartum. *Journal of the American Academy of Audiology* 5: 163-172
- Stuart A, Yang EY, Stenstrom R, Reindorp AG (1993) Auditory brainstem response thresholds to air and bone conducted clicks in neonates and adults. *The American Journal of Otology* 14:176-182
- Valeriotte, H, & Small. S.A. (2015). Comparison of air- and bone-conduction auditory brainstem and multiple 80-Hz auditory steady-state responses in infants with normal hearing and conductive hearing loss. Abstracts of the XXIV Biennial Symposium of the International Evoked Response Audiometry Study Group IERASG (May 10-14, 2015, Busan, Korea.) Retrieved April 4, 2015, from <http://www.ierasg2015.org/>
- Vander Werff, K. R., Prieve, B. A., & Georgantas, L. M. (2009). Infant air and bone conduction tone burst auditory brain stem responses for classification of hearing loss and the relationship to behavioral thresholds. *Ear and Hearing*, 30(3), 350-368.
- Yang EY, Rupert AL, Moushegian G (1987) A developmental study of bone conduction auditory brainstem responses in infants. *Ear and Hearing* 8:244-251

Thank you!