

If you are viewing this course as a recorded course after the live webinar, you can use the scroll bar at the bottom of the player window to pause and navigate the course.

This handout is for reference only. It may not include content identical to the powerpoint. Any links included in the handout are current at the time of the live webinar, but are subject to change and may not be current at a later date.

One Size Does Not Fit All

Selecting and Fitting Devices for Tinnitus Management

Jennifer Martin, AuD
National University of Singapore
Audiology Online July 2016

Learning Outcomes

- Define the role of sound therapy as part of the tinnitus/sound sensitivity management plan.
- Explain when each type of device (hearing aid, sound generator, combination unit) would be used as part of the management plan for patients with tinnitus and sound sensitivity.
- Describe the differences between fitting hearing aids for hearing loss and fitting hearing aids for tinnitus and sound sensitivity management.

6

Recommended Components of Tinnitus Management

- 1) Education & Counseling
- 2) Stress reduction & Relaxation
- 3) Therapeutic sound

Hoare et al. (2014)

7

Sound Therapy

- “Any use of sound where the intention is to alter the tinnitus perception and/or the reactions to tinnitus in a clinically meaningful way” Hoare et al. (2014)
- “using external sounds to provide relief from tinnitus” Folmer & Carroll (2006)
- “sound is used to directly or indirectly shift attention away from the tinnitus”

Henry et al. (2008)

8

Sound Therapy

- ◉ Very broad area
- ◉ Use the general guidelines of sound therapy
- ◉ Focus on the use of ear level devices as part of the sound therapy strategy

9

General Guidelines of Sound Therapy

- ◉ Low level sound that is not bothersome
- ◉ Should not interfere with communication or concentration
- ◉ Reduces the contrast between the tinnitus and the environment

10

Therapeutic Sound

- ◉ Environmental sound
- ◉ Music
- ◉ Speech

Henry et al. (2008)

11

Purposes of Therapeutic Sound

- ◉ Sense of relief from tinnitus-related stress (soothing sound)
- ◉ Passive diversion from tinnitus (background sound)
- ◉ Active diversion from tinnitus (interesting sound)

Henry et al. (2008)

12

Devices for tinnitus management

- ◉ Sound Generators
- ◉ Hearing aids
- ◉ Combination units

13

Devices for tinnitus management

- ◉ Sound Generators
 - › Ear level devices
 - › No amplification
 - › Some have multiple sound options
 - › Some are programmable
 - › Various styles
 - › Less expensive than hearing aids

14

Devices for tinnitus management

- Hearing aids
 - › Amplification only
 - › Wireless devices to stream sounds

15

Benefits of Amplification

- Improved ability to hear soothing sound, background sound and distracting sound to provide relief from tinnitus and tinnitus-related stress
- Reduction in tinnitus loudness
- Less stress associated with straining to hear
- Stimulation of the auditory system

16

Devices for tinnitus management

- ◉ Combination units
 - › Amplification only
 - › Sound generator only
 - › Combination of both
 - › Can also use wireless devices to stream sounds

17

Considerations: Device Type

- ◉ Hearing Aid versus Combination Unit
 - › Starting Point:
 - HA for normal hearing to mild HL in LF
 - CU for greater than mild LF loss
- McNeil et al. (2012)
- ◉ Always be open to trying both to allow patient to decide what works best for them
- ◉ Current technology: order combination unit and use features as needed

18

Considerations: Device Style

- ◉ Open fit
 - › Slim tube BTE
 - › Receiver in the ear
- ◉ If hearing loss is too great for open fit, then use maximum venting possible

Searchfield (2006)

19

Considerations: Monaural vs Binaural

- ◉ Binaural fitting for binaural hearing loss, even if the tinnitus is only present unilaterally
 - › Provides more normal auditory balance

Searchfield (2006)

20

Considerations: Prosthetic vs Therapeutic

- Most hearing aids are fit prosthetically
 - › Fitting to replace the lost hearing
 - › Focus on communication
- Tinnitus often requires fitting therapeutically
 - › Changes are made to the programming to allow the hearing aid to have more effect on the tinnitus perception

21

Considerations: HL vs Tinnitus

- Which is the BIGGER problem?
 - › HL: Fit for the hearing loss and patient's communication needs
 - Tinnitus benefit is secondary goal
 - › Tinnitus: Fit primarily for tinnitus relief
 - Improved communication ability is secondary goal
- Can set two programs if needed
 - › 1 for better hearing ability
 - › 1 for better tinnitus relief

Henry et al. (2008)


22

Considerations: Acoustic Programming

- Use feedback reduction for most open fit
- Disable internal noise reduction (expansion)
- Disable environmental noise reduction
- Low compression knee point
- Omnidirectional microphone setting
- Fitting protocol: DSL I/O v5

Searchfield (2006)
Henry et al. (2010)
Wise (2003)

23




24

Considerations: Acoustic Programming

- Use feedback reduction for most open fit
- Disable internal noise reduction (expansion)
- Disable environmental noise reduction
- Low compression knee point
- Omnidirectional microphone setting
- Fitting protocol: DSL I/O v5

Searchfield (2006)
Henry et al. (2010)
Wise (2003)

25




Considerations: Acoustic Programming



- Use feedback reduction for most open fit
- Disable internal noise reduction (expansion)
- Disable environmental noise reduction
- Low compression knee point
- Omnidirectional microphone setting
- Fitting protocol: DSL I/O v5

Searchfield (2006)
Henry et al. (2010)
Wise (2003)

27



28




continued™

Considerations: Acoustic Programming

- Use feedback reduction for most open fit
- Disable internal noise reduction (expansion)
- Disable environmental noise reduction
- Low compression knee point
- Omnidirectional microphone setting
- Fitting protocol: DSL I/O v5

Searchfield (2006)
Henry et al. (2010)
Wise (2003)

31




32

Considerations: Acoustic Programming

- Use feedback reduction for most open fit
- Disable internal noise reduction (expansion)
- Disable environmental noise reduction
- Low compression knee point
- Omnidirectional microphone setting
- Fitting protocol: DSL I/O v5



Searchfield (2006)
Henry et al. (2010)
Wise (2003)

33



34

continued™




continued™

Considerations: Acoustic Programming



- Use feedback reduction for most open fit
- Disable internal noise reduction (expansion)
- Disable environmental noise reduction
- Low compression knee point
- Omnidirectional microphone setting
- Fitting protocol: DSL v5 versus NAL-NL2

Searchfield (2006)
Henry et al. (2010)
Wise (2003)


37



38



continued™



41

WIDEX COMPASS GPS

Enter patient's ID or full name here

SESSION START SELECTION FITTING FINE TUNING LOG HANDLING

Connect

END SESSION

FITTING

- Overview
- Feedback test
- Sensogram
- Acclimatization

MORE TOOLS

- Audiometry
- Rationale
- Fitting wizard

Fitting rationale

- ☐ Widex rationale
- ☒ NAL-NL2
- ☐ DSL v5.0 Pediatric
- ☐ DSL v5.0 Adult

Gender

- ☐ Male
- ☒ Female
- ☐ Not specified

Hearing aid experience

- ☐ Experienced
- ☒ Inexperienced

Language

- ☒ Non-tonal
- ☐ Tonal

Hearing loss nature

Sensorineural

Assessment of in-situ acoustics (AISA)

Vent compensation off

SoundPlayer

Music_jazz

Play

WIDEX DREAM

42

continued™

Considerations: Sound Sensitivity

- ◉ Sound Generator versus Hearing Aid/Combination Unit
- ◉ Occlusion versus Venting

Searchfield (2006)

43


Considerations: Sound Sensitivity

- ◉ Programming modifications
 - › Low compression knee point + higher than normal compression ratio
 - Use hearing aid as amplifier and limiter
 - As sound sensitivity improves, compression ratios can gradually be reduced


- ◉ Reduce maximum power output (MPO)

Searchfield (2006)

44




45



46

continued™



Special Considerations

- Each tinnitus patient is unique and requires individualized care/empathy
- Tune devices to meet the patient's individual hearing, tinnitus and comfort needs
- Perform probe microphone measurements to verify acoustic fit
- Multi-disciplinary care is best
- Know when to refer to a tinnitus specialist

48

References

- Hoare DJ, Searchfield GD, Refaie AE, Henry JA. (2014) Sound Therapy for Tinnitus Management: Practicable Options. *J Am Acad Audiol* 25: 62-75.
- Henry JA, Zaugg TL, Myers PJ, Schechter MA. (2008) Using Therapeutic Sound With Progressive Audiologic Tinnitus Management. *Trends in Amplification* 12: 188-209.
- Searchfield GD, Kaur M, Martin WH. (2010) Hearing aids as an adjunct to counselling: Tinnitus patients who choose amplification do better than those that don't. *Intl J Aud*
- McNeil C, Tavora-Vieira D, Alnafjan F, Searchfield GD, Welch D. (2012) Tinnitus pitch, masking, and the effectiveness of hearing aids for tinnitus therapy. *Intl j Aud* 12: 914-919.
- Searchfield GD. (2006) Hearing Aids and Tinnitus. In Tyler RS (Editor), *Tinnitus Treatment*
- Shekhawat GS, Searchfield, GD, Kobayashi K, Stinear CM. (2013) Prescription of hearing-aid output for tinnitus relief. *Intl J Aud* 9: 617-625.
- Folmer RL, Carroll, JR. (2006) Long-Term Effectiveness of Ear-Level Devices for Tinnitus. *Oto Head Neck Surg* 134: 132-137.
- Henry JA, Zaugg TL, Myers PJ, Kendall CJ. (2010) *Progressive Tinnitus Management: Clinical Handbook for Audiologists*.

49