



ReSound GN

Quick reference guide

















ReSound LiNX 3D™ and ReSound ENZO²™



DESIGN

Custom Models

RIE/BTE Models

		Non-W	W	Options	LP: Mild-Moderate MP: Mild-Moderate	HP: Moderate - Severe UP: Severe - Profound	Matrices
	IIC	●		10A battery Faceplate color: anthracite, beige tinnitus sound generator, SoundShaper, Acceptance Manager			low power: 115/40
	CIC	●		10A battery VC, push button, tinnitus sound generator, SoundShaper, Acceptance Manager			
	MIH-S	●		10A battery VC, push button, tinnitus sound generator, SoundShaper, Acceptance Manager			
	MIH	●	●	10A, 312 or 13 battery Made for iPhone, Ear-to-Ear communication, push button, VC, telecoil, tinnitus sound generator, SoundShaper, Synchronized Acceptance Manager, 10A (non-wireless only)			low power: 115/40 medium power: 119/50 high power: 121/60 ultra power: 130/70
	ITC	●	●	312 battery Made for iPhone, directional microphones, telecoil, VC, push button, SoundShaper, Acceptance Manager, tinnitus sound generator			
	ITE HS	●	●	312 battery Made for iPhone, Ear-to-Ear communication, directional microphones, telecoil, VC, push button, tinnitus sound generator, SoundShaper, Synchronized Acceptance Manager			
	ITE FS	●	●	312 or 13 battery Made for iPhone, Ear-to-Ear communication, directional microphones, telecoil, VC, push button, tinnitus sound generator, SoundShaper, Synchronized Acceptance Manager			medium power: 119/50 high power: 121/60 ultra power: 130/70

Available Colors















ReSound custom models

ReSound 3D™ RIE and BTE



ReSound ENZO2™ BTE



		Non-W	W	Options	LP: Mild-Moderate MP: Mild-Moderate	HP: Moderate - Severe UP: Severe - Profound	Matrices
	Mini RIE 61		●	312 battery options: Made for iPhone, Ear-to-Ear communication, directional microphones, programmable push button (PB or VC), tinnitus sound generator, SoundShaper, Synchronized Acceptance Manager			low power: 113/52 medium power: 116/58 high power: 122/65 ultra power: 129/75
	RIE 62		●	13 battery Made for iPhone, Ear-to-Ear communication, directional microphones, telecoil, multifunction button (VC, PB, & streaming), DAI, tinnitus sound generator, SoundShaper, Synchronized Acceptance Manager			
	Mini BTE 67		●	312 battery Made for iPhone, Ear-to-Ear communication, directional microphones, earhook or Thin Tube fitting, telecoil, push button, tinnitus sound generator, SoundShaper, Synchronized Acceptance Manager			open: 122/52 closed: 120/56
	Standard BTE 77		●	13 battery Made for iPhone, Ear-to-Ear communication, directional microphones, earhook or Thin Tube fitting, telecoil, VC, push button, DAI, tinnitus sound generator, SoundShaper, Synchronized Acceptance Manager			open: 127/51 closed: 124/57
	Power BTE 88		●	13 battery Made for iPhone, Ear-to-Ear communication, directional microphones, telecoil, VC, push button, DAI, tinnitus sound generator, SoundShaper, Synchronized Acceptance Manager			power: 132/67
	High Power BTE 88		●	13 battery Made for iPhone, Ear-to-Ear communication, directional microphones, telecoil, VC, push button, DAI, tinnitus sound generator, SoundShaper, Synchronized Acceptance Manager			power: 132/73
	Super Power BTE 98		●	675 battery Made for iPhone, Ear-to-Ear communication, directional microphones, telecoil, VC, push button, DAI, tinnitus sound generator, SoundShaper, Synchronized Acceptance Manager			power: 141/83

(#) = Keychain reference



CONNECTIVITY

ReSound wireless connections

Your patients can connect directly to all ReSound wireless accessories for direct audio streaming and discreet control over their listening experience - without using a neckloop. With a full ecosystem of 2.4 GHz wireless connectivity, your patients can hear everything better - from family members to the TV - with less effort.

Easily set up by pressing the pairing button on the accessory, open and close the battery doors of the hearing aids, and the accessory is automatically paired with wireless ReSound hearing aids and ready for wireless communication.

- **ReSound Micro Mic:** Significantly improves signal-to-noise ratio in noisy situations, delivers exceptional speech understanding and extends the hearing range to over 80 feet.
- **ReSound Multi Mic:** With the same hearing range of over 80 feet, doubles as a table microphone, connects with loop and FM systems**, and has a mini-jack input for streaming audio from a computer or mobile device.
- **ReSound Remote Control 2:** Adjust the volume, mute hearing aids, switch programs and see all settings at a glance.
- **ReSound Phone Clip+:** Make crystal clear phone calls, mute background noise while talking, or stream music, podcasts or any other audio from a Bluetooth-enabled phone.
- **ReSound TV 2:** Stream stereo sound from a TV, computer or music system directly to the hearing aids over a range of up to 20 feet. Includes a Dolby® adaptor for newer televisions.

**FM receiver required

© 2017 The GN ReSound Group, all rights reserved. Apple, iPad, iPhone and iPod touch are trademarks of Apple Inc, registered in the U.S. and other countries. App Store is a service mark of Apple Inc. Android and Google Play are registered trademarks of Google Inc. The Android robot is reproduced or modified from work created and shared by Google and used according to terms described in the Creative Commons 3.0 Attribution License. Dolby and the double-D symbol are registered trademarks of Dolby Laboratories.



ReSound Micro Mic



ReSound Multi Mic



ReSound Phone Clip+



ReSound Remote Control 2



ReSound TV Streamer 2

Manufacturer according to FDA:

**ReSound
North America**
8001 Bloomington Freeway
Bloomington, MN 55420
1-800-248-4327
resoundpro.com

Manufacturer according to Health Canada:

**ReSound
Canada**
303 Supertest Road
Toronto, Ontario M3J 2M4
1-888-737-6863
resoundpro.com

ReSound GN



MK604806 Rev A 2017.01