

RECHARGEABLE & CONVENIENT



Stellar RIC Lithium-ion

Rechargeable

Stellar RIC 8C is our lithium-ion rechargeable RIC that provides a full day of wear time, on a 3-hour charge. Stellar RIC is suitable for up to a severe hearing loss and with an IP68 rating, it can stand up to moisture, sweat and dirt.



inoX click and fit

No impression required

inoX 8C is the new premium click and fit that doesn't require an impression. It is perfect for someone who likes discretion without the wait. Our click sleeves allow for a safe and comfortable fit. inoX 8C is remote controllable with the Smart Remote App or Smart Key Remote.

DISCREET & EASY TO USE




Discover the Smart Choice.

Rexton is built on the premise that exceptional hearing technology should be accessible to all. Our products offer superior speech understanding, durability and personalization to deliver unprecedented wearing comfort. Our tenured, passionate people are empowered to provide solution-driven service to meet each customer's individual needs. That is why Rexton is the smart choice.

The *Bluetooth*® word mark and logos are owned by the Bluetooth SIG, Inc., and any use of such marks by Sivantos Inc. is under license. Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. iPhone is a trademark of Apple Inc., registered in the U.S. and other countries. Copyright © 2019 Rexton. All rights reserved.

RX9Q41 Rev 319 22B771 PN RX201902



REXTON

REXTON

PLATFORM OVERVIEW



MyCore 8C Technology:

A perfect combination of natural sound
and direct audio streaming.

Winter/Spring 2019

REXTON



MyCore 8C. Made for Me!

Rexton's MyCore platform combines the best of two worlds: direct audio streaming from an iPhone¹ and advanced binaural signal processing that allows wearers to hear more comfortably and naturally. MyCore 8C technology is available in RIC, BTE, click and fit, and Custom hearing solutions.

Wearers can now enjoy a truly individualized hearing experience.

DIRECT AUDIO STREAMING

Connect MyCore hearing aids directly to iPhones.

Benefits:

Hearing Care Professional:

- Offer wearers connectivity

Wearer:

- Easily connect to iPhones and TV

NATURAL SOUND

Hear more naturally with advanced binaural signal processing.

Benefits:

Hearing Care Professional:

- Higher acceptance at fitting

Wearer:

- Natural sound of their own voice

¹Apple iPhone 5 and later.



Emerald™ XS



Emerald S



Emerald M



Stellar™ RIC

RIC



EASY-TO-USE

Hearing aids can be completely automatic or customized.

Benefits:

Hearing Care Professional:

- No additional programs required

Wearer:

- Completely automatic



BTE

Click and fit

Custom

FEATURE OVERVIEW

Product Details		Benefits	80/60/40
Signal Processing Channels		Provides clean, natural sound	48/32/24
Gain Handles		Audibility and comfort	20/16/12
Wireless Sync ²		Creates an accurate and natural sounding binaural experience	✓
Volume & Program Coupling ²		Ensures optimal hearing instrument system performance	✓
Hearing Programs		Options for the wearer's unique situations	6
Tinnitus Notch Function ⁴		Removes sound in wearer's tinnitus frequency to relax nerve activity and avoid over-stimulation	✓
Features		Benefits	80/60/40
Direct Audio Streaming from iPhone ¹		Stream phone and audio signals from iPhone, directly to hearing aids	✓
Rechargeable ⁵		Rechargeable lithium-ion hearing aid. Provides all day hearing with a single 3-hour charge	✓
Push Button Volume or Program Change		Allows the selection of different programs or volume levels	✓
XPhone ^{2,3}		Used with hand held phones; sound is transferred wirelessly to both hearing aids	✓
Telecoil ³		Picks up signals from an induction loop system via exchangeable battery door	✓
Accessories/Apps		Benefits	80/60/40
Smart Mic		Stream phone calls from Android and other Bluetooth phones.	✓
Smart Transmitter 2.4		Stream sounds from TV or other devices directly to the hearing aids	✓
Smart Key		Multi-function remote with easy-to-use buttons	✓
CROS RIC 8C Transmitter		Compatible with all MyCore 8C devices featuring Direct Streaming (Bluetooth)	✓
CROS inoX CIC 8C Transmitter		Compatible with all devices featuring Wireless Sync but not Direct Streaming	✓
Smart Direct App		Allows hearing aid to be controlled via a smartphone (iOS & Android)	✓
Smart Remote App		Adjusts program, volume, microphone directionality and treble	✓
Program Options		Benefits	80/60/40
Noahlink Wireless Programming		Program without wires	✓
HiPro, NoahLink		Requires battery pill adapters	10/312/13 Adapter
ConnexxAir		Wireless programmer for inoX and IIC only	✓

¹Apple iPhone 5 and later. ²Bilateral fitting required ³Available for Emerald M , Mosaic M and Mosaic P via exchangeable battery door ⁴Availability is country dependent ⁵Available on Stellar RIC