



The Magic of Direct Drive:

- ... What is it that makes the sound quality of direct drive hearing aids so good?

Drew Dundas, Ph.D.



Speaker & Disclosures

- Vice President, Audiology & Product Strategy at Earlens Corporation
- Ph.D., Audiology & Biomedical Engineering, Vanderbilt University
- Research interest: Quantification of listener perceived benefit
- Original research presented was funded in part by Earlens corporation

Learning Objectives

- Identify and define the various approaches to direct drive of the middle ear system.
- Identify advantages and disadvantages of middle ear implants and acoustic hearing aids.

Why bother with Direct Drive?

- Air conduction hearing aids (ACHA) have limitations that adversely impact the listening experience
- Severely limited low frequency response when vented
- Limited audibility and headroom in the high frequencies due to receiver roll off and feedback
- Users of direct drive hearing devices report:
 - Very natural and superior sound quality,
 - Superior performance in challenging listening situations,
 - Superior ability to hear soft sounds
 - Superior comfort for loud sounds


e.g., Kraus et al. 2011, Luetje et al. 2002, Jenkins et al. 2008 , Levy et al. 2016

Contents are confidential and proprietary

"It doesn't sound natural"


- Air Conduction devices cause distortion in the output signal
 - Harmonic Distortion
 - Intermodulation distortion
 - Comb Filtering
 - Spectral Ripples
 - Signal processing purposely distorts the signal to create audibility of high frequency speech cues
- Limited audible bandwidth and dynamic range
- At high levels, the TM creates distortions in the signal at the stapes due to anisotropic properties

Relative Audible Bandwidth – Mild to Mod SNHL




Increased Bandwidth Enhances Sound Quality

- If tolerable, extended high frequency amplification is beneficial
- Rated 'Naturalness'
 - Increases with bandwidth
- Rated Preference
 - Increases with bandwidth




High Frequency Audibility Is Not Enough

- Merely amplifying a narrow region can sound harsh and tinny
- If balanced with additional lows, providing MORE highs doesn't sound tinny, it sounds more natural⁵



Direct Drive Advantage


- Direct coupling to the ossicular chain overcomes the impedance mismatch that limits energy transfer from the air into the ear



ISO 226 Equal-Loudness-Level Contour

However...

- Implantable devices have limited audible bandwidth due to gain limitations, transducer drive efficiency (effective mass) and power budget considerations




Retrieved from: https://s3.medel.com/pdf/VSB_relaunch/28477_10_FactsheetVSBSystem_en_.pdf

Contents are confidential and proprietary

Exception: Earlens

- Consistently produces 100-10,000Hz audible bandwidth
- Maximum output increases smoothly from 500-10KHz, allowing for exceptional audibility even with steeply sloping loss




*data on file at Earlens


High Gain Margin without Trade-offs

With direct vibration instead of acoustic transmission:

- High gain margin with open fit
- Sound quality preserved without use of feedback cancellation at all in some cases
- Broad spectrum audibility without the annoyance of feedback



Relative Audible Bandwidth – Mild to Mod SNHL



Treatment of SNHL


... ACHA vs. Direct Drive Approaches

Approaches




- Fully Implantable
 - e.g., Esteem (Med El)
- Partially Implantable
 - e.g., Soundbridge (Med el)
 - e.g., Maxum (Ototronics)
- Non-implantable
 - e.g., Earlens
- Conventional Acoustic Hearing Aids


How it works - Conventional HA




How it works – Partially Implantable




How it works – Fully Implantable



How it works – Contact Drive



Earlens Contact Hearing Solution



Research Question

- Relative to acoustic stimulation, what is different about the signal transmitted to the stapes footplate via direct drive?

Experiment




Hypotheses

- Direct Drive of the ossicular chain generates a signal at the stapes footplate that is free of distortions to the input signal
- The calibration gain/damping effect of the direct drive mechanism prevents comb filtering


Methods

- In donated cadaver temporal bones:
- Compared stapes footplate responses achieved via:
 - Acoustic stimulation
 - Direct drive with Earlens Contact Hearing Solution
- Measured acceleration of stapes footplate using laser doppler vibrometry
- Derived frequency & phase responses to high level chirp signal


Experimental Setup



Experimental setup – Block Diagram




Measurement Setup



Stimulus

- 50ms Chirp
- 24.4 - 25,000Hz
- Repeated and averaged
160-180x



Signal Capture and Analysis


- Hardware: NI USB-4431 data-acquisition module (National Instruments, Austin, TX) with a maximum sampling rate of 96 KHz.
- Software: LabVIEW based synchronous-averaging measurement software (Gottlieb et al., 2016).
 - Sampling rate: 48KHz,
 - Fast Fourier Transform (FFT) length: 4096
 - Runs averaged/temporal bone: 10
- From the measurements of stapes velocity (VST) and ear canal pressure (PEC), the following quantities are calculated:
 - The baseline sound-driven stapes transfer function without the Tympanic Lens on the TM.
 - The equivalent pressure output of the direct drive system

Results

- At the stapes footplate,
 - Acoustic drive condition demonstrated:
 - Minor comb filtering
 - Spectral Ripples
 - Phase shift with frequency
 - Contact direct drive demonstrated:
 - Smooth, flat spectral response
 - Consistent phase relationship


Results

- Spectral ripples minimized with direct drive
- Superior audible bandwidth to both low and high frequencies



Results

- Phase relationship preserved with direct drive



Conclusions

- Substantial differences in transmitted signal quality are observed between acoustic and direct drive modalities
- Direct drive exhibits:
 - Superior effective bandwidth to both low and high frequencies
 - Smooth spectral shape and no induced ripples
 - Preserved phase relationship between input and output signals
 - Minimal group delay
 - Superior stable gain margin

Next Steps

- Investigation of direct/amplified path interaction effects on stapes signal with open fittings
 - Output SNR
 - Comb filtering effects

Acknowledgements

- Thank you to:
 - Sunil Puria, PhD
 - Suzanne Levy, PhD
 - Brian Moore, PhD
 - Morteza Khalegi, PhD

References

1. Struck CJ, Prusick L. (2017) Comparison of real-world bandwidth in hearing aids vs Earlens light-driven hearing aid system. *Hearing Review*. 24(3):24-29.
2. Moore BCJ, Tan CT (2003). Perceived naturalness of spectrally distorted speech and music. *J Acoust Soc Am*. 114:408-419.
3. Turner CW, Henry BA (2002). Benefits of amplification for speech recognition in background noise. *J Acoust Soc Am*. 112:1675-1680.
4. Hornsby, BW, Ricketts, TA (2003). The effects of hearing loss on the contribution of high- and low-frequency speech information to speech understanding. *J Acoust Soc Am*. 113:1706-1717.
5. Plyler PN, Fleck EL (2006). The effects of high-frequency amplification on the objective and subjective performance of hearing instrument users with varying degrees of high-frequency hearing loss. *J Sp Lang Hear Res*. 49:616-627.
6. Moore BCJ, Füllgrabe C, Stone MA (2010). Effect of spatial separation, extended bandwidth, and compression speed on intelligibility in a competing-speech task. *J Acoust Soc Am*. 128:360-371.
7. Ahlstrom JB, Horwitz AR, Dubno JR. (2014) Spatial separation benefit for unaided and aided listening. *Ear Hear*. 35:72-85.
8. Levy SC, Freed DJ, Nilsson M, Moore BC, Puria S. (2015) Extended high-frequency bandwidth improves speech reception in the presence of spatially separated masking speech. *Ear Hear*. 36:e214-e224.
9. Puria S., (2013), "Middle ear hearing devices", In "The Middle Ear: Science, Otosurgery, and Technology". Edited by Puria, S., Popper, A.N., and Fay, R.R. *Springer Handbook of Auditory Research Series*.
10. Fay JP, Perkins R, Levy SC, Nilsson M, Puria S. (2013) "Preliminary evaluation of a light-based contact hearing device for the hearing impaired." *Otol Neurotol*. Jul;34(5):912-21.
11. Gantz BJ, Perkins R, Murray MT, Levy SC, & Puria S (2017) Light-Driven Contact Hearing Aid for Broad Spectrum Amplification: Safety and Effectiveness Pivotal Study. *Otol Neurotol*. Mar;38(3):352-359
12. McElveen J, Gantz BJ, Perkins R, Murray M, Syms CA, Chen DA, Arriaga M, Mehta R, & Levy SC (2017), High-frequency Gain Is Better Provided By A Light-driven Hearing Aid: A Comparative Study, COSM, San Diego.
13. Khaleghi M, and Puria S (2017). Attenuating the ear canal feedback pressure of a laser-driven hearing aid. *J Acoust Soc Am.*, 141(3), pp.1683-1693.
14. Maier H, Muller M, Salcher R, Lenarz T, (2016) "Improved efficiency in round window stimulation by controlled pretension with the Hannover coupler," IHCON 2016 conference, 2016.
15. Grosse-Himmich M, Busch S, Waldmann B, Lenarz T, Maier H. (2016) "Validating a method to characterize the efficiency of implantable middle ear hearing devices in cadaver ears by intracochlear pressure," ARO midwinter meeting.
16. Kraus EM, Shohet, JA, & Catalano PJ (2011). Envoy Esteem totally implantable hearing system: Phase 2 trial, 1-year hearing results. *Otolaryngology – Head and Neck Surgery*, 145(1), 100–109.
17. Luetje CM, Brackman D, Balkany TJ, Maw J, Baker RS, Kelsall D, Backous D, Miyamoto R, Parisier S, & Arts A (2002). Phase III clinical trial results with the Vibrant Soundbridge implantable middle ear hearing device: A prospective controlled multicenter study. *Otolaryngology – Head and Neck Surgery*, 126(2), 97–107.